



Psychē Systems
EndLIS Possibilities

PSYCHE SYSTEMS

BACKBONE NUCLEOLIS API

REVISION 1.9

10/29/2024

(REQUIRES BACKBONE SERVER V32 OR LATER)

TABLE OF CONTENTS

Overview – getting started.....	6
Theory of operation	6
SECURITY	6
1. Authentication.	7
2. Secure Transport Protocol (HTTP-S).	7
3. SlidinG WINDOW TOKEN EXPIRY	7
Getting Started	7
Step 1 – Verify/CONFIRM Backbone URL Connectivity	7
Step 2 – Logging Into the Backbone NucleoLIS API	8
Step 3- Calling API endpoints within a session.....	10
Step 4- TERMINATE A session and log out.....	11
Accessing Help	12
API Calls Specifications	13
get heartbeat	13
get patient	13
get patients.....	15
get case.....	18
get cases	34
get specimen.....	38
get specimens	40
get test order	53
get test orders.....	61
get physician	65
get physicians	67
set status step	72
set status steps	73
get physician locations.....	75
get report.....	77
get reports	93
get result report	96

get case type definitions	97
get source definitions.....	102
get profile definitions	107
get question definitions	110
set image Document.....	114
set instrument test results.....	116
set instrument Batch test results.....	119
set uploaded text	122
add patient.....	123
update patient	125
add case	128
update case.....	129
add specimen.....	131
update specimen.....	133
add test order.....	135
cancel case	137
cancel specimen.....	138
cancel test order	139
get plate map definitions	140
get plates	144
get plate map.....	146
get charge definitions	149
get charges	151
get charge.....	154
Coding Examples Introduction	157
Logon Coding Examples	159
C – libcurl	159
C# - RestSharp.....	159
CURL.....	159
GO - Native	160
HTTP.....	161

Java - OKhttp	161
Java – UNIREST	161
JAVASCRIPT - FETCH	161
JAVASCRIPT - JQUERY	162
JAVASCRIPT – XHR	162
NODEJS – AXIOS	163
NODEJS – NATIVE	163
NODEJS – REQUEST	164
NODEJS – UNIREST	164
OBJECTIVE-C NSURLSESSION	165
OCAML – COHTTP	166
PHP – CURL	166
PHP – HTTP_Request2	167
PHP – pecl_http	167
POWERSHELL – RESTMETHOD	168
python – http.client	168
python – requests	168
ruby – net::HTTP	169
SHELL – HTTPIE	169
SHELL – WGET	169
SWIFT - URLSESSION	169
VB.NET - RestSHARP	170
Logout Coding Examples	170
C - libcurl	170
C# - RestSharp	171
cURL	171
Go – native	171
HTTP	172
Java – OKHTTP	172
Java – unirest	173
javascript – fetch	173

javaScript – jquery.....	173
javascript – xhr	174
NodeJS – AXIOS.....	174
nodeJS – NATIVE	174
nodejs – request.....	175
nodejs – unirest	176
objective-C – nsurlsession	176
ocaml – cohttp	177
PHP – curl.....	177
php – HTTP_Request2.....	178
php – pecl_http	178
powershell – restmethod.....	179
python – http.client.....	179
python – requests.....	179
ruby – net::HTTP.....	180
shell – httpie.....	180
shell – wget	180
swift – URLSession.....	181
vb.net – restsharp	181
FOR MORE INFORMATION	181

BACKBONE API

OVERVIEW – GETTING STARTED

This is a guide for using the Psyche Systems Backbone application programming interface (API).

This section of the guide covers getting started information for four essential requirements for communicating with the Backbone Server using the Backbone API:

1. Verifying/Confirming the Backbone Server Uniform Resource Locator (URL)
2. Logging in, authenticating, and establishing a session.
3. Calling API endpoints during a session.
4. Terminating a session and log out.

Appendix A provides more details about the available API commands that can be used for interacting with NucleoLIS.

Appendix B provides multiple coding examples illustrating how you can program a client side implementation to call the Backbone API.

THEORY OF OPERATION

The Backbone API is a Representational state transfer (REST) based API that is designed to run on Microsoft Internet Information Server (IIS). Typically, Psyche Systems deploys the server-side of the Backbone on IIS version 7 or later to host the Backbone Server API. The client side implementation of the API can be written in any programming language and run on any operating system as long as the client side implementation adheres to the standards outlined in this document.

To ensure Health Insurance Portability And Accountability Act (HIPAA) privacy compliance and to safeguard Protected Health Information (PHI), all requests to the Backbone Server API must use Hypertext Transfer Protocol Secure (HTTPS) over port 443.

The Backbone Server (being a REST based web server implementation) is considered to be stateless, with the exception of the "session" cookie which is passed back and forth between the client side implementation and the Backbone Server with each request after a successful logon is established.

SECURITY

To further ensure HIPAA compliance and security of PHI, the Backbone API provides session-based security as follows:

1. AUTHENTICATION.

A valid username and password are required to access the Backbone Server API. Once logged in, the API returns an access token that valid for 15 minutes and can be used to make subsequent requests. This access token is stored in a cookie which is transparently passed back and forth with each API request/response. This access token resets every time a successful API call is completed. With this approach, an application can authenticate once, make any number of API calls, and then log out, without needing to pass authentication credentials with every request.

2. SECURE TRANSPORT PROTOCOL (HTTP-S).

All API calls are made using Internet standard Hypertext Transport Protocol Secure (HTTPS) on port 443. The bidirectional [encryption](#) of communications between a client and server protects against [eavesdropping](#) and [tampering](#) of the communication. In practice, this provides a reasonable assurance that one is communicating without interference by attackers with the website that one intended to communicate with, as opposed to an impostor.

This implementation also uses the latest [Transport Layer Security](#) (TLS) encryption. This public-private key encryption overlays the symmetric AES encryption used for encrypting the payload, so essentially, the payload is double encrypted and is extremely secure.

3. SLIDING WINDOW TOKEN EXPIRY

Each time a session token is provided to a requesting client, the token is encoded so that it is only valid for 15 minutes. This provides what is known as “sliding window” expiration. By using this security method, tokens that are copy/pasted can only be used for a limited period of time. Each time a successful request is made to the Backbone Server API, the token is refreshed to provide a new 15 minute window. In the event that no request is made within 15 minutes, the token expires and can no longer be used without logging in again to obtain a new token. This is the equivalent to the “idle timeout” feature that many online systems offer.

GETTING STARTED

For each of the examples shown, this document also provides in Appendix “B”, coding examples in a wide variety of programming languages so the implementers can research options for creating their custom client implementation that communicates with the API.

STEP 1 – VERIFY/CONFIRM BACKBONE URL CONNECTIVITY

When first getting started with a client side implementation, it is helpful to first ensure that you can “reach” the Backbone Server from your client location. To do this, enter the URL of the Backbone Server that you were provided into any web browser. If you are able to connect using that URL, you will see a page that looks like this:

Welcome to Psyche Systems Web API!

Access to this API is granted to Authorized users only.

For more information please contact Psyche Systems at: (508) 478-2047.
Or, via email at: elixasupport@PsycheSystems.com.

If you don't see this page, please contact Psyche Systems support for assistance. If you do see this page, then this is the URL that you will use when making all API requests. This URL will be referred to as "https://www.SomeServer.com/api/" in the rest of this document.

STEP 2 – LOGGING INTO THE BACKBONE NUCLEOLIS API

Logging in, authenticates the requesting client as a valid user, and upon successful login, returns an authorization token that is stored in a session cookie. This cookie must be passed after each subsequent request.

This is accomplished by calling the "login" endpoint with username and password information that is POSTED to the Backbone Server with the following settings:

URI	<API URL>/api/authenticate/login
Host	Provide the hostname of the web server you are connecting with.
Authorization Method	No Auth
Content-Type	<code>application/x-www-form-urlencoded</code>
Content-Length	<p>The length of the body payload in bytes. For example, if you are posting:</p> <pre>username=LAB&password=123456</pre> <p>Then, the Content-Length parameter would be set to 28. (Because there are 28 characters in the payload).</p>
User-Agent	Set this to a unique name for your client application.

Accept */*

Posted Data username=<username>&password=<password>

Here, <username> and <password> are the username and password you were provided to access the API. If you do not have these credentials, please contact Psyche Systems support.

When the called API command succeeds, the response will look like this:

Status 200 OK

Body <?xml version="1.0" encoding="UTF-8" ?><login>Success</login>

Set-Cookie session=session-id=qAuOpF9tIBeBR8SIGLmDVRQGUFJ%2ftxLUbELWr8HtNx%3d&username=346650&epoch=1604331210&timelimit=900; expires=Mon, 02 Nov 2020 15:48:30 GMT; domain=SomeServer.com; path=/; secure; httponly

Content-Type Text/plain; charset=utf-8

Server Microsoft-IIS/8.0

If the username or password is invalid, the response will look like this:

Status 401

Body <?xml version="1.0" encoding="UTF-8" ?><login>username or password invalid</login>

Set-Cookie (No cookie returned)

Content-Type Text/plain; charset=utf-8

Server Microsoft-IIS/8.0

Other error messages can sometimes be shown as well. As of Backbone 1.1.26, errors are returned in the same format as specified in the call (XML vs JSON). If you don't receive the word "Success" within the xml return message, then you should assume that the login failed, and the message returned is the error message.

Note: On initial login, there no cookie specified yet. The cookie gets created by the Backbone Server, and is returned to the client only after a successful logon request is executed. This "session" cookie will contain the authentication token and will look similar to this:

```
Cookie: session=session-  
id=qAuOpF9tIBeBR8SIGLmDVRQGUFJ%2ftxLUbELwr8HtNxo%3d&username=346650&e  
poch=1604331210&timelimit=900; expires=Mon, 02 Nov 2020 15:48:30 GMT;  
domain=SomeServer.com; path=/; secure; httponly
```

This cookie will be updated after each successful request and must be passed back to the server on each subsequent request. The session token contained within the cookie is only valid for 15 minutes. If a subsequent API endpoint request is not made within the session "timeout" limit of 15 minutes, then the session token becomes invalid, and the client will need to login again to get a new session token by calling the "login" API endpoint again.

See Appendix "B" for programming examples that call this endpoint.

STEP 3- CALLING API ENDPOINTS WITHIN A SESSION

Once you are logged in, the next step is to call an API Endpoint using GET. In this example we will show how to call the NucleoLIS api/N/GetHeartbeat API endpoint. This API endpoint returns a heartbeat message so you can determine if the API is able to communicate with NucleoLIS.

URI	<API URL>/api/N/GetHeartbeat
Host	Provide the hostname of the web server you are connecting with.
Authorization Method	No Auth
Cookie	Pass in the cookie you received in the response from the logon API call. It should look similar to this: <pre>session=session- id=H38%2fAImLMVMNssLZsY%2fx7mWkfNsHIM4m1CVO%2bFDc9aU%3 d&username=346650&epoch=1604330631&timelimit=900</pre>
User-Agent	Set this to a unique name for your client application.
Accept	*/*

When the called API command succeeds, the response will look like this:

Status	200 OK
Body	<?xml version="1.0" encoding="UTF-8" ?><Result>OK</Result>
Set-Cookie	session=session-id=xjxAnqMxbu9D%2fbYakV5nRMeKxLlEVnADTiMJqUJycp0%3d&username=346650&epoch=1604330884&timelimit=900; expires=Mon, 02 Nov 2020 15:43:04 GMT; domain=SomeServer.com; path=/; secure; httponly
Content-Type	Text/plain; charset=utf-8
Server	Microsoft-IIS/8.0

See Appendix "B" for programming examples that call this endpoint.

STEP 4- TERMINATE A SESSION AND LOG OUT

When you are done communicating with the Backbone Server API, there are two choices on how to end the session.

1. Do nothing – the cookie containing the latest session token will expire after 15 minutes.
2. Call the Logout API endpoint to destroy the cookie and immediately terminate the session:

URI	<API URL>/api/authenticate/logout
Host	Provide the hostname of the web server you are connecting with.
Authorization Method	No Auth
Cookie	Pass in the cookie you received in the response from the logon API call. It should look similar to this: session=session-id=H38%2fAImLMVMNssLZsY%2fx7mWkfNsHIM4m1CVO%2bFDc9aU%3d&username=346650&epoch=1604330631&timelimit=900

User-Agent Set this to a unique name for your client application.

Accept * / *

When the called API command succeeds, the response will look like this:

Status 200 OK

Body <?xml version="1.0" encoding="UTF-8" ?><logout>Success</logout>

Cookie (The cookie no longer exists)

Content-Type Text/plain; charset=utf-8

ACCESSING HELP

As of Backbone 1.1.24, all available API calls may be viewed online by entering **/help** after the URL you are using. (<https://www.SomeServer.com/help>) Below is a sample of the API Help page.

Web API Help Page

Introduction

This Web API is for authorized users only.
For more information, please contact Psyche Systems Corporation at (508) 473-1500.

Security

API	Description
GET api/authenticate/GetResetPassword?user_id={user_id}&email={email}	api/authenticate/GetResetPassword?user_id={user_id}&email={email} Reset password request - only ONE of the following parameters must be provided: email, user_id, or code. Returns: A Token to be used for facilitating the reset process, else FAIL Remarks: Either the User Id or the Email must be provided. Both parameters are not required.
POST api/authenticate/login	api/authenticate/login Logs in and creates a session Returns: SUCCESS if logon was successful, else an error message. Remarks: Uses post - Account will be locked if there are three failed login attempts in a row.
GET api/authenticate/logout	api/authenticate/logout Logout ends a session. Returns: Success if login credentials are destroyed. Remarks: This is a remark.
GET api/authenticate/GetChangePassword?email={email}&old_password={old_password}&token={token}&new_password={new_password}	api/authenticate/GetChangePassword?email={email}&old_password={old_password}&token={token}&new_password={new_password} Change the users password Returns: SUCCESS if password was changed, else FAIL with message Remarks:
GET api/Security/GetResetPassword?user_id={user_id}&email={email}	api/Security/GetResetPassword?user_id={user_id}&email={email} Reset password request - only ONE of the following parameters must be provided: email, user_id, or code. Returns: A Token to be used for facilitating the reset process, else FAIL Remarks: Either the User Id or the Email must be provided. Both parameters are not required.
POST api/Security/login	api/Security/login Logs in and creates a session Returns: SUCCESS if logon was successful, else an error message. Remarks: Uses post - Account will be locked if there are three failed login attempts in a row.

Logged in to <https://QA-WEB.psyche.prj/WebApiBackboneQA> at 1/26/2023 3:10:51 PM

APPENDIX A – NUCLEOLIS API

ENDPOINTS

API CALLS SPECIFICATIONS

In addition to the logon and logoff API endpoints covered in section 1, the endpoints in this appendix are all of the available API endpoints for NucleoLIS Interaction.

GET HEARTBEAT

Endpoint	api/N/GetHeartbeat
Purpose	Returns a NucleoLIS heartbeat message so you can determine if the API is able to communicate with NucleoLIS system.
Request Type	GET
Parameters	{None}
Example	<code>https://www.SomeServer.com/api/N/GetHeartbeat</code>
Returns	If the heartbeat succeeds, returns: <pre><?xml version="1.0" encoding="UTF8" ?><Result>OK</Result></pre>

GET PATIENT

Endpoint	api/N/GetPatient
Purpose	Returns a single, full NucleoLIS patient record based on a patient unique identifier.
Request Type	GET

Parameters

```
patient_id={patient_id}

user_id={user_id}

return_format={return_format}
```

Where:

patient_id is the objectID of the requested Patient record.

user_id is the user code for the user in the NucleoLIS user table that is requesting the record.

return_format is the format for returning the data.
Where:

0 = XML Format

1 = JSON Format

Example

```
https://www.SomeServer.com/api/N/GetPatient?patient_id=518535440&user_id=245654&return_format=0
```

Returns If the request succeeds, returns a patient record in the following XML format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <PatientTable>
    <Patient_ObjectID>518535440</Patient_ObjectID>
    <Patient_CreationDate>2020-11-02T11:23:31.247-05:00</Patient_CreationDate>
    <Patient_UpdateTime>2020-11-02T11:23:31.247-05:00</Patient_UpdateTime>
    <Patient_Name>Boop, Betty</Patient_Name>
    <Patient_Sex>F</Patient_Sex>
    <Patient_Street />
    <Patient_City />
    <Patient_State />
    <Patient_Zip />
    <Patient_Phone />
    <Patient_Fax />
    <Patient_Code>518535440</Patient_Code>
    <Patient_Soundex>BOPBTY</Patient_Soundex>
    <Patient_Dob>1965-03-23</Patient_Dob>
```

```

    <Patient_Ssn>011-22-3333</Patient_Ssn>
    <Patient_MedicareNumber />
    <Patient_Status>I</Patient_Status>
    <Patient_CREATOR>LAB</Patient_CREATOR>
    <Patient_RACE>White</Patient_RACE>
    <Patient_STREET2 />
    <Patient_USER1 />
    <Patient_USER2 />
    <Patient_USER3 />
    <Patient_USER4 />
    <Patient_MODIFIER>LAB</Patient_MODIFIER>
    <Patient_LastName>Boop</Patient_LastName>
    <Patient_FirstName>Betty</Patient_FirstName>
    <Patient_MiddleName />
    <Case_Status>I</Case_Status>
  </PatientTable>

</NewDataSet>

```

GET PATIENTS

Endpoint	api/N/GetPatients
Purpose	Returns a collection of NucleoLIS patient records based on filter parameters passed.
Request Type	GET
Parameters	last_name = {last_name} (use * for wildcard) first_name = {first_name} (use * for wildcard) activesOnly = {activesOnly} (true, or false) filter_expression = {filter_expression} user_id = {user_id} return_format = {return_format}

Where:

last_name is the last name of the patient(s) to search for. Wildcards are supported. (Example S* will return all patients with a last name beginning with 'S'.)

first_name is the first name of the patient(s) to search for. Wildcards are supported. (Example L* will return all patients with a first name beginning with 'L'.)

activesOnly is a true or false parameter. If true, only active patients will be returned. If false, all patients (both active and inactive) will be returned.

filter_expression is used as a secondary custom filtering criteria applied to the result set . Wildcards are supported. It can be used with any column returned in the result set. Any valid SQL WHERE clause is accepted.

Multiple columns can be used when separated with 'AND'. Example: Patient_State='MA' AND Patient_City='W*' will return all patients that live in MA in a city that begins with 'W'.

user_id is the user code for the user in the NucleoLIS user table that is requesting the record.

return_format is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

Example

`https://www.SomeServer.com/api/N/GetPatients?last_name=S*&first_name=*&activesOnly=true&filter_expression=Patient_Sex='Male'&user_id=245654&return_format=0`

Returns

If the request succeeds, returns a collection of 1..n patient records matching the filtering criteria in the following XML format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <PatientTable>
    <Patient_ObjectID>433016111</Patient_ObjectID>
    <Patient_Name>Claus, Santa</Patient_Name>
    <Patient_Code>433016111</Patient_Code>
```



```

    <Patient_Street />
    <Patient_STREET2 />
    <Patient_City />
    <Patient_State />
    <Patient_Zip />
    <Patient_Sex>Male</Patient_Sex>
    <Patient_Dob>1900-12-25</Patient_Dob>
    <Patient_Ssn />
  </PatientTable>
  <PatientTable>
    <Patient_ObjectID>503609702</Patient_ObjectID>
    <Patient_Name>Claus, Santa</Patient_Name>
    <Patient_Code>503609702</Patient_Code>
    <Patient_Street />
    <Patient_STREET2 />
    <Patient_City />
    <Patient_State />
    <Patient_Zip />
    <Patient_Sex>Male</Patient_Sex>
    <Patient_Dob>1900-12-25</Patient_Dob>
    <Patient_Ssn />
  </PatientTable>
  <PatientTable>
    <Patient_ObjectID>504251809</Patient_ObjectID>
    <Patient_Name>Claus, Santa</Patient_Name>
    <Patient_Code>504251809</Patient_Code>
    <Patient_Street />
    <Patient_STREET2 />
    <Patient_City />
    <Patient_State />
    <Patient_Zip />
    <Patient_Sex>Male</Patient_Sex>
    <Patient_Dob>1900-12-25</Patient_Dob>
    <Patient_Ssn />
  </PatientTable>
  <PatientTable>
    <Patient_ObjectID>4550996</Patient_ObjectID>
    <Patient_Name>Sanding, Sammie</Patient_Name>
    <Patient_Code>4550996</Patient_Code>
    <Patient_Street />
    <Patient_STREET2 />
    <Patient_City />
    <Patient_State />
    <Patient_Zip />
    <Patient_Sex>Male</Patient_Sex>
    <Patient_Dob>1980-04-22</Patient_Dob>
    <Patient_Ssn />
  </PatientTable>

```

```

</PatientTable>
<PatientTable>
  <Patient_ObjectID>507538376</Patient_ObjectID>
  <Patient_Name>Sanding, Sammie</Patient_Name>
  <Patient_Code>507538376</Patient_Code>
  <Patient_Street />
  <Patient_STREET2 />
  <Patient_City />
  <Patient_State />
  <Patient_Zip />
  <Patient_Sex>Male</Patient_Sex>
  <Patient_Dob>1980-04-22</Patient_Dob>
  <Patient_Ssn />
</PatientTable>

</NewDataSet>

```

GET CASE

Endpoint	api/N/GetCase
Purpose	Return a single, full NucleoLIS case record based on case number. The full case record inclusive of all data elements associated with patient, case, specimen, test order and test result records are returned.
Request Type	GET
Parameters	<p>case_number={case_number}</p> <p>user_id={user_id}</p> <p>return_format={return_format}</p> <p>Where:</p> <p>case_number is the number of the case.</p> <p>user_id is the user code for the user in the NucleoLIS user table that is requesting the record.</p> <p>return_format is the format for returning the data.</p> <p>Where:</p> <p>0 = XML Format</p> <p>1 = JSON Format</p>

Example

[https:// www.SomeServer.com
/api/N/GetCase?case_number=M2020-
000032&user_id=245654&return_format=0](https://www.SomeServer.com/api/N/GetCase?case_number=M2020-000032&user_id=245654&return_format=0)

Returns

If the request succeeds, returns the requested case record in the following XML format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <CaseTable>
    <Patient_ObjectID>163307929</Patient_ObjectID>
    <Patient_CreationDate>2016-11-10T16:11:29.84-
05:00</Patient_CreationDate>
    <Patient_UpdateTime>2019-05-22T10:51:15.137-
04:00</Patient_UpdateTime>
    <Patient_Name>O'Hara, Jaime</Patient_Name>
    <Patient_Sex>Female</Patient_Sex>
    <Patient_Street />
    <Patient_City />
    <Patient_State />
    <Patient_Zip />
    <Patient_Phone />
    <Patient_Fax />
    <Patient_Code>163307929</Patient_Code>
    <Patient_Soundex>O'HRJM</Patient_Soundex>
    <Patient_Dob>2018-05-16</Patient_Dob>
    <Patient_Ssn />
    <Patient_MedicareNumber />
    <Patient_Status>I</Patient_Status>
    <Patient_CREATOR>LAB</Patient_CREATOR>
    <Patient_RACE />
    <Patient_STREET2 />
    <Patient_USER1 />
    <Patient_USER2 />
    <Patient_USER3 />
    <Patient_USER4 />
    <Patient_Modifier>Scaron</Patient_Modifier>
    <Patient_LastName>O'Hara</Patient_LastName>
    <Patient_FirstName>Jaime</Patient_FirstName>
    <Visit_ObjectID>498604164</Visit_ObjectID>
    <Visit_CreationDate>2019-05-22T10:51:34.293-
04:00</Visit_CreationDate>
    <Visit_UpdateTime>2019-05-22T10:51:34.293-
04:00</Visit_UpdateTime>

    <Visit_SuperobjectID>163307929</Visit_SuperobjectID>
    <Visit_Patient>163307929</Visit_Patient>
```

```

<Visit_Code>498604164</Visit_Code>
<Visit_Location />
<Visit_Type>OP</Visit_Type>
<Visit_VisitDate>2019-05-22</Visit_VisitDate>
<Visit_Status>I</Visit_Status>
<Visit_DischargeDate />
<Visit_ExtraField1 />
<Visit_ExtraField2 />
<Visit_CREATOR>Scaron</Visit_CREATOR>
<Visit_USER1 />
<Visit_USER2 />
<Visit_USER3 />
<Visit_USER4 />
<Visit_comment />
<Case_ObjectID>506873424</Case_ObjectID>
<Case_CreationDate>2020-03-09T14:50:06.023-
04:00</Case_CreationDate>
<Case_UpdateTime>2020-05-11T14:03:15.34-
04:00</Case_UpdateTime>
<Case_SuperobjectID>498604164</Case_SuperobjectID>
<Case_User1 />
<Case_User2 />
<Case_User3 />
<Case_User4 />
<Case_CREATOR>Scaron</Case_CREATOR>
<Case_MODIFIER>Bishop</Case_MODIFIER>
<Case_Number>M2020-000032</Case_Number>
<Case_ReqNumber />
<Case_ReportedFlag>False</Case_ReportedFlag>

<Case_SubmitterLocation>258012099</Case_SubmitterLocat
ion>
<Case_OrderID />
<Case_Status>I</Case_Status>
<Case_ReportStatus>FINAL</Case_ReportStatus>
<Case_CaseType>14662557</Case_CaseType>
<Case_CorrectedComment />
<Case_AmendedComment />
<Case_ReportedBy>Bishop</Case_ReportedBy>
<Case_CanceledBy />
<Case_CanceledDate />
<Case_ReviewFlag>2020-03-09</Case_ReviewFlag>
<Case_ReviewedBy />
<Case_ReviewedDate />
<Case_ReviewedComment />
<Case_ReviewedTime />
<Case_StatusStep />

```

```

<Case_CanceledFor />
<Case_Signer>370583063</Case_Signer>
<Case_SignerDate>2020-03-09</Case_SignerDate>
<Case_ICDs />
<Case_ClinicalInformation />
<Case_Images><?xml version="1.0" encoding="utf-
16"?>
<ChrysalisImgData
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema"
/></Case_Images>
  <Case_DocumentID />
  <Case_OrderDate>2020-03-09</Case_OrderDate>
  <Case_OrderTime />
  <Case_ExternalID />
  <Case_ExternalPatientID />
  <Case_ReopenReason />
  <Case_ReopenDate />
  <Case_Signatures><?xml version="1.0"
encoding="utf-16"?>
<CaseSignatureHistory
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema">
  <Sigs>
    <CaseSignature>
      <Signer>370583063</Signer>
      <SignerDate>2020-03-09T15:08:36</SignerDate>
      <ReOpenReason />
      <ReOpener>0</ReOpener>
      <ReOpenedDate>0001-01-01T00:00:00</ReOpenedDate>
      <UnSign>false</UnSign>
    </CaseSignature>
    <CaseSignature>
      <Signer>370583063</Signer>
      <SignerDate>2020-03-09T15:10:42</SignerDate>
      <ReOpenReason />
      <ReOpener>0</ReOpener>
      <ReOpenedDate>0001-01-01T00:00:00</ReOpenedDate>
      <UnSign>true</UnSign>
    </CaseSignature>
    <CaseSignature>
      <Signer>370583063</Signer>
      <SignerDate>2020-03-09T15:11:09</SignerDate>
      <ReOpenReason />
      <ReOpener>0</ReOpener>
      <ReOpenedDate>0001-01-01T00:00:00</ReOpenedDate>
      <UnSign>false</UnSign>
    </CaseSignature>
  </Sigs>
</CaseSignatureHistory>
</Case_Signatures>

```

```

        </CaseSignature>
    </Sigs>
</CaseSignatureHistory></Case_Signatures>
    <Case_NeedsCorrectedF>False</Case_NeedsCorrectedF>
    <Case_PrescriptionDrugs><?xml version="1.0"
encoding="utf-16"?>
<PrescribedDrugs
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema"
/></Case_PrescriptionDrugs>
    <Case_DataCollect />
    <Case_ReportIX>1</Case_ReportIX>
    <Case_OInitStatus>I</Case_OInitStatus>

<Case_RequireSignature>True</Case_RequireSignature>
    <Case_SignerTime>15:11</Case_SignerTime>
    <Case_CorrectionRequest />
    <Case_CorrectionRequestDate />
    <Case_CorrectionRequestTime />

<Case_SignerLocation>246584261</Case_SignerLocation>
    <Case_ClearCorrectedReason />
    <Case_QuestionsStash><?xml version="1.0"
encoding="utf-16"?>
<AAOEstasherBase
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema">
    <Responses>
        <AAOEItemBase>
            <QObjectID>302726734</QObjectID>
            <Question>Additional Information</Question>
            <AllowOther>true</AllowOther>
            <SingleSelect>false</SingleSelect>
            <Required>true</Required>
            <FullResponse />
            <GridQuestion>false</GridQuestion>
            <GridQuestions />
            <Answers>
                <string>none</string>
            </Answers>
        </AAOEItemBase>
        <AAOEItemBase>
            <QObjectID>481211124</QObjectID>
            <Question>General Health</Question>
            <AllowOther>true</AllowOther>
            <SingleSelect>true</SingleSelect>
            <Required>false</Required>

```

```

        <FullResponse />
        <GridQuestion>false</GridQuestion>
        <GridQuestions />
        <Answers>
            <string>Good</string>
        </Answers>
    </AAOEItemBase>
    <AAOEItemBase>
        <QObjectID>481215687</QObjectID>
        <Question>Fasting</Question>
        <AllowOther>false</AllowOther>
        <SingleSelect>true</SingleSelect>
        <Required>true</Required>
        <FullResponse />
        <GridQuestion>false</GridQuestion>
        <GridQuestions />
        <Answers>
            <string>Yes</string>
        </Answers>
    </AAOEItemBase>
</Responses>
</AAOEstasherBase></Case_QuestionsStash>
    <Case_Questions>Additional Information:none
    General Health:Good
    Fasting:Yes</Case_Questions>

<Case_QuestionsComplete>True</Case_QuestionsComplete>
    <Specimen_ObjectID>506873435</Specimen_ObjectID>
    <Specimen_CreationDate>2020-03-09T14:50:41.513-
    04:00</Specimen_CreationDate>
    <Specimen_UpdateTime>2020-03-09T14:51:10.153-
    04:00</Specimen_UpdateTime>

<Specimen_SuperobjectID>506873424</Specimen_SuperobjectID>
    <Specimen_User1 />
    <Specimen_User2 />
    <Specimen_User3 />
    <Specimen_User4 />
    <Specimen_CREATOR>Scaron</Specimen_CREATOR>
    <Specimen_MODIFIER>Scaron</Specimen_MODIFIER>
    <Specimen_CollectionDate>2020-03-
    09</Specimen_CollectionDate>

<Specimen_CollectionTime>09:00</Specimen_CollectionTime>
e>

```

```

    <Specimen_ReceivedDate>2020-03-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>14:51:10</Specimen_ReceivedTime
>
    <Specimen_Receiver>Scaron</Specimen_Receiver>
    <Specimen_Number>506873435</Specimen_Number>
    <Specimen_Status>A</Specimen_Status>
    <Specimen_Comment />
    <Specimen_Description />

<Specimen_SourceDefinition>4348467</Specimen_SourceDef
inition>
    <Specimen_Interpretation />
    <Specimen_StatusStep />
    <Specimen_SuperOrders />
    <Specimen_UnInlabDate />
    <Specimen_DataCollect />
    <Specimen_Priority>Routine</Specimen_Priority>
    <Specimen_Site />
    <Specimen_Gross />
    <Specimen_AddOnDate />
    <Specimen_AddOnTime />
    <Specimen_AddOn>False</Specimen_AddOn>
    <Specimen_MultipleCountID />
    <Specimen_MergeComments />
    <Specimen_MergeIntoComments />
    <Specimen_QuestionsStash />
    <Specimen_Questions />

<Specimen_QuestionsComplete>True</Specimen_QuestionsCo
mplete>
    <Specimen_FastingStatus />
    <Profile_TYPE>IMMUNO</Profile_TYPE>
    <Profile_ObjectID>506873438</Profile_ObjectID>
    <Profile_CreationDate>2020-03-09T14:50:49.687-
04:00</Profile_CreationDate>
    <Profile_UpdateTime>2020-03-09T15:07:50.84-
04:00</Profile_UpdateTime>

<Profile_SuperobjectID>506873435</Profile_SuperobjectI
D>
    <Profile_User1 />
    <Profile_User2 />
    <Profile_User3 />
    <Profile_User4 />
    <Profile_CREATOR>Scaron</Profile_CREATOR>

```



```

    <Profile_MODIFIER>Suec</Profile_MODIFIER>
    <Profile_Description />
    <Profile_OrderDate>2020-03-09</Profile_OrderDate>
    <Profile_OrderTime>14:50</Profile_OrderTime>
    <Profile_OrderedBy>Scaron</Profile_OrderedBy>

    <Profile_PanelDefinition>416976655</Profile_PanelDefin
    ition>
        <Profile_R1 />
        <Profile_R2 />
        <Profile_R3 />
        <Profile_R4 />
        <Profile_R5 />
        <Profile_R6 />
        <Profile_Result />
        <Profile_Reportable>False</Profile_Reportable>

    <Profile_DatumEntryDefinition>478882298</Profile_Datum
    EntryDefinition>
        <Profile_Images><?xml version="1.0" encoding="utf-
        16"?>
        <ChrysalisImgDats
        xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
        xmlns:xsd="http://www.w3.org/2001/XMLSchema"
        /></Profile_Images>
        <Profile_Abnormal>False</Profile_Abnormal>
        <Profile_AOV>False</Profile_AOV>
        <Profile_StatusStep />
        <Profile_Range />
        <Profile_RangeComment />
        <Profile_CanceledDate />
        <Profile_CanceledBy />
        <Profile_CanceledFor />
        <Profile_ReflexedBy />
        <Profile_Signer>270861542</Profile_Signer>
        <Profile_SignerDate>2020-03-
    09</Profile_SignerDate>
        <Profile_IOrder />
        <Profile_IOrderDate />
        <Profile_IOrderTime />
        <Profile_Complete>True</Profile_Complete>
        <Profile_DataCollect />
        <Profile_Signatures><?xml version="1.0"
        encoding="utf-16"?>
        <ProfileSignatureHistory
        xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
        xmlns:xsd="http://www.w3.org/2001/XMLSchema">

```

```

<Sigs>
  <ProfileSignature>
    <Signer>270861542</Signer>
    <SignerDate>2020-03-09T15:07:50</SignerDate>
    <UnSign>false</UnSign>
  </ProfileSignature>
</Sigs>
</ProfileSignatureHistory></Profile_Signatures>

<Profile_RequireSignature>True</Profile_RequireSignatu
re>
  <Profile_Approver>162355184</Profile_Approver>
  <Profile_ApproveDate>2020-03-
09</Profile_ApproveDate>
  <Profile_Approvals><?xml version="1.0"
encoding="utf-16"?>
<ProfileApprovalHistory
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema">
  <Sigs>
    <ProfileApproval>
      <Signer>162355184</Signer>
      <SignerDate>2020-03-09T14:56:31</SignerDate>
      <UnSign>false</UnSign>
    </ProfileApproval>
  </Sigs>
</ProfileApprovalHistory></Profile_Approvals>

<Profile_RequireApproval>True</Profile_RequireApproval
>
  <Profile_SignerTime>15:07</Profile_SignerTime>
  <Profile_ApproveTime>14:56</Profile_ApproveTime>
  <Profile_WPDLink />
  <Profile_CompleteDate />
  <Profile_CompleteTime />

<Profile_SignerLocation>246584261</Profile_SignerLocat
ion>

<Profile_PendingReagent>False</Profile_PendingReagent>
  <Profile_ORDLink />
  <Profile_QuestionsStash><?xml version="1.0"
encoding="utf-16"?>
<AAOEstasherBase
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema"
/></Profile_QuestionsStash>

```

```

    <Profile_Questions />

    <Profile_QuestionsComplete>True</Profile_QuestionsComplete>
    <Profile_ROrder>I</Profile_ROrder>
    <Profile_ROrderDate>2020-03-09</Profile_ROrderDate>
    <Profile_ROrderTime>14:50</Profile_ROrderTime>
    <Profile_ROrderLab>111200992</Profile_ROrderLab>
    <Profile_Interpretation />

    <Profile_ProfileDefinitionObjectID>416976655</Profile_ProfileDefinitionObjectID>
    <Profile_ProfileDefinitionCreationDate>2018-10-31T10:16:34.31-04:00</Profile_ProfileDefinitionCreationDate>
    <Profile_ProfileDefinitionUpdateTime>2020-03-09T14:48:48.757-04:00</Profile_ProfileDefinitionUpdateTime>
    <Profile_ProfileDefinitionUser1 />
    <Profile_ProfileDefinitionUser2 />
    <Profile_ProfileDefinitionUser3 />
    <Profile_ProfileDefinitionUser4 />

    <Profile_ProfileDefinitionCREATOR>debbieb</Profile_ProfileDefinitionCREATOR>

    <Profile_ProfileDefinitionMODIFIER>Scaron</Profile_ProfileDefinitionMODIFIER>

    <Profile_ProfileDefinitionCode>AMPH</Profile_ProfileDefinitionCode>

    <Profile_ProfileDefinitionDescription>Amphetamine</Profile_ProfileDefinitionDescription>

    <Profile_ProfileDefinitionRETIRED>False</Profile_ProfileDefinitionRETIRED>

    <Profile_ProfileDefinitionDatumEntryDefinition>478882298</Profile_ProfileDefinitionDatumEntryDefinition>
    <Profile_ProfileDefinitionAlwaysReportMessage />
    <Profile_ProfileDefinitionMethodology />

    <Profile_ProfileDefinitionRequireSignature>True</Profile_ProfileDefinitionRequireSignature>
    <Profile_ProfileDefinitionCodeMap />

```

```

<Profile_ProfileDefinitionLabelCount>1</Profile_ProfileDefinitionLabelCount>

<Profile_ProfileDefinitionAutoComplete>False</Profile_ProfileDefinitionAutoComplete>
  <Profile_ProfileDefinitionConstituentPanels><?xml
version="1.0" encoding="utf-16"?>
<GPathConstituentPanels
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema"
/></Profile_ProfileDefinitionConstituentPanels>

<Profile_ProfileDefinitionAllowDuplicate>False</Profile_ProfileDefinitionAllowDuplicate>
  <Profile_ProfileDefinitionPendingMessage />

<Profile_ProfileDefinitionWPLink>False</Profile_ProfileDefinitionWPLink>

<Profile_ProfileDefinitionWPDiscrete>False</Profile_ProfileDefinitionWPDiscrete>

<Profile_ProfileDefinitionWPCyto>False</Profile_ProfileDefinitionWPCyto>

<Profile_ProfileDefinitionApprovalSignature>True</Profile_ProfileDefinitionApprovalSignature>

<Profile_ProfileDefinitionLabelDefinition>318256</Profile_ProfileDefinitionLabelDefinition>

<Profile_ProfileDefinitionWPDiscreteAdvanceOnly>False</Profile_ProfileDefinitionWPDiscreteAdvanceOnly>

<Profile_ProfileDefinitionRequireLocation>False</Profile_ProfileDefinitionRequireLocation>

<Profile_ProfileDefinitionUseDataCSort>False</Profile_ProfileDefinitionUseDataCSort>

<Profile_ProfileDefinitionAllowBatchResult>True</Profile_ProfileDefinitionAllowBatchResult>

<Profile_ProfileDefinitionIsToxicology>True</Profile_ProfileDefinitionIsToxicology>

```

```

<Profile_ProfileDefinitionIsClinical>False</Profile_ProfileDefinitionIsClinical>

<Profile_ProfileDefinitionSuppressOverlap>False</Profile_ProfileDefinitionSuppressOverlap>

<Profile_ProfileDefinitionOutreachCode>AMPH</Profile_ProfileDefinitionOutreachCode>

<Profile_ProfileDefinitionOutreachDescription>Amphetamine</Profile_ProfileDefinitionOutreachDescription>

<Profile_ProfileDefinitionAllowReferenceLab>True</Profile_ProfileDefinitionAllowReferenceLab>

<SourceDefinition_ObjectID>4348467</SourceDefinition_ObjectID>
  <SourceDefinition_CreationDate>2015-04-23T09:31:01.957-04:00</SourceDefinition_CreationDate>
  <SourceDefinition_UpdateTime>2020-09-18T16:59:13.113-04:00</SourceDefinition_UpdateTime>
  <SourceDefinition_User1 />
  <SourceDefinition_User2 />
  <SourceDefinition_User3 />
  <SourceDefinition_User4 />

<SourceDefinition_CREATOR>deannac</SourceDefinition_CREATOR>

<SourceDefinition_MODIFIER>Scaron</SourceDefinition_MODIFIER>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition_Description>

<SourceDefinition_RETIRED>False</SourceDefinition_RETIRED>

<SourceDefinition_LabelCount>1</SourceDefinition_LabelCount>

<SourceDefinition_WPLink>False</SourceDefinition_WPLink>

```

```

<SourceDefinition_WPRequireProfiles>False</SourceDefin
ition_WPRequireProfiles>

<SourceDefinition_WPAllowNOrdering>False</SourceDefini
tion_WPAllowNOrdering>

<SourceDefinition_ForeignAllowNOrdering>False</SourceD
efinition_ForeignAllowNOrdering>

<SourceDefinition_WPAddDiagnosis>False</SourceDefiniti
on_WPAddDiagnosis>

<SourceDefinition_LabelDefinition>280272780</SourceDef
inition_LabelDefinition>
    <SourceDefinition_Sites><?xml version="1.0"
encoding="utf-16"?>
<XStrings xmlns:xsi="http://www.w3.org/2001/XMLSchema-
instance" xmlns:xsd="http://www.w3.org/2001/XMLSchema"
/></SourceDefinition_Sites>

<SourceDefinition_OverdueDefault>8</SourceDefinition_O
verdueDefault>

<SourceDefinition_Multiples>2</SourceDefinition_Multip
les>

<SourceDefinition_AutoMultiples>False</SourceDefinitio
n_AutoMultiples>

<SourceDefinition_GrossReceived>False</SourceDefinitio
n_GrossReceived>

<SourceDefinition_AAOESet>485355213</SourceDefinition_
AAOESet>

<SourceDefinition_IgnoreLabelDefault>False</SourceDefi
nition_IgnoreLabelDefault>

<CaseTypeDefinition_ObjectID>14662557</CaseTypeDefinit
ion_ObjectID>
    <CaseTypeDefinition_CreationDate>2015-06-
25T11:10:41.087-
04:00</CaseTypeDefinition_CreationDate>
    <CaseTypeDefinition_UpdateTime>2020-10-
27T11:57:35.247-04:00</CaseTypeDefinition_UpdateTime>
    <CaseTypeDefinition_User1 />

```

```

    <CaseTypeDefinition_User2 />
    <CaseTypeDefinition_User3 />
    <CaseTypeDefinition_User4 />

<CaseTypeDefinition_CREATOR>0</CaseTypeDefinition_CREATOR>

<CaseTypeDefinition_MODIFIER>Scaron</CaseTypeDefinition_MODIFIER>

<CaseTypeDefinition_Code>MOLC</CaseTypeDefinition_Code>

<CaseTypeDefinition_Description>Molecular</CaseTypeDefinition_Description>

<CaseTypeDefinition_RETIREDA>False</CaseTypeDefinition_RETIREDA>

<CaseTypeDefinition_ReportDefinition>217397414</CaseTypeDefinition_ReportDefinition>

<CaseTypeDefinition_OverdueDays>1</CaseTypeDefinition_OverdueDays>

<CaseTypeDefinition_RequirePreliminary>False</CaseTypeDefinition_RequirePreliminary>

<CaseTypeDefinition_StepDefinition>245594</CaseTypeDefinition_StepDefinition>

<CaseTypeDefinition_RequireSignature>False</CaseTypeDefinition_RequireSignature>

<CaseTypeDefinition_AccessionLetter>M</CaseTypeDefinition_AccessionLetter>

<CaseTypeDefinition_UseProfileReports>False</CaseTypeDefinition_UseProfileReports>

<CaseTypeDefinition_PrelimReportDefinition>217397414</CaseTypeDefinition_PrelimReportDefinition>

<CaseTypeDefinition_WPLink>False</CaseTypeDefinition_WPLink>

```

```

<CaseTypeDefinition_WPSignOnReport>False</CaseTypeDefinition_WPSignOnReport>

<CaseTypeDefinition_WPLoadReport>False</CaseTypeDefinition_WPLoadReport>

<CaseTypeDefinition_WPAllowSupplemental>False</CaseTypeDefinition_WPAllowSupplemental>

<CaseTypeDefinition_LabelDefinition>348382382</CaseTypeDefinition_LabelDefinition>
  <CaseTypeDefinition_WPSignSuppress><?xml
version="1.0" encoding="utf-16"?>
<XStrings xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance" xmlns:xsd="http://www.w3.org/2001/XMLSchema"/></CaseTypeDefinition_WPSignSuppress>

<CaseTypeDefinition_WPSignSuppressSupplemental><?xml
version="1.0" encoding="utf-16"?>
<XStrings xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance" xmlns:xsd="http://www.w3.org/2001/XMLSchema"/></CaseTypeDefinition_WPSignSuppressSupplemental>

<CaseTypeDefinition_NeverPreliminary>False</CaseTypeDefinition_NeverPreliminary>

<CaseTypeDefinition_AbnormalSpecialFlag>348386126</CaseTypeDefinition_AbnormalSpecialFlag>

<CaseTypeDefinition_WPSetSpecialFlag>True</CaseTypeDefinition_WPSetSpecialFlag>

<CaseTypeDefinition_IgnoreLabelDefault>False</CaseTypeDefinition_IgnoreLabelDefault>
  <Constituent_TYPE>IMMUNO</Constituent_TYPE>

<Constituent_ObjectID>506873439</Constituent_ObjectID>
  <Constituent_CreationDate>2020-03-09T14:50:50.123-04:00</Constituent_CreationDate>
  <Constituent_UpdateTime>2020-03-09T14:52:49.107-04:00</Constituent_UpdateTime>

<Constituent_SuperobjectID>506873438</Constituent_SuperobjectID>
  <Constituent_User1 />
  <Constituent_User2 />

```



```

    <Constituent_User3 />
    <Constituent_User4 />
    <Constituent_CREATOR>Scaron</Constituent_CREATOR>

<Constituent_MODIFIER>Scaron</Constituent_MODIFIER>

<Constituent_Reportable>True</Constituent_Reportable>
    <Constituent_Description />

<Constituent_Definition>414331635</Constituent_Definition>
    <Constituent_Result>Negative</Constituent_Result>
    <Constituent_R1 />
    <Constituent_R2 />
    <Constituent_R3 />
    <Constituent_R4 />
    <Constituent_R5 />
    <Constituent_R6 />
    <Constituent_Abnormal>False</Constituent_Abnormal>
    <Constituent_AOV>False</Constituent_AOV>
    <Constituent_Range />
    <Constituent_RangeComment />

<Constituent_IsConfirm>False</Constituent_IsConfirm>

<Constituent_ConstituentDefinitionObjectID>414331635</Constituent_ConstituentDefinitionObjectID>

<Constituent_ConstituentDefinitionCreationDate>2018-10-26T11:55:05.557-04:00</Constituent_ConstituentDefinitionCreationDate>
    <Constituent_ConstituentDefinitionUpdateTime>2020-11-03T10:33:26.493-05:00</Constituent_ConstituentDefinitionUpdateTime>

<Constituent_ConstituentDefinitionCREATOR>Import</Constituent_ConstituentDefinitionCREATOR>

<Constituent_ConstituentDefinitionMODIFIER>Scaron</Constituent_ConstituentDefinitionMODIFIER>

<Constituent_ConstituentDefinitionCode>AMPH</Constituent_ConstituentDefinitionCode>

<Constituent_ConstituentDefinitionDescription>Amphetamine</Constituent_ConstituentDefinitionDescription>

```

```

<Constituent_ConstituentDefinitionRETIRED>False</Constituent_ConstituentDefinitionRETIRED>
  <Constituent_ConstituentDefinitionCutOffValue>250
ng/dl</Constituent_ConstituentDefinitionCutOffValue>

<Constituent_ConstituentDefinitionDetectedString>DETECTED</Constituent_ConstituentDefinitionDetectedString>

<Constituent_ConstituentDefinitionUndetectedString>NOT
DETECTED</Constituent_ConstituentDefinitionUndetectedString>

<Constituent_ConstituentDefinitionIllicit>False</Constituent_ConstituentDefinitionIllicit>

<Constituent_ConstituentDefinitionDetectionWindow>1-2
days</Constituent_ConstituentDefinitionDetectionWindow>

<Constituent_ConstituentDefinitionOrderConfirmationWhen>Detected</Constituent_ConstituentDefinitionOrderConfirmationWhen>

<Constituent_ConstituentDefinitionConfirmation>4788675
97</Constituent_ConstituentDefinitionConfirmation>

<Constituent_ConstituentDefinitionConfirmationCharge>2
59621897</Constituent_ConstituentDefinitionConfirmationCharge>

<Constituent_ConstituentDefinitionIsTextual>False</Constituent_ConstituentDefinitionIsTextual>
  </CaseTable>

</NewDataSet>

```

GET CASES

Endpoint	api/N/GetCases
Purpose	Returns a collection of NucleoLIS case records based on filter parameters passed.

Request Type	GET
Parameters	<p>case_number={case_number}</p> <p>activesOnly={activesOnly}</p> <p>order_date_from={order_date_from}</p> <p>order_date_to={order_date_to}</p> <p>submitting_physicians={submitting_physicians}</p> <p>submitting_groups={submitting_groups}</p> <p>patient_objectid</p> <p>filter_expression={filter_expression}</p> <p>user_ID={user_ID}</p> <p>return_format={return_format}</p> <p>Where:</p> <p>case_number is the requested NucleoLIS case_number. (Wildcards are supported)</p> <p>activesOnly is a true or false value. If true, only active cases will be returned. If false, all cases will be returned.</p> <p>order_date_from is the order from date for filtering cases. Dates must be in the format YYYY-MM-DD.</p> <p>order_date_to is the order to date for filtering cases. Dates must be in the format YYYY-MM-DD.</p> <p>submitting_physicians is a list of submitting physicians.</p> <p>submitting_groups is a list of submitting groups.</p> <p>patient_objectid is the ObjectID of the patient record or patient code.</p> <p>filter_expression is used as a secondary custom filtering criteria applied to the result set . Wildcards are supported. Any valid SQL WHERE clause is accepted. It can be used with any column returned in the result set. Multiple columns can</p>

be used when separated with 'AND'. Filtering criteria supported include: Patient_ObjectID, Patient_Name, Patient_Code, Patient_Dob, Patient_Sex, Case_Number, Case_OrderDate, Case_CreationDate, Case_ReportStatus, Case_Status, Case_ObjectID.

user_id is the user code for the user in the NucleoLIS user table that is requesting the record.

return_format is the format for returning the data. Where:

0 = XML Format

1 = JSON Format

Example

```
https://SomeServer.com/api/N/GetCases?case_number=CL2020*&activesOnly=true&return_format=0&user_id=245654&order_date_from=2020-09-01&order_date_to=2020-11-02&filter_expression=case_status='A' AND patient_sex='Female'&submitting_physicians=*&submitting_groups=*
```

Returns

If the request succeeds, returns a collection of 1..n case records matching the filtering criteria in the following XML format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
<CaseTable>
  <Patient_ObjectID>506408213</Patient_ObjectID>
  <Patient_Name>Boop, Betty</Patient_Name>
  <Patient_Code>506408213</Patient_Code>
  <Patient_Dob>1965-03-23</Patient_Dob>
  <Patient_Sex>Female</Patient_Sex>
  <Case_Number>CL2020-000080</Case_Number>
  <Case_OrderDate>2020-09-17</Case_OrderDate>
  <Case_CreationDate>2020-09-17T17:27:27.473-04:00</Case_CreationDate>
  <Case_ReportStatus>PENDING</Case_ReportStatus>
  <Case_Status>A</Case_Status>
  <Case_ObjectID>514638163</Case_ObjectID>
  <Submitter_PhysicianCode>PYNE</Submitter_PhysicianCode>
  <Submitter_PhysicianName>Pyne, Warren</Submitter_PhysicianName>
</CaseTable>
<CaseTable>
  <Patient_ObjectID>503799913</Patient_ObjectID>
  <Patient_Name>Bell, Tinker</Patient_Name>
```

```

<Patient_Code>503799913</Patient_Code>
<Patient_Dob>1950-01-01</Patient_Dob>
<Patient_Sex>Female</Patient_Sex>
<Case_Number>CL2020-000082</Case_Number>
<Case_OrderDate>2020-09-18</Case_OrderDate>
<Case_CreationDate>2020-09-18T14:53:08.89-04:00</Case_CreationDate>
<Case_ReportStatus>PENDING</Case_ReportStatus>
<Case_Status>A</Case_Status>
<Case_ObjectID>514704625</Case_ObjectID>
<Submitter_GroupCode>PC</Submitter_GroupCode>
<Submitter_PhysicianCode>PHIL</Submitter_PhysicianCode>
<Submitter_GroupName>PrimaCare's, P.C.</Submitter_GroupName>
<Submitter_PhysicianName>Philip, Linsey</Submitter_PhysicianName>
</CaseTable>
<CaseTable>
  <Patient_ObjectID>503799913</Patient_ObjectID>
  <Patient_Name>Bell, Tinker</Patient_Name>
  <Patient_Code>503799913</Patient_Code>
  <Patient_Dob>1950-01-01</Patient_Dob>
  <Patient_Sex>Female</Patient_Sex>
  <Case_Number>CL2020-000084</Case_Number>
  <Case_OrderDate>2020-09-28</Case_OrderDate>
  <Case_CreationDate>2020-09-28T13:43:54.22-04:00</Case_CreationDate>
  <Case_ReportStatus>PENDING</Case_ReportStatus>
  <Case_Status>A</Case_Status>
  <Case_ObjectID>515443274</Case_ObjectID>
  <Submitter_PhysicianCode>PHIL</Submitter_PhysicianCode>
  <Submitter_PhysicianName>Philip, Linsey</Submitter_PhysicianName>
</CaseTable>
<CaseTable>
  <Patient_ObjectID>503625913</Patient_ObjectID>
  <Patient_Name>Accession, Annie</Patient_Name>
  <Patient_Code>503625913</Patient_Code>
  <Patient_Dob>2018-12-11</Patient_Dob>
  <Patient_Sex>Female</Patient_Sex>
  <Case_Number>CL2020-000089</Case_Number>
  <Case_OrderDate>2020-09-30</Case_OrderDate>
  <Case_CreationDate>2020-09-30T13:32:44.8-04:00</Case_CreationDate>
  <Case_ReportStatus>PENDING</Case_ReportStatus>
  <Case_Status>A</Case_Status>
  <Case_ObjectID>515593040</Case_ObjectID>
  <Submitter_PhysicianCode>ROBINS</Submitter_PhysicianCode>
  <Submitter_PhysicianName>Robins, Tim</Submitter_PhysicianName>

```

```

</CaseTable>
<CaseTable>
  <Patient_ObjectID>433509601</Patient_ObjectID>
  <Patient_Name>Friday, Girl</Patient_Name>
  <Patient_Code>433509601</Patient_Code>
  <Patient_Dob>2009-03-13</Patient_Dob>
  <Patient_Sex>Female</Patient_Sex>
  <Case_Number>CL2020-000091</Case_Number>
  <Case_OrderDate>2020-10-05</Case_OrderDate>
  <Case_CreationDate>2020-10-05T15:56:57.673-04:00</Case_CreationDate>
  <Case_ReportStatus>PENDING</Case_ReportStatus>
  <Case_Status>A</Case_Status>
  <Case_ObjectID>515975948</Case_ObjectID>
  <Submitter_PhysicianCode>PYNE</Submitter_PhysicianCode>
  <Submitter_PhysicianName>Pyne, Warren</Submitter_PhysicianName>
</CaseTable>
<CaseTable>
  <Patient_ObjectID>516131995</Patient_ObjectID>
  <Patient_Name>Boop, Betty</Patient_Name>
  <Patient_Code>516131995</Patient_Code>
  <Patient_Dob>1965-03-23</Patient_Dob>
  <Patient_Sex>Female</Patient_Sex>
  <Case_Number>CL2020-000105</Case_Number>
  <Case_OrderDate>2020-10-07</Case_OrderDate>
  <Case_CreationDate>2020-10-07T14:49:44.863-04:00</Case_CreationDate>
  <Case_ReportStatus>PENDING</Case_ReportStatus>
  <Case_Status>A</Case_Status>
  <Case_ObjectID>516132019</Case_ObjectID>
  <Submitter_PhysicianCode>PYNE</Submitter_PhysicianCode>
  <Submitter_PhysicianName>Pyne, Warren</Submitter_PhysicianName>
</CaseTable>
</NewDataSet>

```

GET SPECIMEN

Endpoint api/N/GetSpecimen

Purpose Return a single, full NucleoLIS specimen record based on specimen unique identifier. Note: The full specimen record is inclusive of all specimen level data elements.

Request Type	GET
Parameters	<p>specimen_ID={specimen_ID}</p> <p>user_ID={user_ID}</p> <p>return_format={return_format}</p> <p>Where:</p> <p>specimen_id is the specimen id number of the specimen to return.</p> <p>user_id is the user code for the user in the NucleoLIS user table that is requesting the record.</p> <p>return_format is the format for returning the data.</p> <p>Where:</p> <p>0 = XML Format</p> <p>1 = JSON Format</p>
Example	https://www.SomeServer.com/api/N/GetSpecimen?specimen_id=467291&user_id=245654&return_format=0
Returns	<p>If the request succeeds, returns a specimen record matching the filtering criteria in the following XML format:</p> <pre> <?xml version="1.0" encoding="UTF-8" ?><NewDataSet> <SpecimenTable> <Patient_ObjectID>347631</Patient_ObjectID> <Patient_Name>Farrell, Ivy</Patient_Name> <Patient_Code>347631</Patient_Code> <Patient_Sex>F</Patient_Sex> <Patient_Dob>1970-01-01</Patient_Dob> <Case_Number>2</Case_Number> <Case_Status>I</Case_Status> <Specimen_ObjectID>467291</Specimen_ObjectID> <Specimen_CreationDate>2015-03-27T15:56:44.83-04:00</Specimen_CreationDate> <Specimen_UpdateTime>2015-03-30T07:07:02.71-04:00</Specimen_UpdateTime> <Specimen_SuperobjectID>467290</Specimen_SuperobjectID> </pre>

```

        <Specimen_CREATOR>HIS</Specimen_CREATOR>
        <Specimen_MODIFIER>0</Specimen_MODIFIER>
        <Specimen_CollectionDate>2015-03-
13</Specimen_CollectionDate>

        <Specimen_CollectionTime>16:14</Specimen_CollectionTime>
        <Specimen_ReceivedDate>2015-03-
30</Specimen_ReceivedDate>

        <Specimen_ReceivedTime>07:07:02</Specimen_ReceivedTime>
        <Specimen_Number>467291</Specimen_Number>
        <Specimen_Status>A</Specimen_Status>
        <Specimen_Receiver>0</Specimen_Receiver>
    </SpecimenTable>
</NewDataSet>

```

GET SPECIMENS

Endpoint	api/N/GetSpecimens
Purpose	Returns a collection of NucleoLIS Patient, Case and Specimen level data.
Request Type	GET
Parameters	<pre> received_date_from={received_date_from} received_date_to={received_date_to} activesOnly={activesOnly} specimen_status_steps={specimen_status_steps} specimen_sources={specimen_sources} filter_expression={filter_expression} user_ID={user_ID} return_format={return_format} </pre>

Where:

received_date_from is the "from" date for received specimens. Dates must be in YYYY-MM-DD format.

received_date_to is the "to" date for received specimens. Dates must be in YYYY-MM-DD format.

activesOnly is a true or false value. If true, only active cases are returned. If false, all cases are returned.

specimen_status_steps is a pipe delimited list of specimen status steps.

specimen_sources is a pipe delimited list of specimen sources

filter_expression is used as a secondary custom filtering criteria applied to the result set. Wildcards are supported. It can be used with any column returned in the result set. Any valid SQL WHERE clause is accepted. Multiple columns can be used when separated with 'AND'. Filtering criteria supported include: Patient_ObjectID, Patient_Name, Patient_Code, Patient_Dob, Patient_Sex, Case_Number, Specimen_ObjectID, Specimen_Received date, Specimen_Received time, Specimen_CollectionDate, Specimen_CollectionTime, Specimen_StatusStep, SourceDefinition_Code, and SourceDefinition_Description.

user_id is the user code for the user in the NucleoLIS user table that is requesting the record.

return_format is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

Example

https://www.SomeServer.com/api/N/GetSpecimens?received_date_from=2020-10-01&received_date_to=2020-11-02&activesOnly=true&specimen_status_steps=&specimen_sources=Blood|Serum&filter_expression=Patient_Sex='Male'&user_id=245654&return_format=0

Returns

If the request succeeds, returns a collection of 1..n specimen records matching the filtering criteria in the following XML format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <SpecimenTable>
    <Patient_ObjectID>507179564</Patient_ObjectID>
    <Patient_Name>Bishop, Tom</Patient_Name>
    <Patient_Code>507179564</Patient_Code>
    <Patient_Dob>1963-08-02</Patient_Dob>
    <Patient_Sex>Male</Patient_Sex>
    <Case_Number>507179566OR</Case_Number>
    <Specimen_ReceivedDate>2020-10-
13</Specimen_ReceivedDate>

    <Specimen_ReceivedTime>14:07:31</Specimen_ReceivedTim
e>
    <Specimen_CollectionDate>2020-10-
13</Specimen_CollectionDate>
    <Specimen_CollectionTime />
    <Specimen_StatusStep />

    <SourceDefinition_Code>Blood</SourceDefinition_Code>

    <SourceDefinition_Description>Blood</SourceDefinition
_Description>
    <Specimen_ObjectID>516628021</Specimen_ObjectID>
  </SpecimenTable>
  <SpecimenTable>
    <Patient_ObjectID>514700814</Patient_ObjectID>
    <Patient_Name>Barksalot, Bruno</Patient_Name>
    <Patient_Code>514700814</Patient_Code>
    <Patient_Dob>2016-03-30</Patient_Dob>
    <Patient_Sex>Male</Patient_Sex>
    <Case_Number>CL2020-000086</Case_Number>
    <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

    <Specimen_ReceivedTime>18:15:58</Specimen_ReceivedTim
e>
    <Specimen_CollectionDate>2020-09-
18</Specimen_CollectionDate>

    <Specimen_CollectionTime>13:38</Specimen_CollectionTi
me>

    <SourceDefinition_Code>Blood</SourceDefinition_Code>
```

```

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
  <Specimen_ObjectID>515221988</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
  <Patient_ObjectID>503782367</Patient_ObjectID>
  <Patient_Name>Bishop, Tom</Patient_Name>
  <Patient_Code>503782367</Patient_Code>
  <Patient_Dob>1963-08-02</Patient_Dob>
  <Patient_Sex>Male</Patient_Sex>
  <Case_Number>CL2020-000093</Case_Number>
  <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:15:34</Specimen_ReceivedTim
e>
  <Specimen_CollectionDate>2020-10-
06</Specimen_CollectionDate>
  <Specimen_CollectionTime />
  <Specimen_StatusStep />

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
  <Specimen_ObjectID>516056679</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
  <Patient_ObjectID>511955302</Patient_ObjectID>
  <Patient_Name>Cricket, Jimimy</Patient_Name>
  <Patient_Code>511955302</Patient_Code>
  <Patient_Dob>2018-11-13</Patient_Dob>
  <Patient_Sex>Male</Patient_Sex>
  <Case_Number>CL2020-000110</Case_Number>
  <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:16:15</Specimen_ReceivedTim
e>
  <Specimen_CollectionDate>2020-07-
21</Specimen_CollectionDate>

<Specimen_CollectionTime>17:28</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

```

```

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
  <Specimen_ObjectID>511968657</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
  <Patient_ObjectID>510788319</Patient_ObjectID>
  <Patient_Name>Caron, Jack</Patient_Name>
  <Patient_Code>510788319</Patient_Code>
  <Patient_Dob>2019-06-15</Patient_Dob>
  <Patient_Sex>Male</Patient_Sex>
  <Case_Number>CL2020-000111</Case_Number>
  <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:16:18</Specimen_ReceivedTim
e>
  <Specimen_CollectionDate>2020-06-
08</Specimen_CollectionDate>

<Specimen_CollectionTime>16:20</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
  <Specimen_ObjectID>510788325</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
  <Patient_ObjectID>506987816</Patient_ObjectID>
  <Patient_Name>Duck, Daffy</Patient_Name>
  <Patient_Code>506987816</Patient_Code>
  <Patient_Dob>1950-03-03</Patient_Dob>
  <Patient_Sex>Male</Patient_Sex>
  <Case_Number>CL2020-000113</Case_Number>
  <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:16:20</Specimen_ReceivedTim
e>
  <Specimen_CollectionDate>2020-03-
12</Specimen_CollectionDate>

<Specimen_CollectionTime>14:40</Specimen_CollectionTi
me>

```

```

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
  <Specimen_ObjectID>506987822</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
  <Patient_ObjectID>506407748</Patient_ObjectID>
  <Patient_Name>Cricket, Jimimy</Patient_Name>
  <Patient_Code>506407748</Patient_Code>
  <Patient_Dob>2018-11-13</Patient_Dob>
  <Patient_Sex>Male</Patient_Sex>
  <Case_Number>CL2020-000115</Case_Number>
  <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:16:24</Specimen_ReceivedTim
e>
  <Specimen_CollectionDate>2020-02-
28</Specimen_CollectionDate>

<Specimen_CollectionTime>14:57</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
  <Specimen_ObjectID>506407755</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
  <Patient_ObjectID>298771370</Patient_ObjectID>
  <Patient_Name>Jones, Tom</Patient_Name>
  <Patient_Code>298771370</Patient_Code>
  <Patient_Dob>1944-01-01</Patient_Dob>
  <Patient_Sex>Male</Patient_Sex>
  <Case_Number>CL2020-000155</Case_Number>
  <Specimen_ReceivedDate>2020-10-
21</Specimen_ReceivedDate>

<Specimen_ReceivedTime>17:10:05</Specimen_ReceivedTim
e>
  <Specimen_CollectionDate>2020-10-
21</Specimen_CollectionDate>
  <Specimen_CollectionTime />
  <Specimen_StatusStep />

```

```

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
  <Specimen_ObjectID>517413942</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
  <Patient_ObjectID>517962619</Patient_ObjectID>
  <Patient_Name>Caron, Jack</Patient_Name>
  <Patient_Code>517962619</Patient_Code>
  <Patient_Dob>2019-06-15</Patient_Dob>
  <Patient_Sex>Male</Patient_Sex>
  <Case_Number>CL2020-000156</Case_Number>
  <Specimen_ReceivedDate>2020-10-
27</Specimen_ReceivedDate>

<Specimen_ReceivedTime>10:50:07</Specimen_ReceivedTim
e>
  <Specimen_CollectionDate>2020-10-
27</Specimen_CollectionDate>

<Specimen_CollectionTime>10:46</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
  <Specimen_ObjectID>517962628</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
  <Patient_ObjectID>518533883</Patient_ObjectID>
  <Patient_Name>Caron, Jack</Patient_Name>
  <Patient_Code>518533883</Patient_Code>
  <Patient_Dob>2019-06-15</Patient_Dob>
  <Patient_Sex>Male</Patient_Sex>
  <Case_Number>CL2020-000160</Case_Number>
  <Specimen_ReceivedDate>2020-11-
02</Specimen_ReceivedDate>

<Specimen_ReceivedTime>11:05:39</Specimen_ReceivedTim
e>
  <Specimen_CollectionDate>2020-11-
02</Specimen_CollectionDate>

```

<Specimen_CollectionTime>10:59</Specimen_CollectionTime>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition_Description>

<Specimen_ObjectID>518533891</Specimen_ObjectID>

</SpecimenTable>

<SpecimenTable>

<Patient_ObjectID>516213756</Patient_ObjectID>

<Patient_Name>Caron, Jack</Patient_Name>

<Patient_Code>516213756</Patient_Code>

<Patient_Dob>2019-06-15</Patient_Dob>

<Patient_Sex>Male</Patient_Sex>

<Case_Number>CLI2020-000002</Case_Number>

<Specimen_ReceivedDate>2020-10-08</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:15:11</Specimen_ReceivedTime>

<Specimen_CollectionDate>2020-10-08</Specimen_CollectionDate>

<Specimen_CollectionTime>18:10</Specimen_CollectionTime>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition_Description>

<Specimen_ObjectID>516213763</Specimen_ObjectID>

</SpecimenTable>

<SpecimenTable>

<Patient_ObjectID>516287111</Patient_ObjectID>

<Patient_Name>Caron, Jack</Patient_Name>

<Patient_Code>516287111</Patient_Code>

<Patient_Dob>2019-06-15</Patient_Dob>

<Patient_Sex>Male</Patient_Sex>

<Case_Number>CLI2020-000005</Case_Number>

<Specimen_ReceivedDate>2020-10-09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>14:45:44</Specimen_ReceivedTime>

```

        <Specimen_CollectionDate>2020-10-
09</Specimen_CollectionDate>

<Specimen_CollectionTime>14:43</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
    <Specimen_ObjectID>516287120</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
    <Patient_ObjectID>516287344</Patient_ObjectID>
    <Patient_Name>Cricket, Jimimy</Patient_Name>
    <Patient_Code>516287344</Patient_Code>
    <Patient_Dob>2018-11-13</Patient_Dob>
    <Patient_Sex>Male</Patient_Sex>
    <Case_Number>CLI2020-000006</Case_Number>
    <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>15:44:18</Specimen_ReceivedTim
e>
    <Specimen_CollectionDate>2020-10-
09</Specimen_CollectionDate>

<Specimen_CollectionTime>14:46</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
    <Specimen_ObjectID>516287379</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
    <Patient_ObjectID>516291416</Patient_ObjectID>
    <Patient_Name>Duck, Daffy</Patient_Name>
    <Patient_Code>516291416</Patient_Code>
    <Patient_Dob>1950-03-03</Patient_Dob>
    <Patient_Sex>Male</Patient_Sex>
    <Case_Number>CLI2020-000009</Case_Number>
    <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

```



```

<Specimen_ReceivedTime>18:11:49</Specimen_ReceivedTime>
    <Specimen_CollectionDate>2020-10-
09</Specimen_CollectionDate>

<Specimen_CollectionTime>15:54</Specimen_CollectionTime>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
    <Specimen_ObjectID>516291421</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
    <Patient_ObjectID>516276481</Patient_ObjectID>
    <Patient_Name>Duck, Daffy</Patient_Name>
    <Patient_Code>516276481</Patient_Code>
    <Patient_Dob>1950-03-03</Patient_Dob>
    <Patient_Sex>Male</Patient_Sex>
    <Case_Number>CLI2020-000010</Case_Number>
    <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:12:16</Specimen_ReceivedTime>
    <Specimen_CollectionDate>2020-10-
09</Specimen_CollectionDate>

<Specimen_CollectionTime>11:46</Specimen_CollectionTime>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
    <Specimen_ObjectID>516276511</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
    <Patient_ObjectID>516298160</Patient_ObjectID>
    <Patient_Name>QA, Final</Patient_Name>
    <Patient_Code>516298160</Patient_Code>
    <Patient_Dob>2019-10-23</Patient_Dob>
    <Patient_Sex>Male</Patient_Sex>
    <Case_Number>CLI2020-000013</Case_Number>

```

```

        <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:13:24</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-10-
09</Specimen_CollectionDate>

<Specimen_CollectionTime>17:45</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>516298169</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
        <Patient_ObjectID>518078029</Patient_ObjectID>
        <Patient_Name>QA, Final</Patient_Name>
        <Patient_Code>518078029</Patient_Code>
        <Patient_Dob>2019-10-23</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>CLI2020-000016</Case_Number>
        <Specimen_ReceivedDate>2020-10-
28</Specimen_ReceivedDate>

<Specimen_ReceivedTime>15:48:59</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-10-
28</Specimen_CollectionDate>

<Specimen_CollectionTime>15:44</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>518078227</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
        <Patient_ObjectID>510676219</Patient_ObjectID>
        <Patient_Name>Cricket, Jiminy</Patient_Name>
        <Patient_Code>510676219</Patient_Code>
        <Patient_Dob>2011-07-01</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>

```

```

        <Case_Number>FC2020-000002</Case_Number>
        <Specimen_ReceivedDate>2020-10-
13</Specimen_ReceivedDate>

<Specimen_ReceivedTime>08:30:47</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-06-
04</Specimen_CollectionDate>

<Specimen_CollectionTime>15:26</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>510676225</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
        <Patient_ObjectID>516298310</Patient_ObjectID>
        <Patient_Name>QA, Final</Patient_Name>
        <Patient_Code>516298310</Patient_Code>
        <Patient_Dob>2019-10-23</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>HE2020-000003</Case_Number>
        <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:14:16</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-10-
09</Specimen_CollectionDate>

<Specimen_CollectionTime>17:49</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>516298330</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
        <Patient_ObjectID>506647493</Patient_ObjectID>
        <Patient_Name>Mouse, Mickey</Patient_Name>
        <Patient_Code>506647493</Patient_Code>
        <Patient_Dob>2000-07-08</Patient_Dob>

```

```

        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>M2020-000030</Case_Number>
        <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:16:21</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-03-
04</Specimen_CollectionDate>

<Specimen_CollectionTime>16:09</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>506647513</Specimen_ObjectID>
    </SpecimenTable>
    <SpecimenTable>
        <Patient_ObjectID>510588219</Patient_ObjectID>
        <Patient_Name>Covid, Charlie</Patient_Name>
        <Patient_Code>510588219</Patient_Code>
        <Patient_Dob>2019-12-15</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>M2020-000059</Case_Number>
        <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:36:05</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-06-
01</Specimen_CollectionDate>

<Specimen_CollectionTime>10:00</Specimen_CollectionTi
me>
        <Specimen_StatusStep />

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>510589737</Specimen_ObjectID>
    </SpecimenTable>
</NewDataSet>

```

GET TEST ORDER

Endpoint `api/N/GetTestOrder`

Purpose Returns a single, full test record and available results based on test order unique identifier.

Request Type GET

Parameters `test_order_ID={test_order_ID}`
`user_ID={user_ID}`
`return_format={return_format}`

Where:

test_order_id is the id of the test order to return.

user_id is the user code for the user in the NucleoLIS user table that is requesting the record.

return_format is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

Example `https://www.SomeServer.com/api/N/GetTestOrder?test_order_id=8606534&user_id=245654&return_format=0`

Returns If the request succeeds, it returns the requested full profile test order record inclusive of all test order and available test result data in the following XML format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <ProfileTable>
    <Patient_ObjectID>4180414</Patient_ObjectID>
    <Patient_Name>Wednesday, Wendy</Patient_Name>
    <Patient_Code>4180414</Patient_Code>
    <Patient_Sex />
    <Patient_Dob />
    <Case_Number>PSY2015-000042</Case_Number>
```

```

        <Case_Status>I</Case_Status>
        <Specimen_ObjectID>8606398</Specimen_ObjectID>
        <Specimen_CollectionDate />
        <Specimen_CollectionTime />
        <Specimen_ReceivedDate>2015-05-
19</Specimen_ReceivedDate>

<Specimen_ReceivedTime>15:37:45</Specimen_ReceivedTime>

        <Specimen_StatusStep />

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Profile_TYPE>FISH</Profile_TYPE>
        <Profile_ObjectID>8606534</Profile_ObjectID>
        <Profile_CreationDate>2015-05-19T15:37:18.65-
04:00</Profile_CreationDate>
        <Profile_UpdateTime>2015-05-19T15:42:05.01-
04:00</Profile_UpdateTime>

<Profile_SuperobjectID>8606398</Profile_SuperobjectID
>
        <Profile_User1 />
        <Profile_User2 />
        <Profile_User3 />
        <Profile_User4 />
        <Profile_CREATOR>deannac</Profile_CREATOR>
        <Profile_MODIFIER>deannac</Profile_MODIFIER>
        <Profile_Description />
        <Profile_OrderDate>2015-05-19</Profile_OrderDate>
        <Profile_OrderTime>15:37</Profile_OrderTime>
        <Profile_OrderedBy>deannac</Profile_OrderedBy>

<Profile_PanelDefinition>249815</Profile_PanelDefinit
ion>
        <Profile_R1>Clinical Significance: ROS1 gene
rearrangements are found in 1-2% of non-small cell
lung carcinoma (NSCLC). Pre-clinical and early
clinical evidence suggests ROS1-rearranged tumors may
be sensitive to the dual ALK/MET inhibitor
crizotinib.</Profile_R1>
        <Profile_R2 />
        <Profile_R3 />
        <Profile_R4 />
        <Profile_R5 />

```

```

    <Profile_R6 />
    <Profile_Result>Present</Profile_Result>
    <Profile_Reportable>True</Profile_Reportable>
    <Profile_DatumDef>245550</Profile_DatumDef>

<Profile_DatumEntryDefinition>4380952</Profile_DatumE
ntryDefinition>
    <Profile_Images><?xml version="1.0"
encoding="utf-16"?>
    <ChrysalisImgDats
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema">
    <ImgDat />
</ChrysalisImgDats></Profile_Images>
    <Profile_Abnormal>False</Profile_Abnormal>
    <Profile_AOV>False</Profile_AOV>
    <Profile_StatusStep />
    <Profile_Range />
    <Profile_RangeComment />
    <Profile_CanceledDate />
    <Profile_CanceledBy />
    <Profile_CanceledFor />
    <Profile_ReflexedBy />
    <Profile_IOrder />
    <Profile_IOrderDate />
    <Profile_IOrderTime />
    <Profile_Complete>True</Profile_Complete>

<Profile_RequireSignature>>false</Profile_RequireSigna
ture>

<Profile_RequireApproval>>false</Profile_RequireApprov
al>

<Profile_ProfileDefinitionObjectID>249815</Profile_Pr
ofileDefinitionObjectID>
    <Profile_ProfileDefinitionCreationDate>2014-12-
31T10:55:34.78-
05:00</Profile_ProfileDefinitionCreationDate>
    <Profile_ProfileDefinitionUpdateTime>2015-10-
01T15:08:20.637-
04:00</Profile_ProfileDefinitionUpdateTime>
    <Profile_ProfileDefinitionUser1 />
    <Profile_ProfileDefinitionUser2 />
    <Profile_ProfileDefinitionUser3 />
    <Profile_ProfileDefinitionUser4 />

```

<Profile_ProfileDefinitionCREATOR>0</Profile_ProfileDefinitionCREATOR>

<Profile_ProfileDefinitionMODIFIER>deannac</Profile_ProfileDefinitionMODIFIER>

<Profile_ProfileDefinitionCode>NSmCCA</Profile_ProfileDefinitionCode>

<Profile_ProfileDefinitionDescription>Non-Small
Cell Carcinoma,
Lung</Profile_ProfileDefinitionDescription>

<Profile_ProfileDefinitionRETIRED>False</Profile_ProfileDefinitionRETIRED>

<Profile_ProfileDefinitionDatumLink>245550</Profile_ProfileDefinitionDatumLink>

<Profile_ProfileDefinitionConstituentDatumLink>342477
12</Profile_ProfileDefinitionConstituentDatumLink>

<Profile_ProfileDefinitionDatumEntryDefinition>347149
41</Profile_ProfileDefinitionDatumEntryDefinition>

<Profile_ProfileDefinitionAlwaysReportMessage>ROS1
FISH uses a dual-color break-apart probe which
produces single green, red, and yellow or fusion
signals in rearranged cells.

References:

Kwak EL et al. New Engl J Med. 2010;363(18):1693-1703.

Vysis ALK Break Apart FISH Probe Kit package insert. Vysis ALK United States website. Accessed 12-31-2012.

Bergethon K, Shaw AT, Ignatius Ou S, et al. J Clin Oncol. 2012;30(8):863-870.

Takeuchi K, Soda M, Togashi Y, et al. Nature Med. 2012;18(3):378-381.

Davies KD, Le AT, Theodoro MF, et al. Clin Cancer Res. 2012;18(17):4570-4579.

Shaw AT, Camidge DR, Engelman JA, et al. J Clin Oncol. 2012;30 (suppl): abstract 7508.

Bang YJ. Arch Pathol Lab Med. 2012;136:1201-1204.

NCCN Guidelines, Non-Small Cell Lung Cancer, version 1.2013. National Comprehensive Cancer Network website. Accessed 12-31-2012.

ALK page. Atlas of Genetics and Cytogenetics in Oncology and Haematology. Accessed 12-31-2012; last updated 2-2010.

Stumpfova M, Janne PA. Clin Cancer Res. 2012;18(16):4222-4224.

Sasaki T et al. Eur J Cancer. 2010;46:1773-1780.
</Profile_ProfileDefinitionAlwaysReportMessage>

<Profile_ProfileDefinitionMethodology>FISH</Profile_ProfileDefinitionMethodology>

<Profile_ProfileDefinitionCodeMap />

<Profile_ProfileDefinitionLabelCount>1</Profile_ProfileDefinitionLabelCount>

<Profile_ProfileDefinitionAutoComplete>False</Profile_ProfileDefinitionAutoComplete>

<Profile_ProfileDefinitionWPDDiscreteAdvanceOnly>>false</Profile_ProfileDefinitionWPDDiscreteAdvanceOnly>

<Profile_DatumDefinitionObjectID>245550</Profile_DatumDefinitionObjectID>

<Profile_DatumDefinitionCreationDate>2014-12-31T10:54:49.03-05:00</Profile_DatumDefinitionCreationDate>

<Profile_DatumDefinitionUpdateTime>2019-04-09T15:17:02.137-04:00</Profile_DatumDefinitionUpdateTime>

<Profile_DatumDefinitionUser1 />

<Profile_DatumDefinitionUser2 />

<Profile_DatumDefinitionUser3 />

<Profile_DatumDefinitionUser4 />

<Profile_DatumDefinitionCREATOR>0</Profile_DatumDefinitionCREATOR>

<Profile_DatumDefinitionMODIFIER>0</Profile_DatumDefinitionMODIFIER>

<Profile_DatumDefinitionCode>416980</Profile_DatumDefinitionCode>

<Profile_DatumDefinitionDescription />

```
<Profile_DatumDefinitionRETIRED>False</Profile_DatumDefinitionRETIRED>
```

```
<Profile_DatumDefinitionR1>Comment</Profile_DatumDefinitionR1>
```

```
<Profile_DatumDefinitionR2 />
```

```
<Profile_DatumDefinitionR3 />
```

```
<Profile_DatumDefinitionR4 />
```

```
<Profile_DatumDefinitionR5 />
```

```
<Profile_DatumDefinitionR6 />
```

```
<Profile_DatumDefinitionHideResult>False</Profile_DatumDefinitionHideResult>
```

```
<Profile_DatumDefinitionRange />
```

```
<Profile_DatumDefinitionRangeComment />
```

```
<Profile_DatumDefinitionExtensions><?xml  
version="1.0" encoding="utf-16"?>
```

```
<DatumExtensionDefinition
```

```
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
```

```
xmlns:xsd="http://www.w3.org/2001/XMLSchema">
```

```
<Items>
```

```
<DatumExtensionItemDefinition>
```

```
<FieldName>Result</FieldName>
```

```
<IsReadOnly>>false</IsReadOnly>
```

```
</DatumExtensionItemDefinition>
```

```
<DatumExtensionItemDefinition>
```

```
<FieldName>R1</FieldName>
```

```
<IsReadOnly>>false</IsReadOnly>
```

```
</DatumExtensionItemDefinition>
```

```
<DatumExtensionItemDefinition>
```

```
<FieldName>R2</FieldName>
```

```
<IsReadOnly>>false</IsReadOnly>
```

```
</DatumExtensionItemDefinition>
```

```
<DatumExtensionItemDefinition>
```

```
<FieldName>R3</FieldName>
```

```
<IsReadOnly>>false</IsReadOnly>
```

```
</DatumExtensionItemDefinition>
```

```
<DatumExtensionItemDefinition>
```

```
<FieldName>R4</FieldName>
```

```
<IsReadOnly>>false</IsReadOnly>
```

```
</DatumExtensionItemDefinition>
```

```
<DatumExtensionItemDefinition>
```

```
<FieldName>R5</FieldName>
```

```
<IsReadOnly>>false</IsReadOnly>
```

```
</DatumExtensionItemDefinition>
```

```
<DatumExtensionItemDefinition>
```

```

        <FieldName>R6</FieldName>
        <IsReadOnly>>false</IsReadOnly>
    </DatumExtensionItemDefinition>
    <DatumExtensionItemDefinition>
        <FieldName>Range</FieldName>
        <IsReadOnly>>true</IsReadOnly>
    </DatumExtensionItemDefinition>
    <DatumExtensionItemDefinition>
        <FieldName>RangeComment</FieldName>
        <IsReadOnly>>true</IsReadOnly>
    </DatumExtensionItemDefinition>
</Items>
</DatumExtensionDefinition></Profile_DatumDefinitionE
xtensions>
    <Profile_DatumDefinitionTableDisplayDefinition />
    <Profile_DatumDefinitionResultFieldMap />
    <Profile_DatumDefinitionInstrumentAutoOrder />
    <Constituent_TYPE>FISH</Constituent_TYPE>

<Constituent_ObjectID>8606535</Constituent_ObjectID>
    <Constituent_CreationDate>2015-05-
19T15:37:18.843-04:00</Constituent_CreationDate>
    <Constituent_UpdateTime>2015-05-19T15:39:04.96-
04:00</Constituent_UpdateTime>

<Constituent_SuperobjectID>8606534</Constituent_Super
objectID>
    <Constituent_User1 />
    <Constituent_User2 />
    <Constituent_User3 />
    <Constituent_User4 />

<Constituent_CREATOR>deannac</Constituent_CREATOR>

<Constituent_MODIFIER>deannac</Constituent_MODIFIER>

<Constituent_Reportable>True</Constituent_Reportable>
    <Constituent_Description />

<Constituent_Definition>249789</Constituent_Definitio
n>
    <Constituent_Result>Present</Constituent_Result>
    <Constituent_R1 />
    <Constituent_R2 />
    <Constituent_R3 />
    <Constituent_R4 />
    <Constituent_R5 />

```

```

        <Constituent_R6 />

    <Constituent_Abnormal>False</Constituent_Abnormal>
        <Constituent_AOV>False</Constituent_AOV>
        <Constituent_Range />
        <Constituent_RangeComment />

    <Constituent_ConstituentDefinitionObjectID>249789</Co
nstituent_ConstituentDefinitionObjectID>

    <Constituent_ConstituentDefinitionCreationDate>2014-
12-31T10:55:34.527-
05:00</Constituent_ConstituentDefinitionCreationDate>

    <Constituent_ConstituentDefinitionUpdateTime>2014-12-
31T10:55:34.527-
05:00</Constituent_ConstituentDefinitionUpdateTime>
        <Constituent_ConstituentDefinitionUser1 />
        <Constituent_ConstituentDefinitionUser2 />
        <Constituent_ConstituentDefinitionUser3 />
        <Constituent_ConstituentDefinitionUser4 />

    <Constituent_ConstituentDefinitionCREATOR>0</Constitu
ent_ConstituentDefinitionCREATOR>

    <Constituent_ConstituentDefinitionMODIFIER>0</Constit
uent_ConstituentDefinitionMODIFIER>

    <Constituent_ConstituentDefinitionCode>ROS1</Constitu
ent_ConstituentDefinitionCode>

    <Constituent_ConstituentDefinitionDescription>break-
apart probe at
(6q22)</Constituent_ConstituentDefinitionDescription>

    <Constituent_ConstituentDefinitionRETIRED>False</Cons
tituent_ConstituentDefinitionRETIRED>
        <Constituent_ConstituentDefinitionCodeMap />
        </ProfileTable>

</NewDataSet>

```

GET TEST ORDERS

Endpoint	api/N/GetTestOrders
Purpose	Returns a collection of Patient, Case, Specimen and Test level data including Patient_ObjectID, Patient_Name, Patient_Code, Patient_Dob, Patient_Sex, Case_Number, Specimen_ObjectID, Specimen_Received date, Specimen_Received time, Specimen_CollectionDate, Specimen_CollectionTime, Specimen_StatusStep, SourceDefinition_Code, SourceDefinition_Description, Test order date, status, completion status, test code and unique identifier of all ordered tests records meeting filtering requirements.

.

Request Type	GET
---------------------	-----

Parameters	<code>order_date_from={order_date_from}</code> <code>order_date_to={order_date_to}</code> <code>activesOnly={activesOnly}</code> <code>order_status_steps={order_status_steps}</code> <code>order_codes={order_codes}</code> <code>filter_expression={filter_expression}</code> <code>user_ID={user_ID}</code> <code>return_format={return_format}</code>
-------------------	--

Where:

order_date_from is the "from" order date. Dates must be in YYYY-MM-DD format.

order_date_to is the "to" order date. Dates must be in YYYY-MM-DD format.

activesOnly is a true or false value. If set to true, only active orders records are returned if false, all order records are returned.

order_status_steps is pipe delimited list of order status steps

order_codes is a pipe delimited list of order codes

filter_expression is used as a secondary custom filtering criteria applied to the result set . Wildcards are supported. It can be used with any column returned in the result set. Any valid SQL WHERE clause is accepted. Multiple columns can be used when separated with 'AND'. Filtering criteria supported include: Patient_ObjectID, Patient_Name, Patient_Code, Patient_DOB, Patient_Sex, Case_Number, Specimen_ObjectID, Specimen_RetreivalDate, Specimen_RetreivalTime, Specimen_CollectionDate, Specimen_CollectionTime, Specimen_StatusStep, SourceDefinition_Code, SourceDefinition_Description, Profile_ObjectID, Profile_OrderDate, Profile_OrderTime, Profile_ProfileDefinitionCode, Profile_ProfileDefinitionDescription, Profile_Complete, Profile_CompleteDate, Profile_SignerDate, and Profile_StatusStep.

user_id is the user code for the user in the NucleoLIS user table that is requesting the record.

return_format is the format for returning the data. Where:

0 = XML Format

1 = JSON Format

Example

```
https://www.SomeServer.com/api/N/GetTestOrders?order_date_from=2020-09-17&order_date_to=2020-09-18&activesOnly=true&order_status_steps=*&order_codes=CDI&filter_expression=Patient_DOB like '1960%*&user_ID=245654&return_format=0
```

Returns

If the request succeeds, returns a collection of 1..n test order profile records matching the filtering criteria in the following XML format:

```
<?xml version="1.0" encoding="UTF-8"?><NewDataSet>
<ProfileTable>
  <Patient_ObjectID>506646654</Patient_ObjectID>
  <Patient_Name>Buckley, Betty Sue Joe</Patient_Name>
```

```

<Patient_Code>506646654</Patient_Code>
<Patient_Dob>1960-07-14</Patient_Dob>
<Patient_Sex>Female</Patient_Sex>
<Case_Number>CL2020-000013</Case_Number>
<Specimen_ObjectID>506686900</Specimen_ObjectID>
<Specimen_ReceivedDate>2020-03-04</Specimen_ReceivedDate>
<Specimen_ReceivedTime>15:49:27</Specimen_ReceivedTime>
<Specimen_CollectionDate>2020-03-04</Specimen_CollectionDate>
<Specimen_CollectionTime>15:46</Specimen_CollectionTime>
<Specimen_StatusStep />
<SourceDefinition_Code>Blood</SourceDefinition_Code>
<SourceDefinition_Description>Blood</SourceDefinition_Description>
<Profile_ObjectID>514633950</Profile_ObjectID>
<Profile_OrderDate>2020-09-17</Profile_OrderDate>
<Profile_OrderTime>16:05</Profile_OrderTime>
<Profile_ProfileDefinitionCode>CDCI</Profile_ProfileDefinitionCode>

<Profile_ProfileDefinitionDescription>CDCI</Profile_ProfileDefinitionDescription>
</ProfileTable>
<ProfileTable>
  <Patient_ObjectID>291217280</Patient_ObjectID>
  <Patient_Name>Ashbey, Clywd</Patient_Name>
  <Patient_Code>291217280</Patient_Code>
  <Patient_Dob>1960-07-07</Patient_Dob>
  <Patient_Sex>Male</Patient_Sex>
  <Case_Number>CL2020-000068</Case_Number>
  <Specimen_ObjectID>511969438</Specimen_ObjectID>
  <Specimen_ReceivedDate>2020-08-06</Specimen_ReceivedDate>
  <Specimen_ReceivedTime>08:22:09</Specimen_ReceivedTime>
  <Specimen_CollectionDate>2020-08-06</Specimen_CollectionDate>
  <Specimen_CollectionTime />
  <Specimen_StatusStep />
  <SourceDefinition_Code>Blood</SourceDefinition_Code>
  <SourceDefinition_Description>Blood</SourceDefinition_Description>
  <Profile_ObjectID>514715461</Profile_ObjectID>
  <Profile_OrderDate>2020-09-18</Profile_OrderDate>
  <Profile_OrderTime>18:21</Profile_OrderTime>
  <Profile_ProfileDefinitionCode>CDCI</Profile_ProfileDefinitionCode>

```

<Profile_ProfileDefinitionDescription>CDCl</Profile_ProfileDefinitionDescription>

<Profile_Complete>False</Profile_Complete>

<Profile_CompleteDate />

<Profile_SignerDate />

<Profile_StatusStep />

</ProfileTable>

<ProfileTable>

<Patient_ObjectID>132099125</Patient_ObjectID>

<Patient_Name>Miller, Leonard</Patient_Name>

<Patient_Code>132099125</Patient_Code>

<Patient_Dob>1960-12-12</Patient_Dob>

<Patient_Sex>Male</Patient_Sex>

<Case_Number>CL2020-000079</Case_Number>

<Specimen_ObjectID>514704804</Specimen_ObjectID>

<Specimen_ReceivedDate>2020-09-18</Specimen_ReceivedDate>

<Specimen_ReceivedTime>16:57:21</Specimen_ReceivedTime>

<Specimen_CollectionDate>2020-09-18</Specimen_CollectionDate>

<Specimen_CollectionTime />

<Specimen_StatusStep />

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition_Description>

<Profile_ObjectID>514704809</Profile_ObjectID>

<Profile_OrderDate>2020-09-18</Profile_OrderDate>

<Profile_OrderTime>14:56</Profile_OrderTime>

<Profile_ProfileDefinitionCode>CDCl</Profile_ProfileDefinitionCode>

<Profile_ProfileDefinitionDescription>CDCl</Profile_ProfileDefinitionDescription>

<Profile_Complete>False</Profile_Complete>

<Profile_CompleteDate />

<Profile_SignerDate />

<Profile_StatusStep />

</ProfileTable>

<ProfileTable>

<Patient_ObjectID>376500</Patient_ObjectID>

<Patient_Name>Farrell, Tiffany</Patient_Name>

<Patient_Code>376500</Patient_Code>

<Patient_Dob>1960-07-01</Patient_Dob>

<Patient_Sex>Female</Patient_Sex>

<Case_Number>PSY2015-000029</Case_Number>

<Specimen_ObjectID>4350748</Specimen_ObjectID>


```

<Specimen_ReceivedDate>2015-06-02</Specimen_ReceivedDate>
<Specimen_ReceivedTime>07:40:06</Specimen_ReceivedTime>
<Specimen_CollectionDate>2015-04-22</Specimen_CollectionDate>
<Specimen_CollectionTime>10:00</Specimen_CollectionTime>
<Specimen_StatusStep />
<SourceDefinition_Code>Blood</SourceDefinition_Code>
<SourceDefinition_Description>Blood</SourceDefinition_Description>
<Profile_ObjectID>514703397</Profile_ObjectID>
<Profile_OrderDate>2020-09-18</Profile_OrderDate>
<Profile_OrderTime>14:30</Profile_OrderTime>
<Profile_ProfileDefinitionCode>CDCl</Profile_ProfileDefinitionCode>

<Profile_ProfileDefinitionDescription>CDCl</Profile_ProfileDefinitionDescription>
</ProfileTable>
</NewDataSet>

```

GET PHYSICIAN

Endpoint	api/N/GetPhysician
Purpose	Returns a single, full physician record based on a unique code.
Request Type	GET
Parameters	<pre> physician_code={physician_code} user_ID={user_ID} return_format={return_format} </pre>

Where:

physician_code is the unique identifier for the requested physician.

user_id is the user code for the user in the NucleoLIS user table that is requesting the record.

return_format is the format for returning the data.
Where:

0 = XML Format

1 = JSON Format

Example

`https://www.SomeServer.com/
api/N/GetPhysician?physician_code=381596&user_ID=24565
4&return_format=0`

Returns

If the request succeeds, returns the requested physician record in the following format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <PhysicianTable>
    <Physician_ObjectID>230022</Physician_ObjectID>
    <Physician_CreationDate>2004-07-30T12:14:57.687-
04:00</Physician_CreationDate>
    <Physician_UpdateTime>2020-01-21T11:25:23.293-
05:00</Physician_UpdateTime>
    <Physician_Name>test, physician</Physician_Name>
    <Physician_Sex />
    <Physician_Street />
    <Physician_City />
    <Physician_State />
    <Physician_Zip />
    <Physician_Phone />
    <Physician_Code>TEST1</Physician_Code>
    <Physician_Soundex>TESTPHYSCN</Physician_Soundex>
    <Physician_UPIN />
    <Physician_STOP />
    <Physician_STREET2 />
    <Physician_RETIRED>False</Physician_RETIRED>
    <Physician_GROUP>False</Physician_GROUP>
    <Physician_STREET3 />
    <Physician_MODIFIER>0</Physician_MODIFIER>
    <Physician_NPI />
    <Physician_WebPathStartDate />

    <Physician_WithholdGroupResults>False</Physician_With
holdGroupResults>
    <Physician_WebPathReturnAs />
    <Physician_WebCopyTo>False</Physician_WebCopyTo>
```

```

        <Physician_Instructions />

<Physician_WebPathSupP>False</Physician_WebPathSupP>

<Physician_WebPathALL>False</Physician_WebPathALL>
    <Physician_WebPathReturnCC><ALL MY
GROUPS></Physician_WebPathReturnCC>

<Physician_UsePreliminaryDistribution>False</Physicia
n_UsePreliminaryDistribution>
    <Physician_PreliminaryReportDistribution><?xml
version="1.0" encoding="utf-16"?>
<DistributionDefinition
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema">
    <AsSubmtterDistribution>
        <PrintEnabled>>false</PrintEnabled>
        <FaxEnabled>>false</FaxEnabled>
        <SecureMailEnabled>>false</SecureMailEnabled>
    </AsSubmtterDistribution>
    <AsOtherDistribution>
        <PrintEnabled>>false</PrintEnabled>
        <FaxEnabled>>false</FaxEnabled>
        <SecureMailEnabled>>false</SecureMailEnabled>
    </AsOtherDistribution>
</DistributionDefinition></Physician_PreliminaryRepor
tDistribution>
    </PhysicianTable>

</NewDataSet>

```

GET PHYSICIANS

Endpoint	api/N/GetPhysicians
Purpose	Returns a collection of 1..n physician records based on filter parameters passed.
Request Type	GET
Parameters	last_name={last_name}

```
first_name={first_name}

filter_expression={filter_expression}

user_id={user_id}

return_format={return_format}

with_retired={true/false}
```

Where:

last_name is the last name of the Physician. Wildcards are supported.

first_name is the first name of the physician. Wildcards are supported.

filter_expression is is used as a secondary custom filtering criteria applied to the result set . Wildcards are supported. Any valid SQL WHERE clause is accepted. It can be used with any column returned in the result set. Multiple columns can be used when separated with 'AND'.

user_id is the user code for the user in the NucleoLIS user table that is requesting the record.

return_format is the format for returning the data.
Where:

0 = XML Format

1 = JSON Format

with_retired is a true or false value. If false, only non-retired physicians are returned. If true, all physicians are returned.

Example

```
https://www.SomeServer.com/
api/N/GetPhysicians?last_name=*&first_name=T*&filter_e
xpression=Physician_GROUP='false'&user_id=245654&retur
n_format=0&with_retired=false
```

Returns

If the request succeeds, returns the requested physician records in the following format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
```

```

<PhysicianTable>
  <Physician_ObjectID>15460945</Physician_ObjectID>
  <Physician_Name>Jones, Tim</Physician_Name>
  <Physician_Street>This physician has a single
"none" location called Office</Physician_Street>
  <Physician_STREET2 />
  <Physician_STREET3 />
  <Physician_City />
  <Physician_State />
  <Physician_Zip />
  <Physician_Code>JONES</Physician_Code>
  <Physician_GROUP>False</Physician_GROUP>
</PhysicianTable>
<PhysicianTable>

<Physician_ObjectID>503790496</Physician_ObjectID>
  <Physician_Name>New physician, testing
roles</Physician_Name>
  <Physician_Street />
  <Physician_STREET2 />
  <Physician_STREET3 />
  <Physician_City />
  <Physician_State />
  <Physician_Zip />
  <Physician_Code>NP1</Physician_Code>
  <Physician_GROUP>False</Physician_GROUP>
</PhysicianTable>
<PhysicianTable>

<Physician_ObjectID>503612491</Physician_ObjectID>
  <Physician_Name>Sig, Text Only</Physician_Name>
  <Physician_Street />
  <Physician_STREET2 />
  <Physician_STREET3 />
  <Physician_City />
  <Physician_State />
  <Physician_Zip />
  <Physician_Code>Phy01</Physician_Code>
  <Physician_GROUP>False</Physician_GROUP>
</PhysicianTable>
<PhysicianTable>
  <Physician_ObjectID>390657</Physician_ObjectID>
  <Physician_Name>Physician, Test</Physician_Name>
  <Physician_Street />
  <Physician_STREET2 />
  <Physician_STREET3 />
  <Physician_City />

```

```

        <Physician_State />
        <Physician_Zip />
        <Physician_Code>PHY2</Physician_Code>
        <Physician_GROUP>False</Physician_GROUP>
    </PhysicianTable>
    <PhysicianTable>
        <Physician_ObjectID>14661001</Physician_ObjectID>
        <Physician_Name>Physician, Test</Physician_Name>
        <Physician_Street />
        <Physician_STREET2 />
        <Physician_STREET3 />
        <Physician_City />
        <Physician_State />
        <Physician_Zip />
        <Physician_Code>PHY3</Physician_Code>
        <Physician_GROUP>False</Physician_GROUP>
    </PhysicianTable>
    <PhysicianTable>
        <Physician_ObjectID>16932034</Physician_ObjectID>
        <Physician_Name>Physician,
TestAgain</Physician_Name>
        <Physician_Street />
        <Physician_STREET2 />
        <Physician_STREET3 />
        <Physician_City />
        <Physician_State />
        <Physician_Zip />
        <Physician_Code>PHY5</Physician_Code>
        <Physician_GROUP>False</Physician_GROUP>
    </PhysicianTable>
    <PhysicianTable>
        <Physician_ObjectID>135677272</Physician_ObjectID>
        <Physician_Name>Robins, Tim</Physician_Name>
        <Physician_Street />
        <Physician_STREET2 />
        <Physician_STREET3 />
        <Physician_City />
        <Physician_State />
        <Physician_Zip />
        <Physician_Code>ROBINS</Physician_Code>
        <Physician_GROUP>False</Physician_GROUP>
    </PhysicianTable>
    <PhysicianTable>
        <Physician_ObjectID>232073067</Physician_ObjectID>
        <Physician_Name>Robins, Tim Psy</Physician_Name>

```

```

        <Physician_Street />
        <Physician_STREET2 />
        <Physician_STREET3 />
        <Physician_City />
        <Physician_State />
        <Physician_Zip />
        <Physician_Code>ROBINS/PSY</Physician_Code>
        <Physician_GROUP>False</Physician_GROUP>
    </PhysicianTable>
    <PhysicianTable>

<Physician_ObjectID>135676995</Physician_ObjectID>
    <Physician_Name>Smith, Tim</Physician_Name>
    <Physician_Street />
    <Physician_STREET2 />
    <Physician_STREET3 />
    <Physician_City />
    <Physician_State />
    <Physician_Zip />
    <Physician_Code>SMITH</Physician_Code>
    <Physician_GROUP>False</Physician_GROUP>
    </PhysicianTable>
    <PhysicianTable>

<Physician_ObjectID>232073138</Physician_ObjectID>
    <Physician_Name>Smith, Tim Psy</Physician_Name>
    <Physician_Street />
    <Physician_STREET2 />
    <Physician_STREET3 />
    <Physician_City />
    <Physician_State />
    <Physician_Zip />
    <Physician_Code>SMITH/PSY</Physician_Code>
    <Physician_GROUP>False</Physician_GROUP>
    </PhysicianTable>
    <PhysicianTable>

<Physician_ObjectID>119430985</Physician_ObjectID>
    <Physician_Name>Doctor, Thursday</Physician_Name>
    <Physician_Street />
    <Physician_STREET2 />
    <Physician_STREET3 />
    <Physician_City />
    <Physician_State />
    <Physician_Zip />
    <Physician_Code>THURS</Physician_Code>
    <Physician_GROUP>False</Physician_GROUP>

```

```

        </PhysicianTable>
    </NewDataSet>

```

SET STATUS STEP

Endpoint	api/N/ SetStatusStep
Purpose	Set a processing status of a single case, specimen or test order record.
Request Type	GET
Parameters	<pre> object_id={object_id} status_set={status_set} status_advance={status_advance} user_ID={user_ID} return_format={return_format} </pre>

Where:

object_id is the objectId of the record to update.

status_set is the new status setting.

status_advance is a true or false value. If true, the status of the record to update will be advanced to the next defined step. When false, the status_advance argument is ignored.

user_id is the user code for the user in the NucleoLIS user table that is requesting the status change.

return_format is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

Example `https://www.SomeServer.com/api/N/SetStatusStep?object_id=503766463&status_set=&status_advance=true&user_ID=200001&return_format=0`

Returns If the request succeeds, returns the results of processing the status update:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <StatusStepUpdateTable>
    <object_id>503766463</object_id>
    <object_type>PCR</object_type>
    <status_set>Purification</status_set>
    <original_status>Accessioned</original_status>
    <message>success</message>
  </StatusStepUpdateTable>
</NewDataSet>
```

SET STATUS STEPS

Endpoint `api/N/SetStatusSteps`

Purpose Set a processing status of a multiple case, specimen or test order records.

Request Type GET

Parameters

```
object_ids={object_ids}
status_set={status_set}
status_advance={status_advance}
user_ID={user_ID}
return_format={return_format}
```

Where:

object_ids is a pipe delimited list of object ids that specifies which cases should be updated.

status_set is the new processing status used to update the specified cases.

status_advance is a true or false value. If true, the specified cases will be advanced. If false the specified cases will not be advanced.

user_id is the user code for the user in the NucleoLIS user table that is requesting the status change(s).

return_format is the format for returning the data. Where:

0 = XML Format

1 = JSON Format

Example

`https://www.SomeServer.com/api/N/SetStatusSteps?object_ids=504252001|504165951&status_set=Amplify&status_advance=false&user_ID=200001&return_format=0`

Returns

If the request succeeds, returns the results of processing the status updates:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <StatusStepUpdateTable>
    <object_id>504252001</object_id>
    <object_type>PCR</object_type>
    <status_set>Amplify</status_set>
    <original_status>Accessioned</original_status>
    <message>success</message>
  </StatusStepUpdateTable>
  <StatusStepUpdateTable>
    <object_id>504165951</object_id>
    <object_type>PCR</object_type>
    <status_set>Amplify</status_set>
    <original_status>Accessioned</original_status>
    <message>success</message>
```

```

        </StatusStepUpdateTable>
    </NewDataSet>

```

GET PHYSICIAN LOCATIONS

Endpoint `api/N/GetPhysicianLocations`

Purpose Returns all physician location records for the specified physician, based on a unique physician code.

Request Type GET

Parameters

```

physician_code={physician_code}

user_ID={user_ID}

return_format={return_format}

```

Where:

Physician_code is the code of the physician record to return.

user_id is the user code for the user in the NucleoLIS user table that is uploading the text.

return_format is the format for returning the data.
Where:

0 = XML Format

1 = JSON Format

Example

```

https://www.SomeServer.com/api/N/
GetPhysicianLocations?physician_code=00002&user_ID=154
92125&return_format=0

```

Returns If the request succeeds, returns the status of the text that was parsed and processed.:

```

<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <PhysicianLocationTable>
    <Location_ObjectID>2291844</Location_ObjectID>

```

```

        <Location_CreationDate>2012-10-08T09:58:17.373-
04:00</Location_CreationDate>
        <Location_UpdateTime>2016-08-02T08:47:04.33-
04:00</Location_UpdateTime>
        <Location_Location>191919-
GROUP</Location_Location>
        <Location_Phone />
        <Location_Fax />
        <Location_CODE>00002-2291844</Location_CODE>
        <Location_USER1 />
        <Location_USER2 />
        <Location_USER3 />
        <Location_USER4 />
        <Location_EMAIL />
        <Location_Retired>True</Location_Retired>
        <Group_ObjectID>1709111</Group_ObjectID>
        <Group_CreationDate>2007-05-29T16:17:29-
04:00</Group_CreationDate>
        <Group_UpdateTime>2016-11-30T07:37:15.85-
05:00</Group_UpdateTime>
        <Group_Name>Physicians Preferred</Group_Name>
        <Group_Street>127 Plain Street</Group_Street>
        <Group_City>Milford</Group_City>
        <Group_State>MA</Group_State>
        <Group_Zip>01568</Group_Zip>
        <Group_Code>PPM</Group_Code>
        <Group_STREET2>APT 2</Group_STREET2>
        <Group_RETIREDD>False</Group_RETIREDD>
    </PhysicianLocationTable>
    <PhysicianLocationTable>
        <Location_ObjectID>1933150</Location_ObjectID>
        <Location_CreationDate>2010-05-27T15:11:59.63-
04:00</Location_CreationDate>
        <Location_UpdateTime>2016-08-02T08:47:02.77-
04:00</Location_UpdateTime>
        <Location_Location>OFFICE</Location_Location>
        <Location_Phone />
        <Location_CODE>1933150-OFFICE</Location_CODE>
        <Location_Retired>False</Location_Retired>
    </PhysicianLocationTable>
</NewDataSet>

```

GET REPORT

Endpoint `api/N/GetReport`

Purpose Return a single, full case reported result record based on case number and report Id.

Request Type GET

Parameters `case_number={case_number}`
`report_id={report_id}`
`user_ID={user_ID}`
`return_format={return_format}`

Where:

case_number is the case number of the report record to return.

report_id is the report object id of the report record to return.

user_id is the user code for the user in the NucleoLIS user table that is requesting the report.

return_format is the format for returning the data.
Where:

0 = XML Format

1 = JSON Format

Example [https://www.SomeServer.com/api/N/GetReport?](https://www.SomeServer.com/api/N/GetReport?case_number=184&report_id=15386635&user_ID=15492125&return_format=0)
`case_number=184&report_id=15386635&user_ID=15492125&re`
`turn_format=0`

Returns If the request succeeds, returns the full case record inclusive of all data elements associated with patient, case, specimen, test order and test result records:

```

<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>

  <ResultReportTable>

    <Location_AdmitterCity />
    <Location_AdmitterCode />
    <Location_AdmitterName />
    <Location_AdmitterNPI />
    <Location_AdmitterPhone />
    <Location_AdmitterState />
    <Location_AdmitterStreet />
    <Location_AdmitterStreet2 />
    <Location_AdmitterStreet3 />
    <Location_AdmitterZip />
    <Case_Abnormal>False</Case_Abnormal>
    <Case_AmendedComment />
    <Case_CanceledBy />
    <Case_CanceledDate />
    <Case_CaseType>2327730</Case_CaseType>
    <Case_CaseTypeCode>Cytogenetics</Case_CaseTypeCode>

    <Case_CaseTypeDescription>Cytogenetics</Case_CaseTypeD
escription>

    <Case_ChargeString>123 x 2, BCCHARGE x
2</Case_ChargeString>

    <Case_ClinicalInformation />
    <Case_Comment />
    <Case_CorrectedComment />
    <Case_CreationDate>2018-05-07
09:42:53</Case_CreationDate>

    <Case_CREATOR>johnb</Case_CREATOR>

```

```

<Case_ExternalID />
<Case_ExternalPatientID />
<Case_ICDCodesDescriptions />
<Case_ICDDescriptions />
<Case_ICDS />
<Case_Images />
<Case_MODIFIER>0</Case_MODIFIER>
<Case_Number>184</Case_Number>
<Case_ObjectID>13767203</Case_ObjectID>
<Case_OrderDate />
<Case_OrderID />
<Case_OtherPhysician1GreetingDisplay />
<Case_OtherPhysician1Name />
<Case_OtherPhysician2GreetingDisplay />
<Case_OtherPhysician2Name />
<Case_OtherPhysician3GreetingDisplay />
<Case_OtherPhysician3Name />
<Case_OtherPhysician4GreetingDisplay />
<Case_OtherPhysician4Name />
<Case_OtherPhysician5GreetingDisplay />
<Case_OtherPhysician5Name />
<Case_OtherPhysicianString />
<Case_PDFs />
<Case_PendingResultMessage />
<Case_PrescribedMedications />
<Case_PriorReportDate>2019-09-
24</Case_PriorReportDate>

```

<Case_PriorReportTime>11:26</Case_PriorReportTime>
<Case_PriorReportType>FINAL</Case_PriorReportType>
<Case_Questions />
<Case_ReportDate>2019-09-24</Case_ReportDate>
<Case_ReportDefinition />
<Case_ReportedBy>0</Case_ReportedBy>
<Case_ReportedByTitle />
<Case_ReportedFlag>False</Case_ReportedFlag>
<Case_ReportIndex>2</Case_ReportIndex>
<Case_ReportingComment />
<Case_ReportStatus>FINAL</Case_ReportStatus>
<Case_ReportTime>11:25</Case_ReportTime>
<Case_ReportTitleModifier />
<Case_ReportType>PRELIMINARY</Case_ReportType>
<Case_ReqNumber />
<Case_ReviewedBy />
<Case_ReviewedComment />
<Case_ReviewedDate />
<Case_ReviewedTime />
<Case_ReviewFlag>2019-09-24</Case_ReviewFlag>
<Case_Signer />
<Case_SignerDate />
<Case_SignerLocationAddress />
<Case_SignerLocationAddress2 />
<Case_SignerLocationCity />
<Case_SignerLocationDescription />
<Case_SignerLocationName />


```

<Case_SignerLocationObjectID>0</Case_SignerLocationObjectID>

    <Case_SignerLocationState />
    <Case_SignerLocationZip />
    <Case_SignerObjectID>0</Case_SignerObjectID>
    <Case_SignerTime />
    <Case_SignerTitle />
    <Case_Status>I</Case_Status>
    <Case_StatusStep />

<Case_SubmitterLocation>1437553</Case_SubmitterLocation>

    <Case_SuperobjectID>13767202</Case_SuperobjectID>
    <Case_UpdateTime>2019-09-25
05:06:39</Case_UpdateTime>
    <Case_User1 />
    <Case_User2 />
    <Case_User3 />
    <Case_User4 />
    <Constituent_Abnormal>False</Constituent_Abnormal>
    <Constituent_AlwaysReport />

<Constituent_CategorySequence>2</Constituent_CategorySequence>

    <Constituent_Definition>OFBUPRENORPHINE-
C</Constituent_Definition>
    <Constituent_Description>Buprenorphine
</Constituent_Description>
    <Constituent_ImageName1 />
    <Constituent_ImageName2 />

```

```

    <Constituent_IsolateLabel />
    <Constituent_IsolateOrganism />
    <Constituent_IsolateOrganismDescription />
    <Constituent_IsolateOrganismModifier />
    <Constituent_KB />
    <Constituent_MIC />

    <Constituent_ObjectID>15382509</Constituent_ObjectID>

    <Constituent_PriorAbnormal>False</Constituent_PriorAbnormal>

        <Constituent_PriorR1>NEG</Constituent_PriorR1>

    <Constituent_PriorR2>CONSISTENT</Constituent_PriorR2>

        <Constituent_PriorR3>1.5
    ng/ml</Constituent_PriorR3>

        <Constituent_PriorR4>6 to 18
    Hours</Constituent_PriorR4>

        <Constituent_PriorR5 />

        <Constituent_PriorR6 />

        <Constituent_PriorRange />

        <Constituent_PriorRangeComment />

    <Constituent_PriorResult>0.1</Constituent_PriorResult>

        <Constituent_PriorResultDescription />

        <Constituent_PType>Immuno</Constituent_PType>

        <Constituent_R1>NEG</Constituent_R1>

        <Constituent_R1Desc>Detected</Constituent_R1Desc>

        <Constituent_R2>CONSISTENT</Constituent_R2>

```

```

<Constituent_R2Desc>Consistent</Constituent_R2Desc>
    <Constituent_R3>1.5 ng/ml</Constituent_R3>
    <Constituent_R3Desc>Cut Off</Constituent_R3Desc>
    <Constituent_R4>6 to 18 Hours</Constituent_R4>
    <Constituent_R4Desc>Detection
Window</Constituent_R4Desc>
    <Constituent_R5 />
    <Constituent_R5Desc>Result
Text</Constituent_R5Desc>
    <Constituent_R6 />
    <Constituent_R6Desc />
    <Constituent_Range />
    <Constituent_RangeComment />
    <Constituent_RangeCommentDesc />
    <Constituent_RangeDesc />
    <Constituent_ReportCategory />

<Constituent_ReportCategorySequence>283</Constituent_R
eportCategorySequence>
    <Constituent_Result>0.1</Constituent_Result>
    <Constituent_ResultDate />
    <Constituent_ResultDescription />
    <Constituent_ResultTime />
    <Constituent_Sensitivity />
    <Location_ConsultorCity />
    <Location_ConsultorCode />
    <Location_ConsultorName />
    <Location_ConsultorNPI />

```

```

<Location_ConsultorPhone />
<Location_ConsultorState />
<Location_ConsultorStreet />
<Location_ConsultorStreet2 />
<Location_ConsultorStreet3 />
<Location_ConsultorZip />
<Location_GroupCity />
<Location_GroupCode>HARTM</Location_GroupCode>
<Location_GroupCounty />
<Location_GroupFax />
  <Location_GroupName>Last,
First</Location_GroupName>
  <Location_GroupPhone />
  <Location_GroupState />
  <Location_GroupStreet />
  <Location_GroupStreet2 />
  <Location_GroupStreet3 />
  <Location_GroupUser1 />
  <Location_GroupZip />
  <Patient_City />
  <Patient_Code>13767197</Patient_Code>
  <Patient_County />
  <Patient_CreationDate>2018-05-07
12:42:08</Patient_CreationDate>
  <Patient_CREATOR>johnb</Patient_CREATOR>
  <Patient_Dob>1970-01-01</Patient_Dob>
  <Patient_Ethnicity />
  <Patient_Fax />

```

```
<Patient_FirstName>J</Patient_FirstName>
<Patient_LastName>Smith</Patient_LastName>
<Patient_MedicareNumber />
<Patient_MiddleName />
<Patient_MODIFIER>johnb</Patient_MODIFIER>
<Patient_Name>Smith, J</Patient_Name>
<Patient_ObjectID>13767197</Patient_ObjectID>
<Patient_Phone />
<Patient_Race />
<Patient_Sex />
<Patient_Soundex>SMTHJ</Patient_Soundex>
<Patient_Ssn>111-11-1111</Patient_Ssn>
<Patient_State />
<Patient_Street />
<Patient_Street2 />
<Patient_UpdateTime>2018-05-07
09:42:08</Patient_UpdateTime>
<Patient_USER1 />
<Patient_USER2 />
<Patient_USER3 />
<Patient_USER4 />
<Patient_Zip />
<Profile_Abnormal>False</Profile_Abnormal>
<Profile_AlwaysReport />
<Profile_ApproveDate />
<Profile_Approver />

<Profile_ApproverObjectID>0</Profile_ApproverObjectID>
```

```

    <Profile_ApproverTitle />

    <Profile_CompleteDate>2019-09-
24</Profile_CompleteDate>

    <Profile_CompletedBy>Vanger,
John</Profile_CompletedBy>

<Profile_CompleterObjectID>2325137</Profile_CompleterO
bjectID>

    <Profile_CompleteTime />

    <Profile_Definition>OFBASIC</Profile_Definition>

    <Profile_DefinitionUser1 />

    <Profile_DefinitionUser2 />

    <Profile_DefinitionUser3 />

    <Profile_Description>OFBASIC (ORAL
FLUID)</Profile_Description>

    <Profile_Designation />

    <Profile_DNASequenceInterpretation />

    <Profile_ImageName1 />

    <Profile_ImageName2 />

    <Profile_ImageName3 />

    <Profile_ImageName4 />

    <Profile_ImageName5 />

    <Profile_ImageName6 />

    <Profile_Images />

    <Profile_Interpretation />

    <Profile_KaryotypeInterpretation />

    <Profile_Methodology />

    <Profile_ObjectID>15382508</Profile_ObjectID>

    <Profile_OrderDate />

```

```

    <Profile_OrderTime />
    <Profile_PDFs />

    <Profile_PriorAbnormal>False</Profile_PriorAbnormal>
    <Profile_PriorDNASequenceInterpretation />
    <Profile_PriorInterpretation />
    <Profile_PriorR1 />
    <Profile_PriorR2 />
    <Profile_PriorR3 />
    <Profile_PriorR4 />
    <Profile_PriorR5 />
    <Profile_PriorR6 />
    <Profile_PriorRange />
    <Profile_PriorRangeComment />
    <Profile_PriorResult />
    <Profile_PriorResultDescription />
    <Profile_PType>Immuno</Profile_PType>
    <Profile_Questions />
    <Profile_R1 />
    <Profile_R1Desc />
    <Profile_R2 />
    <Profile_R2Desc />
    <Profile_R3 />
    <Profile_R3Desc />
    <Profile_R4 />
    <Profile_R4Desc />
    <Profile_R5 />

```

```

<Profile_R5Desc />
<Profile_R6 />
<Profile_R6Desc />
<Profile_Range />
<Profile_RangeComment />
<Profile_RangeCommentDesc />
<Profile_RangeDesc />
<Profile_ReportCategory />
<Profile_Result />
<Profile_ResultDate />
<Profile_ResultDescription />
<Profile_ResultTime />
<Profile_Signer />
<Profile_SignerDate />
<Profile_SignerLocationAddress />
<Profile_SignerLocationAddress2 />
<Profile_SignerLocationCity />
<Profile_SignerLocationDescription />
<Profile_SignerLocationName />

<Profile_SignerLocationObjectID>0</Profile_SignerLocationObjectID>

<Profile_SignerLocationState />
<Profile_SignerLocationZip />
<Profile_SignerObjectID>0</Profile_SignerObjectID>
<Profile_SignerTitle />
<Profile_StartDate />
<Profile_StartTime />

```



```

    <Profile_TestingLocationAddress />
    <Profile_TestingLocationAddress2 />
    <Profile_TestingLocationCity />
    <Profile_TestingLocationDescription />
    <Profile_TestingLocationName />

    <Profile_TestingLocationObjectID>0</Profile_TestingLoc
ationObjectID>

    <Profile_TestingLocationState />
    <Profile_TestingLocationZip />
    <Specimen_CollectionDate>2018-06-
29</Specimen_CollectionDate>
    <Specimen_CollectionTime />
    <Specimen_Comment />
    <Specimen_Description />
    <Specimen_Designation />
    <Specimen_DiagnosisFull />

    <Specimen_DiagnosisMalignant>False</Specimen_Diagnosis
Malignant>

    <Specimen_FastingStatus />
    <SpecimenFrozenSection_CompleteDate />
    <SpecimenFrozenSection_CompletedBy />
    <SpecimenFrozenSection_CompleteTime />
    <SpecimenFrozenSection_Diagnosis />
    <SpecimenFrozenSection_Gross />
    <SpecimenFrozenSection_ImageName1 />
    <SpecimenFrozenSection_ImageName2 />
    <SpecimenFrozenSection_Images />

```

<SpecimenFrozenSection_ObjectID>0</SpecimenFrozenSection_ObjectID>

<SpecimenFrozenSection_PDFs />

<SpecimenFrozenSection_Signer />

<SpecimenFrozenSection_SignerDate />

<SpecimenFrozenSection_SignerTime />

<SpecimenFrozenSection_SignerTitle />

<SpecimenFrozenSection_Specimen>0</SpecimenFrozenSection_Specimen>

<Specimen_Gross />

<Specimen_ImageName1 />

<Specimen_ImageName2 />

<Specimen_Images />

<Specimen_IndexNumber>1</Specimen_IndexNumber>

<Specimen_Interpretation />

<Specimen_Microscopic />

<Specimen_Number>13770892</Specimen_Number>

<Specimen_ObjectID>13770892</Specimen_ObjectID>

<Specimen_OID>13770892</Specimen_OID>

<Specimen_PatientAge>48</Specimen_PatientAge>

<Specimen_PDFs />

<Specimen_PendingResultMessage />

<Specimen_Priority>ROUTINEx</Specimen_Priority>

<Specimen_Questions />

<Specimen_ReceivedDate>2018-06-29</Specimen_ReceivedDate>

<Specimen_ReceivedTime>05:55:13</Specimen_ReceivedTime>

<Specimen_Site />

<Specimen_SourceCode>APPENDIX</Specimen_SourceCode>

<Specimen_SourceDescription>APPENDIX</Specimen_SourceDescription>

<Specimen_User1 />

<Specimen_User2 />

<Specimen_User3 />

<Specimen_User4 />

<Location_PhysicianCity />

<Location_PhysicianCode>T124</Location_PhysicianCode>

<Location_PhysicianCounty />

<Location_PhysicianDegree />

<Location_PhysicianFax />

<Location_PhysicianGreetingDisplay>Earl Gentry
M.D.</Location_PhysicianGreetingDisplay>

<Location_PhysicianLocationCity />

<Location_PhysicianLocationLocation>None</Location_PhysicianLocationLocation>

<Location_PhysicianLocationPhone />

<Location_PhysicianLocationState />

<Location_PhysicianLocationStreet />

<Location_PhysicianLocationStreet2 />

<Location_PhysicianLocationStreet3 />

<Location_PhysicianLocationZip />

```

    <Location_PhysicianName>Gentry M.D.,
Earl</Location_PhysicianName>

    <Location_PhysicianNPI />

    <Location_PhysicianPhone />

    <Location_PhysicianState />

    <Location_PhysicianStreet />

    <Location_PhysicianStreet2 />

    <Location_PhysicianStreet3 />

    <Location_PhysicianUPIN />

    <Location_PhysicianUser1 />

    <Location_PhysicianZip />

    <Visit_ADMITTERFAX />

    <Visit_Code>13767202</Visit_Code>

    <Visit_Comment />

    <Visit_CONSULTORFAX />

    <Visit_CreationDate>2018-05-07
12:42:52</Visit_CreationDate>

    <Visit_DischargeDate />

    <Visit_Location />

    <Visit_MODIFIER>johnb</Visit_MODIFIER>

    <Visit_ObjectID>13767202</Visit_ObjectID>

<Visit_SuperObjectID>13767197</Visit_SuperObjectID>

    <Visit_Type />

    <Visit_UpdateTime>2018-05-07
09:42:52</Visit_UpdateTime>

    <Visit_USER1 />

    <Visit_USER2 />

    <Visit_USER3 />

```

```

    <Visit_USER4 />

    <Visit_VisitDate>2018-05-07</Visit_VisitDate>

  </ResultReportTable>

</NewDataSet>

```

GET REPORTS

Endpoint	api/N/GetReports
Purpose	Returns a collection of report records based on filter parameters passed.
Request Type	GET
Parameters	<pre> case_number={case_number} actives_only = {actives_only} report_date_from = {report_date_from} report_date_to = {report_date_to} submitting_physicians = {submitting_physicians} submitting_groups = {submitting_groups} filter_expression = {filter_expression} user_ID={user_ID} return_format={return_format} </pre>

Where:

case_number is a record identifier which has the specified value in the case number. Wildcards can be embedded in the argument as "*".

actives_only if set to true, only returns active cases. If set to false, returns all cases.

report_date_from is the oldest inclusive report date in the format: (MMDDYYYY OR YYYY-MM-DD).

report_date_to is the newest inclusive report date in the format: (MMDDYYYY OR YYYY-MM-DD).

submitting_physicians is a pipe delimited string specifying the physician identifying codes.

submitting_groups is a pipe delimited string specifying the physician facility identifying codes.

filter_expression is a dynamic filter expression which may pertain to any data element returned in the dataset.

user_id is the user code for the user in the NucleoLIS user table that is requesting the reports.

return_format is the format for returning the data. Where:

0 = XML Format

1 = JSON Format

Example

https://www.SomeServer.com/api/N/GetReports?case_number=18*&activesOnly=false&report_date_from=2000-01-01&report_date_to=2021-02-21&submitting_physicians=&submitting_groups=&filter_expression=&user_ID=15492125&return_format=0

Returns

If the request succeeds, returns a collection of Patient, Case and Report level data based on filtering parameters passed including Patient_ObjectID, Patient_Name, Patient_Code, Patient_Dob, Patient_Sex, Case_Number, Case_OrderDate, Case_CreationDate, Case_ReportStatus, Case_Status, Case_ObjectID, Report_ObjectID, Report_Type, Report_Date, Report_Time and Report_TitleModifier:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <ResultReportTable>
    <Patient_ObjectID>13767197</Patient_ObjectID>
    <Patient_Name>Smith, J</Patient_Name>
    <Patient_Code>13767197</Patient_Code>
    <Patient_Dob>1970-01-01</Patient_Dob>
```

```

    <Patient_Sex />
    <Case_Number>183</Case_Number>
    <Case_OrderDate />
    <Case_CreationDate>2018-05-07T09:42:40.91-
04:00</Case_CreationDate>
    <Case_ReportStatus>FINAL</Case_ReportStatus>
    <Case_Status>I</Case_Status>
    <Case_ObjectID>13767201</Case_ObjectID>
    <Report_ObjectID>13781492</Report_ObjectID>
    <Report_Type>FINAL</Report_Type>
    <Report_Date>2018-09-24</Report_Date>
    <Report_Time>05:57</Report_Time>
</ResultReportTable>
<ResultReportTable>
    <Patient_ObjectID>13767197</Patient_ObjectID>
    <Patient_Name>Smith, J</Patient_Name>
    <Patient_Code>13767197</Patient_Code>
    <Patient_Dob>1970-01-01</Patient_Dob>
    <Patient_Sex />
    <Case_Number>184</Case_Number>
    <Case_OrderDate />
    <Case_CreationDate>2018-05-07T09:42:53.527-
04:00</Case_CreationDate>
    <Case_ReportStatus>FINAL</Case_ReportStatus>
    <Case_Status>I</Case_Status>
    <Case_ObjectID>13767203</Case_ObjectID>
    <Report_ObjectID>15386635</Report_ObjectID>
    <Report_Type>PRELIMINARY</Report_Type>

```

```

    <Report_Date>2019-09-24</Report_Date>

    <Report_Time>11:25</Report_Time>

  </ResultReportTable>

</NewDataSet>

```

GET RESULT REPORT

Endpoint	api/N/GetResultReport
Purpose	Return a PDF result report file based on case number and report Id.
Request Type	GET
Parameters	case_number={case_number} report_id={report_id} user_ID={user_ID} return_format={return_format}

Where:

case_number is the case number of the report record to return.

report_id is the report object id of the report record to return.

user_id is the user code for the user in the NucleoLIS user table that is requesting the result report.

return_format is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

Example [https://www.SomeServer.com/api/N/GetResultReport?](https://www.SomeServer.com/api/N/GetResultReport?case_number=FSH2020-)
[case_number=FSH2020-](https://www.SomeServer.com/api/N/GetResultReport?case_number=FSH2020-)

000002&report_id=15498801&user_ID=15492125&return_format=0

Returns

If the request succeeds, returns a base64 encoded copy of the requested PDF result report:

```
<?xml version="1.0" encoding="UTF-8" ?><ReportRecords>

  <Reports>

    <report_id>15498801</report_id>

    <case_number>FSH2020-000002</case_number>

    <base_64>JVBERi0xLjQNCiWxsrO0DQolQ3Jl....o=</base_64>

    <message />

  </Reports>

</ReportRecords>
```

GET CASE TYPE DEFINITIONS

Endpoint	api/N/GetCaseTypeDefinitions
Purpose	Returns a collection of case type definitions in the system build based on filter parameters passed.
Request Type	GET
Parameters	user_ID={user_ID} return_format={return_format} with_retired = {with_retired}

Where:

user_id is the user code for the user in the NucleoLIS user table that is requesting the case type definitions.

return_format is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

with_retired If false, will exclude retired entries in the system build. If true, will return all entries in the system build.

Example

https://www.SomeServer.com/api/N/GetCaseTypeDefinitions?user_ID=15492125&return_format=0&with_retired=false

Returns

If the request succeeds, returns the system unique identifier, code, description and retired state for each item defined in the system build:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>

  <CaseTypeDefinitionTable>

    <CaseTypeDefinition_ObjectID>2327075</CaseTypeDefinition_ObjectID>

    <CaseTypeDefinition_Code>GEN</CaseTypeDefinition_Code>

    <CaseTypeDefinition_Description>GENERAL</CaseTypeDefinition_Description>

    <CaseTypeDefinition_Retired>False</CaseTypeDefinition_Retired>

  </CaseTypeDefinitionTable>

</CaseTypeDefinitionTable>
```

<CaseTypeDefinition_ObjectID>2327641</CaseTypeDefinition_ObjectID>

<CaseTypeDefinition_Code>Prenatal</CaseTypeDefinition_Code>

<CaseTypeDefinition_Description>Prenatal Testing</CaseTypeDefinition_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_Retired>

</CaseTypeDefinitionTable>

<CaseTypeDefinitionTable>

<CaseTypeDefinition_ObjectID>2327727</CaseTypeDefinition_ObjectID>

<CaseTypeDefinition_Code>Enzyme</CaseTypeDefinition_Code>

<CaseTypeDefinition_Description>Enzyme</CaseTypeDefinition_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_Retired>

</CaseTypeDefinitionTable>

<CaseTypeDefinitionTable>

<CaseTypeDefinition_ObjectID>2327728</CaseTypeDefinition_ObjectID>

<CaseTypeDefinition_Code>DNA Seq</CaseTypeDefinition_Code>

<CaseTypeDefinition_Description>DNA Sequence</CaseTypeDefinition_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_Retired>

</CaseTypeDefinitionTable>

<CaseTypeDefinitionTable>

<CaseTypeDefinition_ObjectID>2327729</CaseTypeDefinition_ObjectID>

<CaseTypeDefinition_Code>Disorder</CaseTypeDefinition_Code>

<CaseTypeDefinition_Description>Genetic Disorders</CaseTypeDefinition_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_Retired>

</CaseTypeDefinitionTable>

<CaseTypeDefinitionTable>

<CaseTypeDefinition_ObjectID>2327730</CaseTypeDefinition_ObjectID>

<CaseTypeDefinition_Code>Cytogenetics</CaseTypeDefinition_Code>

<CaseTypeDefinition_Description>Cytogenetics</CaseTypeDefinition_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_Retired>

</CaseTypeDefinitionTable>

<CaseTypeDefinitionTable>

<CaseTypeDefinition_ObjectID>2327732</CaseTypeDefinition_ObjectID>

<CaseTypeDefinition_Code>Cytogenetics2</CaseTypeDefinition_Code>

<CaseTypeDefinition_Description>FISH</CaseTypeDefinition_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_Retired>

</CaseTypeDefinitionTable>

<CaseTypeDefinitionTable>

<CaseTypeDefinition_ObjectID>9758103</CaseTypeDefinition_ObjectID>

<CaseTypeDefinition_Code>Cytogenetics_1</CaseTypeDefinition_Code>

<CaseTypeDefinition_Description>Karyotype_1</CaseTypeDefinition_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_Retired>

</CaseTypeDefinitionTable>

<CaseTypeDefinitionTable>

<CaseTypeDefinition_ObjectID>10866436</CaseTypeDefinition_ObjectID>

<CaseTypeDefinition_Code>pGX</CaseTypeDefinition_Code>

<CaseTypeDefinition_Description>TSI
Integration</CaseTypeDefinition_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_Retired>

```

</CaseTypeDefinitionTable>

<CaseTypeDefinitionTable>

<CaseTypeDefinition_ObjectID>11363917</CaseTypeDefinition_ObjectID>

<CaseTypeDefinition_Code>SURG</CaseTypeDefinition_Code>

<CaseTypeDefinition_Description>SURGICAL</CaseTypeDefinition_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_Retired>

</CaseTypeDefinitionTable>

<CaseTypeDefinitionTable>

<CaseTypeDefinition_ObjectID>13792557</CaseTypeDefinition_ObjectID>

<CaseTypeDefinition_Code>Gastro</CaseTypeDefinition_Code>

<CaseTypeDefinition_Description>Gastro</CaseTypeDefinition_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_Retired>

</CaseTypeDefinitionTable>

</NewDataSet>

```

GET SOURCE DEFINITIONS

Endpoint `api/N/GetSourceDefinitions`

Purpose	Returns a collection of source definitions in the system build based on filter parameters passed.
Request Type	GET
Parameters	<p>user_ID={user_ID}</p> <p>return_format={return_format}</p> <p>with_retired = {with_retired}</p> <p>Where:</p> <p>user_id is the user code for the user in the NucleoLIS user table that is requesting the source definitions.</p> <p>return_format is the format for returning the data.</p> <p>Where:</p> <p>0 = XML Format</p> <p>1 = JSON Format</p> <p>with_retired If false, will exclude retired entries in the system build. If true, will return all entries in the system build.</p>
Example	<pre>https://www.SomeServer.com/api/N/GetSourceDefinitions? user_ID=15492125&return_format=0&with_retired=false</pre>
Returns	<p>If the request succeeds, returns the system unique identifier, code, description and retired state for each item defined in the system build:</p> <pre><?xml version="1.0" encoding="UTF-8" ?><NewDataSet> <SourceDefinitionTable> <SourceDefinition_ObjectID>2326047</SourceDefinition_ObjectID></pre>

<SourceDefinition_Code>APPENDIX</SourceDefinition_Code>

<SourceDefinition_Description>APPENDIX</SourceDefinition_Description>

<SourceDefinition_Retired>False</SourceDefinition_Retired>

</SourceDefinitionTable>

<SourceDefinitionTable>

<SourceDefinition_ObjectID>2326116</SourceDefinition_ObjectID>

<SourceDefinition_Code>BLOOD</SourceDefinition_Code>

<SourceDefinition_Description>BLOOD</SourceDefinition_Description>

<SourceDefinition_Retired>False</SourceDefinition_Retired>

</SourceDefinitionTable>

<SourceDefinitionTable>

<SourceDefinition_ObjectID>2326606</SourceDefinition_ObjectID>

<SourceDefinition_Code>DBlood</SourceDefinition_Code>

<SourceDefinition_Description>Dried Blood
</SourceDefinition_Description>

<SourceDefinition_Retired>False</SourceDefinition_Retired>

</SourceDefinitionTable>


```

<SourceDefinitionTable>

<SourceDefinition_ObjectID>2326607</SourceDefinition_ObjectID>

<SourceDefinition_Code>Serum</SourceDefinition_Code>

<SourceDefinition_Description>Serum</SourceDefinition_Description>

<SourceDefinition_Retired>False</SourceDefinition_Retired>

</SourceDefinitionTable>

<SourceDefinitionTable>

<SourceDefinition_ObjectID>2326608</SourceDefinition_ObjectID>

<SourceDefinition_Code>Urine</SourceDefinition_Code>

<SourceDefinition_Description>Urine</SourceDefinition_Description>

<SourceDefinition_Retired>False</SourceDefinition_Retired>

</SourceDefinitionTable>

<SourceDefinitionTable>

<SourceDefinition_ObjectID>2326609</SourceDefinition_ObjectID>

<SourceDefinition_Code>Plasma</SourceDefinition_Code>

<SourceDefinition_Description>Plasma</SourceDefinition_Description>

```

```

<SourceDefinition_Retired>False</SourceDefinition_Retired>

  </SourceDefinitionTable>

  <SourceDefinitionTable>

    <SourceDefinition_ObjectID>2326610</SourceDefinition_ObjectID>

    <SourceDefinition_Code>Leukocytes</SourceDefinition_Code>

    <SourceDefinition_Description>Leukocytes</SourceDefinition_Description>

    <SourceDefinition_Retired>False</SourceDefinition_Retired>

    </SourceDefinitionTable>

    <SourceDefinitionTable>

      <SourceDefinition_ObjectID>2326611</SourceDefinition_ObjectID>

      <SourceDefinition_Code>Fibroblasts</SourceDefinition_Code>

      <SourceDefinition_Description>Fibroblasts</SourceDefinition_Description>

      <SourceDefinition_Retired>False</SourceDefinition_Retired>

      </SourceDefinitionTable>

    </NewDataSet>

```

GET PROFILE DEFINITIONS

Endpoint	api/N/GetProfileDefinitions
Purpose	Returns a collection of profile definitions in the system build based on filter parameters passed.
Request Type	GET
Parameters	<code>user_ID={user_ID}</code> <code>return_format={return_format}</code> <code>with_retired = {with_retired}</code>

Where:

user_id is the user code for the user in the NucleoLIS user table that is requesting the profile definitions.

return_format is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

with_retired If false, will exclude retired entries in the system build. If true, will return all entries in the system build.

Example	<code>https://www.SomeServer.com/api/N/GetProfileDefinitions?user_ID=15492125&return_format=0&with_retired=false</code>
----------------	---

Returns	If the request succeeds, returns the system unique identifier, code, description and retired state for each item defined in the system build.
----------------	---

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <ProfileDefinitionTable>
```

<ProfileDefinition_ObjectID>2325220</ProfileDefinition_ObjectID>

<ProfileDefinition_Code>BAESO</ProfileDefinition_Code>

<ProfileDefinition_Description>Barrett's Esophagus
Progression: check for Esophageal
Adenocarcinoma</ProfileDefinition_Description>

<ProfileDefinition_Retired>False</ProfileDefinition_Retired>

<ProfileDefinition_LIBRARYTYPE>FISH</ProfileDefinition_LIBRARYTYPE>

</ProfileDefinitionTable>

<ProfileDefinitionTable>

<ProfileDefinition_ObjectID>2335695</ProfileDefinition_ObjectID>

<ProfileDefinition_Code>BCL-
6</ProfileDefinition_Code>

<ProfileDefinition_Description>Break Apart
Probe</ProfileDefinition_Description>

<ProfileDefinition_Retired>False</ProfileDefinition_Retired>

<ProfileDefinition_LIBRARYTYPE>FISH</ProfileDefinition_LIBRARYTYPE>

</ProfileDefinitionTable>

<ProfileDefinitionTable>

<ProfileDefinition_ObjectID>2326451</ProfileDefinition_ObjectID>

```

<ProfileDefinition_Code>OVCA</ProfileDefinition_Code>

    <ProfileDefinition_Description>Ovarian Cancer
    Kit</ProfileDefinition_Description>

<ProfileDefinition_Retired>False</ProfileDefinition_Re
tired>

<ProfileDefinition_LIBRARYTYPE>PCR</ProfileDefinition_
LIBRARYTYPE>

    </ProfileDefinitionTable>

    <ProfileDefinitionTable>

<ProfileDefinition_ObjectID>2326532</ProfileDefinition
_ObjectID>

<ProfileDefinition_Code>RESPVIRUS</ProfileDefinition_C
ode>

    <ProfileDefinition_Description>Respiratory Virus
    Panel</ProfileDefinition_Description>

<ProfileDefinition_Retired>False</ProfileDefinition_Re
tired>

<ProfileDefinition_LIBRARYTYPE>PCR</ProfileDefinition_
LIBRARYTYPE>

    </ProfileDefinitionTable>

    <ProfileDefinitionTable>

<ProfileDefinition_ObjectID>2326565</ProfileDefinition
_ObjectID>

<ProfileDefinition_Code>BL</ProfileDefinition_Code>

```

```

        <ProfileDefinition_Description>Beta-Lactamase
Detection Comp Gram Neg
Panel</ProfileDefinition_Description>

<ProfileDefinition_Retired>False</ProfileDefinition_Re
tired>

<ProfileDefinition_LIBRARYTYPE>PCR</ProfileDefinition_
LIBRARYTYPE>

    </ProfileDefinitionTable>

    <ProfileDefinitionTable>

<ProfileDefinition_ObjectID>2326570</ProfileDefinition
_ObjectID>

<ProfileDefinition_Code>HPV</ProfileDefinition_Code>
    <ProfileDefinition_Description>HPV Genotyping-
TMA</ProfileDefinition_Description>

<ProfileDefinition_Retired>False</ProfileDefinition_Re
tired>

<ProfileDefinition_LIBRARYTYPE>PCR</ProfileDefinition_
LIBRARYTYPE>

    </ProfileDefinitionTable>

</NewDataSet>

```

GET QUESTION DEFINITIONS

Endpoint	api/N/GetQuestionDefinitions
Purpose	Returns a collection of profile definitions in the system build based on filter parameters passed.
Request Type	GET

Parameters

```
user_ID={user_ID}

return_format={return_format}

with_retired={with_retired}
```

Where:

user_id is the user code for the user in the NucleoLIS user table that is requesting the question definitions.

return_format is the format for returning the data.
Where:

0 = XML Format

1 = JSON Format

with_retired If false, will exclude retired entries in the system build. If true, will return all entries in the system build.

Example

```
https://www.SomeServer.com/api/N/GetQuestionDefinitions?user_ID=15492125&return_format=0&with_retired=false
```

Returns

If the request succeeds, returns the system unique identifier, code, description and retired state for each item defined in the system build for each question and response defined:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
```

```
<QuestionDefinitionTable>
```

```
<AAOEQuestionDefinition_ObjectID>1722188</AAOEQuestionDefinition_ObjectID>
```

```
<AAOEQuestionDefinition_Code>CATCH</AAOEQuestionDefinition_Code>
```

```
<AAOEQuestionDefinition_Description>Enter catch method</AAOEQuestionDefinition_Description>
```

<AAOEQuestionDefinition_RETIREDD>False</AAOEQuestionDef
inition_RETIREDD>

<AAOEQuestionDefinition_AllowOther>False</AAOEQuestion
Definition_AllowOther>

<AAOEQuestionDefinition_SingleSelect>False</AAOEQuesti
onDefinition_SingleSelect>

<AAOEQuestionDefinition_Required>False</AAOEQuestionDe
finition_Required>

<AAOEResponseDefinition_Code>CLEAN</AAOEResponseDefini
tion_Code>

<AAOEResponseDefinition_Description>Clean
Catch</AAOEResponseDefinition_Description>

<AAOEResponseDefinition_RETIREDD>False</AAOEResponseDef
inition_RETIREDD>

</QuestionDefinitionTable>

<QuestionDefinitionTable>

<AAOEQuestionDefinition_ObjectID>1722188</AAOEQuestion
Definition_ObjectID>

<AAOEQuestionDefinition_Code>CATCH</AAOEQuestionDefini
tion_Code>

<AAOEQuestionDefinition_Description>Enter catch
method</AAOEQuestionDefinition_Description>

<AAOEQuestionDefinition_RETIREDD>False</AAOEQuestionDef
inition_RETIREDD>

<AAOEQuestionDefinition_AllowOther>False</AAOEQuestion
Definition_AllowOther>

<AAOEQuestionDefinition_SingleSelect>False</AAOEQuestionDefinition_SingleSelect>

<AAOEQuestionDefinition_Required>False</AAOEQuestionDefinition_Required>

<AAOEResponseDefinition_Code>RANDOM</AAOEResponseDefinition_Code>

<AAOEResponseDefinition_Description>Random</AAOEResponseDefinition_Description>

<AAOEResponseDefinition_RETIREDD>False</AAOEResponseDefinition_RETIREDD>

</QuestionDefinitionTable>

<QuestionDefinitionTable>

<AAOEQuestionDefinition_ObjectID>1722191</AAOEQuestionDefinition_ObjectID>

<AAOEQuestionDefinition_Code>TEXT</AAOEQuestionDefinition_Code>

<AAOEQuestionDefinition_Description>Enter Text</AAOEQuestionDefinition_Description>

<AAOEQuestionDefinition_RETIREDD>False</AAOEQuestionDefinition_RETIREDD>

</QuestionDefinitionTable>

<QuestionDefinitionTable>

<AAOEQuestionDefinition_ObjectID>1722239</AAOEQuestionDefinition_ObjectID>

<AAOEQuestionDefinition_Code>HOW</AAOEQuestionDefinition_Code>

```

        <AAOEQuestionDefinition_Description>How are you
today?</AAOEQuestionDefinition_Description>

<AAOEQuestionDefinition_RETIREDD>False</AAOEQuestionDef
inition_RETIREDD>

<AAOEQuestionDefinition_AllowOther>False</AAOEQuestion
Definition_AllowOther>

<AAOEQuestionDefinition_SingleSelect>True</AAOEQuestio
nDefinition_SingleSelect>

<AAOEQuestionDefinition_Required>True</AAOEQuestionDef
inition_Required>

<AAOEQuestionDefinition_DefaultResponse>Good</AAOEQues
tionDefinition_DefaultResponse>

<AAOEResponseDefinition_Code>GOOD</AAOEResponseDefinit
ion_Code>

<AAOEResponseDefinition_Description>Good</AAOEResponse
Definition_Description>

<AAOEResponseDefinition_RETIREDD>False</AAOEResponseDef
inition_RETIREDD>

    </QuestionDefinitionTable>

</NewDataSet>

```

SET IMAGE DOCUMENT

Endpoint	api/N/SetImageDocument
Purpose	Adds an image to case, specimen, or profile
Request Type	POST

Parameters

```
object_id={object_id}
id_type={case,specimen,profile}
file_format={ pdf,bmp,jpg,gif,png,tif}
content={POSTCONTENT:location and path of image file}
reportable={true,false}
friendly_name={name_value}
notes={notes_value}
user_ID={user_ID}
return_format={return_format}
```

Where:

object_ids is a (mandatory) String, containing a unique identifier of case, specimen or profile record.

id_type is a (mandatory) String - case, specimen or profile.

file_format is a specific image file type being uploaded - pdf,bmp,jpg,gif,png,tif

content is the full path and filename of the image to attach.

reportable is to set the image as reportable or not. (true,false).

friendly_name is the name of the file displayed in NucleoLIS.

notes is the notes displayed in NucleoLIS.

user_id is a (mandatory) String, the authorized user

return_format is the desired return format (mandatory)
- Integer, 0 = XML, 1 = JSON

Example

```
https://www.someserver.com/api/N/SetImageDocument?object_id=15631506&id_type=Case&file_format=jpg&content=PO
```

```
ST_CONTENT:C:\Users\johnb\Desktop\Capture.JPG&content_
bytes=&reportable=false&friendly_name=my name&notes=my
notes&user_id=561544911&return_format=0
```

Returns

If the request succeeds, returns the success status:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>

  <ImageUploadTable>

    <object_type>ImageUploadTable</object_type>

    <message>Image_Attached</message>

    <status>Success</status>

  </ImageUploadTable>

</NewDataSet>  </ImageUploadTable>

</NewDataSet>
```

SET INSTRUMENT TEST RESULTS

Endpoint	api/N/SetInstrumentTestResults
Purpose	Sets result values for test orders.
Request Type	GET
Parameters	<pre>object_ids={object_ids} at_profile={at_profile} result_value={result_value} r1_value={r1_value} r2_value={r2_value} r3_value={r3_value} r4_value={r4_value}</pre>

```
r5_value={r5_value}

r6_value={r6_value}range_value is a (mandatory)
String, the a range value to set


rangeComment_value={rangeComment_value}

abnormal_value={abnormal_value}

sent_from={sent_from}

batch_id={batch_id}

user_ID={user_ID}

return_format={return_format}
```

Where:

object_ids is a (mandatory) String, containing a pipe delimited list of unique identifiers of profile records.

test_code is a (mandatory) String, containing profile or constituent code to result.

at_profile is a (mandatory) Boolean, true for profile level result, false for constituent result.

result_value is a (mandatory) String, a result value to set.

r1_value is a (mandatory) String, the r1 result value to set.

r2_value is a (mandatory) String, the r2 result value to set.

r3_value is a (mandatory) String, the r3 result value to set.

r4_value is a (mandatory) String, the r4 result value to set.

r5_value is a (mandatory) String, the r5 result value to set.

r6_value is a (mandatory) String, the r6 result value to set.

range_value is a (mandatory) String, the a range value to set.

rangeComment_value is a (mandatory) String, the result comment value to set.

abnormal_value is a (mandatory) Boolean, true or false signifying if this value is abnormal.

sent_from is (mandatory), the instrument or system in which the data originated

batch_id is a (mandatory) String, the id of the batch associated with the results (if in use)

user_id is a (mandatory) String, the authorized user

return_format is the desired return format (mandatory)
- Integer, 0 = XML, 1 = JSON

Example

```
https://www.SomeServer.com/api/N/SetInstrumentTestResults?object_ids=15403525|15494579|15494594&test_code=Escherichia coli  
0157&at_profile=false&result_value=Negative&r1_value=&r2_value=&r3_value=&r4_value=&r5_value=&r6_value=&range_value=&rangeComment_value=&abnormal_value=False&sent_from=Luminex&batch_id=Q123456&user_id=15492125&return_format=0
```

Returns

If the request succeeds, returns the results of processing the result updates:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>  
  <InstrumentResultUpdate>  
    <object_id>15403525</object_id>  
    <object_type>PCR</object_type>  
    <message>success</message>  
  </InstrumentResultUpdate>  
  <InstrumentResultUpdate>
```

```

        <object_id>15494579</object_id>

        <object_type>PCR</object_type>

        <message>success</message>

    </InstrumentResultUpdate>

    <InstrumentResultUpdate>

        <object_id>15494594</object_id>

        <object_type>PCR</object_type>

        <message>success</message>

    </InstrumentResultUpdate>

</NewDataSet>

```

SET INSTRUMENT BATCH TEST RESULTS

Endpoint	api/N/SetInstrumentBatchTestResults
Purpose	Sets result values for instrument orders.
Request Type	POST
Parameters (payload)	order_type=(order_type) sent_from=(sent_from) at_profile=(at_profile) user_ID=(user_ID) return_format=(return_format) file_text=object_id,test_code,result_value,r1_value,r2_value,r3_value,r4_value,r5_value,r6_value,range_value,rangeComment_value,abnormal_value,batch_id object_id,test_code,result_value,r1_value,r2_value,r3_value,r4_value,r5_value,r6_value,range_value,rangeComment_value,abnormal_value,batch_id object_id,test_code,result_value,r1_value,r2_value,r3_value,r4_value,r5_value,r6_value,range_value,rangeComment_value,abnormal_value,batch_id

object_id,test_code,result_value,r1_value,r2_value,r3_value,r4_value,r5_value,r6_value,range_value,rangeComment_value,abnormal_value,batch_id

Where:

object_ids is a (mandatory) String, containing a pipe delimited list of unique identifiers of profile records.

test_code is a (mandatory) String, containing profile or constituent code to result.

at_profile is a (mandatory) Boolean, true for profile level result, false for constituent result.

result_value is a (mandatory) String, a result value to set.

r1_value is a (mandatory) String, the r1 result value to set.

r2_value is a (mandatory) String, the r2 result value to set.

r3_value is a (mandatory) String, the r3 result value to set.

r4_value is a (mandatory) String, the r4 result value to set.

r5_value is a (mandatory) String, the r5 result value to set.

r6_value is a (mandatory) String, the r6 result value to set.

range_value is a (mandatory) String, the a range value to set.

rangeComment_value is a (mandatory) String, the result comment value to set.

abnormal_value is a (mandatory) Boolean, true or false signifying if this value is abnormal.

sent_from is (mandatory), the instrument or system in which the data originated

batch_id is a (mandatory) String, the id of the batch associated with the results (if in use)

user_id is a (mandatory) String, the authorized user
return_format is the desired return format (mandatory)
- Integer, 0 = XML, 1 = JSON

Example

```
https://www.SomeServer.com/api/N/SetInstrumentBatchTestResults?order_type=PCR
```

```
sent_from=Luxinda
```

```
at_profile=false
```

```
user_ID=3088751
```

```
return_format=0
```

```
file_text=object_id,test_code,result_value,r1_value,r2_value,r3_value,r4_value,r5_value,r6_value,range_value,rangeComment_value,abnormal_value,batch_id
```

```
object_id,test_code,result_value,r1_value,r2_value,r3_value,r4_value,r5_value,r6_value,range_value,rangeComment_value,abnormal_value,batch_id
```

```
object_id,test_code,result_value,r1_value,r2_value,r3_value,r4_value,r5_value,r6_value,range_value,rangeComment_value,abnormal_value,batch_id
```

```
object_id,test_code,result_value,r1_value,r2_value,r3_value,r4_value,r5_value,r6_value,range_value,rangeComment_value,abnormal_value,batch_id
```

Returns

If the request succeeds, returns the results of processing the result updates:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>  
  <InstrumentResultUpdate>  
    <object_id>15403525</object_id>  
    <object_type>PCR</object_type>  
    <message>success</message>  
  </InstrumentResultUpdate>  
  <InstrumentResultUpdate>
```

```

        <object_id>15494579</object_id>

        <object_type>PCR</object_type>

        <message>success</message>

    </InstrumentResultUpdate>

    <InstrumentResultUpdate>

        <object_id>15494594</object_id>

        <object_type>PCR</object_type>

        <message>success</message>

    </InstrumentResultUpdate>

</NewDataSet>

```

SET UPLOADED TEXT

Endpoint	api/N/SetUploadedText
Purpose	Sends inbound text and test processor identification for customized parsing and processing.
Request Type	GET
Parameters	processor_id={processor_id} file_text={file_text} user_ID={user_ID} return_format={return_format}

Where:

processor_id is the (mandatory) the id of the text processor to use.

file_text is the (mandatory) text to parse and process.

user_id is the user code for the user in the NucleoLIS user table that is uploading the text.

return_format is the format for returning the data.
Where:

0 = XML Format

1 = JSON Format

Example

`https://www.SomeServer.com/api/N/SetUploadedText?processor_id=MyProcessor&file_text=MSH|&user_ID=15492125&return_format=0`

Returns

If the request succeeds, returns the status of the text that was parsed and processed.:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <CustomTable>
    <object_type>CustomTable</object_type>
    <message>File Processed</message>
    <status>Success</status>
  </PatientTable>
</NewDataSet>
```

ADD PATIENT

Endpoint `api/N/Add Patient`

Purpose Adds a patient record

Request Type GET

Parameters

```
patient_lastname={patient_lastname}
patient_firstname={patient_firstname}
patient_middlename={patient_middlename}
patient_street={patient_street}
```

```
patient_street2={patient_street2}
patient_city={patient_city}
patient_state={patient_state}
patient_zip={patient_zip}
patient_sex={patient_sex}
patient_dob={patient_dob}
patient_ssn={patient_ssn}
additional_fields={field1^value1|field2^value2}
user_ID={user_ID}
return_format={return_format}
```

Where:

patient_lastname is the last name of the patient.

patient_firstname is the first name of the patient.

patient_middleName is the middle name of the patient.

patient_street is the street of the patient.

patient_street2 is the street of the patient.

patient_city is the city of the patient.

patient_state is the state of the patient.

patient_zip is the zip code of the patient.

patient_sex is the sex of the patient.

patient_dob is the date of birth of the patient in the
format: yyyy-mm-dd.

patient_ssn is the Social Security Number of the
patient.

additional_fields is Pipe delimited fields to update
when in use. FieldName^FieldValue|

user_id is the user code for the user in the NucleoLIS user table that is adding the patient.

return_format is the format for returning the data.
Where:

0 = XML Format

1 = JSON Format

Example

```
https://www.SomeServer.com/api/N/
Addpatient?patient_lastname=Smith&patient_firstname=Fr
ancisco&patient_middleName=John&patient_street=126
Maple Street&patient_street2=Apt
9&patient_city=Milford&patient_state=MA&patient_zip=01
757&patient_Sex=M&patient_dob=1970-01-
01&patient_ssn=123-456-
7890&additional_fields=&user_ID=15492125&return_format
=0
```

Returns

If the request succeeds, returns the status of the patient add request:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <PatientTable>
    <object_type>PatientTable</object_type>
    <message>Patient Added</message>
    <status>Success</status>
  </PatientTable>
</NewDataSet>
```

UPDATE PATIENT

Endpoint	api/N/UpdatePatient
Purpose	Updates an existing patient record.
Request Type	GET
Parameters	patient_id={patient_id} patient_lastname={patient_lastname}

```
patient_firstname={patient_firstname}
patient_middlename={patient_middlename}
patient_street={patient_street}
patient_street2={patient_street2}
patient_city={patient_city}
patient_state={patient_state}
patient_zip={patient_zip}
patient_sex={patient_sex}
patient_dob={patient_dob}
patient_ssn={patient_ssn}
additional_fields={field1^value1|field2^value2}
user_ID={user_ID}
return_format={return_format}
```

Where:

patient_id is the unique system identifier of the patient to update.

patient_lastname is the last name of the patient.

patient_firstname is the first name of the patient.

patient_middleName is the middle name of the patient.

patient_street is the street of the patient.

patient_street2 is the street of the patient.

patient_city is the city of the patient.

patient_state is the state of the patient.

patient_zip is the zip code of the patient.

patient_sex is the sex of the patient.

patient_dob is the date of birth of the patient in the format: yyyy-mm-dd.

patient_ssn is the Social Security Number of the patient.

additional_fields is Pipe delimited fields to update when in use. FieldName^FieldValue|

user_id is the user code for the user in the NucleoLIS user table that is adding the patient.

return_format is the format for returning the data. Where:

0 = XML Format

1 = JSON Format

Example

```
https://www.SomeServer.com/api/N/
UpdatePatient?patient_id=15498920&patient_lastname=Smith&patient_firstname=John&patient_middleName=J&patient_street=1 Pine Street&patient_street2=&patient_city=Milford&patient_state=MA&patient_zip=01757&patient_Sex=M&patient_dob=1976-07-01&patient_ssn=123-456-0000&additional_fields=&user_ID=15492125&return_format=0
```

Returns

If the request succeeds, returns the updated patient identifier and status of processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <PatientTable>
    <object_type>PatientTable</object_type>
    <message>Patient Updated</message>
    <status>Success</status>
  </PatientTable>
</NewDataSet>
```

ADD CASE

Endpoint	api/N/AddCase
Purpose	Adds a case record.
Request Type	GET
Parameters	<pre>patient_id={patient_id} caseType_code={caseType_code} submittingPhysician_code={submittingPhysician_code} icds={icds} clinical_information={clinical_information} aoe_questions={aoe_questions} additional_fields={field1^value1 field2^value2} user_ID={user_ID} return_format={return_format}</pre>

Where:

patient_id is The unique system identifier of the patient in which to create a case.

caseType_code is The unique code of the case type definition.

submittingPhysician_code is The unique code of the submitting physician.

submittingLocation is The location of the submitting physician.

icds is Pipe Delimited icd codes.

clinical_information is Clinical Information.

aoe_questions is Not yet supported - Pipe Delimited AOE questions and answers Question^Response.

additional_fields is Pipe delimited fields to update when in use. FieldName^FieldValue|

user_id is the user code for the user in the NucleoLIS user table that is adding the case.

return_format is the format for returning the data.
Where:

0 = XML Format

1 = JSON Format

Example

```
https://www.SomeServer.com/api/N/
AddCase?patient_id=15498920&caseType_code=GEN&submittingPhysician_code=000000&submittingLocation=OFFICE&icds=001.0|001.1|001.9|002.0&clinical_information=assorted ailments&aoe_questions=&additional_fields=&user_ID=15492125&return_format=0
```

Returns

If the request succeeds, returns the system generated case identifier and status of processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <CaseTable>
    <object_type>CaseTable</object_type>
    <message>Case Added</message>
    <status>Success</status>
  </CaseTable>
</NewDataSet>
```

UPDATE CASE

Endpoint	api/N/UpdateCase
Purpose	Updates a case record.
Request Type	GET
Parameters	patient_id={patient_id}

```

case_id={case_id}

caseType_code={caseType_code}

submittingPhysician_code={submittingPhysician_code}

icds={icds}

clinical_information={clinical_information}

aoe_questions={aoe_questions}

additional_fields={field1^value1|field2^value2}

user_ID={user_ID}

return_format={return_format}

```

Where:

patient_id is the unique system identifier of the patient in which to create a case.

case_id is the unique system identifier of the case to update.

caseType_code is the unique code of the case type definition.

submittingPhysician_code is the unique code of the submitting physician.

submittingLocation is the location of the submitting physician.

icds is a list of pipe delimited icd codes.

clinical_information is Clinical Information.

aoe_questions is not yet supported - Pipe Delimited AOE questions and answers Question^Response.

additional_fields is Pipe delimited fields to update when in use. FieldName^FieldValue|

user_id is the user code for the user in the NucleoLIS user table that is updating the case.

return_format is the format for returning the data.
Where:

0 = XML Format
1 = JSON Format

Example `https://www.SomeServer.com/api/N/UpdateCase?patient_id=15498920&case_id=15498922&caseType_code=GEN&submittingPhysician_code=000000&submittingLocation=000000&icds=001.0&clinical_information=&aoe_questions=&additional_fields=&user_ID=15492125&return_format=0`

Returns If the request succeeds, returns the system generated case identifier and status of processing. Note: Case updates may be disallowed depending on the state of case processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <CaseTable>
    <object_type>CaseTable</object_type>
    <message>Case Updated</message>
    <status>Success</status>
  </CaseTable>
</NewDataSet>
```

ADD SPECIMEN

Endpoint	api/N/AddSpecimen
Purpose	Add a specimen record.
Request Type	GET
Parameters	patient_id={patient_id} case_id={case_id} sourceDefinition_code={sourceDefinition_code} collection_date={collection_date} collection_time={collection_time}

```

test_orders={test_orders}

test_order={test_order}

constituentDefinition_codes={constituent_code|constitu
ent_code}

aoe_questions={aoe_questions}

additional_fields={field1^value1|field2^value2}

user_ID={user_ID}

return_format={return_format}

```

Where:

patient_id is the unique system identifier of the patient in which to add a specimen.

case_id is the unique system identifier of the case in which to add a specimen.

sourceDefinition_code is the unique system code of the source of the specimen.

collection_date is the collection date in yyyy-mm-dd format.

collection_time is the collection time in hh:mm format.

test_orders is a list of pipe delimited Test Codes to order on the specimen.

test_order is a single Test Code to order on the specimen and used with constituentDefinition_codes.

constituentDefinition_codes is a list of pipe delimited Constituent Codes to order on the specimen. Must use test_order, not test_orders.

aoe_questions is Pipe Delimited AOE questions and answers Question^Response.

additional_fields is Pipe delimited fields to update when in use. FieldName^FieldValue|

user_id is the user code for the user in the NucleoLIS user table that is updating the case.

return_format is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

Example

```
https://www.SomeServer.com/api/N/
AddSpecimen?patient_id=15498920&case_id=15498922&sourceDefinition_code=BLOOD&collection_date=2021-04-
19&collection_time=11:25&test_Orders=ZIKA|XTAG|RPR&aoe
_questions=&additional_fields=&user_ID=15492125&return
_format=0
```

Returns

If the request succeeds, returns the system generated case identifier and status of processing. Note: Case updates may be disallowed depending on the state of processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <SpecimenTable>
    <object_type>SpecimenTable</object_type>
    <message>Specimen Added</message>
    <status>Success</status>
  </SpecimenTable>
</NewDataSet>
```

UPDATE SPECIMEN

Endpoint	api/N/UpdateSpecimen
Purpose	Updates a specimen record.
Request Type	GET
Parameters	patient_id={patient_id} case_id={case_id} specimen_id = {specimen_id}

```

sourceDefinition_code={sourceDefinition_code}

collection_date={collection_date}

collection_time={collection_time}

aoe_questions={aoe_questions}

additional_fields={field1^value1|field2^value2}

user_ID={user_ID}

return_format={return_format}

```

Where:

patient_id is the unique system identifier of the patient in which to add a specimen.

case_id is the unique system identifier of the case in which to add a specimen.

specimen_id is the unique system identifier of the specimen in which to update.

sourceDefinition_code is the unique system code of the source of the specimen.

collection_date is the collection date in yyyy-mm-dd format.

collection_time is the collection time in hh:mm format.

aoe_questions is Pipe Delimited AOE questions and answers Question^Response.

additional_fields is Pipe delimited fields to update when in use. FieldName^FieldValue|

user_id is the user code for the user in the NucleoLIS user table that is updating the specimen.

return_format is the format for returning the data.
Where:

0 = XML Format

1 = JSON Format

Example `https://www.SomeServer.com/api/N/UpdateSpecimen?patient_id=15498920&case_id=15498922&specimen_id=15498927&sourceDefinition_code=BLOOD&collection_date=2021-04-20&collection_time=11:00&aoe_questions=Last Dose^15mg|Last Does Date^2021-03-19&additional_fields=&user_ID=15492125&return_format=0`

Returns If the request succeeds, returns the system generated specimen identifier and status of processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <SpecimenTable>
    <object_type>SpecimenTable</object_type>
    <message>Specimen Updated</message>
    <status>Success</status>
  </SpecimenTable>
</NewDataSet>
```

ADD TEST ORDER

Endpoint	api/N/AddTestOrder
Purpose	Adds a test order.
Request Type	GET
Parameters	<pre>patient_id={patient_id} case_id={case_id} specimen_id = {specimen_id} profileDefinition_code={profileDefinition_code} constituentDefinition_codes={constituent_code constituent_code} user_ID={user_ID} return_format={return_format}</pre>

Where:

patient_id is the unique system identifier of the patient in which to add a specimen.

case_id is the unique system identifier of the case in which to add a test.

specimen_id is the unique system identifier of the specimen in which to add to a test.

profileDefinition_code is the unique system code of the profile to order.

constituentDefinition_codes is a list of pipe delimited Constituent Codes to order on the specimen.

user_id is the user code for the user in the NucleoLIS user table that is adding the test order.

return_format is the format for returning the data.
Where:

0 = XML Format

1 = JSON Format

Example

```
https://www.SomeServer.com/api/N/
AddTestOrder?patient_id=15498920&case_id=15498948&specimen_id=15498949&profileDefinition_code=CLL&user_ID=15492125&return_format=0
```

Returns

If the request succeeds, returns the status of processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <SpecimenTable>
    <object_type>SpecimenTable</object_type>
    <message>Test Added</message>
    <status>Success</status>
  </SpecimenTable>
</NewDataSet>
```


CANCEL CASE

Endpoint api/N/CancelCase

Purpose Cancels a case record.

Request Type GET

Parameters

```
patient_id={patient_id}
case_id={case_id}
user_ID={user_ID}
return_format={return_format}
```

Where:

patient_id is the unique system identifier of the patient in which to cancel the case.

case_id is the unique system identifier of the case to cancel.

user_id is the user code for the user in the NucleoLIS user table that is canceling the case.

return_format is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

Example

```
https://www.SomeServer.com/api/N/
CancelCase?patient_id=15498977&case_id=15498979&user_I
D=15492125&return_format=0
```

Returns If the request succeeds, returns the status of processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <CaseTable>
    <object_type>CaseTable</object_type>
    <message>Case Canceled</message>
    <status>Success</status>
  </CaseTable>
</NewDataSet>
```

CANCEL SPECIMEN

Endpoint `api/N/CancelSpecimen`

Purpose Cancels a specimen record.

Request Type GET

Parameters `patient_id={patient_id}`
`case_id={case_id}`
`specimen_id={specimen_id}`
`user_ID={user_ID}`
`return_format={return_format}`

Where:

patient_id is the unique system identifier of the patient in which to cancel a specimen.

case_id is the unique system identifier of the case of the specimen to cancel.

specimen_id is the unique system identifier of the specimen to cancel.

user_id is the user code for the user in the NucleoLIS user table that is canceling the specimen.

return_format is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

Example `https://www.SomeServer.com/api/N/CancelSpecimen?patient_id=15498977&case_id=15498979&specimen_id=15498980&user_ID=15492125&return_format=0`

Returns If the request succeeds, returns the status of processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <SpecimenTable>
    <object_type>SpecimenTable</object_type>
    <message>Specimen Canceled</message>
    <status>Success</status>
  </SpecimenTable>
</NewDataSet>
```

CANCEL TEST ORDER

Endpoint api/N/CancelTestOrder

Purpose Cancels a test order.

Request Type GET

Parameters patient_id={patient_id}
 case_id={case_id}
 specimen_id={specimen_id}
 profileDefinition_Code={ profileDefinition_Code}
 user_ID={user_ID}
 return_format={return_format}

Where:

patient_id is the unique system identifier of the patient in which to cancel a test.

case_id is the unique system identifier of the case in which to cancel a test.

specimen_id is the unique system identifier of the specimen to cancel a test.

profileDefinition_Code is the unique system code of the profile to cancel.

user_id is the user code for the user in the NucleoLIS user table that is canceling the specimen.

return_format is the format for returning the data.
Where:

0 = XML Format

1 = JSON Format

Example

```
https://www.SomeServer.com/api/N/
CancelTestOrder?patient_id=15498977&case_id=15498979&specimen_id=15498980&profileDefinition_Code=ZIKA&user_ID=15492125&return_format=0
```

Returns

If the request succeeds, returns the status of processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <SpecimenTable>
    <object_type>SpecimenTable</object_type>
    <message>Test Canceled</message>
    <status>Success</status>
  </SpecimenTable>
</NewDataSet>
```

GET PLATE MAP DEFINITIONS

Endpoint

api/N/GetPlateMapDefinitions

Purpose

Returns a collection of worksheet templates defined as plate maps.

Request Type

GET

Parameters

user_ID={user_ID}

return_format={return_format}

with_retired={with_retired}

Where:

user_id is the user code for the user in the NucleoLIS user table that is requesting the plate map definitions.

return_format is the format for returning the data.
Where:

0 = XML Format

1 = JSON Format

with_retired if false will exclude retired entries in the system build. If true, all entries will be returned.

Example

```
https://www.SomeServer.com/api/N/
GetPlateMapDefinitions?user_ID=15492125&return_format=
0&with_retired=false
```

Returns

If the request succeeds, returns the system unique identifier, and other descriptive information:

```
<PlateMapDefinitionTable>
```

```
<WorksheetTemplate_ObjectID>15394236</WorksheetTemplate_ObjectID>
```

```
<WorksheetTemplate_CreationDate>2019-12-
03T06:05:40.123-
```

```
05:00</WorksheetTemplate_CreationDate>
```

```
<WorksheetTemplate_UpdateTime>2021-04-
```

```
23T07:06:39.687-04:00</WorksheetTemplate_UpdateTime>
```

```
<WorksheetTemplate_User1 />
```

```
<WorksheetTemplate_User2 />
```

```
<WorksheetTemplate_User3 />
```

```
<WorksheetTemplate_User4 />
```

```
<WorksheetTemplate_CREATOR>johnb</WorksheetTemplate_CREATOR>
```

```
<WorksheetTemplate_MODIFIER>johnb</WorksheetTemplate_MODIFIER>
```

```
<WorksheetTemplate_Code>Extraction</WorksheetTemplate_Code>
```

```

<WorksheetTemplate_Description>Extraction</WorksheetTemplate_Description>

<WorksheetTemplate_RETIREDD>False</WorksheetTemplate_RETIREDD>

<WorksheetTemplate_CurrentVersion>8</WorksheetTemplate_CurrentVersion>

<WorksheetTemplate_Locked>False</WorksheetTemplate_Locked>
    <WorksheetTemplate_LockedBy />
    <WorksheetTemplate_LockedDate />
    <WorksheetTemplate_LockedTime />

<WorksheetTemplate_IsVariant>False</WorksheetTemplate_IsVariant>
    <WorksheetTemplate_VariantMap />

<WorksheetTemplate_HasPlateMap>True</WorksheetTemplate_HasPlateMap>

<WorksheetTemplate_PlateMapTabID>Extraction</WorksheetTemplate_PlateMapTabID>

<WorksheetTemplate_PlateMapStartRow>2</WorksheetTemplate_PlateMapStartRow>

<WorksheetTemplate_PlateMapEndRow>25</WorksheetTemplate_PlateMapEndRow>

<WorksheetTemplate_PlateMapPositionColumn>A</WorksheetTemplate_PlateMapPositionColumn>

<WorksheetTemplate_PlateMapIDColumn>B</WorksheetTemplate_PlateMapIDColumn>

<WorksheetTemplate_PlateMapTypeColumn>C</WorksheetTemplate_PlateMapTypeColumn>
    </PlateMapDefinitionTable>
    <PlateMapDefinitionTable>

<WorksheetTemplate_ObjectID>15398875</WorksheetTemplate_ObjectID>

```

```

        <WorksheetTemplate_CreationDate>2020-05-
06T10:07:53.493-
04:00</WorksheetTemplate_CreationDate>
        <WorksheetTemplate_UpdateTime>2021-04-
23T07:15:00.657-04:00</WorksheetTemplate_UpdateTime>
        <WorksheetTemplate_User1 />
        <WorksheetTemplate_User2 />
        <WorksheetTemplate_User3 />
        <WorksheetTemplate_User4 />

<WorksheetTemplate_CREATOR>johnb</WorksheetTemplate_C
REATOR>

<WorksheetTemplate_MODIFIER>johnb</WorksheetTemplate_
MODIFIER>

<WorksheetTemplate_Code>PCRFILL</WorksheetTemplate_Co
de>

<WorksheetTemplate_Description>PCRFILL</WorksheetTemp
late_Description>

<WorksheetTemplate_RETIRED>False</WorksheetTemplate_R
ETIRED>

<WorksheetTemplate_CurrentVersion>3</WorksheetTemplat
e_CurrentVersion>

<WorksheetTemplate_Locked>False</WorksheetTemplate_Lo
cked>
        <WorksheetTemplate_LockedBy />
        <WorksheetTemplate_LockedDate />
        <WorksheetTemplate_LockedTime />

<WorksheetTemplate_IsVariant>False</WorksheetTemplate
_IsVariant>
        <WorksheetTemplate_VariantMap />

<WorksheetTemplate_HasPlateMap>True</WorksheetTemplat
e_HasPlateMap>

<WorksheetTemplate_PlateMapTabID>Sheet1</WorksheetTem
plate_PlateMapTabID>

<WorksheetTemplate_PlateMapStartRow>2</WorksheetTempl
ate_PlateMapStartRow>

```

```
<WorksheetTemplate_PlateMapEndRow>97</WorksheetTemplate_PlateMapEndRow>
```

```
<WorksheetTemplate_PlateMapPositionColumn>F</WorksheetTemplate_PlateMapPositionColumn>
```

```
<WorksheetTemplate_PlateMapIDColumn>C</WorksheetTemplate_PlateMapIDColumn>
```

```
<WorksheetTemplate_PlateMapTypeColumn>D</WorksheetTemplate_PlateMapTypeColumn>  
</PlateMapDefinitionTable>  
</NewDataSet>
```

GET PLATES

Endpoint	api/N/GetPlates
Purpose	Returns a collection of committed plate maps.
Request Type	GET
Parameters	<pre>template_code={template_code} min_date={min_date} max_date={max_date} user_ID={user_ID} return_format={return_format}</pre>

Where:

template_code is the system unique code value for the plate map query.

min_date is the earliest commit date to include in the result set, inclusive in the format YYYY-MM-DD.

max_date is the latest commit date to include in the result set, inclusive in the format YYYY-MM-DD.

profileDefinition_Code is the unique system code of the profile to cancel.

user_id is the user code for the user in the NucleoLIS user table that is making the plate map request.

return_format is the format for returning the data.
Where:

0 = XML Format

1 = JSON Format

Example

```
https://www.SomeServer.com/api/N/
GetPlates?template_code=PCRFILL&min_date=2021-04-
23&max_date=2021-04-
26&user_ID=15492125&return_format=0
```

Returns

If the request succeeds, returns the system unique plate identifier, committed date and time, notes and template definition code for plate matching the criteria:

```
<PlateMapTable>
  <Worksheet_ObjectID>15499389</Worksheet_ObjectID>

  <Worksheet_WorksheetName>15499389</Worksheet_WorksheetName>

  <Worksheet_CommittedBy>johnb</Worksheet_CommittedBy>
    <Worksheet_CommittedDate>2021-04-
23</Worksheet_CommittedDate>

  <Worksheet_CommittedTime>07:24</Worksheet_CommittedTime>
    <Worksheet_Notes />

  <WorksheetTemplate_Code>PCRFILL</WorksheetTemplate_Code>

  <WorksheetTemplate_ObjectID>15398875</WorksheetTemplate_ObjectID>
</PlateMapTable>
</NewDataSet>
```

GET PLATE MAP

Endpoint	api/N/GetPlateMap
Purpose	Returns the map for the unique identifier passed.
Request Type	GET
Parameters	<code>platemap_id={platemap_id}</code> <code>user_ID={user_ID}</code> <code>return_format={return_format}</code>

Where:

platemap_id is the unique identifier of the plate.

user_id is the user code for the user in the NucleoLIS user table that is requesting the specified plate map.

return_format is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

Example	<code>https://www.SomeServer.com/api/N/GetPlateMap?platemad_id=15499389&user_ID=15492125&return_format=0</code>
----------------	---

Returns	If the request succeeds, returns the plate map including plate position, specimen identifier, and additional configured plate map related information:
----------------	--

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <PlateMapTable>
    <plate_map_position>A1</plate_map_position>
    <plate_map_id>GEN2020-1000185</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>A2</plate_map_position>
    <plate_map_id>GEN2020-1000185</plate_map_id>
```

```

    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>A3</plate_map_position>
    <plate_map_id>GEN2020-1000188</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>A4</plate_map_position>
    <plate_map_id>GEN2020-1000188</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>A5</plate_map_position>
    <plate_map_id>GEN2020-1000187</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>A6</plate_map_position>
    <plate_map_id>GEN2020-1000187</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>A7</plate_map_position>
    <plate_map_id>GEN2020-1000191</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>A8</plate_map_position>
    <plate_map_id>GEN2020-1000191</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>B1</plate_map_position>
    <plate_map_id>GEN2020-1000200</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>B2</plate_map_position>
    <plate_map_id>GEN2020-1000200</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>B3</plate_map_position>
    <plate_map_id>GEN2020-1000195</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>

```

```

<PlateMapTable>
  <plate_map_position>B4</plate_map_position>
  <plate_map_id>GEN2020-1000195</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>B5</plate_map_position>
  <plate_map_id>GEN2020-1000201</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>B6</plate_map_position>
  <plate_map_id>GEN2020-1000201</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>B7</plate_map_position>
  <plate_map_id>GEN2020-1000202</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>B8</plate_map_position>
  <plate_map_id>GEN2020-1000202</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>C1</plate_map_position>
  <plate_map_id>GEN2020-1000193</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>C2</plate_map_position>
  <plate_map_id>GEN2020-1000193</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>C3</plate_map_position>
  <plate_map_id>GEN2020-1000203</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>C4</plate_map_position>
  <plate_map_id>GEN2020-1000203</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>C5</plate_map_position>

```

```

    <plate_map_id>GEN2020-1000196</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>C6</plate_map_position>
    <plate_map_id>GEN2020-1000196</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>C7</plate_map_position>
    <plate_map_id>GEN2020-1000190</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>C8</plate_map_position>
    <plate_map_id>GEN2020-1000190</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
</NewDataSet>

```

GET CHARGE DEFINITIONS

Endpoint	api/N/GetChargeDefinitions
Purpose	Gets data associated with the charge code library (general call access)
Request Type	GET
Parameters	user_ID={user_ID} return_format={return_format}

Where:

user_id is the user code for the user in the NucleoLIS user table that is making the charge definitions request.

return_format is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

Example

https://www.SomeServer.com/api/N/GetChargeDefinitions?user_ID=15492125&return_format=0&with_retired=false

Returns

If the request succeeds, returns the system charge definitions defined:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <ChargeDefinitionTable>

    <ChargeDefinition_ObjectID>22711758</ChargeDefinition_ObjectID>

    <ChargeDefinition_Code>88305</ChargeDefinition_Code>
      <ChargeDefinition_Description>Level IV - Surgical pathology, gross and microscopic examination</ChargeDefinition_Description>

    <ChargeDefinition_Retired>False</ChargeDefinition_Retired>
  </ChargeDefinitionTable>
  <ChargeDefinitionTable>

    <ChargeDefinition_ObjectID>22868612</ChargeDefinition_ObjectID>

    <ChargeDefinition_Code>88300</ChargeDefinition_Code>
      <ChargeDefinition_Description>Level I - Surgical pathology, gross examination only</ChargeDefinition_Description>

    <ChargeDefinition_ChargeType>Examination</ChargeDefinition_ChargeType>

    <ChargeDefinition_Retired>False</ChargeDefinition_Retired>
  </ChargeDefinitionTable>
  <ChargeDefinitionTable>

    <ChargeDefinition_ObjectID>22868641</ChargeDefinition_ObjectID>

    <ChargeDefinition_Code>88302</ChargeDefinition_Code>
```

```

    <ChargeDefinition_Description>Level II - Surgical
    pathology, gross and microscopic
    examination</ChargeDefinition_Description>

    <ChargeDefinition_Retired>False</ChargeDefinition_Ret
    ired>
    </ChargeDefinitionTable>
</NewDataSet>

```

GET CHARGES

Endpoint	api/N/GetCharges
Purpose	Gets summary of charge data (requires billing module access)
Request Type	GET
Parameters	charge_date_from={min_date} charge_date_to={max_date} actives_only={true/false} filter_expression={filter expression} user_ID={user_ID} return_format={return_format}

Where:

charge_date_from is the earliest charge date to include in the result set, inclusive in the format YYYY-MM-DD.

charge_date_to is the latest charge date to include in the result set, inclusive in the format YYYY-MM-DD.

activesOnly is whether the case is active or not

filter_expression is used as a secondary custom filtering criteria applied to the result set. Wildcards are supported. It can be used with any

column returned in the result set. Any valid SQL WHERE clause is accepted.

user_id is the user code for the user in the NucleoLIS user table that is making the charges request.

return_format is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

Example

https://www.SomeServer.com/api/N/GetCharges?charge_date_from=2021-04-23&charge_date_to=2021-04-26&actives_only=false&filter_expression=&user_ID=15492125&return_format=0

Returns

If the request succeeds, returns the system charges matching the criteria:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <ChargeSlipTable>
    <Patient_ObjectID>549052794</Patient_ObjectID>
    <Patient_Name>Caron, Jack</Patient_Name>
    <Patient_Code>549052794</Patient_Code>
    <Patient_Sex>Male</Patient_Sex>
    <Patient_Dob>2019-06-15</Patient_Dob>
    <Case_Number>CL2022-000002</Case_Number>

    <ChargeSlip_ObjectID>572543969</ChargeSlip_ObjectID>
    <ChargeSlip_Parent>572543311</ChargeSlip_Parent>

    <ChargeSlip_ParentType>Specimen</ChargeSlip_ParentType>
  </ChargeSlipTable>

  <ChargeSlip_Code>80333</ChargeSlip_Code>
  <ChargeSlip_ChargeDate>2022-08-23</ChargeSlip_ChargeDate>
  <ChargeSlip_Credit>False</ChargeSlip_Credit>

  <ChargeSlip_ChargeType>Staining</ChargeSlip_ChargeType>
  <ChargeSlip_Description>Immuno/Clinical-HH</ChargeSlip_Description>
</NewDataSet>
```



```

<ChargeSlipTable>
  <Patient_ObjectID>536351343</Patient_ObjectID>
  <Patient_Name>Caron, Bobby</Patient_Name>
  <Patient_Code>536351343</Patient_Code>
  <Patient_Sex>Male</Patient_Sex>
  <Patient_Dob>2005-11-05</Patient_Dob>
  <Case_Number>CL2022-000069</Case_Number>

<ChargeSlip_ObjectID>573231128</ChargeSlip_ObjectID>
  <ChargeSlip_Parent>536351358</ChargeSlip_Parent>

<ChargeSlip_ParentType>Specimen</ChargeSlip_ParentType>
e>
  <ChargeSlip_Code>80333</ChargeSlip_Code>
  <ChargeSlip_ChargeDate>2022-08-
24</ChargeSlip_ChargeDate>
  <ChargeSlip_Credit>False</ChargeSlip_Credit>

<ChargeSlip_ChargeType>Staining</ChargeSlip_ChargeType>
e>
  <ChargeSlip_Description>Immuno/Clinical-
HH</ChargeSlip_Description>
</ChargeSlipTable>
<ChargeSlipTable>
  <Patient_ObjectID>574218209</Patient_ObjectID>
  <Patient_Name>Caron, Jack</Patient_Name>
  <Patient_Code>574218209</Patient_Code>
  <Patient_Sex>Male</Patient_Sex>
  <Patient_Dob>2019-06-15</Patient_Dob>
  <Case_Number>CL2022-000074</Case_Number>

<ChargeSlip_ObjectID>574218878</ChargeSlip_ObjectID>
  <ChargeSlip_Parent>574218842</ChargeSlip_Parent>

<ChargeSlip_ParentType>Specimen</ChargeSlip_ParentType>
e>
  <ChargeSlip_Code>80333</ChargeSlip_Code>
  <ChargeSlip_ChargeDate>2022-08-
25</ChargeSlip_ChargeDate>
  <ChargeSlip_Credit>False</ChargeSlip_Credit>

<ChargeSlip_ChargeType>Staining</ChargeSlip_ChargeType>
e>
  <ChargeSlip_Description>Immuno/Clinical-
HH</ChargeSlip_Description>
</ChargeSlipTable>
<ChargeSlipTable>

```

```

    <Patient_ObjectID>574341267</Patient_ObjectID>
    <Patient_Name>Cricket, Jiminy</Patient_Name>
    <Patient_Code>574341267</Patient_Code>
    <Patient_Sex>M</Patient_Sex>
    <Patient_Dob>2018-11-13</Patient_Dob>
    <Case_Number>CL2022-000075</Case_Number>

    <ChargeSlip_ObjectID>574352641</ChargeSlip_ObjectID>
    <ChargeSlip_Parent>574341274</ChargeSlip_Parent>

    <ChargeSlip_ParentType>Specimen</ChargeSlip_ParentType>
    <ChargeSlip_Code>80333</ChargeSlip_Code>
    <ChargeSlip_ChargeDate>2022-08-25</ChargeSlip_ChargeDate>
    <ChargeSlip_Credit>False</ChargeSlip_Credit>

    <ChargeSlip_ChargeType>Staining</ChargeSlip_ChargeType>
    <ChargeSlip_Description>Immuno/Clinical-HH</ChargeSlip_Description>
  </ChargeSlipTable>
</NewDataSet>

```

GET CHARGE

Endpoint	api/N/GetCharge
Purpose	Gets all data associated with a charge (requires billing module access)
Request Type	GET
Parameters	<pre> charge_id={charge slip objectID} user_ID={user_ID} return_format={return_format} </pre>

Where:

charge_id is the charge slip ObjectID to include in the result set.

user_id is the user code for the user in the NucleoLIS user table that is making the charge request.

return_format is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

Example

https://www.SomeServer.com/api/N/GetCharge?charge_id=11371897&user_ID=15492125&return_format=0

Returns

If the request succeeds, returns the system unique charge identifier matching the criteria:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <ChargeSlipTable>
    <Patient_ObjectID>175277649</Patient_ObjectID>
    <Patient_Name>Knowles, Beyonce</Patient_Name>
    <Patient_Code>175277649</Patient_Code>
    <Patient_Sex>Female</Patient_Sex>
    <Patient_Dob>1981-09-04</Patient_Dob>
    <Case_Number>M2022-000070</Case_Number>
    <Case_Status>A</Case_Status>

    <ChargeSlip_ObjectID>581034845</ChargeSlip_ObjectID>
    <ChargeSlip_CreationDate>2022-09-02T16:56:05.1-04:00</ChargeSlip_CreationDate>
    <ChargeSlip_UpdateTime>2022-09-02T16:56:05.1-04:00</ChargeSlip_UpdateTime>

    <ChargeSlip_SuperobjectID>559484547</ChargeSlip_SuperobjectID>
    <ChargeSlip_User1 />
    <ChargeSlip_User2 />
    <ChargeSlip_User3 />
    <ChargeSlip_User4 />
    <ChargeSlip_CREATOR>Scaron</ChargeSlip_CREATOR>
    <ChargeSlip_MODIFIER>Scaron</ChargeSlip_MODIFIER>
    <ChargeSlip_Parent>559484549</ChargeSlip_Parent>

    <ChargeSlip_ParentType>Specimen</ChargeSlip_ParentType>
  </ChargeSlipTable>
</NewDataSet>
```

```

    <ChargeSlip_Code>87633</ChargeSlip_Code>
    <ChargeSlip_ChargeDate>2022-09-
02</ChargeSlip_ChargeDate>
    <ChargeSlip_Billed>none</ChargeSlip_Billed>
    <ChargeSlip_Credit>False</ChargeSlip_Credit>
    <ChargeSlip_ChargeType />

<ChargeSlip_ChargeObject>559484549</ChargeSlip_Charge
Object>
    <ChargeSlip_Description>Specimen-
NP</ChargeSlip_Description>
    </ChargeSlipTable>
</NewDataSet>

```

APPENDIX B – CODING EXAMPLES

CODING EXAMPLES INTRODUCTION

Note: The coding examples in this appendix assume that the URL being connected to is <https://www.SomeServer.com>, and the username = LAB, and the password = 123456. These examples show how to connect to the Backbone API using a variety of programming languages and techniques. These examples are solely to help you get started communicating with the API, and are not intended to be used as production code.

For more information on the various frameworks mentioned in the code examples, please refer to the links below:

Language	Framework
C	https://curl.haxx.se/libcurl/c/
C#	http://restsharp.org/
cURL	https://curl.haxx.se/
Go	https://golang.org/pkg/net/http/
HTTP	None (Raw HTTP request)
Java OkHTTP	https://github.com/square/okhttp
Java UniREST	https://github.com/Kong/unirest-java
JavaScript Fetch	https://developer.mozilla.org/en-US/docs/Web/API/Fetch_API
JavaScript jQuery	http://api.jquery.com/jquery.ajax/
JavaScript XHR	https://developer.mozilla.org/en-US/docs/Web/API/XMLHttpRequest
NodeJS Axios	https://github.com/axios/axios
NodeJS Http	https://nodejs.org/api/http.html
NodeJS Request	https://github.com/request/request

NodeJS UniRest	https://github.com/Kong/unirest-nodejs
Objective-C	https://developer.apple.com/library/ios/documentation/Foundation/Reference/NSURLSession_class/
OCaml	https://github.com/mirage/ocaml-cohttp
PHP CURl	http://php.net/manual/en/ref.curl.php
PHP HttpRequest	https://www.php.net/manual/en/reserved.variables.request.php
PHP pecl_http	https://mdref.m6w6.name/http
PowerShell RestMethod	https://docs.microsoft.com/en-us/powershell/module/microsoft.powershell.utility/invoke-restmethod?view=powershell-7
Python Built-in http.client (Python 3)	https://docs.python.org/3/library/http.client.html
Python Requests	https://requests.readthedocs.io/en/master/
Ruby Built-in NET::Http	http://docs.ruby-lang.org/en/2.0.0/Net/HTTP.html
Shell HTTPie	https://github.com/jkbrzt/httpie
Shell wget	https://www.gnu.org/software/wget/
Shell cURL	https://curl.haxx.se/
Swift Built-in NSURLSession	https://developer.apple.com/documentation/foundation/urlsession
VB.NET RestSharp	http://restsharp.org/

LOGON CODING EXAMPLES

Note: The Logon endpoint is a POST command. Any API call that uses POST would follow a similar pattern.

C – LIBCURL

```
CURL *curl;
CURLcode res;
curl = curl_easy_init();
if(curl) {
    curl_easy_setopt(curl, CURLOPT_CUSTOMREQUEST, "POST");
    curl_easy_setopt(curl, CURLOPT_URL, "https://www.SomeServer.com/api/authenticate/logon");
    curl_easy_setopt(curl, CURLOPT_FOLLOWLOCATION, 1L);
    curl_easy_setopt(curl, CURLOPT_DEFAULT_PROTOCOL, "https");
    struct curl_slist *headers = NULL;
    headers = curl_slist_append(headers, "Content-Type: application/x-www-form-urlencoded");
    curl_easy_setopt(curl, CURLOPT_HTTPHEADER, headers);
    const char *data = "username=LAB&password=123456";
    curl_easy_setopt(curl, CURLOPT_POSTFIELDS, data);
    res = curl_easy_perform(curl);
}
curl_easy_cleanup(curl);
```

C# - RESTSHARP

```
var client = new RestClient("https://www.SomeServer.com/api/authenticate/logon");
client.Timeout = -1;
var request = new RestRequest(Method.POST);
request.AddHeader("Content-Type", "application/x-www-form-urlencoded");
request.AddParameter("username", "LAB");
request.AddParameter("password", "123456");
IRestResponse response = client.Execute(request);
Console.WriteLine(response.Content);
```

CURL

```
curl --location --request POST 'https://www.SomeServer.com/api/authenticate/logon' \
--header 'Content-Type: application/x-www-form-urlencoded' \
--data-urlencode 'username=LAB' \
--data-urlencode 'password=123456'
```

GO - NATIVE

```
package main

import (
    "fmt"
    "strings"
    "net/http"
    "io/ioutil"
)

func main() {

    url := "https://www.SomeServer.com/api/authenticate/logon"
    method := "POST"

    payload := strings.NewReader("username=LAB&password=123456")

    client := &http.Client {
    }
    req, err := http.NewRequest(method, url, payload)

    if err != nil {
        fmt.Println(err)
        return
    }
    req.Header.Add("Content-Type", "application/x-www-form-urlencoded")

    res, err := client.Do(req)
    if err != nil {
        fmt.Println(err)
        return
    }
    defer res.Body.Close()

    body, err := ioutil.ReadAll(res.Body)
    if err != nil {
        fmt.Println(err)
        return
    }
    fmt.Println(string(body))
}
```


HTTP

POST /WebApiBackboneQA/api/authenticate/logon? HTTP/1.1

Host: SomeServer.com

Content-Type: application/x-www-form-urlencoded

Content-Length: 28

username=LAB&password=123456

JAVA - OKHTTP

```
OkHttpClient client = new OkHttpClient().newBuilder()
    .build();
MediaType mediaType = MediaType.parse("application/x-www-form-urlencoded");
RequestBody body = RequestBody.create(mediaType, "username=LAB&password=123456");
Request request = new Request.Builder()
    .url("https://www.SomeServer.com/api/authenticate/logon")
    .method("POST", body)
    .addHeader("Content-Type", "application/x-www-form-urlencoded")
    .build();
Response response = client.newCall(request).execute();
```

JAVA – UNIREST

```
Unirest.setTimeouts(0, 0);
HttpResponse<String> response = Unirest.post("https://www.SomeServer.com/api/authenticate/logo
n")
    .header("Content-Type", "application/x-www-form-urlencoded")
    .field("username", "LAB")
    .field("password", "123456")
    .asString();
```

JAVASCRIPT - FETCH

```
var myHeaders = new Headers();
myHeaders.append("Content-Type", "application/x-www-form-urlencoded");

var urlencoded = new URLSearchParams();
urlencoded.append("username", "LAB");
urlencoded.append("password", "123456");

var requestOptions = {
```

```

    method: 'POST',
    headers: myHeaders,
    body: urlencoded,
    redirect: 'follow'
  };

  fetch("https://www.SomeServer.com/api/authenticate/logon", requestOptions)
    .then(response => response.text())
    .then(result => console.log(result))
    .catch(error => console.log('error', error));

```

JAVASCRIPT - JQUERY

```

var settings = {
  "url": "https://www.SomeServer.com/api/authenticate/logon",
  "method": "POST",
  "timeout": 0,
  "headers": {
    "Content-Type": "application/x-www-form-urlencoded",
  },
  "data": {
    "username": "LAB",
    "password": "123456"
  }
};

$.ajax(settings).done(function (response) {
  console.log(response);
});

```

JAVASCRIPT – XHR

```

var data = "username=LAB&password=123456";
var xhr = new XMLHttpRequest();
xhr.withCredentials = true;

xhr.addEventListener("readystatechange", function() {
  if(this.readyState === 4) {
    console.log(this.responseText);
  }
});

xhr.open("POST", "https://www.SomeServer.com/api/authenticate/logon");
xhr.setRequestHeader("Content-Type", "application/x-www-form-urlencoded");

xhr.send(data);

```

NODEJS – AXIOS

```
var axios = require('axios');
var qs = require('qs');
var data = qs.stringify({
  'username': 'LAB',
  'password': '123456'
});
var config = {
  method: 'post',
  url: 'https://www.SomeServer.com/api/authenticate/logon',
  headers: {
    'Content-Type': 'application/x-www-form-urlencoded',
  },
  data : data
};

axios(config)
  .then(function (response) {
    console.log(JSON.stringify(response.data));
  })
  .catch(function (error) {
    console.log(error);
  });
```

NODEJS – NATIVE

```
var https = require('follow-redirects').https;
var fs = require('fs');

var qs = require('querystring');

var options = {
  'method': 'POST',
  'hostname': 'SomeServer.com',
  'path': '/WebApiBackboneQA/api/authenticate/logon?',
  'headers': {
    'Content-Type': 'application/x-www-form-urlencoded',
  },
  'maxRedirects': 20
};

var req = https.request(options, function (res) {
  var chunks = [];

  res.on("data", function (chunk) {
```

```

        chunks.push(chunk);
    });

    res.on("end", function (chunk) {
        var body = Buffer.concat(chunks);
        console.log(body.toString());
    });

    res.on("error", function (error) {
        console.error(error);
    });
});

var postData = qs.stringify({
    'username': 'LAB',
    'password': '123456'
});

req.write(postData);

req.end();

```

NODEJS – REQUEST

```

var request = require('request');
var options = {
    'method': 'POST',
    'url': 'https://www.SomeServer.com/api/authenticate/logon',
    'headers': {
        'Content-Type': 'application/x-www-form-urlencoded',
    },
    form: {
        'username': 'LAB',
        'password': '123456'
    }
};
request(options, function (error, response) {
    if (error) throw new Error(error);
    console.log(response.body);
});

```

NODEJS – UNIREST

```

var unirest = require('unirest');
var req = unirest('POST', 'https://www.SomeServer.com/api/authenticate/logon')
    .headers({
        'Content-Type': 'application/x-www-form-urlencoded',
    })

```

```

    })
    .send('username=LAB')
    .send('password=123456')
    .end(function (res) {
        if (res.error) throw new Error(res.error);
        console.log(res.raw_body);
    });

```

OBJECTIVE-C NSURLSESSION

```
#import <Foundation/Foundation.h>
```

```
dispatch_semaphore_t sema = dispatch_semaphore_create(0);
```

```

NSMutableURLRequest *request = [NSMutableURLRequest requestWithURL:[NSURL URLWithString:@"http
s://www.SomeServer.com/api/authenticate/login"]
    cachePolicy:NSURLRequestUseProtocolCachePolicy
    timeoutInterval:10.0];
NSDictionary *headers = @{
    @"Content-Type": @"application/x-www-form-urlencoded",
};

```

```

[request setAllHTTPHeaderFields:headers];
NSMutableData *postData = [[NSMutableData alloc] initWithData:@"username=LAB" dataUsingEncoding:
:NSUTF8StringEncoding];
[postData appendData:@"&password=123456" dataUsingEncoding:NSUTF8StringEncoding];
[request setHTTPBody:postData];

```

```
[request setHTTPMethod:@"POST"];
```

```

NSURLSession *session = [NSURLSession sharedSession];
NSURLSessionDataTask *dataTask = [session dataTaskWithRequest:request
    completionHandler:^(NSData *data, NSURLResponse *response, NSError *error) {
    if (error) {
        NSLog(@"%@", error);
    } else {
        NSHTTPURLResponse *httpResponse = (NSHTTPURLResponse *) response;
        NSError *parseError = nil;
        NSDictionary *responseDictionary = [NSJSONSerialization JSONObjectWithData:data options:0
            error:&parseError];
        NSLog(@"%@", responseDictionary);
        dispatch_semaphore_signal(sema);
    }
}];
[dataTask resume];
dispatch_semaphore_wait(sema, DISPATCH_TIME_FOREVER);

```

OCAML – COHTTP

```
open Lwt
open Cohttp
open Cohttp_lwt_unix

let postData = ref "username=LAB&password=123456";;

let reqBody =
  let uri = Uri.of_string "https://www.SomeServer.com/api/authenticate/logon" in
  let headers = Header.init ()
  |> fun h -> Header.add h "Content-Type" "application/x-www-form-urlencoded"
  in
  let body = Cohttp_lwt.Body.of_string !postData in

  Client.call ~headers ~body `POST uri >>= fun (_resp, body) ->
  body |> Cohttp_lwt.Body.to_string >|= fun body -> body

let () =
  let respBody = Lwt_main.run reqBody in
  print_endline (respBody)
```

PHP – CURL

<?php

```
$curl = curl_init();

curl_setopt_array($curl, array(
  CURLOPT_URL => "https://www.SomeServer.com/api/authenticate/logon",
  CURLOPT_RETURNTRANSFER => true,
  CURLOPT_ENCODING => "",
  CURLOPT_MAXREDIRS => 10,
  CURLOPT_TIMEOUT => 0,
  CURLOPT_FOLLOWLOCATION => true,
  CURLOPT_HTTP_VERSION => CURL_HTTP_VERSION_1_1,
  CURLOPT_CUSTOMREQUEST => "POST",
  CURLOPT_POSTFIELDS => "username=LAB&password=123456",
  CURLOPT_HTTPHEADER => array(
    "Content-Type: application/x-www-form-urlencoded",
  ),
));

$response = curl_exec($curl);

curl_close($curl);
echo $response;
```

PHP – HTTP_REQUEST2

```
<?php
require_once 'HTTP/Request2.php';
$request = new HTTP_Request2();
$request->setUrl('https://www.SomeServer.com/api/authenticate/logon');
$request->setMethod(HTTP_Request2::METHOD_POST);
$request->setConfig(array(
    'follow_redirects' => TRUE
));
$request->setHeader(array(
    'Content-Type' => 'application/x-www-form-urlencoded',
));
$request->addPostParameter(array(
    'username' => 'LAB',
    'password' => '123456'
));
try {
    $response = $request->send();
    if ($response->getStatus() == 200) {
        echo $response->getBody();
    }
    else {
        echo 'Unexpected HTTP status: ' . $response->getStatus() . ' ' .
            $response->getReasonPhrase();
    }
}
catch(HTTP_Request2_Exception $e) {
    echo 'Error: ' . $e->getMessage();
}
```

PHP – PECL_HTTP

```
<?php
$client = new http\Client;
$request = new http\Client\Request;
$request->setRequestUrl('https://www.SomeServer.com/api/authenticate/logon');
$request->setRequestMethod('POST');
$body = new http\Message\Body;
$body->append(new http\QueryString(array(
    'username' => 'LAB',
    'password' => '123456')));
$request->setBody($body);
$request->setOptions(array());
$request->setHeaders(array(
    'Content-Type' => 'application/x-www-form-urlencoded',
));
$client->enqueue($request)->send();
```

```
$response = $client->getResponse();  
echo $response->getBody();
```

POWERSHELL – RESTMETHOD

```
$headers = New-Object "System.Collections.Generic.Dictionary[[String],[String]]"  
$headers.Add("Content-Type", "application/x-www-form-urlencoded")  
  
$body = "username=LAB&password=123456"  
  
$response = Invoke-RestMethod 'https://www.SomeServer.com/api/authenticate/logon' -  
Method 'POST' -Headers $headers -Body $body
```

PYTHON – HTTP.CLIENT

```
import http.client  
import mimetypes  
conn = http.client.HTTPSConnection("SomeServer.com")  
payload = 'username=LAB&password=123456'  
headers = {  
    'Content-Type': 'application/x-www-form-urlencoded',  
}  
conn.request("POST", "/WebApiBackboneQA/api/authenticate/logon", payload, headers)  
res = conn.getresponse()  
data = res.read()  
print(data.decode("utf-8"))
```

PYTHON – REQUESTS

```
import requests  
  
url = "https://www.SomeServer.com/api/authenticate/logon"  
  
payload='username=LAB&password=123456'  
headers = {  
    'Content-Type': 'application/x-www-form-urlencoded',  
}  
  
response = requests.request("POST", url, headers=headers, data=payload)  
  
print(response.text)
```

*Note: For an example Python integration with the API,
please go to the following link: <https://github.com/marcuslehr/NucleoLIS-API>*

RUBY – NET::HTTP

```
require "uri"
require "net/http"

url = URI("https://www.someserver.com/api/authenticate/logon")

https = Net::HTTP.new(url.host, url.port)
https.use_ssl = true

request = Net::HTTP::Post.new(url)
request["Content-Type"] = "application/x-www-form-urlencoded"
request.body = "username=LAB&password=123456"

response = https.request(request)
puts response.read_body
```

SHELL – HTTPIE

```
http --ignore-stdin --form --follow --
timeout 3600 POST https://www.someserver.com/api/authenticate/logon \
'username='LAB' \
'password='123456' \
Content-Type:'application/x-www-form-urlencoded' \
```

SHELL – WGET

```
wget --no-check-certificate --quiet \
--method POST \
--timeout=0 \
--header 'Content-Type: application/x-www-form-urlencoded' \
--body-data 'username=LAB&password=123456' \
'https://www.someserver.com/api/authenticate/logon'
```

SWIFT - URLSESSION

```
import Foundation

var semaphore = DispatchSemaphore(value: 0)

let parameters = "username=LAB&password=123456"
let postData = parameters.data(using: .utf8)

var request = URLRequest(url: URL(string: "https://www.someserver.com/api/authenticate/logon")
!,timeoutInterval: Double.infinity)
request.addValue("application/x-www-form-urlencoded", forHTTPHeaderField: "Content-Type")
```

```

request.httpMethod = "POST"
request.httpBody = postData

let task = URLSession.shared.dataTask(with: request) { data, response, error in
    guard let data = data else {
        print(String(describing: error))
        semaphore.signal()
        return
    }
    print(String(data: data, encoding: .utf8)!)
    semaphore.signal()
}

task.resume()
semaphore.wait()

```

VB.NET - RESTSHARP

```

Dim client = New
RestClient("https://Www.SomeServer.com/api/authenticate/logon")
client.Timeout = -1
Dim request = New RestRequest(Method.POST)
request.AddHeader("Content-Type", "application/x-www-form-urlencoded")
request.AddParameter("username", "LAB")
request.AddParameter("password", "123456")
Dim response As IRestResponse = client.Execute(request)
Console.WriteLine(response.Content)

```

LOGOUT CODING EXAMPLES

Note: The Logout endpoint is a GET command. Any API call that uses GET would follow a similar pattern.

C - LIBCURL

```

CURL *curl;
CURLcode res;
curl = curl_easy_init();
if(curl) {
    curl_easy_setopt(curl, CURLOPT_CUSTOMREQUEST, "GET");
    curl_easy_setopt(curl, CURLOPT_URL, "https://Www.SomeServer.com/api/authenticate/logout");
    curl_easy_setopt(curl, CURLOPT_FOLLOWLOCATION, 1L);
    curl_easy_setopt(curl, CURLOPT_DEFAULT_PROTOCOL, "https");
    struct curl_slist *headers = NULL;
    headers = curl_slist_append(headers, "Cookie: session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900");
    curl_easy_setopt(curl, CURLOPT_HTTPHEADER, headers);

```

```

    res = curl_easy_perform(curl);
}
curl_easy_cleanup(curl);

```

C# - RESTSHARP

```

var client = new RestClient("https://www.SomeServer.com/api/authenticate/logout");
client.Timeout = -1;
var request = new RestRequest(Method.GET);
request.AddHeader("Cookie", "session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900");
IRestResponse response = client.Execute(request);
Console.WriteLine(response.Content);

```

CURL

```

curl --location --request GET 'https://www.SomeServer.com/api/authenticate/logout' \
--header 'Cookie: session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'

```

GO – NATIVE

```

package main

import (
    "fmt"
    "net/http"
    "io/ioutil"
)

func main() {

    url := "https://www.SomeServer.com/api/authenticate/logout"
    method := "GET"

    client := &http.Client {
    }
    req, err := http.NewRequest(method, url, nil)

    if err != nil {
        fmt.Println(err)
        return
    }
}

```

```

req.Header.Add("Cookie", "session=session-
id=kUn8axg5aDzfBOuVqs2QwF2pZiaI%2fv%2fQSsQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900")

res, err := client.Do(req)
if err != nil {
    fmt.Println(err)
    return
}
defer res.Body.Close()

body, err := ioutil.ReadAll(res.Body)
if err != nil {
    fmt.Println(err)
    return
}
fmt.Println(string(body))
}

```

HTTP

GET /WebApiBackboneQA/api/authenticate/logout HTTP/1.1

Host: SomeServer.com

Cookie: session=session-

id=kUn8axg5aDzfBOuVqs2QwF2pZiaI%2fv%2fQSsQt2mFq000%3d&username=346650&epoch=1604350665&timelimit=900

JAVA – OKHTTP

```

OkHttpClient client = new OkHttpClient().newBuilder()
    .build();
Request request = new Request.Builder()
    .url("https://www.SomeServer.com/api/authenticate/logout")
    .method("GET", null)
    .addHeader("Cookie", "session=session-
id=kUn8axg5aDzfBOuVqs2QwF2pZiaI%2fv%2fQSsQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900")
    .build();
Response response = client.newCall(request).execute();

```

JAVA – UNIREST

```
Unirest.setTimeouts(0, 0);
HttpResponse<String> response = Unirest.get("https://www.SomeServer.com/api/authenticate/logout")
    .header("Cookie", "session=session-id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelimit=900")
    .asString();
```

JAVASCRIPT –FETCH

```
var myHeaders = new Headers();
myHeaders.append("Cookie", "session=session-id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelimit=900");

var requestOptions = {
  method: 'GET',
  headers: myHeaders,
  redirect: 'follow'
};

fetch("https://www.SomeServer.com/api/authenticate/logout", requestOptions)
  .then(response => response.text())
  .then(result => console.log(result))
  .catch(error => console.log('error', error));
```

JAVASCRIPT – JQUERY

```
var settings = {
  "url": "https://www.SomeServer.com/api/authenticate/logout",
  "method": "GET",
  "timeout": 0,
  "headers": {
    "Cookie": "session=session-id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelimit=900"
  },
};

$.ajax(settings).done(function (response) {
  console.log(response);
});
```

JAVASCRIPT – XHR

```
var xhr = new XMLHttpRequest();
xhr.withCredentials = true;

xhr.addEventListener("readystatechange", function() {
  if(this.readyState === 4) {
    console.log(this.responseText);
  }
});

xhr.open("GET", "https://www.SomeServer.com/api/authenticate/logout");
xhr.setRequestHeader("Cookie", "session=session-id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelimit=900");

xhr.send();
```

NODEJS – AXIOS

```
var axios = require('axios');

var config = {
  method: 'get',
  url: 'https://www.SomeServer.com/api/authenticate/logout',
  headers: {
    'Cookie': 'session=session-id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelimit=900'
  }
};

axios(config)
  .then(function (response) {
    console.log(JSON.stringify(response.data));
  })
  .catch(function (error) {
    console.log(error);
  });
```

NODEJS – NATIVE

```
var https = require('follow-redirects').https;
var fs = require('fs');
```

```

var options = {
  'method': 'GET',
  'hostname': 'SomeServer.com',
  'path': '/WebApiBackboneQA/api/authenticate/logout',
  'headers': {
    'Cookie': 'session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQsQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'
  },
  'maxRedirects': 20
};

var req = https.request(options, function (res) {
  var chunks = [];

  res.on("data", function (chunk) {
    chunks.push(chunk);
  });

  res.on("end", function (chunk) {
    var body = Buffer.concat(chunks);
    console.log(body.toString());
  });

  res.on("error", function (error) {
    console.error(error);
  });
});

req.end();

```

NODEJS – REQUEST

```

var request = require('request');
var options = {
  'method': 'GET',
  'url': 'https://www.SomeServer.com/api/authenticate/logout',
  'headers': {
    'Cookie': 'session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQsQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'
  }
};
request(options, function (error, response) {
  if (error) throw new Error(error);
  console.log(response.body);
});

```

NODEJS – UNIREST

```
var unirest = require('unirest');
var req = unirest('GET', 'https://www.SomeServer.com/api/authenticate/logout')
    .headers({
        'Cookie': 'session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSqt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'
    })
    .end(function (res) {
        if (res.error) throw new Error(res.error);
        console.log(res.raw_body);
    });
```

OBJECTIVE-C – NSURLSESSION

```
#import <Foundation/Foundation.h>

dispatch_semaphore_t sema = dispatch_semaphore_create(0);

NSMutableURLRequest *request = [NSMutableURLRequest requestWithURL:[NSURL URLWithString:@"http
s://www.SomeServer.com/api/authenticate/logout"]
    cachePolicy:NSURLRequestUseProtocolCachePolicy
    timeoutInterval:10.0];
NSDictionary *headers = @{
    @"Cookie": @"session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSqt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900"
};

[request setAllHTTPHeaderFields:headers];

[request setHTTPMethod:@"GET"];

NSURLSession *session = [NSURLSession sharedSession];
NSURLSessionDataTask *dataTask = [session dataTaskWithRequest:request
    completionHandler:^(NSData *data, NSURLResponse *response, NSError *error) {
    if (error) {
        NSLog(@"%@", error);
    } else {
        NSHTTPURLResponse *httpResponse = (NSHTTPURLResponse *) response;
        NSError *parseError = nil;
        NSDictionary *responseDictionary = [NSJSONSerialization JSONObjectWithData:data options:0
            error:&parseError];
        NSLog(@"%@", responseDictionary);
        dispatch_semaphore_signal(sema);
    }
}];
```



```

    }
  }];
[dataTask resume];
dispatch_semaphore_wait(sema, DISPATCH_TIME_FOREVER);

```

OCAML – COHTTP

```

open Lwt
open Cohttp
open Cohttp_lwt_unix

let reqBody =
  let uri = Uri.of_string "https://www.SomeServer.com/api/authenticate/logout" in
  let headers = Header.init ()
  |> fun h -> Header.add h "Cookie" "session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQsQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900"
  in
  Client.call ~headers `GET uri >= fun (_resp, body) ->
  body |> Cohttp_lwt.Body.to_string >|= fun body -> body

let () =
  let respBody = Lwt_main.run reqBody in
  print_endline (respBody)

```

PHP – CURL

```

<?php

$curl = curl_init();

curl_setopt_array($curl, array(
  CURLOPT_URL => "https://www.SomeServer.com/api/authenticate/logout",
  CURLOPT_RETURNTRANSFER => true,
  CURLOPT_ENCODING => "",
  CURLOPT_MAXREDIRS => 10,
  CURLOPT_TIMEOUT => 0,
  CURLOPT_FOLLOWLOCATION => true,
  CURLOPT_HTTP_VERSION => CURL_HTTP_VERSION_1_1,
  CURLOPT_CUSTOMREQUEST => "GET",
  CURLOPT_HTTPHEADER => array(
    "Cookie: session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQsQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900"
  ),
));

```

```
$response = curl_exec($curl);
```

```
curl_close($curl);
```

```
echo $response;
```

PHP – HTTP_REQUEST2

```
<?php
require_once 'HTTP/Request2.php';
$request = new HTTP_Request2();
$request->setUrl('https://www.SomeServer.com/api/authenticate/logout');
$request->setMethod(HTTP_Request2::METHOD_GET);
$request->setConfig(array(
    'follow_redirects' => TRUE
));
$request->setHeader(array(
    'Cookie' => 'session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSqt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'
));
try {
    $response = $request->send();
    if ($response->getStatus() == 200) {
        echo $response->getBody();
    }
    else {
        echo 'Unexpected HTTP status: ' . $response->getStatus() . ' ' .
        $response->getReasonPhrase();
    }
}
catch(HTTP_Request2_Exception $e) {
    echo 'Error: ' . $e->getMessage();
}
```

PHP – PECL_HTTP

```
<?php
$client = new http\Client;
$request = new http\Client\Request;
$request->setRequestUrl('https://www.SomeServer.com/api/authenticate/logout');
$request->setRequestMethod('GET');
$request->setOptions(array());
$request->setHeaders(array(
```

```

    'Cookie' => 'session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'
));
$client->enqueue($request)->send();
$response = $client->getResponse();
echo $response->getBody();

```

POWERSHELL – RESTMETHOD

```

$headers = New-Object "System.Collections.Generic.Dictionary[[String],[String]]"
$headers.Add("Cookie", "session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900")

$response = Invoke-RestMethod 'https://www.SomeServer.com/api/authenticate/logout' -
Method 'GET' -Headers $headers -Body $body

```

PYTHON – HTTP.CLIENT

```

import http.client
import mimetypes
conn = http.client.HTTPSConnection("SomeServer.com")
payload = ''
headers = {
    'Cookie': 'session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'
}
conn.request("GET", "/WebApiBackboneQA/api/authenticate/logout", payload, headers)
res = conn.getresponse()
data = res.read()
print(data.decode("utf-8"))

```

PYTHON – REQUESTS

```

import requests

url = "https://www.SomeServer.com/api/authenticate/logout"

payload={}
headers = {
    'Cookie': 'session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'
}

```

```
response = requests.request("GET", url, headers=headers, params=payload)
```

```
print(response.text)
```

RUBY – NET::HTTP

```
require "uri"
```

```
require "net/http"
```

```
url = URI("https://www.SomeServer.com/api/authenticate/logout")
```

```
https = Net::HTTP.new(url.host, url.port)
```

```
https.use_ssl = true
```

```
request = Net::HTTP::Get.new(url)
```

```
request["Cookie"] = "session=session-
```

```
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim  
it=900"
```

```
response = https.request(request)
```

```
puts response.read_body
```

SHELL – HTTPIE

```
http --follow --timeout 3600 GET https://www.SomeServer.com/api/authenticate/logout \
```

```
Cookie:'session=session-
```

```
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim  
it=900'
```

SHELL – WGET

```
wget --no-check-certificate --quiet \
```

```
--method GET \
```

```
--timeout=0 \
```

```
--header 'Cookie: session=session-
```

```
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim  
it=900' \
```

```
'https://www.SomeServer.com/api/authenticate/logout'
```

SWIFT – URLSESSION

```
import Foundation

var semaphore = DispatchSemaphore (value: 0)

var request = URLRequest(url: URL(string: "https://www.SomeServer.com/api/authenticate/logout"
)!, timeoutInterval: Double.infinity)
request.addValue("session=session-
id=kUn8axg5aDzfBOuVqs2QwF2pZiaI%2fv%2fQSsQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900", forHTTPHeaderField: "Cookie")

request.httpMethod = "GET"

let task = URLSession.shared.dataTask(with: request) { data, response, error in
    guard let data = data else {
        print(String(describing: error))
        semaphore.signal()
        return
    }
    print(String(data: data, encoding: .utf8)!)
    semaphore.signal()
}

task.resume()
semaphore.wait()
```

VB.NET – RESTSHARP

```
Dim client = New
RestClient("https://www.SomeServer.com/api/authenticate/logout")
client.Timeout = -1
Dim request = New RestRequest(Method.[GET])
request.AddHeader("Cookie", "session=session-
id=kUn8axg5aDzfBOuVqs2QwF2pZiaI%2fv%2fQSsQt2mFq000%3d&username=346650&epoch=1
604350665&timelimit=900")
Dim response As IRestResponse = client.Execute(request)
Console.WriteLine(response.Content)
```

FOR MORE INFORMATION

Please contact Psyche Systems technical support.