



PSYCHE SYSTEMS

# BACKBONE NUCLEOLIS API

REVISION 1.3

04/26/2021

(REQUIRES BACKBONE SERVER V19 OR LATER)

## TABLE OF CONTENTS

Overview – getting started.....	6
Theory of operation .....	6
SECURITY .....	6
1. Authentication. ....	7
2. Secure Transport Protocol (HTTP-S). ....	7
3. Sliding WINDOW TOKEN EXPIRY .....	7
Getting Started .....	7
Step 1 – Verify/CONFIRM Backbone URL Connectivity .....	7
Step 2 – Logging Into the Backbone NucleoLIS API .....	8
Step 3- Calling API endpoints within a session.....	10
Step 4- TERMINATE A session and log out.....	11
API Calls Specifications .....	12
get heartbeat .....	13
get patient .....	13
get patients.....	15
get case.....	18
get cases .....	35
get specimen.....	39
get specimens .....	41
get test order .....	53
get test orders.....	61
get physician .....	66
get physicians .....	68
set status step .....	73
set status steps .....	74
get physician locations.....	76
get report.....	78
get reports .....	95

get result report .....	98
get case type definitions .....	99
get source definitions.....	105
get profile definitions .....	109
get question definitions .....	113
set instrument test results.....	117
set uploaded text .....	120
add patient.....	121
update patient .....	123
add case .....	126
update case.....	127
add specimen.....	129
update specimen.....	131
add test order.....	133
cancel case .....	135
cancel specimen.....	136
cancel test order .....	137
get plate map definitions .....	139
get plates .....	143
get plate map.....	144
Coding Examples Introduction .....	148
Logon Coding Examples .....	150
C – libcurl .....	150
C# - RestSharp.....	150
CURL.....	151
GO - Native .....	151
HTTP.....	152
Java - OKhttp .....	152
Java – UNIREST.....	152
JAVASCRIPT - FETCH .....	152
JAVASCRIPT - JQUERY .....	153

JAVASCRIPT – XHR .....	153
NODEJS – AXIOS .....	154
NODEJS – NATIVE .....	154
NODEJS – REQUEST .....	155
NODEJS – UNIREST .....	156
OBJECTIVE-C NSURLSESSION .....	156
OCAML – COHTTP .....	157
PHP – CURL .....	157
PHP – HTTP_Request2 .....	158
PHP – pecl_http .....	158
POWERSHELL – RESTMETHOD .....	159
python – http.client .....	159
python – requests .....	159
ruby – net::HTTP .....	160
SHELL – HTTPIE .....	160
SHELL – WGET .....	160
SWIFT - URLSESSION .....	161
VB.NET - RestSHARP .....	161
Logout Coding Examples .....	161
C - libcurl .....	161
C# - RestSharp .....	162
cURL .....	162
Go – native .....	162
HTTP .....	163
Java – OKHTTP .....	163
Java – unirest .....	164
javascript – fetch .....	164
javaScript – jquery .....	164
javascript – xhr .....	165
NodeJS – AXIOS .....	165
nodeJS – NATIVE .....	166

nodejs – request .....	166
nodejs – unirest .....	167
objective-C – NSURLSession .....	167
ocaml – cohttp .....	168
PHP – curl .....	168
php – HTTP_Request2 .....	169
php – pecl_http .....	170
powershell – restmethod .....	170
python – http.client .....	170
python – requests .....	171
ruby – net::HTTP .....	171
shell – httpie .....	171
shell – wget .....	171
swift – URLSession .....	172
vb.net – restsharp .....	172
FOR MORE INFORMATION .....	173

# BACKBONE API

## OVERVIEW – GETTING STARTED

This is a guide for using the Psyche Systems Backbone application programming interface (API).

This section of the guide covers getting started information for four essential requirements for communicating with the Backbone Server using the Backbone API:

1. Verifying/Confirming the Backbone Server Uniform Resource Locator (URL)
2. Logging in, authenticating, and establishing a session.
3. Calling API endpoints during a session.
4. Terminating a session and log out.

Appendix A provides more details about the available API commands that can be used for interacting with NucleoLIS.

Appendix B provides multiple coding examples illustrating how you can program a client side implementation to call the Backbone API.

## THEORY OF OPERATION

The Backbone API is a Representational state transfer (REST) based API that is designed to run on Microsoft Internet Information Server (IIS). Typically, Psyche Systems deploys the server-side of the Backbone on IIS version 7 or later to host the Backbone Server API. The client side implementation of the API can be written in any programming language and run on any operating system as long as the client side implementation adheres to the standards outlined in this document.

To ensure Health Insurance Portability And Accountability Act (HIPAA) privacy compliance and to safeguard Protected Health Information (PHI), all requests to the Backbone Server API must use Hypertext Transfer Protocol Secure (HTTPS) over port 443.

The Backbone Server (being a REST based web server implementation) is considered to be stateless, with the exception of the "session" cookie which is passed back and forth between the client side implementation and the Backbone Server with each request after a successful logon is established.

## SECURITY

To further ensure HIPAA compliance and security of PHI, the Backbone API provides session-based security as follows:

## 1. AUTHENTICATION.

A valid username and password is required to access the Backbone Server API. Once logged in, the API returns an access token that valid for 15 minutes and can be used to make subsequent requests. This access token is stored in a cookie which is transparently passed back and forth with each API request/response. This access token resets every time a successful API call is completed. With this approach, an application can authenticate once, make any number of API calls, and then log out, without needing to pass authentication credentials with every request.

## 2. SECURE TRANSPORT PROTOCOL (HTTP-S).

All API calls are made using Internet standard Hypertext Transport Protocol Secure (HTTPS) on port 443. The bidirectional [encryption](#) of communications between a client and server protects against [eavesdropping](#) and [tampering](#) of the communication. In practice, this provides a reasonable assurance that one is communicating without interference by attackers with the website that one intended to communicate with, as opposed to an impostor.

This implementation also uses the latest [Transport Layer Security](#) (TLS) encryption. This public-private key encryption overlays the symmetric AES encryption used for encrypting the payload, so essentially, the payload is double encrypted and is extremely secure.

## 3. SLIDING WINDOW TOKEN EXPIRY

Each time a session token is provided to a requesting client, the token is encoded so that it is only valid for 15 minutes. This provides what is known as “sliding window” expiration. By using this security method, tokens that are copy/pasted can only be used for a limited period of time. Each time a successful request is made to the Backbone Server API, the token is refreshed to provide a new 15 minute window. In the event that no request is made within 15 minutes, the token expires and can no longer be used without logging in again to obtain a new token. This is the equivalent to the “idle timeout” feature that many online systems offer.

## GETTING STARTED

For each of the examples shown, this document also provides in Appendix “B”, coding examples in a wide variety of programming languages so the implementers can research options for creating their custom client implementation that communicates with the API.

### STEP 1 – VERIFY/CONFIRM BACKBONE URL CONNECTIVITY

When first getting started with a client side implementation, it is helpful to first ensure that you can “reach” the Backbone Server from your client location. To do this, enter the URL of the Backbone Server that you were provided into any web browser. If you are able to connect using that URL, you will see a page that looks like this:

## Welcome to Psyche Systems Web API!

### Access to this API is granted to Authorized users only.

For more information please contact Psyche Systems at: (508) 478-2047.  
Or, via email at: [elixasupport@PsycheSystems.com](mailto:elixasupport@PsycheSystems.com).

If you don't see this page, please contact Psyche Systems support for assistance. If you do see this page, then this is the URL that you will use when making all API requests. This URL will be referred to as "https://www.SomeServer.com/api/" in the rest of this document.

#### STEP 2 – LOGGING INTO THE BACKBONE NUCLEOLIS API

Logging in, authenticates the requesting client as a valid user, and upon successful login, returns an authorization token that is stored in a session cookie. This cookie must be passed after each subsequent request.

This is accomplished by calling the "logon" endpoint with username and password information that is POSTED to the Backbone Server with the following settings:

URI	<a href="https://www.SomeServer.com/api/authenticate/logon">&lt;API URL&gt;/api/authenticate/logon</a>
Host	Provide the hostname of the web server you are connecting with.
Authorization Method	No Auth
Content-Type	application/x-www-form-urlencoded
Content-Length	The length of the body payload in bytes. For example, if you are posting:  <code>username=LAB&amp;password=123456</code>  Then, the Content-Length parameter would be set to 28. (Because there are 28 characters in the payload).



User-Agent	Set this to a unique name for your client application.
Accept	*/*
Posted Data	username=<username>&password=<password>  Here, <username> and <password> are the username and password you were provided to access the API. If you do not have these credentials, please contact Psyche Systems support.

When the called API command succeeds, the response will look like this:

Status	200 OK
Body	<?xml version="1.0" encoding="UTF-8" ?><login>Success</login>
Set-Cookie	session=session-id=qAuOpF9tIBeBR8SIGLmDVRQGUFJ%2ftxLUbELwr8HtNx0%3d&username=346650&epoch=1604331210&timelimit=900; expires=Mon, 02 Nov 2020 15:48:30 GMT; domain=SomeServer.com; path=/; secure; httponly
Content-Type	Text/plain; charset=utf-8
Server	Microsoft-IIS/8.0

If the username or password is invalid, the response will look like this:

Status	401
Body	<?xml version="1.0" encoding="UTF-8" ?><login>username or password invalid</login>
Set-Cookie	(No cookie returned)
Content-Type	Text/plain; charset=utf-8

Other error messages can sometimes be shown as well. If you don't receive the word "Success" within the xml return message, then you should assume that the login failed, and the message returned is the error message.

Note: On initial login, there no cookie specified yet. The cookie gets created by the Backbone Server, and is returned to the client only after a successful logon request is executed. This "session" cookie will contain the authentication token and will look similar to this:

```
Cookie: session=session-  
id=qAuOpF9tIBeBR8SIGLmDVRQGUFJ%2ftxLUbELwr8HtNxo%3d&username=346650&e  
poch=1604331210&timelimit=900; expires=Mon, 02 Nov 2020 15:48:30 GMT;  
domain=SomeServer.com; path=/; secure; httponly
```

This cookie will be updated after each successful request, and must be passed back to the server on each subsequent request. The session token contained within the cookie is only valid for 15 minutes. If a subsequent API endpoint request is not made within the session "timeout" limit of 15 minutes, then the session token becomes invalid, and the client will need to login again to get a new session token by calling the "logon" API endpoint again.

See Appendix "B" for programming examples that call this endpoint.

### STEP 3- CALLING API ENDPOINTS WITHIN A SESSION

Once you are logged in, the next step is to call an API Endpoint using GET. In this example we will show how to call the NucleoLIS api/N/GetHeartbeat API endpoint. This API endpoint returns a heartbeat message so you can determine if the API is able to communicate with NucleoLIS.

URI	<API URL>/ api/N/GetHeartbeat
Host	Provide the hostname of the web server you are connecting with.
Authorization Method	No Auth
Cookie	Pass in the cookie you received in the response from the logon API call. It should look similar to this:

```
session=session-  
id=H38%2fAImLMVMNssLZsY%2fx7mWkfNsHIM4m1CVO%2bFDc9aU%3  
d&username=346650&epoch=1604330631&timelimit=900
```

User-Agent      Set this to a unique name for your client application.

Accept            \*/\*

When the called API command succeeds, the response will look like this:

Status	200 OK
Body	<?xml version="1.0" encoding="UTF-8" ?><Result>OK</Result>
Set-Cookie	session=session- id=xjxAnqMxbu9D%2fbYakV5nRMeKxLlEVnADTiMJqUJycp0%3 d&username=346650&epoch=1604330884&timelimit=900; expires=Mon, 02 Nov 2020 15:43:04 GMT; domain=SomeServer.com; path=/; secure; httponly
Content-Type	Text/plain; charset=utf-8
Server	Microsoft-IIS/8.0

See Appendix "B" for programming examples that call this endpoint.

#### STEP 4- TERMINATE A SESSION AND LOG OUT

When you are done communicating with the Backbone Server API, there are two choices on how to end the session.

1. Do nothing – the cookie containing the latest session token will expire after 15 minutes.
2. Call the Logout API endpoint to destroy the cookie and immediately terminate the session:

URI                <API URL>/api/authenticate/logout

Host               Provide the hostname of the web server you are connecting with.

Authorization Method	No Auth
Cookie	<p>Pass in the cookie you received in the response from the logon API call. It should look similar to this:</p> <pre>session=session-id=H38%2fAImLMVMNssLZsY%2fx7mWkfNsHIM4m1CVO%2bFDc9aU%3d&amp;username=346650&amp;epoch=1604330631&amp;timelimit=900</pre>
User-Agent	Set this to a unique name for your client application.
Accept	*/*

When the called API command succeeds, the response will look like this:

Status	200 OK
Body	<?xml version="1.0" encoding="UTF-8" ?><logout>Success</logout>
Cookie	(The cookie no longer exists)
Content-Type	Text/plain; charset=utf-8

## APPENDIX A – NUCLEOLIS API ENDPOINTS

### API CALLS SPECIFICATIONS

In addition to the logon and logoff API endpoints covered in section 1, the endpoints in this appendix are all of the available API endpoints for NucleoLIS Interaction.

## GET HEARTBEAT

<b>Endpoint</b>	api/N/GetHeartbeat
<b>Purpose</b>	Returns a NucleoLIS heartbeat message so you can determine if the API is able to communicate with NucleoLIS system.
<b>Request Type</b>	GET
<b>Parameters</b>	{None}
<b>Example</b>	<code>https://www.SomeServer.com/api/N/GetHeartbeat</code>
<b>Returns</b>	<p>If the heartbeat succeeds, returns:</p> <pre>&lt;?xml version="1.0" encoding="UTF8" ?&gt;&lt;Result&gt;OK&lt;/Result&gt;</pre>

## GET PATIENT

<b>Endpoint</b>	api/N/GetPatient
<b>Purpose</b>	Returns a single, full NucleoLIS patient record based on a patient unique identifier.
<b>Request Type</b>	GET
<b>Parameters</b>	<pre>patient_id={patient_id} user_id={user_id} return_format={return_format}</pre>

Where:

**patient\_id** is the objectID of the requested Patient record.

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the record.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

**Example** `https://www.SomeServer.com/api/N/GetPatient?patient_id=518535440&user_id=245654&return_format=0`

**Returns** If the request succeeds, returns a patient record in the following XML format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <PatientTable>
    <Patient_ObjectID>518535440</Patient_ObjectID>
    <Patient_CreationDate>2020-11-02T11:23:31.247-
05:00</Patient_CreationDate>
    <Patient_UpdateTime>2020-11-02T11:23:31.247-
05:00</Patient_UpdateTime>
    <Patient_Name>Boop, Betty</Patient_Name>
    <Patient_Sex>F</Patient_Sex>
    <Patient_Street />
    <Patient_City />
    <Patient_State />
    <Patient_Zip />
    <Patient_Phone />
    <Patient_Fax />
    <Patient_Code>518535440</Patient_Code>
    <Patient_Soundex>BOPBTY</Patient_Soundex>
    <Patient_Dob>1965-03-23</Patient_Dob>
    <Patient_Ssn>011-22-3333</Patient_Ssn>
    <Patient_MedicareNumber />
    <Patient_Status>I</Patient_Status>
    <Patient_CREATOR>LAB</Patient_CREATOR>
    <Patient_RACE>White</Patient_RACE>
    <Patient_STREET2 />
```

```

    <Patient_USER1 />
    <Patient_USER2 />
    <Patient_USER3 />
    <Patient_USER4 />
    <Patient_MODIFIER>LAB</Patient_MODIFIER>
    <Patient_LastName>Boop</Patient_LastName>
    <Patient_FirstName>Betty</Patient_FirstName>
    <Patient_MiddleName />
    <Case_Status>I</Case_Status>
  </PatientTable>

</NewDataSet>

```

## GET PATIENTS

<b>Endpoint</b>	api/N/GetPatients
<b>Purpose</b>	Returns a collection of NucleoLIS patient records based on filter parameters passed.
<b>Request Type</b>	GET
<b>Parameters</b>	last_name = {last_name} (use * for wildcard) first_name = {first_name} (use * for wildcard) activesOnly = {activesOnly} (true, or false) filter_expression = {filter_expression} user_id = {user_id} return_format = {return_format}

Where:

**last\_name** is the last name of the patient(s) to search for. Wildcards are supported. (Example S\* will return all patients with a last name beginning with 'S'.)

**first\_name** is the first name of the patient(s) to search for. Wildcards are supported. (Example L\* will return all patients with a first name beginning with 'L'.)

**activesOnly** is a true or false parameter. If true, only active patients will be returned. If false, all patients (both active and inactive) will be returned.

**filter\_expression** is used as a secondary custom filtering criteria applied to the result set . Wildcards are supported. It can be used with any column returned in the result set. Any valid SQL WHERE clause is accepted.

Multiple columns can be used when separated with 'AND'. Example: Patient\_State='MA' AND Patient\_City='W\*' will return all patients that live in MA in a city that begins with 'W'.

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the record.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

#### Example

```
https://www.SomeServer.com/api/N/GetPatients?last_name=S
*&first_name=*&activesOnly=true&filter_expression=Patien
t_Sex='Male'&user_id=245654&return_format=0
```

#### Returns

If the request succeeds, returns a collection of 1..n patient records matching the filtering criteria in the following XML format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <PatientTable>
    <Patient_ObjectID>433016111</Patient_ObjectID>
    <Patient_Name>Claus, Santa</Patient_Name>
    <Patient_Code>433016111</Patient_Code>
```



```

    <Patient_Street />
    <Patient_STREET2 />
    <Patient_City />
    <Patient_State />
    <Patient_Zip />
    <Patient_Sex>Male</Patient_Sex>
    <Patient_Dob>1900-12-25</Patient_Dob>
    <Patient_Ssn />
  </PatientTable>
  <PatientTable>
    <Patient_ObjectID>503609702</Patient_ObjectID>
    <Patient_Name>Claus, Santa</Patient_Name>
    <Patient_Code>503609702</Patient_Code>
    <Patient_Street />
    <Patient_STREET2 />
    <Patient_City />
    <Patient_State />
    <Patient_Zip />
    <Patient_Sex>Male</Patient_Sex>
    <Patient_Dob>1900-12-25</Patient_Dob>
    <Patient_Ssn />
  </PatientTable>
  <PatientTable>
    <Patient_ObjectID>504251809</Patient_ObjectID>
    <Patient_Name>Claus, Santa</Patient_Name>
    <Patient_Code>504251809</Patient_Code>
    <Patient_Street />
    <Patient_STREET2 />
    <Patient_City />
    <Patient_State />
    <Patient_Zip />
    <Patient_Sex>Male</Patient_Sex>
    <Patient_Dob>1900-12-25</Patient_Dob>
    <Patient_Ssn />
  </PatientTable>
  <PatientTable>
    <Patient_ObjectID>4550996</Patient_ObjectID>
    <Patient_Name>Sanding, Sammie</Patient_Name>
    <Patient_Code>4550996</Patient_Code>
    <Patient_Street />
    <Patient_STREET2 />
    <Patient_City />
    <Patient_State />
    <Patient_Zip />
    <Patient_Sex>Male</Patient_Sex>
    <Patient_Dob>1980-04-22</Patient_Dob>
    <Patient_Ssn />
  </PatientTable>

```

```

</PatientTable>
<PatientTable>
  <Patient_ObjectID>507538376</Patient_ObjectID>
  <Patient_Name>Sanding, Sammie</Patient_Name>
  <Patient_Code>507538376</Patient_Code>
  <Patient_Street />
  <Patient_STREET2 />
  <Patient_City />
  <Patient_State />
  <Patient_Zip />
  <Patient_Sex>Male</Patient_Sex>
  <Patient_Dob>1980-04-22</Patient_Dob>
  <Patient_Ssn />
</PatientTable>

</NewDataSet>

```

## GET CASE

<b>Endpoint</b>	api/N/GetCase
<b>Purpose</b>	Return a single, full NucleoLIS case record based on case number. The full case record inclusive of all data elements associated with patient, case, specimen, test order and test result records are returned.
<b>Request Type</b>	GET
<b>Parameters</b>	case_number={case_number} user_id={user_id} return_format={return_format}

Where:

**case\_number** is the number of the case.

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the record.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

### Example

```
https:// www.SomeServer.com  
/api/N/GetCase?case_number=M2020-  
000032&user_id=245654&return_format=0
```

### Returns

If the request succeeds, returns the requested case record in the following XML format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>  
  <CaseTable>  
    <Patient_ObjectID>163307929</Patient_ObjectID>  
    <Patient_CreationDate>2016-11-10T16:11:29.84-  
05:00</Patient_CreationDate>  
    <Patient_UpdateTime>2019-05-22T10:51:15.137-  
04:00</Patient_UpdateTime>  
    <Patient_Name>O'Hara, Jaime</Patient_Name>  
    <Patient_Sex>Female</Patient_Sex>  
    <Patient_Street />  
    <Patient_City />  
    <Patient_State />  
    <Patient_Zip />  
    <Patient_Phone />  
    <Patient_Fax />  
    <Patient_Code>163307929</Patient_Code>  
    <Patient_Soundex>O'HRJM</Patient_Soundex>  
    <Patient_Dob>2018-05-16</Patient_Dob>  
    <Patient_Ssn />  
    <Patient_MedicareNumber />  
    <Patient_Status>I</Patient_Status>  
    <Patient_CREATOR>LAB</Patient_CREATOR>  
    <Patient_RACE />  
    <Patient_STREET2 />  
    <Patient_USER1 />  
    <Patient_USER2 />  
    <Patient_USER3 />  
    <Patient_USER4 />  
    <Patient_Modifier>Scaron</Patient_Modifier>
```

```

    <Patient_LastName>O'Hara</Patient_LastName>
    <Patient_FirstName>Jaime</Patient_FirstName>
    <Visit_ObjectID>498604164</Visit_ObjectID>
    <Visit_CreationDate>2019-05-22T10:51:34.293-
04:00</Visit_CreationDate>
    <Visit_UpdateTime>2019-05-22T10:51:34.293-
04:00</Visit_UpdateTime>

<Visit_SuperobjectID>163307929</Visit_SuperobjectID>
    <Visit_Patient>163307929</Visit_Patient>
    <Visit_Code>498604164</Visit_Code>
    <Visit_Location />
    <Visit_Type>OP</Visit_Type>
    <Visit_VisitDate>2019-05-22</Visit_VisitDate>
    <Visit_Status>I</Visit_Status>
    <Visit_DischargeDate />
    <Visit_ExtraField1 />
    <Visit_ExtraField2 />
    <Visit_CREATOR>Scaron</Visit_CREATOR>
    <Visit_USER1 />
    <Visit_USER2 />
    <Visit_USER3 />
    <Visit_USER4 />
    <Visit_comment />
    <Case_ObjectID>506873424</Case_ObjectID>
    <Case_CreationDate>2020-03-09T14:50:06.023-
04:00</Case_CreationDate>
    <Case_UpdateTime>2020-05-11T14:03:15.34-
04:00</Case_UpdateTime>
    <Case_SuperobjectID>498604164</Case_SuperobjectID>
    <Case_User1 />
    <Case_User2 />
    <Case_User3 />
    <Case_User4 />
    <Case_CREATOR>Scaron</Case_CREATOR>
    <Case_MODIFIER>Bishop</Case_MODIFIER>
    <Case_Number>M2020-000032</Case_Number>
    <Case_ReqNumber />
    <Case_ReportedFlag>False</Case_ReportedFlag>

<Case_SubmitterLocation>258012099</Case_SubmitterLocat
ion>
    <Case_OrderID />
    <Case_Status>I</Case_Status>
    <Case_ReportStatus>FINAL</Case_ReportStatus>
    <Case_CaseType>14662557</Case_CaseType>
    <Case_CorrectedComment />

```

```

<Case_AmendedComment />
<Case_ReportedBy>Bishop</Case_ReportedBy>
<Case_CanceledBy />
<Case_CanceledDate />
<Case_ReviewFlag>2020-03-09</Case_ReviewFlag>
<Case_ReviewedBy />
<Case_ReviewedDate />
<Case_ReviewedComment />
<Case_ReviewedTime />
<Case_StatusStep />
<Case_CanceledFor />
<Case_Signer>370583063</Case_Signer>
<Case_SignerDate>2020-03-09</Case_SignerDate>
<Case_ICDs />
<Case_ClinicalInformation />
<Case_Images><?xml version="1.0" encoding="utf-
16"?>
<ChrysalisImgDats
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema"
/></Case_Images>
  <Case_DocumentID />
  <Case_OrderDate>2020-03-09</Case_OrderDate>
  <Case_OrderTime />
  <Case_ExternalID />
  <Case_ExternalPatientID />
  <Case_ReopenReason />
  <Case_ReopenDate />
  <Case_Signatures><?xml version="1.0"
encoding="utf-16"?>
<CaseSignatureHistory
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema">
  <Sigs>
    <CaseSignature>
      <Signer>370583063</Signer>
      <SignerDate>2020-03-09T15:08:36</SignerDate>
      <ReOpenReason />
      <ReOpener>0</ReOpener>
      <ReOpenedDate>0001-01-01T00:00:00</ReOpenedDate>
      <UnSign>false</UnSign>
    </CaseSignature>
    <CaseSignature>
      <Signer>370583063</Signer>
      <SignerDate>2020-03-09T15:10:42</SignerDate>
      <ReOpenReason />
      <ReOpener>0</ReOpener>

```

```

        <ReOpenedDate>0001-01-01T00:00:00</ReOpenedDate>
        <UnSign>true</UnSign>
    </CaseSignature>
    <CaseSignature>
        <Signer>370583063</Signer>
        <SignerDate>2020-03-09T15:11:09</SignerDate>
        <ReOpenReason />
        <ReOpener>0</ReOpener>
        <ReOpenedDate>0001-01-01T00:00:00</ReOpenedDate>
        <UnSign>false</UnSign>
    </CaseSignature>
</Sigs>
</CaseSignatureHistory></Case_Signatures>
    <Case_NeedsCorrectedF>False</Case_NeedsCorrectedF>
    <Case_PrescriptionDrugs><?xml version="1.0"
encoding="utf-16"?>
<PrescribedDrugs
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema"
/></Case_PrescriptionDrugs>
    <Case_DataCollect />
    <Case_ReportIX>1</Case_ReportIX>
    <Case_OInitStatus>I</Case_OInitStatus>

<Case_RequireSignature>True</Case_RequireSignature>
    <Case_SignerTime>15:11</Case_SignerTime>
    <Case_CorrectionRequest />
    <Case_CorrectionRequestDate />
    <Case_CorrectionRequestTime />

<Case_SignerLocation>246584261</Case_SignerLocation>
    <Case_ClearCorrectedReason />
    <Case_QuestionsStash><?xml version="1.0"
encoding="utf-16"?>
<AAOEstasherBase
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema">
    <Responses>
        <AAOEItemBase>
            <QObjectID>302726734</QObjectID>
            <Question>Additional Information</Question>
            <AllowOther>true</AllowOther>
            <SingleSelect>false</SingleSelect>
            <Required>true</Required>
            <FullResponse />
            <GridQuestion>false</GridQuestion>
            <GridQuestions />

```

```

    <Answers>
      <string>none</string>
    </Answers>
  </AAOEItemBase>
  <AAOEItemBase>
    <QObjectID>481211124</QObjectID>
    <Question>General Health</Question>
    <AllowOther>true</AllowOther>
    <SingleSelect>true</SingleSelect>
    <Required>false</Required>
    <FullResponse />
    <GridQuestion>false</GridQuestion>
    <GridQuestions />
    <Answers>
      <string>Good</string>
    </Answers>
  </AAOEItemBase>
  <AAOEItemBase>
    <QObjectID>481215687</QObjectID>
    <Question>Fasting</Question>
    <AllowOther>false</AllowOther>
    <SingleSelect>true</SingleSelect>
    <Required>true</Required>
    <FullResponse />
    <GridQuestion>false</GridQuestion>
    <GridQuestions />
    <Answers>
      <string>Yes</string>
    </Answers>
  </AAOEItemBase>
</Responses>
</AAOEstasherBase></Case_QuestionsStash>
  <Case_Questions>Additional Information:none
  General Health:Good
  Fasting:Yes</Case_Questions>

<Case_QuestionsComplete>True</Case_QuestionsComplete>
  <Specimen_ObjectID>506873435</Specimen_ObjectID>
  <Specimen_CreationDate>2020-03-09T14:50:41.513-
04:00</Specimen_CreationDate>
  <Specimen_UpdateTime>2020-03-09T14:51:10.153-
04:00</Specimen_UpdateTime>

<Specimen_SuperobjectID>506873424</Specimen_Superobjec
tID>
  <Specimen_User1 />
  <Specimen_User2 />

```

```

    <Specimen_User3 />
    <Specimen_User4 />
    <Specimen_CREATOR>Scaron</Specimen_CREATOR>
    <Specimen_MODIFIER>Scaron</Specimen_MODIFIER>
    <Specimen_CollectionDate>2020-03-
09</Specimen_CollectionDate>

    <Specimen_CollectionTime>09:00</Specimen_CollectionTime>

    <Specimen_ReceivedDate>2020-03-
09</Specimen_ReceivedDate>

    <Specimen_ReceivedTime>14:51:10</Specimen_ReceivedTime>
    >
        <Specimen_Receiver>Scaron</Specimen_Receiver>
        <Specimen_Number>506873435</Specimen_Number>
        <Specimen_Status>A</Specimen_Status>
        <Specimen_Comment />
        <Specimen_Description />

    <Specimen_SourceDefinition>4348467</Specimen_SourceDefinition>
    <Specimen_Interpretation />
    <Specimen_StatusStep />
    <Specimen_SuperOrders />
    <Specimen_UnInlabDate />
    <Specimen_DataCollect />
    <Specimen_Priority>Routine</Specimen_Priority>
    <Specimen_Site />
    <Specimen_Gross />
    <Specimen_AddOnDate />
    <Specimen_AddOnTime />
    <Specimen_AddOn>False</Specimen_AddOn>
    <Specimen_MultipleCountID />
    <Specimen_MergeComments />
    <Specimen_MergeIntoComments />
    <Specimen_QuestionsStash />
    <Specimen_Questions />

    <Specimen_QuestionsComplete>True</Specimen_QuestionsComplete>
    <Specimen_FastingStatus />
    <Profile_TYPE>IMMUNO</Profile_TYPE>
    <Profile_ObjectID>506873438</Profile_ObjectID>
    <Profile_CreationDate>2020-03-09T14:50:49.687-
04:00</Profile_CreationDate>

```



```

    <Profile_UpdateTime>2020-03-09T15:07:50.84-
04:00</Profile_UpdateTime>

<Profile_SuperobjectID>506873435</Profile_SuperobjectI
D>
    <Profile_User1 />
    <Profile_User2 />
    <Profile_User3 />
    <Profile_User4 />
    <Profile_CREATOR>Scaron</Profile_CREATOR>
    <Profile_MODIFIER>Suec</Profile_MODIFIER>
    <Profile_Description />
    <Profile_OrderDate>2020-03-09</Profile_OrderDate>
    <Profile_OrderTime>14:50</Profile_OrderTime>
    <Profile_OrderedBy>Scaron</Profile_OrderedBy>

<Profile_PanelDefinition>416976655</Profile_PanelDefin
ition>
    <Profile_R1 />
    <Profile_R2 />
    <Profile_R3 />
    <Profile_R4 />
    <Profile_R5 />
    <Profile_R6 />
    <Profile_Result />
    <Profile_Reportable>False</Profile_Reportable>

<Profile_DatumEntryDefinition>478882298</Profile_Datum
EntryDefinition>
    <Profile_Images><?xml version="1.0" encoding="utf-
16"?>
<ChrysalisImgDats
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema"
/></Profile_Images>
    <Profile_Abnormal>False</Profile_Abnormal>
    <Profile_AOV>False</Profile_AOV>
    <Profile_StatusStep />
    <Profile_Range />
    <Profile_RangeComment />
    <Profile_CanceledDate />
    <Profile_CanceledBy />
    <Profile_CanceledFor />
    <Profile_ReflexedBy />
    <Profile_Signer>270861542</Profile_Signer>
    <Profile_SignerDate>2020-03-
09</Profile_SignerDate>

```

```

        <Profile_IOrder />
        <Profile_IOrderDate />
        <Profile_IOrderTime />
        <Profile_Complete>True</Profile_Complete>
        <Profile_DataCollect />
        <Profile_Signatures><?xml version="1.0"
encoding="utf-16"?>
<ProfileSignatureHistory
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema">
    <Sigs>
        <ProfileSignature>
            <Signer>270861542</Signer>
            <SignerDate>2020-03-09T15:07:50</SignerDate>
            <UnSign>false</UnSign>
        </ProfileSignature>
    </Sigs>
</ProfileSignatureHistory></Profile_Signatures>

<Profile_RequireSignature>True</Profile_RequireSignatu
re>
    <Profile_Approver>162355184</Profile_Approver>
    <Profile_ApproveDate>2020-03-
09</Profile_ApproveDate>
    <Profile_Approvals><?xml version="1.0"
encoding="utf-16"?>
<ProfileApprovalHistory
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema">
    <Sigs>
        <ProfileApproval>
            <Signer>162355184</Signer>
            <SignerDate>2020-03-09T14:56:31</SignerDate>
            <UnSign>false</UnSign>
        </ProfileApproval>
    </Sigs>
</ProfileApprovalHistory></Profile_Approvals>

<Profile_RequireApproval>True</Profile_RequireApproval
>
    <Profile_SignerTime>15:07</Profile_SignerTime>
    <Profile_ApproveTime>14:56</Profile_ApproveTime>
    <Profile_WPDLink />
    <Profile_CompleteDate />
    <Profile_CompleteTime />

```

```

<Profile_SignerLocation>246584261</Profile_SignerLocation>

<Profile_PendingReagent>False</Profile_PendingReagent>
  <Profile_ORDLink />
  <Profile_QuestionsStash><?xml version="1.0"
encoding="utf-16"?>
<AAOEstasherBase
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema"
/></Profile_QuestionsStash>
  <Profile_Questions />

<Profile_QuestionsComplete>True</Profile_QuestionsComplete>
  <Profile_ROrder>I</Profile_ROrder>
  <Profile_ROrderDate>2020-03-
09</Profile_ROrderDate>
  <Profile_ROrderTime>14:50</Profile_ROrderTime>
  <Profile_ROrderLab>111200992</Profile_ROrderLab>
  <Profile_Interpretation />

<Profile_ProfileDefinitionObjectID>416976655</Profile_ProfileDefinitionObjectID>
  <Profile_ProfileDefinitionCreationDate>2018-10-
31T10:16:34.31-
04:00</Profile_ProfileDefinitionCreationDate>
  <Profile_ProfileDefinitionUpdateTime>2020-03-
09T14:48:48.757-
04:00</Profile_ProfileDefinitionUpdateTime>
  <Profile_ProfileDefinitionUser1 />
  <Profile_ProfileDefinitionUser2 />
  <Profile_ProfileDefinitionUser3 />
  <Profile_ProfileDefinitionUser4 />

<Profile_ProfileDefinitionCREATOR>debbieb</Profile_ProfileDefinitionCREATOR>

<Profile_ProfileDefinitionMODIFIER>Scaron</Profile_ProfileDefinitionMODIFIER>

<Profile_ProfileDefinitionCode>AMPH</Profile_ProfileDefinitionCode>

<Profile_ProfileDefinitionDescription>Amphetamine</Profile_ProfileDefinitionDescription>

```

```

<Profile_ProfileDefinitionRETIRED>False</Profile_ProfileDefinitionRETIRED>

<Profile_ProfileDefinitionDatumEntryDefinition>4788822
98</Profile_ProfileDefinitionDatumEntryDefinition>
  <Profile_ProfileDefinitionAlwaysReportMessage />
  <Profile_ProfileDefinitionMethodology />

<Profile_ProfileDefinitionRequireSignature>True</Profile_ProfileDefinitionRequireSignature>
  <Profile_ProfileDefinitionCodeMap />

<Profile_ProfileDefinitionLabelCount>1</Profile_ProfileDefinitionLabelCount>

<Profile_ProfileDefinitionAutoComplete>False</Profile_ProfileDefinitionAutoComplete>
  <Profile_ProfileDefinitionConstituentPanels><?xml
version="1.0" encoding="utf-16"?>
<GPathConstituentPanels
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema"
/></Profile_ProfileDefinitionConstituentPanels>

<Profile_ProfileDefinitionAllowDuplicate>False</Profile_ProfileDefinitionAllowDuplicate>
  <Profile_ProfileDefinitionPendingMessage />

<Profile_ProfileDefinitionWPLink>False</Profile_ProfileDefinitionWPLink>

<Profile_ProfileDefinitionWPDiscrete>False</Profile_ProfileDefinitionWPDiscrete>

<Profile_ProfileDefinitionWPCyto>False</Profile_ProfileDefinitionWPCyto>

<Profile_ProfileDefinitionApprovalSignature>True</Profile_ProfileDefinitionApprovalSignature>

<Profile_ProfileDefinitionLabelDefinition>318256</Profile_ProfileDefinitionLabelDefinition>

<Profile_ProfileDefinitionWPDiscreteAdvanceOnly>False</Profile_ProfileDefinitionWPDiscreteAdvanceOnly>

```

```

<Profile_ProfileDefinitionRequireLocation>False</Profile_ProfileDefinitionRequireLocation>

<Profile_ProfileDefinitionUseDataCSort>False</Profile_ProfileDefinitionUseDataCSort>

<Profile_ProfileDefinitionAllowBatchResult>True</Profile_ProfileDefinitionAllowBatchResult>

<Profile_ProfileDefinitionIsToxicology>True</Profile_ProfileDefinitionIsToxicology>

<Profile_ProfileDefinitionIsClinical>False</Profile_ProfileDefinitionIsClinical>

<Profile_ProfileDefinitionSuppressOverlap>False</Profile_ProfileDefinitionSuppressOverlap>

<Profile_ProfileDefinitionOutreachCode>AMPH</Profile_ProfileDefinitionOutreachCode>

<Profile_ProfileDefinitionOutreachDescription>Amphetamine</Profile_ProfileDefinitionOutreachDescription>

<Profile_ProfileDefinitionAllowReferenceLab>True</Profile_ProfileDefinitionAllowReferenceLab>

<SourceDefinition_ObjectID>4348467</SourceDefinition_ObjectID>
  <SourceDefinition_CreationDate>2015-04-23T09:31:01.957-04:00</SourceDefinition_CreationDate>
  <SourceDefinition_UpdateTime>2020-09-18T16:59:13.113-04:00</SourceDefinition_UpdateTime>
  <SourceDefinition_User1 />
  <SourceDefinition_User2 />
  <SourceDefinition_User3 />
  <SourceDefinition_User4 />

<SourceDefinition_CREATOR>deannac</SourceDefinition_CREATOR>

<SourceDefinition_MODIFIER>Scaron</SourceDefinition_MODIFIER>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

```

<SourceDefinition\_Description>Blood</SourceDefinition\_Description>

<SourceDefinition\_RETIRED>False</SourceDefinition\_RETIRED>

<SourceDefinition\_LabelCount>1</SourceDefinition\_LabelCount>

<SourceDefinition\_WPLink>False</SourceDefinition\_WPLink>

<SourceDefinition\_WPRequireProfiles>False</SourceDefinition\_WPRequireProfiles>

<SourceDefinition\_WPAllowNOrdering>False</SourceDefinition\_WPAllowNOrdering>

<SourceDefinition\_ForeignAllowNOrdering>False</SourceDefinition\_ForeignAllowNOrdering>

<SourceDefinition\_WPAddDiagnosis>False</SourceDefinition\_WPAddDiagnosis>

<SourceDefinition\_LabelDefinition>280272780</SourceDefinition\_LabelDefinition>

<SourceDefinition\_Sites><?xml version="1.0" encoding="utf-16"?>  
<XStrings xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance" xmlns:xsd="http://www.w3.org/2001/XMLSchema" /></SourceDefinition\_Sites>

<SourceDefinition\_OverdueDefault>8</SourceDefinition\_OverdueDefault>

<SourceDefinition\_Multiples>2</SourceDefinition\_Multiples>

<SourceDefinition\_AutoMultiples>False</SourceDefinition\_AutoMultiples>

<SourceDefinition\_GrossReceived>False</SourceDefinition\_GrossReceived>

<SourceDefinition\_AAOESet>485355213</SourceDefinition\_AAOESet>

```

<SourceDefinition_IgnoreLabelDefault>False</SourceDefinition_IgnoreLabelDefault>

<CaseTypeDefinition_ObjectID>14662557</CaseTypeDefinition_ObjectID>
  <CaseTypeDefinition_CreationDate>2015-06-25T11:10:41.087-04:00</CaseTypeDefinition_CreationDate>
  <CaseTypeDefinition_UpdateTime>2020-10-27T11:57:35.247-04:00</CaseTypeDefinition_UpdateTime>
  <CaseTypeDefinition_User1 />
  <CaseTypeDefinition_User2 />
  <CaseTypeDefinition_User3 />
  <CaseTypeDefinition_User4 />

<CaseTypeDefinition_CREATOR>0</CaseTypeDefinition_CREATOR>

<CaseTypeDefinition_MODIFIER>Scaron</CaseTypeDefinition_MODIFIER>

<CaseTypeDefinition_Code>MOLC</CaseTypeDefinition_Code>

<CaseTypeDefinition_Description>Molecular</CaseTypeDefinition_Description>

<CaseTypeDefinition_RETIRED>False</CaseTypeDefinition_RETIRED>

<CaseTypeDefinition_ReportDefinition>217397414</CaseTypeDefinition_ReportDefinition>

<CaseTypeDefinition_OverdueDays>1</CaseTypeDefinition_OverdueDays>

<CaseTypeDefinition_RequirePreliminary>False</CaseTypeDefinition_RequirePreliminary>

<CaseTypeDefinition_StepDefinition>245594</CaseTypeDefinition_StepDefinition>

<CaseTypeDefinition_RequireSignature>False</CaseTypeDefinition_RequireSignature>

```

```

<CaseTypeDefinition_AccessionLetter>M</CaseTypeDefinition_AccessionLetter>

<CaseTypeDefinition_UseProfileReports>False</CaseTypeDefinition_UseProfileReports>

<CaseTypeDefinition_PrelimReportDefinition>217397414</CaseTypeDefinition_PrelimReportDefinition>

<CaseTypeDefinition_WPLink>False</CaseTypeDefinition_WPLink>

<CaseTypeDefinition_WPSignOnReport>False</CaseTypeDefinition_WPSignOnReport>

<CaseTypeDefinition_WPLoadReport>False</CaseTypeDefinition_WPLoadReport>

<CaseTypeDefinition_WPAllowSupplemental>False</CaseTypeDefinition_WPAllowSupplemental>

<CaseTypeDefinition_LabelDefinition>348382382</CaseTypeDefinition_LabelDefinition>
  <CaseTypeDefinition_WPSignSuppress><?xml
version="1.0" encoding="utf-16"?>
<XStrings xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance" xmlns:xsd="http://www.w3.org/2001/XMLSchema"
/></CaseTypeDefinition_WPSignSuppress>

<CaseTypeDefinition_WPSignSuppressSupplemental><?xml
version="1.0" encoding="utf-16"?>
<XStrings xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance" xmlns:xsd="http://www.w3.org/2001/XMLSchema"
/></CaseTypeDefinition_WPSignSuppressSupplemental>

<CaseTypeDefinition_NeverPreliminary>False</CaseTypeDefinition_NeverPreliminary>

<CaseTypeDefinition_AbnormalSpecialFlag>348386126</CaseTypeDefinition_AbnormalSpecialFlag>

<CaseTypeDefinition_WPSetSpecialFlag>True</CaseTypeDefinition_WPSetSpecialFlag>

<CaseTypeDefinition_IgnoreLabelDefault>False</CaseTypeDefinition_IgnoreLabelDefault>

```



```

    <Constituent_TYPE>IMMUNO</Constituent_TYPE>

    <Constituent_ObjectID>506873439</Constituent_ObjectID>
    <Constituent_CreationDate>2020-03-09T14:50:50.123-
04:00</Constituent_CreationDate>
    <Constituent_UpdateTime>2020-03-09T14:52:49.107-
04:00</Constituent_UpdateTime>

    <Constituent_SuperobjectID>506873438</Constituent_Supe
robjectID>
    <Constituent_User1 />
    <Constituent_User2 />
    <Constituent_User3 />
    <Constituent_User4 />
    <Constituent_CREATOR>Scaron</Constituent_CREATOR>

    <Constituent_MODIFIER>Scaron</Constituent_MODIFIER>

    <Constituent_Reportable>True</Constituent_Reportable>
    <Constituent_Description />

    <Constituent_Definition>414331635</Constituent_Definit
ion>
    <Constituent_Result>Negative</Constituent_Result>
    <Constituent_R1 />
    <Constituent_R2 />
    <Constituent_R3 />
    <Constituent_R4 />
    <Constituent_R5 />
    <Constituent_R6 />
    <Constituent_Abnormal>False</Constituent_Abnormal>
    <Constituent_AOV>False</Constituent_AOV>
    <Constituent_Range />
    <Constituent_RangeComment />

    <Constituent_IsConfirm>False</Constituent_IsConfirm>

    <Constituent_ConstituentDefinitionObjectID>414331635</
Constituent_ConstituentDefinitionObjectID>

    <Constituent_ConstituentDefinitionCreationDate>2018-
10-26T11:55:05.557-
04:00</Constituent_ConstituentDefinitionCreationDate>
    <Constituent_ConstituentDefinitionUpdateTime>2020-
11-03T10:33:26.493-
05:00</Constituent_ConstituentDefinitionUpdateTime>

```

```

<Constituent_ConstituentDefinitionCREATOR>Import</Constituent_ConstituentDefinitionCREATOR>

<Constituent_ConstituentDefinitionMODIFIER>Scaron</Constituent_ConstituentDefinitionMODIFIER>

<Constituent_ConstituentDefinitionCode>AMPH</Constituent_ConstituentDefinitionCode>

<Constituent_ConstituentDefinitionDescription>Amphetamine</Constituent_ConstituentDefinitionDescription>

<Constituent_ConstituentDefinitionRETIRED>False</Constituent_ConstituentDefinitionRETIRED>
    <Constituent_ConstituentDefinitionCutOffValue>250 ng/dl</Constituent_ConstituentDefinitionCutOffValue>

<Constituent_ConstituentDefinitionDetectedString>DETECTED</Constituent_ConstituentDefinitionDetectedString>

<Constituent_ConstituentDefinitionUndetectedString>NOT DETECTED</Constituent_ConstituentDefinitionUndetectedString>

<Constituent_ConstituentDefinitionIllicit>False</Constituent_ConstituentDefinitionIllicit>

<Constituent_ConstituentDefinitionDetectionWindow>1-2 days</Constituent_ConstituentDefinitionDetectionWindow>

<Constituent_ConstituentDefinitionOrderConfirmationWhen>Detected</Constituent_ConstituentDefinitionOrderConfirmationWhen>

<Constituent_ConstituentDefinitionConfirmation>478867597</Constituent_ConstituentDefinitionConfirmation>

<Constituent_ConstituentDefinitionConfirmationCharge>259621897</Constituent_ConstituentDefinitionConfirmationCharge>

<Constituent_ConstituentDefinitionIsTextual>False</Constituent_ConstituentDefinitionIsTextual>
</CaseTable>

```

</NewDataSet>

## GET CASES

<b>Endpoint</b>	api/N/GetCases
<b>Purpose</b>	Returns a collection of NucleoLIS case records based on filter parameters passed.
<b>Request Type</b>	GET
<b>Parameters</b>	<p>case_number={case_number}</p> <p>activesOnly={activesOnly}</p> <p>order_date_from={order_date_from}</p> <p>order_date_to={order_date_to}</p> <p>submitting_physicians={submitting_physicians}</p> <p>submitting_groups={submitting_groups}</p> <p>filter_expression={filter_expression}</p> <p>user_ID={user_ID}</p> <p>return_format={return_format}</p> <p>Where:</p> <p><b>case_number</b> is the requested NucleoLIS case_number. (Wildcards are supported)</p> <p><b>activesOnly</b> is a true or false value. If true, only active cases will be returned. If false, all cases will be returned.</p> <p><b>order_date_from</b> is the order from date for filtering cases. Dates must be in the format YYYY-MM-DD.</p>

**order\_date\_to** is the order to date for filtering cases. Dates must be in the format YYYY-MM-DD.

**submitting\_physicians** is a list of submitting physicians.

**submitting\_groups** is a list of submitting groups.

**filter\_expression** is used as a secondary custom filtering criteria applied to the result set . Wildcards are supported. Any valid SQL WHERE clause is accepted. It can be used with any column returned in the result set. Multiple columns can be used when separated with 'AND'. Filtering criteria supported include: Patient\_ObjectID, Patient\_Name, Patient\_Code, Patient\_Dob, Patient\_Sex, Case\_Number, Case\_OrderDate, Case\_CreationDate, Case\_ReportStatus, Case\_Status, Case\_ObjectID.

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the record.

**return\_format** is the format for returning the data. Where:

0 = XML Format

1 = JSON Format

#### Example

```
https://SomeServer.com/api/N/GetCases?case_number=CL2020*&activesOnly=true&return_format=0&user_id=245654&order_date_from=2020-09-01&order_date_to=2020-11-02&filter_expression=case_status='A' AND patient_sex='Female'&submitting_physicians=*&submitting_groups=*
```

#### Returns

If the request succeeds, returns a collection of 1..n case records matching the filtering criteria in the following XML format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
<CaseTable>
  <Patient_ObjectID>506408213</Patient_ObjectID>
  <Patient_Name>Boop, Betty</Patient_Name>
  <Patient_Code>506408213</Patient_Code>
  <Patient_Dob>1965-03-23</Patient_Dob>
```

```

<Patient_Sex>Female</Patient_Sex>
<Case_Number>CL2020-000080</Case_Number>
<Case_OrderDate>2020-09-17</Case_OrderDate>
<Case_CreationDate>2020-09-17T17:27:27.473-04:00</Case_CreationDate>
<Case_ReportStatus>PENDING</Case_ReportStatus>
<Case_Status>A</Case_Status>
<Case_ObjectID>514638163</Case_ObjectID>
</CaseTable>
<CaseTable>
  <Patient_ObjectID>503799913</Patient_ObjectID>
  <Patient_Name>Bell, Tinker</Patient_Name>
  <Patient_Code>503799913</Patient_Code>
  <Patient_Dob>1950-01-01</Patient_Dob>
  <Patient_Sex>Female</Patient_Sex>
  <Case_Number>CL2020-000082</Case_Number>
  <Case_OrderDate>2020-09-18</Case_OrderDate>
  <Case_CreationDate>2020-09-18T14:53:08.89-04:00</Case_CreationDate>
  <Case_ReportStatus>PENDING</Case_ReportStatus>
  <Case_Status>A</Case_Status>
  <Case_ObjectID>514704625</Case_ObjectID>
</CaseTable>
<CaseTable>
  <Patient_ObjectID>503799913</Patient_ObjectID>
  <Patient_Name>Bell, Tinker</Patient_Name>
  <Patient_Code>503799913</Patient_Code>
  <Patient_Dob>1950-01-01</Patient_Dob>
  <Patient_Sex>Female</Patient_Sex>
  <Case_Number>CL2020-000084</Case_Number>
  <Case_OrderDate>2020-09-28</Case_OrderDate>
  <Case_CreationDate>2020-09-28T13:43:54.22-04:00</Case_CreationDate>
  <Case_ReportStatus>PENDING</Case_ReportStatus>
  <Case_Status>A</Case_Status>
  <Case_ObjectID>515443274</Case_ObjectID>
</CaseTable>
<CaseTable>
  <Patient_ObjectID>503625913</Patient_ObjectID>
  <Patient_Name>Accession, Annie</Patient_Name>
  <Patient_Code>503625913</Patient_Code>
  <Patient_Dob>2018-12-11</Patient_Dob>
  <Patient_Sex>Female</Patient_Sex>
  <Case_Number>CL2020-000089</Case_Number>
  <Case_OrderDate>2020-09-30</Case_OrderDate>
  <Case_CreationDate>2020-09-30T13:32:44.8-04:00</Case_CreationDate>

```

```

<Case_ReportStatus>PENDING</Case_ReportStatus>
<Case_Status>A</Case_Status>
<Case_ObjectID>515593040</Case_ObjectID>
</CaseTable>
<CaseTable>
  <Patient_ObjectID>433509601</Patient_ObjectID>
  <Patient_Name>Friday, Girl</Patient_Name>
  <Patient_Code>433509601</Patient_Code>
  <Patient_Dob>2009-03-13</Patient_Dob>
  <Patient_Sex>Female</Patient_Sex>
  <Case_Number>CL2020-000091</Case_Number>
  <Case_OrderDate>2020-10-05</Case_OrderDate>
  <Case_CreationDate>2020-10-05T15:56:57.673-04:00</Case_CreationDate>
  <Case_ReportStatus>PENDING</Case_ReportStatus>
  <Case_Status>A</Case_Status>
  <Case_ObjectID>515975948</Case_ObjectID>
</CaseTable>
<CaseTable>
  <Patient_ObjectID>516131995</Patient_ObjectID>
  <Patient_Name>Boop, Betty</Patient_Name>
  <Patient_Code>516131995</Patient_Code>
  <Patient_Dob>1965-03-23</Patient_Dob>
  <Patient_Sex>Female</Patient_Sex>
  <Case_Number>CL2020-000105</Case_Number>
  <Case_OrderDate>2020-10-07</Case_OrderDate>
  <Case_CreationDate>2020-10-07T14:49:44.863-04:00</Case_CreationDate>
  <Case_ReportStatus>PENDING</Case_ReportStatus>
  <Case_Status>A</Case_Status>
  <Case_ObjectID>516132019</Case_ObjectID>
</CaseTable>
<CaseTable>
  <Patient_ObjectID>516207223</Patient_ObjectID>
  <Patient_Name>Bell, Tinker</Patient_Name>
  <Patient_Code>516207223</Patient_Code>
  <Patient_Dob>1950-01-01</Patient_Dob>
  <Patient_Sex>Female</Patient_Sex>
  <Case_Number>CL2020-000107</Case_Number>
  <Case_OrderDate>2020-10-08</Case_OrderDate>
  <Case_CreationDate>2020-10-08T13:26:04.687-04:00</Case_CreationDate>
  <Case_ReportStatus>PRELIMINARY</Case_ReportStatus>
  <Case_Status>A</Case_Status>
  <Case_ObjectID>516207240</Case_ObjectID>
</CaseTable>

```

```

<CaseTable>
  <Patient_ObjectID>514701722</Patient_ObjectID>
  <Patient_Name>Boop, Betty</Patient_Name>
  <Patient_Code>514701722</Patient_Code>
  <Patient_Dob>1965-03-23</Patient_Dob>
  <Patient_Sex>Female</Patient_Sex>
  <Case_Number>CL2020-000109</Case_Number>
  <Case_OrderDate>2020-09-18</Case_OrderDate>
  <Case_CreationDate>2020-09-25T11:56:06.06-04:00</Case_CreationDate>
  <Case_ReportStatus>PRELIMINARY</Case_ReportStatus>
  <Case_Status>A</Case_Status>
  <Case_ObjectID>515221990</Case_ObjectID>
</CaseTable>

</NewDataSet>

```

## GET SPECIMEN

<b>Endpoint</b>	api/N/GetSpecimen
<b>Purpose</b>	Return a single, full NucleoLIS specimen record based on specimen unique identifier. Note: The full specimen record is inclusive of all specimen level data elements.
<b>Request Type</b>	GET
<b>Parameters</b>	specimen_ID={specimen_ID} user_ID={user_ID} return_format={return_format}

Where:

**specimen\_id** is the specimen id number of the specimen to return.

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the record.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

**Example**

[https://www.SomeServer.com/api/N/GetSpecimen?specimen\\_id=467291&user\\_id=245654&return\\_format=0](https://www.SomeServer.com/api/N/GetSpecimen?specimen_id=467291&user_id=245654&return_format=0)

**Returns**

If the request succeeds, returns a specimen record matching the filtering criteria in the following XML format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <SpecimenTable>
    <Patient_ObjectID>347631</Patient_ObjectID>
    <Patient_Name>Farrell, Ivy</Patient_Name>
    <Patient_Code>347631</Patient_Code>
    <Patient_Sex>F</Patient_Sex>
    <Patient_Dob>1970-01-01</Patient_Dob>
    <Case_Number>2</Case_Number>
    <Case_Status>I</Case_Status>
    <Specimen_ObjectID>467291</Specimen_ObjectID>
    <Specimen_CreationDate>2015-03-27T15:56:44.83-
04:00</Specimen_CreationDate>
    <Specimen_UpdateTime>2015-03-30T07:07:02.71-
04:00</Specimen_UpdateTime>

    <Specimen_SuperobjectID>467290</Specimen_SuperobjectI
D>
    <Specimen_CREATOR>HIS</Specimen_CREATOR>
    <Specimen_MODIFIER>0</Specimen_MODIFIER>
    <Specimen_CollectionDate>2015-03-
13</Specimen_CollectionDate>

    <Specimen_CollectionTime>16:14</Specimen_CollectionTi
me>
    <Specimen_ReceivedDate>2015-03-
30</Specimen_ReceivedDate>

    <Specimen_ReceivedTime>07:07:02</Specimen_ReceivedTim
e>
    <Specimen_Number>467291</Specimen_Number>
    <Specimen_Status>A</Specimen_Status>
    <Specimen_Receiver>0</Specimen_Receiver>
```



```
</SpecimenTable>

</NewDataSet>
```

## GET SPECIMENS

<b>Endpoint</b>	api/N/GetSpecimens
<b>Purpose</b>	Returns a collection of NucleoLIS Patient, Case and Specimen level data.
<b>Request Type</b>	GET
<b>Parameters</b>	<pre>received_date_from={received_date_from} received_date_to={received_date_to} activesOnly={activesOnly} specimen_status_steps={specimen_status_steps} specimen_sources={specimen_sources} filter_expression={filter_expression} user_ID={user_ID} return_format={return_format}</pre>

Where:

**received\_date\_from** is the "from" date for received specimens. Dates must be in YYYY-MM-DD format.

**received\_date\_to** is the "to" date for received specimens. Dates must be in YYYY-MM-DD format.

**activesOnly** is a true or false value. If true, only active cases are returned. If false, all cases are returned.

**specimen\_status\_steps** is a pipe delimited list of specimen status steps.

**specimen\_sources** is a pipe delimited list of specimen sources

**filter\_expression** is used as a secondary custom filtering criteria applied to the result set . Wildcards are supported. It can be used with any column returned in the result set. Any valid SQL WHERE clause is accepted. Multiple columns can be used when separated with 'AND'. Filtering criteria supported include: Patient\_ObjectID, Patient\_Name, Patient\_Code, Patient\_Dob, Patient\_Sex, Case\_Number, Specimen\_ObjectID, Specimen\_Received date, Specimen\_Received time, Specimen\_CollectionDate, Specimen\_CollectionTime, Specimen\_StatusStep, SourceDefinition\_Code, and SourceDefinition\_Description.

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the record.

**return\_format** is the format for returning the data. Where:

0 = XML Format

1 = JSON Format

**Example** `https://www.SomeServer.com/api/N/GetSpecimens?received_date_from=2020-10-01&received_date_to=2020-11-02&activesOnly=true&specimen_status_steps=&specimen_sources=Blood|Serum&filter_expression=Patient_Sex='Male'&user_id=245654&return_format=0`

**Returns** If the request succeeds, returns a collection of 1..n specimen records matching the filtering criteria in the following XML format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <SpecimenTable>
    <Patient_ObjectID>507179564</Patient_ObjectID>
    <Patient_Name>Bishop, Tom</Patient_Name>
    <Patient_Code>507179564</Patient_Code>
```

```

        <Patient_Dob>1963-08-02</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>5071795660R</Case_Number>
        <Specimen_ReceivedDate>2020-10-
13</Specimen_ReceivedDate>

<Specimen_ReceivedTime>14:07:31</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-10-
13</Specimen_CollectionDate>
        <Specimen_CollectionTime />
        <Specimen_StatusStep />

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>516628021</Specimen_ObjectID>
    </SpecimenTable>
    <SpecimenTable>
        <Patient_ObjectID>514700814</Patient_ObjectID>
        <Patient_Name>Barksalot, Bruno</Patient_Name>
        <Patient_Code>514700814</Patient_Code>
        <Patient_Dob>2016-03-30</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>CL2020-000086</Case_Number>
        <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:15:58</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-09-
18</Specimen_CollectionDate>

<Specimen_CollectionTime>13:38</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>515221988</Specimen_ObjectID>
    </SpecimenTable>
    <SpecimenTable>
        <Patient_ObjectID>503782367</Patient_ObjectID>
        <Patient_Name>Bishop, Tom</Patient_Name>
        <Patient_Code>503782367</Patient_Code>

```

```

        <Patient_Dob>1963-08-02</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>CL2020-000093</Case_Number>
        <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:15:34</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-10-
06</Specimen_CollectionDate>
        <Specimen_CollectionTime />
        <Specimen_StatusStep />

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>516056679</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
        <Patient_ObjectID>511955302</Patient_ObjectID>
        <Patient_Name>Cricket, Jimimy</Patient_Name>
        <Patient_Code>511955302</Patient_Code>
        <Patient_Dob>2018-11-13</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>CL2020-000110</Case_Number>
        <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:16:15</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-07-
21</Specimen_CollectionDate>

<Specimen_CollectionTime>17:28</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>511968657</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
        <Patient_ObjectID>510788319</Patient_ObjectID>
        <Patient_Name>Caron, Jack</Patient_Name>
        <Patient_Code>510788319</Patient_Code>

```

```

        <Patient_Dob>2019-06-15</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>CL2020-000111</Case_Number>
        <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:16:18</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-06-
08</Specimen_CollectionDate>

<Specimen_CollectionTime>16:20</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>510788325</Specimen_ObjectID>
    </SpecimenTable>
    <SpecimenTable>
        <Patient_ObjectID>506987816</Patient_ObjectID>
        <Patient_Name>Duck, Daffy</Patient_Name>
        <Patient_Code>506987816</Patient_Code>
        <Patient_Dob>1950-03-03</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>CL2020-000113</Case_Number>
        <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:16:20</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-03-
12</Specimen_CollectionDate>

<Specimen_CollectionTime>14:40</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>506987822</Specimen_ObjectID>
    </SpecimenTable>
    <SpecimenTable>
        <Patient_ObjectID>506407748</Patient_ObjectID>
        <Patient_Name>Cricket, Jiminy</Patient_Name>

```

```

        <Patient_Code>506407748</Patient_Code>
        <Patient_Dob>2018-11-13</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>CL2020-000115</Case_Number>
        <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:16:24</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-02-
28</Specimen_CollectionDate>

<Specimen_CollectionTime>14:57</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>506407755</Specimen_ObjectID>
    </SpecimenTable>
    <SpecimenTable>
        <Patient_ObjectID>298771370</Patient_ObjectID>
        <Patient_Name>Jones, Tom</Patient_Name>
        <Patient_Code>298771370</Patient_Code>
        <Patient_Dob>1944-01-01</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>CL2020-000155</Case_Number>
        <Specimen_ReceivedDate>2020-10-
21</Specimen_ReceivedDate>

<Specimen_ReceivedTime>17:10:05</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-10-
21</Specimen_CollectionDate>
        <Specimen_CollectionTime />
        <Specimen_StatusStep />

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>517413942</Specimen_ObjectID>
    </SpecimenTable>
    <SpecimenTable>
        <Patient_ObjectID>517962619</Patient_ObjectID>
        <Patient_Name>Caron, Jack</Patient_Name>

```

```

        <Patient_Code>517962619</Patient_Code>
        <Patient_Dob>2019-06-15</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>CL2020-000156</Case_Number>
        <Specimen_ReceivedDate>2020-10-
27</Specimen_ReceivedDate>

<Specimen_ReceivedTime>10:50:07</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-10-
27</Specimen_CollectionDate>

<Specimen_CollectionTime>10:46</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>517962628</Specimen_ObjectID>
    </SpecimenTable>
    <SpecimenTable>
        <Patient_ObjectID>518533883</Patient_ObjectID>
        <Patient_Name>Caron, Jack</Patient_Name>
        <Patient_Code>518533883</Patient_Code>
        <Patient_Dob>2019-06-15</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>CL2020-000160</Case_Number>
        <Specimen_ReceivedDate>2020-11-
02</Specimen_ReceivedDate>

<Specimen_ReceivedTime>11:05:39</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-11-
02</Specimen_CollectionDate>

<Specimen_CollectionTime>10:59</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>518533891</Specimen_ObjectID>
    </SpecimenTable>
    <SpecimenTable>
        <Patient_ObjectID>516213756</Patient_ObjectID>

```

```

        <Patient_Name>Caron, Jack</Patient_Name>
        <Patient_Code>516213756</Patient_Code>
        <Patient_Dob>2019-06-15</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>CLI2020-000002</Case_Number>
        <Specimen_ReceivedDate>2020-10-
08</Specimen_ReceivedDate>

<Specimen_ReceivedTime>18:15:11</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-10-
08</Specimen_CollectionDate>

<Specimen_CollectionTime>18:10</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>516213763</Specimen_ObjectID>
    </SpecimenTable>
    <SpecimenTable>
        <Patient_ObjectID>516287111</Patient_ObjectID>
        <Patient_Name>Caron, Jack</Patient_Name>
        <Patient_Code>516287111</Patient_Code>
        <Patient_Dob>2019-06-15</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>CLI2020-000005</Case_Number>
        <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

<Specimen_ReceivedTime>14:45:44</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-10-
09</Specimen_CollectionDate>

<Specimen_CollectionTime>14:43</Specimen_CollectionTi
me>

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>516287120</Specimen_ObjectID>
    </SpecimenTable>
    <SpecimenTable>

```



```

        <Patient_ObjectID>516287344</Patient_ObjectID>
        <Patient_Name>Cricket, Jiminy</Patient_Name>
        <Patient_Code>516287344</Patient_Code>
        <Patient_Dob>2018-11-13</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>CLI2020-000006</Case_Number>
        <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

    <Specimen_ReceivedTime>15:44:18</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-10-
09</Specimen_CollectionDate>

    <Specimen_CollectionTime>14:46</Specimen_CollectionTi
me>

    <SourceDefinition_Code>Blood</SourceDefinition_Code>

    <SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>516287379</Specimen_ObjectID>
    </SpecimenTable>
    <SpecimenTable>
        <Patient_ObjectID>516291416</Patient_ObjectID>
        <Patient_Name>Duck, Daffy</Patient_Name>
        <Patient_Code>516291416</Patient_Code>
        <Patient_Dob>1950-03-03</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>CLI2020-000009</Case_Number>
        <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

    <Specimen_ReceivedTime>18:11:49</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-10-
09</Specimen_CollectionDate>

    <Specimen_CollectionTime>15:54</Specimen_CollectionTi
me>

    <SourceDefinition_Code>Blood</SourceDefinition_Code>

    <SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>516291421</Specimen_ObjectID>
    </SpecimenTable>

```

```

    <SpecimenTable>
      <Patient_ObjectID>516276481</Patient_ObjectID>
      <Patient_Name>Duck, Daffy</Patient_Name>
      <Patient_Code>516276481</Patient_Code>
      <Patient_Dob>1950-03-03</Patient_Dob>
      <Patient_Sex>Male</Patient_Sex>
      <Case_Number>CLI2020-000010</Case_Number>
      <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

    <Specimen_ReceivedTime>18:12:16</Specimen_ReceivedTim
e>
      <Specimen_CollectionDate>2020-10-
09</Specimen_CollectionDate>

    <Specimen_CollectionTime>11:46</Specimen_CollectionTi
me>

    <SourceDefinition_Code>Blood</SourceDefinition_Code>

    <SourceDefinition_Description>Blood</SourceDefinition
_Description>
      <Specimen_ObjectID>516276511</Specimen_ObjectID>
    </SpecimenTable>
    <SpecimenTable>
      <Patient_ObjectID>516298160</Patient_ObjectID>
      <Patient_Name>QA, Final</Patient_Name>
      <Patient_Code>516298160</Patient_Code>
      <Patient_Dob>2019-10-23</Patient_Dob>
      <Patient_Sex>Male</Patient_Sex>
      <Case_Number>CLI2020-000013</Case_Number>
      <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

    <Specimen_ReceivedTime>18:13:24</Specimen_ReceivedTim
e>
      <Specimen_CollectionDate>2020-10-
09</Specimen_CollectionDate>

    <Specimen_CollectionTime>17:45</Specimen_CollectionTi
me>

    <SourceDefinition_Code>Blood</SourceDefinition_Code>

    <SourceDefinition_Description>Blood</SourceDefinition
_Description>
      <Specimen_ObjectID>516298169</Specimen_ObjectID>

```

```

</SpecimenTable>
<SpecimenTable>
  <Patient_ObjectID>518078029</Patient_ObjectID>
  <Patient_Name>QA, Final</Patient_Name>
  <Patient_Code>518078029</Patient_Code>
  <Patient_Dob>2019-10-23</Patient_Dob>
  <Patient_Sex>Male</Patient_Sex>
  <Case_Number>CLI2020-000016</Case_Number>
  <Specimen_ReceivedDate>2020-10-
28</Specimen_ReceivedDate>

  <Specimen_ReceivedTime>15:48:59</Specimen_ReceivedTim
e>
  <Specimen_CollectionDate>2020-10-
28</Specimen_CollectionDate>

  <Specimen_CollectionTime>15:44</Specimen_CollectionTi
me>

  <SourceDefinition_Code>Blood</SourceDefinition_Code>

  <SourceDefinition_Description>Blood</SourceDefinition
_Description>
    <Specimen_ObjectID>518078227</Specimen_ObjectID>
  </SpecimenTable>
  <SpecimenTable>
    <Patient_ObjectID>510676219</Patient_ObjectID>
    <Patient_Name>Cricket, Jiminy</Patient_Name>
    <Patient_Code>510676219</Patient_Code>
    <Patient_Dob>2011-07-01</Patient_Dob>
    <Patient_Sex>Male</Patient_Sex>
    <Case_Number>FC2020-000002</Case_Number>
    <Specimen_ReceivedDate>2020-10-
13</Specimen_ReceivedDate>

    <Specimen_ReceivedTime>08:30:47</Specimen_ReceivedTim
e>
    <Specimen_CollectionDate>2020-06-
04</Specimen_CollectionDate>

    <Specimen_CollectionTime>15:26</Specimen_CollectionTi
me>

    <SourceDefinition_Code>Blood</SourceDefinition_Code>

    <SourceDefinition_Description>Blood</SourceDefinition
_Description>

```

```

        <Specimen_ObjectID>510676225</Specimen_ObjectID>
    </SpecimenTable>
    <SpecimenTable>
        <Patient_ObjectID>516298310</Patient_ObjectID>
        <Patient_Name>QA, Final</Patient_Name>
        <Patient_Code>516298310</Patient_Code>
        <Patient_Dob>2019-10-23</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>HE2020-000003</Case_Number>
        <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

    <Specimen_ReceivedTime>18:14:16</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-10-
09</Specimen_CollectionDate>

    <Specimen_CollectionTime>17:49</Specimen_CollectionTi
me>

    <SourceDefinition_Code>Blood</SourceDefinition_Code>

    <SourceDefinition_Description>Blood</SourceDefinition
_Description>
        <Specimen_ObjectID>516298330</Specimen_ObjectID>
    </SpecimenTable>
    <SpecimenTable>
        <Patient_ObjectID>506647493</Patient_ObjectID>
        <Patient_Name>Mouse, Mickey</Patient_Name>
        <Patient_Code>506647493</Patient_Code>
        <Patient_Dob>2000-07-08</Patient_Dob>
        <Patient_Sex>Male</Patient_Sex>
        <Case_Number>M2020-000030</Case_Number>
        <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

    <Specimen_ReceivedTime>18:16:21</Specimen_ReceivedTim
e>
        <Specimen_CollectionDate>2020-03-
04</Specimen_CollectionDate>

    <Specimen_CollectionTime>16:09</Specimen_CollectionTi
me>

    <SourceDefinition_Code>Blood</SourceDefinition_Code>

```

```

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
  <Specimen_ObjectID>506647513</Specimen_ObjectID>
</SpecimenTable>
<SpecimenTable>
  <Patient_ObjectID>510588219</Patient_ObjectID>
  <Patient_Name>Covid, Charlie</Patient_Name>
  <Patient_Code>510588219</Patient_Code>
  <Patient_Dob>2019-12-15</Patient_Dob>
  <Patient_Sex>Male</Patient_Sex>
  <Case_Number>M2020-000059</Case_Number>
  <Specimen_ReceivedDate>2020-10-
09</Specimen_ReceivedDate>

  <Specimen_ReceivedTime>18:36:05</Specimen_ReceivedTim
e>
  <Specimen_CollectionDate>2020-06-
01</Specimen_CollectionDate>

  <Specimen_CollectionTime>10:00</Specimen_CollectionTi
me>
  <Specimen_StatusStep />

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
  <Specimen_ObjectID>510589737</Specimen_ObjectID>
</SpecimenTable>
</NewDataSet>

```

## GET TEST ORDER

<b>Endpoint</b>	api/N/GetTestOrder
<b>Purpose</b>	Returns a single, full test record and available results based on test order unique identifier.

**Request Type** GET

**Parameters**

```
test_order_ID={test_order_ID}

user_ID={user_ID}

return_format={return_format}
```

Where:

**test\_order\_id** is the id of the test order to return.

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the record.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

**Example**

```
https://www.SomeServer.com/api/N/GetTestOrder?test_order_id=8606534&user_id=245654&return_format=0
```

**Returns**

If the request succeeds, it returns the requested full profile test order record inclusive of all test order and available test result data in the following XML format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <ProfileTable>
    <Patient_ObjectID>4180414</Patient_ObjectID>
    <Patient_Name>Wednesday, Wendy</Patient_Name>
    <Patient_Code>4180414</Patient_Code>
    <Patient_Sex />
    <Patient_Dob />
    <Case_Number>PSY2015-000042</Case_Number>
    <Case_Status>I</Case_Status>
    <Specimen_ObjectID>8606398</Specimen_ObjectID>
```

```

        <Specimen_CollectionDate />
        <Specimen_CollectionTime />
        <Specimen_ReceivedDate>2015-05-
19</Specimen_ReceivedDate>

<Specimen_ReceivedTime>15:37:45</Specimen_ReceivedTim
e>
    <Specimen_StatusStep />

<SourceDefinition_Code>Blood</SourceDefinition_Code>

<SourceDefinition_Description>Blood</SourceDefinition
_Description>
    <Profile_TYPE>FISH</Profile_TYPE>
    <Profile_ObjectID>8606534</Profile_ObjectID>
    <Profile_CreationDate>2015-05-19T15:37:18.65-
04:00</Profile_CreationDate>
    <Profile_UpdateTime>2015-05-19T15:42:05.01-
04:00</Profile_UpdateTime>

<Profile_SuperobjectID>8606398</Profile_SuperobjectID
>
    <Profile_User1 />
    <Profile_User2 />
    <Profile_User3 />
    <Profile_User4 />
    <Profile_CREATOR>deannac</Profile_CREATOR>
    <Profile_MODIFIER>deannac</Profile_MODIFIER>
    <Profile_Description />
    <Profile_OrderDate>2015-05-19</Profile_OrderDate>
    <Profile_OrderTime>15:37</Profile_OrderTime>
    <Profile_OrderedBy>deannac</Profile_OrderedBy>

<Profile_PanelDefinition>249815</Profile_PanelDefinit
ion>
    <Profile_R1>Clinical Significance: ROS1 gene
rearrangements are found in 1-2% of non-small cell
lung carcinoma (NSCLC). Pre-clinical and early
clinical evidence suggests ROS1-rearranged tumors may
be sensitive to the dual ALK/MET inhibitor
crizotinib.</Profile_R1>
    <Profile_R2 />
    <Profile_R3 />
    <Profile_R4 />
    <Profile_R5 />
    <Profile_R6 />
    <Profile_Result>Present</Profile_Result>

```

```

        <Profile_Reportable>True</Profile_Reportable>
        <Profile_DatumDef>245550</Profile_DatumDef>

<Profile_DatumEntryDefinition>4380952</Profile_DatumEntryDefinition>
    <Profile_Images><?xml version="1.0"
encoding="utf-16"?>
<ChrysalisImgDats
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema">
    <ImgDat />
</ChrysalisImgDats></Profile_Images>
    <Profile_Abnormal>False</Profile_Abnormal>
    <Profile_AOV>False</Profile_AOV>
    <Profile_StatusStep />
    <Profile_Range />
    <Profile_RangeComment />
    <Profile_CanceledDate />
    <Profile_CanceledBy />
    <Profile_CanceledFor />
    <Profile_ReflexedBy />
    <Profile_IOrder />
    <Profile_IOrderDate />
    <Profile_IOrderTime />
    <Profile_Complete>True</Profile_Complete>

<Profile_RequireSignature>>false</Profile_RequireSignature>

<Profile_RequireApproval>>false</Profile_RequireApproval>

<Profile_ProfileDefinitionObjectID>249815</Profile_ProfileDefinitionObjectID>
    <Profile_ProfileDefinitionCreationDate>2014-12-31T10:55:34.78-05:00</Profile_ProfileDefinitionCreationDate>
    <Profile_ProfileDefinitionUpdateTime>2015-10-01T15:08:20.637-04:00</Profile_ProfileDefinitionUpdateTime>
    <Profile_ProfileDefinitionUser1 />
    <Profile_ProfileDefinitionUser2 />
    <Profile_ProfileDefinitionUser3 />
    <Profile_ProfileDefinitionUser4 />

<Profile_ProfileDefinitionCREATOR>0</Profile_ProfileDefinitionCREATOR>

```



<Profile\_ProfileDefinitionMODIFIER>deannac</Profile\_ProfileDefinitionMODIFIER>

<Profile\_ProfileDefinitionCode>NSmCCA</Profile\_ProfileDefinitionCode>

<Profile\_ProfileDefinitionDescription>Non-Small Cell Carcinoma, Lung</Profile\_ProfileDefinitionDescription>

<Profile\_ProfileDefinitionRETIRED>False</Profile\_ProfileDefinitionRETIRED>

<Profile\_ProfileDefinitionDatumLink>245550</Profile\_ProfileDefinitionDatumLink>

<Profile\_ProfileDefinitionConstituentDatumLink>34247712</Profile\_ProfileDefinitionConstituentDatumLink>

<Profile\_ProfileDefinitionDatumEntryDefinition>34714941</Profile\_ProfileDefinitionDatumEntryDefinition>

<Profile\_ProfileDefinitionAlwaysReportMessage>ROS1 FISH uses a dual-color break-apart probe which produces single green, red, and yellow or fusion signals in rearranged cells.

#### References:

Kwak EL et al. New Engl J Med. 2010;363(18):1693-1703.

Vysis ALK Break Apart FISH Probe Kit package insert. Vysis ALK United States website. Accessed 12-31-2012.

Bergethon K, Shaw AT, Ignatius Ou S, et al. J Clin Oncol. 2012;30(8):863-870.

Takeuchi K, Soda M, Togashi Y, et al. Nature Med. 2012;18(3):378-381.

Davies KD, Le AT, Theodoro MF, et al. Clin Cancer Res. 2012;18(17):4570-4579.

Shaw AT, Camidge DR, Engelman JA, et al. J Clin Oncol. 2012;30 (suppl): abstract 7508.

Bang YJ. Arch Pathol Lab Med. 2012;136:1201-1204.

NCCN Guidelines, Non-Small Cell Lung Cancer, version 1.2013. National Comprehensive Cancer Network website. Accessed 12-31-2012.

ALK page. Atlas of Genetics and Cytogenetics in Oncology and Haematology. Accessed 12-31-2012; last updated 2-2010.

Stumpfova M, Janne PA. Clin Cancer Res. 2012;18(16):4222-4224.

Sasaki T et al. Eur J Cancer. 2010;46:1773-1780.  
</Profile\_ProfileDefinitionAlwaysReportMessage>

<Profile\_ProfileDefinitionMethodology>FISH</Profile\_ProfileDefinitionMethodology>

<Profile\_ProfileDefinitionCodeMap />

<Profile\_ProfileDefinitionLabelCount>1</Profile\_ProfileDefinitionLabelCount>

<Profile\_ProfileDefinitionAutoComplete>False</Profile\_ProfileDefinitionAutoComplete>

<Profile\_ProfileDefinitionWPDiscreteAdvanceOnly>>false</Profile\_ProfileDefinitionWPDiscreteAdvanceOnly>

<Profile\_DatumDefinitionObjectID>245550</Profile\_DatumDefinitionObjectID>

<Profile\_DatumDefinitionCreationDate>2014-12-31T10:54:49.03-05:00</Profile\_DatumDefinitionCreationDate>

<Profile\_DatumDefinitionUpdateTime>2019-04-09T15:17:02.137-04:00</Profile\_DatumDefinitionUpdateTime>

<Profile\_DatumDefinitionUser1 />

<Profile\_DatumDefinitionUser2 />

<Profile\_DatumDefinitionUser3 />

<Profile\_DatumDefinitionUser4 />

<Profile\_DatumDefinitionCREATOR>0</Profile\_DatumDefinitionCREATOR>

<Profile\_DatumDefinitionMODIFIER>0</Profile\_DatumDefinitionMODIFIER>

<Profile\_DatumDefinitionCode>416980</Profile\_DatumDefinitionCode>

<Profile\_DatumDefinitionDescription />

<Profile\_DatumDefinitionRETIRED>False</Profile\_DatumDefinitionRETIRED>

```

<Profile_DatumDefinitionR1>Comment</Profile_DatumDefinitionR1>
    <Profile_DatumDefinitionR2 />
    <Profile_DatumDefinitionR3 />
    <Profile_DatumDefinitionR4 />
    <Profile_DatumDefinitionR5 />
    <Profile_DatumDefinitionR6 />

<Profile_DatumDefinitionHideResult>False</Profile_DatumDefinitionHideResult>
    <Profile_DatumDefinitionRange />
    <Profile_DatumDefinitionRangeComment />
    <Profile_DatumDefinitionExtensions><?xml
version="1.0" encoding="utf-16"?>
<DatumExtensionDefinition
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema">
    <Items>
        <DatumExtensionItemDefinition>
            <FieldName>Result</FieldName>
            <IsReadOnly>>false</IsReadOnly>
        </DatumExtensionItemDefinition>
        <DatumExtensionItemDefinition>
            <FieldName>R1</FieldName>
            <IsReadOnly>>false</IsReadOnly>
        </DatumExtensionItemDefinition>
        <DatumExtensionItemDefinition>
            <FieldName>R2</FieldName>
            <IsReadOnly>>false</IsReadOnly>
        </DatumExtensionItemDefinition>
        <DatumExtensionItemDefinition>
            <FieldName>R3</FieldName>
            <IsReadOnly>>false</IsReadOnly>
        </DatumExtensionItemDefinition>
        <DatumExtensionItemDefinition>
            <FieldName>R4</FieldName>
            <IsReadOnly>>false</IsReadOnly>
        </DatumExtensionItemDefinition>
        <DatumExtensionItemDefinition>
            <FieldName>R5</FieldName>
            <IsReadOnly>>false</IsReadOnly>
        </DatumExtensionItemDefinition>
        <DatumExtensionItemDefinition>
            <FieldName>R6</FieldName>
            <IsReadOnly>>false</IsReadOnly>
        </DatumExtensionItemDefinition>
    </Items>
</DatumExtensionDefinition>

```

```

        <DatumExtensionItemDefinition>
            <FieldName>Range</FieldName>
            <IsReadOnly>true</IsReadOnly>
        </DatumExtensionItemDefinition>
        <DatumExtensionItemDefinition>
            <FieldName>RangeComment</FieldName>
            <IsReadOnly>true</IsReadOnly>
        </DatumExtensionItemDefinition>
    </Items>
</DatumExtensionDefinition></Profile_DatumDefinitionE
xtensions>
    <Profile_DatumDefinitionTableDisplayDefinition />
    <Profile_DatumDefinitionResultFieldMap />
    <Profile_DatumDefinitionInstrumentAutoOrder />
    <Constituent_TYPE>FISH</Constituent_TYPE>

<Constituent_ObjectID>8606535</Constituent_ObjectID>
    <Constituent_CreationDate>2015-05-
19T15:37:18.843-04:00</Constituent_CreationDate>
    <Constituent_UpdateTime>2015-05-19T15:39:04.96-
04:00</Constituent_UpdateTime>

<Constituent_SuperobjectID>8606534</Constituent_Super
objectID>
    <Constituent_User1 />
    <Constituent_User2 />
    <Constituent_User3 />
    <Constituent_User4 />

<Constituent_CREATOR>deannac</Constituent_CREATOR>

<Constituent_MODIFIER>deannac</Constituent_MODIFIER>

<Constituent_Reportable>True</Constituent_Reportable>
    <Constituent_Description />

<Constituent_Definition>249789</Constituent_Definitio
n>
    <Constituent_Result>Present</Constituent_Result>
    <Constituent_R1 />
    <Constituent_R2 />
    <Constituent_R3 />
    <Constituent_R4 />
    <Constituent_R5 />
    <Constituent_R6 />

<Constituent_Abnormal>False</Constituent_Abnormal>

```

```

    <Constituent_AOV>False</Constituent_AOV>
    <Constituent_Range />
    <Constituent_RangeComment />

    <Constituent_ConstituentDefinitionObjectID>249789</Co
nstituent_ConstituentDefinitionObjectID>

    <Constituent_ConstituentDefinitionCreationDate>2014-
12-31T10:55:34.527-
05:00</Constituent_ConstituentDefinitionCreationDate>

    <Constituent_ConstituentDefinitionUpdateTime>2014-12-
31T10:55:34.527-
05:00</Constituent_ConstituentDefinitionUpdateTime>
    <Constituent_ConstituentDefinitionUser1 />
    <Constituent_ConstituentDefinitionUser2 />
    <Constituent_ConstituentDefinitionUser3 />
    <Constituent_ConstituentDefinitionUser4 />

    <Constituent_ConstituentDefinitionCREATOR>0</Constitu
ent_ConstituentDefinitionCREATOR>

    <Constituent_ConstituentDefinitionMODIFIER>0</Constit
uent_ConstituentDefinitionMODIFIER>

    <Constituent_ConstituentDefinitionCode>ROS1</Constitu
ent_ConstituentDefinitionCode>

    <Constituent_ConstituentDefinitionDescription>break-
apart probe at
(6q22)</Constituent_ConstituentDefinitionDescription>

    <Constituent_ConstituentDefinitionRETIRED>False</Cons
tituent_ConstituentDefinitionRETIRED>
    <Constituent_ConstituentDefinitionCodeMap />
    </ProfileTable>

</NewDataSet>

```

## GET TEST ORDERS

**Endpoint**      `api/N/GetTestOrders`

**Purpose** Returns a collection of Patient, Case, Specimen and Test level data including Patient\_ObjectID, Patient\_Name, Patient\_Code, Patient\_Dob, Patient\_Sex, Case\_Number, Specimen\_ObjectID, Specimen\_Received date, Specimen\_Received time, Specimen\_CollectionDate, Specimen\_CollectionTime, Specimen\_StatusStep, SourceDefinition\_Code, SourceDefinition\_Description, Test order date, status, completion status, test code and unique identifier of all ordered tests records meeting filtering requirements.

.

**Request Type** GET

**Parameters**

order\_date\_from={order\_date\_from}

order\_date\_to={order\_date\_to}

activesOnly={activesOnly}

order\_status\_steps={order\_status\_steps}

order\_codes={order\_codes}

filter\_expression={filter\_expression}

user\_ID={user\_ID}

return\_format={return\_format}

Where:

**order\_date\_from** is the "from" order date. Dates must be in YYYY-MM-DD format.

**order\_date\_to** is the "to" order date. Dates must be in YYYY-MM-DD format.

**activesOnly** is a true or false value. If set to true, only active orders records are returned if false, all order records are returned.

**order\_status\_steps** is pipe delimited list of order status steps

**order\_codes** is a pipe delimited list of order codes

**filter\_expression** is used as a secondary custom filtering criteria applied to the result set . Wildcards are supported. It can be used with any column returned in the result set. Any valid SQL WHERE clause is accepted. Multiple columns can be used when separated with 'AND'. Filtering criteria supported include: Patient\_ObjectID, Patient\_Name, Patient\_Code, Patient\_DOB, Patient\_Sex, Case\_Number, Specimen\_ObjectID, Specimen\_RetreivalDate, Specimen\_RetreivalTime, Specimen\_CollectionDate, Specimen\_CollectionTime, Specimen\_StatusStep, SourceDefinition\_Code, SourceDefinition\_Description, Profile\_ObjectID, Profile\_OrderDate, Profile\_OrderTime, Profile\_ProfileDefinitionCode, Profile\_ProfileDefinitionDescription, Profile\_Complete, Profile\_CompleteDate, Profile\_SignerDate, and Profile\_StatusStep.

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the record.

**return\_format** is the format for returning the data. Where:

0 = XML Format

1 = JSON Format

**Example** `https://www.SomeServer.com/api/N/GetTestOrders?order_date_from=2020-09-17&order_date_to=2020-09-18&activesOnly=true&order_status_steps=*&order_codes=CDI&filter_expression=Patient_DOB like '1960% '&user_ID=245654&return_format=0`

**Returns** If the request succeeds, returns a collection of 1..n test order profile records matching the filtering criteria in the following XML format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
<ProfileTable>
  <Patient_ObjectID>506646654</Patient_ObjectID>
  <Patient_Name>Buckley, Betty Sue Joe</Patient_Name>
  <Patient_Code>506646654</Patient_Code>
```

```

<Patient_Dob>1960-07-14</Patient_Dob>
<Patient_Sex>Female</Patient_Sex>
<Case_Number>CL2020-000013</Case_Number>
<Specimen_ObjectID>506686900</Specimen_ObjectID>
<Specimen_ReceivedDate>2020-03-04</Specimen_ReceivedDate>
<Specimen_ReceivedTime>15:49:27</Specimen_ReceivedTime>
<Specimen_CollectionDate>2020-03-04</Specimen_CollectionDate>
<Specimen_CollectionTime>15:46</Specimen_CollectionTime>
<Specimen_StatusStep />
<SourceDefinition_Code>Blood</SourceDefinition_Code>
<SourceDefinition_Description>Blood</SourceDefinition_Description>
<Profile_ObjectID>514633950</Profile_ObjectID>
<Profile_OrderDate>2020-09-17</Profile_OrderDate>
<Profile_OrderTime>16:05</Profile_OrderTime>
<Profile_ProfileDefinitionCode>CDCI</Profile_ProfileDefinitionCode>

<Profile_ProfileDefinitionDescription>CDCI</Profile_ProfileDefinitionDescription>
ion>
  <Profile_Complete>False</Profile_Complete>
  <Profile_CompleteDate />
  <Profile_SignerDate />
  <Profile_StatusStep />
</ProfileTable>
<ProfileTable>
  <Patient_ObjectID>291217280</Patient_ObjectID>
  <Patient_Name>Ashbey, Clywd</Patient_Name>
  <Patient_Code>291217280</Patient_Code>
  <Patient_Dob>1960-07-07</Patient_Dob>
  <Patient_Sex>Male</Patient_Sex>
  <Case_Number>CL2020-000068</Case_Number>
  <Specimen_ObjectID>511969438</Specimen_ObjectID>
  <Specimen_ReceivedDate>2020-08-06</Specimen_ReceivedDate>
  <Specimen_ReceivedTime>08:22:09</Specimen_ReceivedTime>
  <Specimen_CollectionDate>2020-08-06</Specimen_CollectionDate>
  <Specimen_CollectionTime />
  <Specimen_StatusStep />
  <SourceDefinition_Code>Blood</SourceDefinition_Code>
  <SourceDefinition_Description>Blood</SourceDefinition_Description>
  <Profile_ObjectID>514715461</Profile_ObjectID>
  <Profile_OrderDate>2020-09-18</Profile_OrderDate>
  <Profile_OrderTime>18:21</Profile_OrderTime>
  <Profile_ProfileDefinitionCode>CDCI</Profile_ProfileDefinitionCode>

```



<Profile\_ProfileDefinitionDescription>CDCl</Profile\_ProfileDefinitionDescription>

<Profile\_Complete>False</Profile\_Complete>

<Profile\_CompleteDate />

<Profile\_SignerDate />

<Profile\_StatusStep />

</ProfileTable>

<ProfileTable>

<Patient\_ObjectID>132099125</Patient\_ObjectID>

<Patient\_Name>Miller, Leonard</Patient\_Name>

<Patient\_Code>132099125</Patient\_Code>

<Patient\_Dob>1960-12-12</Patient\_Dob>

<Patient\_Sex>Male</Patient\_Sex>

<Case\_Number>CL2020-000079</Case\_Number>

<Specimen\_ObjectID>514704804</Specimen\_ObjectID>

<Specimen\_ReceivedDate>2020-09-18</Specimen\_ReceivedDate>

<Specimen\_ReceivedTime>16:57:21</Specimen\_ReceivedTime>

<Specimen\_CollectionDate>2020-09-18</Specimen\_CollectionDate>

<Specimen\_CollectionTime />

<Specimen\_StatusStep />

<SourceDefinition\_Code>Blood</SourceDefinition\_Code>

<SourceDefinition\_Description>Blood</SourceDefinition\_Description>

<Profile\_ObjectID>514704809</Profile\_ObjectID>

<Profile\_OrderDate>2020-09-18</Profile\_OrderDate>

<Profile\_OrderTime>14:56</Profile\_OrderTime>

<Profile\_ProfileDefinitionCode>CDCl</Profile\_ProfileDefinitionCode>

<Profile\_ProfileDefinitionDescription>CDCl</Profile\_ProfileDefinitionDescription>

<Profile\_Complete>False</Profile\_Complete>

<Profile\_CompleteDate />

<Profile\_SignerDate />

<Profile\_StatusStep />

</ProfileTable>

<ProfileTable>

<Patient\_ObjectID>376500</Patient\_ObjectID>

<Patient\_Name>Farrell, Tiffany</Patient\_Name>

<Patient\_Code>376500</Patient\_Code>

<Patient\_Dob>1960-07-01</Patient\_Dob>

<Patient\_Sex>Female</Patient\_Sex>

<Case\_Number>PSY2015-000029</Case\_Number>

<Specimen\_ObjectID>4350748</Specimen\_ObjectID>

```

<Specimen_ReceivedDate>2015-06-02</Specimen_ReceivedDate>
<Specimen_ReceivedTime>07:40:06</Specimen_ReceivedTime>
<Specimen_CollectionDate>2015-04-22</Specimen_CollectionDate>
<Specimen_CollectionTime>10:00</Specimen_CollectionTime>
<Specimen_StatusStep />
<SourceDefinition_Code>Blood</SourceDefinition_Code>
<SourceDefinition_Description>Blood</SourceDefinition_Description>
<Profile_ObjectID>514703397</Profile_ObjectID>
<Profile_OrderDate>2020-09-18</Profile_OrderDate>
<Profile_OrderTime>14:30</Profile_OrderTime>
<Profile_ProfileDefinitionCode>CDCl</Profile_ProfileDefinitionCode>

<Profile_ProfileDefinitionDescription>CDCl</Profile_ProfileDefinitionDescription>
</ProfileTable>
</NewDataSet>

```

## GET PHYSICIAN

<b>Endpoint</b>	api/N/GetPhysician
<b>Purpose</b>	Returns a single, full physician record based on a unique code.
<b>Request Type</b>	GET
<b>Parameters</b>	<pre> physician_code={physician_code}  user_ID={user_ID}  return_format={return_format} </pre>

Where:

**physician\_code** is the unique identifier for the requested physician.

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the record.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

#### Example

```
https://www.SomeServer.com/  
api/N/GetPhysician?physician_code=381596&user_ID=24565  
4&return_format=0
```

#### Returns

If the request succeeds, returns the requested physician record in the following format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>  
  <PhysicianTable>  
    <Physician_ObjectID>230022</Physician_ObjectID>  
    <Physician_CreationDate>2004-07-30T12:14:57.687-  
04:00</Physician_CreationDate>  
    <Physician_UpdateTime>2020-01-21T11:25:23.293-  
05:00</Physician_UpdateTime>  
    <Physician_Name>test, physician</Physician_Name>  
    <Physician_Sex />  
    <Physician_Street />  
    <Physician_City />  
    <Physician_State />  
    <Physician_Zip />  
    <Physician_Phone />  
    <Physician_Code>TEST1</Physician_Code>  
    <Physician_Soundex>TESTPHYSCN</Physician_Soundex>  
    <Physician_UPIN />  
    <Physician_STOP />  
    <Physician_STREET2 />  
    <Physician_RETIRED>False</Physician_RETIRED>  
    <Physician_GROUP>False</Physician_GROUP>  
    <Physician_STREET3 />
```

```

        <Physician_MODIFIER>0</Physician_MODIFIER>
        <Physician_NPI />
        <Physician_WebPathStartDate />

    <Physician_WithholdGroupResults>False</Physician_With
holdGroupResults>
        <Physician_WebPathReturnAs />
        <Physician_WebCopyTo>False</Physician_WebCopyTo>
        <Physician_Instructions />

    <Physician_WebPathSupP>False</Physician_WebPathSupP>

    <Physician_WebPathALL>False</Physician_WebPathALL>
        <Physician_WebPathReturnCC><ALL MY
GROUPS></Physician_WebPathReturnCC>

    <Physician_UsePreliminaryDistribution>False</Physicia
n_UsePreliminaryDistribution>
        <Physician_PreliminaryReportDistribution><?xml
version="1.0" encoding="utf-16"?>
        <DistributionDefinition
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema">
            <AsSubmtterDistribution>
                <PrintEnabled>>false</PrintEnabled>
                <FaxEnabled>>false</FaxEnabled>
                <SecureMailEnabled>>false</SecureMailEnabled>
            </AsSubmtterDistribution>
            <AsOtherDistribution>
                <PrintEnabled>>false</PrintEnabled>
                <FaxEnabled>>false</FaxEnabled>
                <SecureMailEnabled>>false</SecureMailEnabled>
            </AsOtherDistribution>
        </DistributionDefinition></Physician_PreliminaryRepor
tDistribution>
        </PhysicianTable>

</NewDataSet>

```

## GET PHYSICIANS

**Endpoint**      `api/N/GetPhysicians`

<b>Purpose</b>	Returns a collection of 1..n physician records based on filter parameters passed.
<b>Request Type</b>	GET
<b>Parameters</b>	<p> last_name={last_name}  first_name={first_name}  filter_expression={filter_expression}  user_id={user_id}  return_format={return_format} </p> <p>Where:</p> <p><b>last_name</b> is the last name of the Physician. Wildcards are supported.</p> <p><b>first_name</b> is the first name of the physician. Wildcards are supported.</p> <p><b>filter_expression</b> is is used as a secondary custom filtering criteria applied to the result set . Wildcards are supported. Any valid SQL WHERE clause is accepted. It can be used with any column returned in the result set. Multiple columns can be used when separated with 'AND'.</p> <p><b>user_id</b> is the user code for the user in the NucleoLIS user table that is requesting the record.</p> <p><b>return_format</b> is the format for returning the data.</p> <p>Where:</p> <p>0 = XML Format</p> <p>1 = JSON Format</p>

**Example**

```
https://www.SomeServer.com/  
api/N/GetPhysicians?last_name=*&first_name=T*&filter_e  
xpression=Physician_GROUP='false'&user_id=245654&retur  
n_format=0
```

**Returns**

If the request succeeds, returns the requested physician records in the following format:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>  
  <PhysicianTable>  
    <Physician_ObjectID>15460945</Physician_ObjectID>  
    <Physician_Name>Jones, Tim</Physician_Name>  
    <Physician_Street>This physician has a single  
"none" location called Office</Physician_Street>  
    <Physician_STREET2 />  
    <Physician_STREET3 />  
    <Physician_City />  
    <Physician_State />  
    <Physician_Zip />  
    <Physician_Code>JONES</Physician_Code>  
    <Physician_GROUP>False</Physician_GROUP>  
  </PhysicianTable>  
  <PhysicianTable>  
    <Physician_ObjectID>503790496</Physician_ObjectID>  
    <Physician_Name>New physician, testing  
roles</Physician_Name>  
    <Physician_Street />  
    <Physician_STREET2 />  
    <Physician_STREET3 />  
    <Physician_City />  
    <Physician_State />  
    <Physician_Zip />  
    <Physician_Code>NP1</Physician_Code>  
    <Physician_GROUP>False</Physician_GROUP>  
  </PhysicianTable>  
  <PhysicianTable>  
    <Physician_ObjectID>503612491</Physician_ObjectID>  
    <Physician_Name>Sig, Text Only</Physician_Name>  
    <Physician_Street />  
    <Physician_STREET2 />  
    <Physician_STREET3 />  
    <Physician_City />  
    <Physician_State />  
    <Physician_Zip />
```

```

        <Physician_Code>Phy01</Physician_Code>
        <Physician_GROUP>False</Physician_GROUP>
</PhysicianTable>
<PhysicianTable>
    <Physician_ObjectID>390657</Physician_ObjectID>
    <Physician_Name>Physician, Test</Physician_Name>
    <Physician_Street />
    <Physician_STREET2 />
    <Physician_STREET3 />
    <Physician_City />
    <Physician_State />
    <Physician_Zip />
    <Physician_Code>PHY2</Physician_Code>
    <Physician_GROUP>False</Physician_GROUP>
</PhysicianTable>
<PhysicianTable>
    <Physician_ObjectID>14661001</Physician_ObjectID>
    <Physician_Name>Physician, Test</Physician_Name>
    <Physician_Street />
    <Physician_STREET2 />
    <Physician_STREET3 />
    <Physician_City />
    <Physician_State />
    <Physician_Zip />
    <Physician_Code>PHY3</Physician_Code>
    <Physician_GROUP>False</Physician_GROUP>
</PhysicianTable>
<PhysicianTable>
    <Physician_ObjectID>16932034</Physician_ObjectID>
    <Physician_Name>Physician,
TestAgain</Physician_Name>
    <Physician_Street />
    <Physician_STREET2 />
    <Physician_STREET3 />
    <Physician_City />
    <Physician_State />
    <Physician_Zip />
    <Physician_Code>PHY5</Physician_Code>
    <Physician_GROUP>False</Physician_GROUP>
</PhysicianTable>
<PhysicianTable>
    <Physician_ObjectID>135677272</Physician_ObjectID>
    <Physician_Name>Robins, Tim</Physician_Name>
    <Physician_Street />
    <Physician_STREET2 />
    <Physician_STREET3 />

```

```

        <Physician_City />
        <Physician_State />
        <Physician_Zip />
        <Physician_Code>ROBINS</Physician_Code>
        <Physician_GROUP>False</Physician_GROUP>
    </PhysicianTable>
<PhysicianTable>

    <Physician_ObjectID>232073067</Physician_ObjectID>
        <Physician_Name>Robins, Tim Psy</Physician_Name>
        <Physician_Street />
        <Physician_STREET2 />
        <Physician_STREET3 />
        <Physician_City />
        <Physician_State />
        <Physician_Zip />
        <Physician_Code>ROBINS/PSY</Physician_Code>
        <Physician_GROUP>False</Physician_GROUP>
    </PhysicianTable>
<PhysicianTable>

    <Physician_ObjectID>135676995</Physician_ObjectID>
        <Physician_Name>Smith, Tim</Physician_Name>
        <Physician_Street />
        <Physician_STREET2 />
        <Physician_STREET3 />
        <Physician_City />
        <Physician_State />
        <Physician_Zip />
        <Physician_Code>SMITH</Physician_Code>
        <Physician_GROUP>False</Physician_GROUP>
    </PhysicianTable>
<PhysicianTable>

    <Physician_ObjectID>232073138</Physician_ObjectID>
        <Physician_Name>Smith, Tim Psy</Physician_Name>
        <Physician_Street />
        <Physician_STREET2 />
        <Physician_STREET3 />
        <Physician_City />
        <Physician_State />
        <Physician_Zip />
        <Physician_Code>SMITH/PSY</Physician_Code>
        <Physician_GROUP>False</Physician_GROUP>
    </PhysicianTable>
<PhysicianTable>

```



```

<Physician_ObjectID>119430985</Physician_ObjectID>
  <Physician_Name>Doctor, Thursday</Physician_Name>
  <Physician_Street />
  <Physician_STREET2 />
  <Physician_STREET3 />
  <Physician_City />
  <Physician_State />
  <Physician_Zip />
  <Physician_Code>THURS</Physician_Code>
  <Physician_GROUP>False</Physician_GROUP>
</PhysicianTable>

</NewDataSet>

```

## SET STATUS STEP

<b>Endpoint</b>	api/N/ SetStatusStep
<b>Purpose</b>	Set a processing status of a single case, specimen or test order record.
<b>Request Type</b>	GET
<b>Parameters</b>	<pre> object_id={object_id}  status_set={status_set}  status_advance={status_advance}  user_ID={user_ID}  return_format={return_format} </pre>

Where:

**object\_id** is the objectId of the record to update.

**status\_set** is the new status setting.

**status\_advance** is a true or false value. If true, the status of the record to update will be advanced to the next defined step. When false, the status\_advance argument is ignored.

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the status change.

**return\_format** is the format for returning the data. Where:

0 = XML Format

1 = JSON Format

**Example** `https://www.SomeServer.com/api/N/SetStatusStep?object_id=503766463&status_set=&status_advance=true&user_ID=200001&return_format=0`

**Returns** If the request succeeds, returns the results of processing the status update:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <StatusStepUpdateTable>
    <object_id>503766463</object_id>
    <object_type>PCR</object_type>
    <status_set>Purification</status_set>
    <original_status>Accessioned</original_status>
    <message>success</message>
  </StatusStepUpdateTable>
</NewDataSet>
```

## SET STATUS STEPS

**Endpoint** `api/N/SetStatusSteps`

**Purpose** Set a processing status of a multiple case, specimen or test order records.

**Request Type** GET

## Parameters

```
object_ids={object_ids}

status_set={status_set}

status_advance={status_advance}

user_ID={user_ID}

return_format={return_format}
```

Where:

**object\_ids** is a pipe delimited list of object ids that specifies which cases should be updated.

**status\_set** is the new processing status used to update the specified cases.

**status\_advance** is a true or false value. If true, the specified cases will be advanced. If false the specified cases will not be advanced.

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the status change(s).

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

## Example

```
https://www.SomeServer.com/api/N/SetStatusSteps?object_ids=504252001|504165951&status_set=Amplify&status_advance=false&user_ID=200001&return_format=0
```

## Returns

If the request succeeds, returns the results of processing the status updates:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
```

```

<StatusStepUpdateTable>
  <object_id>504252001</object_id>
  <object_type>PCR</object_type>
  <status_set>Amplify</status_set>
  <original_status>Accessioned</original_status>
  <message>success</message>
</StatusStepUpdateTable>
<StatusStepUpdateTable>
  <object_id>504165951</object_id>
  <object_type>PCR</object_type>
  <status_set>Amplify</status_set>
  <original_status>Accessioned</original_status>
  <message>success</message>
</StatusStepUpdateTable>
</NewDataSet>

```

## GET PHYSICIAN LOCATIONS

<b>Endpoint</b>	api/N/GetPhysicianLocations
<b>Purpose</b>	Returns all physician location records for the specified physician, based on a unique physician code.
<b>Request Type</b>	GET
<b>Parameters</b>	physician_code={physician_code} user_ID={user_ID} return_format={return_format}

Where:

**Physician\_code** is the code of the physician record to return.

**user\_id** is the user code for the user in the NucleoLIS user table that is uploading the text.

**return\_format** is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

#### Example

```
https://www.SomeServer.com/api/N/
GetPhysicianLocations?physician_code=00002&user_ID=154
92125&return_format=0
```

#### Returns

If the request succeeds, returns the status of the text that was parsed and processed.:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <PhysicianLocationTable>
    <Location_ObjectID>2291844</Location_ObjectID>
    <Location_CreationDate>2012-10-08T09:58:17.373-
04:00</Location_CreationDate>
    <Location_UpdateTime>2016-08-02T08:47:04.33-
04:00</Location_UpdateTime>
    <Location_Location>191919-
GROUP</Location_Location>
    <Location_Phone />
    <Location_Fax />
    <Location_CODE>00002-2291844</Location_CODE>
    <Location_USER1 />
    <Location_USER2 />
    <Location_USER3 />
    <Location_USER4 />
    <Location_EMAIL />
    <Location_Retired>True</Location_Retired>
    <Group_ObjectID>1709111</Group_ObjectID>
    <Group_CreationDate>2007-05-29T16:17:29-
04:00</Group_CreationDate>
```

```

    <Group_UpdateTime>2016-11-30T07:37:15.85-
05:00</Group_UpdateTime>
    <Group_Name>Physicians Preferred</Group_Name>
    <Group_Street>127 Plain Street</Group_Street>
    <Group_City>Milford</Group_City>
    <Group_State>MA</Group_State>
    <Group_Zip>01568</Group_Zip>
    <Group_Code>PPM</Group_Code>
    <Group_STREET2>APT 2</Group_STREET2>
    <Group_RETIREDD>False</Group_RETIREDD>
  </PhysicianLocationTable>
  <PhysicianLocationTable>
    <Location_ObjectID>1933150</Location_ObjectID>
    <Location_CreationDate>2010-05-27T15:11:59.63-
04:00</Location_CreationDate>
    <Location_UpdateTime>2016-08-02T08:47:02.77-
04:00</Location_UpdateTime>
    <Location_Location>OFFICE</Location_Location>
    <Location_Phone />
    <Location_CODE>1933150-OFFICE</Location_CODE>
    <Location_Retired>False</Location_Retired>
  </PhysicianLocationTable>
</NewDataSet>

```

## GET REPORT

<b>Endpoint</b>	api/N/GetReport
<b>Purpose</b>	Return a single, full case reported result record based on case number and report Id.
<b>Request Type</b>	GET
<b>Parameters</b>	case_number={case_number} report_id={report_id} user_ID={user_ID}

```
return_format={return_format}
```

Where:

**case\_number** is the case number of the report record to return.

**report\_id** is the report object id of the report record to return.

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the report.

**return\_format** is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

#### Example

[https://www.SomeServer.com/api/N/GetReport?](https://www.SomeServer.com/api/N/GetReport?case_number=184&report_id=15386635&user_ID=15492125&return_format=0)  
case\_number=184&report\_id=15386635&user\_ID=15492125&re  
turn\_format=0

#### Returns

If the request succeeds, returns the full case record inclusive of all data elements associated with patient, case, specimen, test order and test result records:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
```

```
<ResultReportTable>
```

```
<Location_AdmitterCity />
```

```
<Location_AdmitterCode />
```

```
<Location_AdmitterName />
```

```
<Location_AdmitterNPI />
```

```
<Location_AdmitterPhone />
```

```

    <Location_AdmitterState />
    <Location_AdmitterStreet />
    <Location_AdmitterStreet2 />
    <Location_AdmitterStreet3 />
    <Location_AdmitterZip />
    <Case_Abnormal>False</Case_Abnormal>
    <Case_AmendedComment />
    <Case_CanceledBy />
    <Case_CanceledDate />
    <Case_CaseType>2327730</Case_CaseType>
    <Case_CaseTypeCode>Cytogenetics</Case_CaseTypeCode>

    <Case_CaseTypeDescription>Cytogenetics</Case_CaseTypeD
escription>

    <Case_ChargeString>123 x 2, BCCHARGE x
2</Case_ChargeString>

    <Case_ClinicalInformation />
    <Case_Comment />
    <Case_CorrectedComment />
    <Case_CreationDate>2018-05-07
09:42:53</Case_CreationDate>
    <Case_CREATOR>johnb</Case_CREATOR>
    <Case_ExternalID />
    <Case_ExternalPatientID />
    <Case_ICDCodesDescriptions />
    <Case_ICDDescriptions />
    <Case_ICDS />
    <Case_Images />
    <Case_MODIFIER>0</Case_MODIFIER>

```



```

<Case_Number>184</Case_Number>
<Case_ObjectID>13767203</Case_ObjectID>
<Case_OrderDate />
<Case_OrderID />
<Case_OtherPhysician1GreetingDisplay />
<Case_OtherPhysician1Name />
<Case_OtherPhysician2GreetingDisplay />
<Case_OtherPhysician2Name />
<Case_OtherPhysician3GreetingDisplay />
<Case_OtherPhysician3Name />
<Case_OtherPhysician4GreetingDisplay />
<Case_OtherPhysician4Name />
<Case_OtherPhysician5GreetingDisplay />
<Case_OtherPhysician5Name />
<Case_OtherPhysicianString />
<Case_PDFs />
<Case_PendingResultMessage />
<Case_PrescribedMedications />
<Case_PriorReportDate>2019-09-
24</Case_PriorReportDate>
<Case_PriorReportTime>11:26</Case_PriorReportTime>
<Case_PriorReportType>FINAL</Case_PriorReportType>
<Case_Questions />
<Case_ReportDate>2019-09-24</Case_ReportDate>
<Case_ReportDefinition />
<Case_ReportedBy>0</Case_ReportedBy>
<Case_ReportedByTitle />

```

```

<Case_ReportedFlag>False</Case_ReportedFlag>
<Case_ReportIndex>2</Case_ReportIndex>
<Case_ReportingComment />
<Case_ReportStatus>FINAL</Case_ReportStatus>
<Case_ReportTime>11:25</Case_ReportTime>
<Case_ReportTitleModifier />
<Case_ReportType>PRELIMINARY</Case_ReportType>
<Case_ReqNumber />
<Case_ReviewedBy />
<Case_ReviewedComment />
<Case_ReviewedDate />
<Case_ReviewedTime />
<Case_ReviewFlag>2019-09-24</Case_ReviewFlag>
<Case_Signer />
<Case_SignerDate />
<Case_SignerLocationAddress />
<Case_SignerLocationAddress2 />
<Case_SignerLocationCity />
<Case_SignerLocationDescription />
<Case_SignerLocationName />

<Case_SignerLocationObjectID>0</Case_SignerLocationObjectID>

<Case_SignerLocationState />
<Case_SignerLocationZip />
<Case_SignerObjectID>0</Case_SignerObjectID>
<Case_SignerTime />
<Case_SignerTitle />

```

```

    <Case_Status>I</Case_Status>

    <Case_StatusStep />

<Case_SubmitterLocation>1437553</Case_SubmitterLocation>

    <Case_SuperobjectID>13767202</Case_SuperobjectID>

    <Case_UpdateTime>2019-09-25
05:06:39</Case_UpdateTime>

    <Case_User1 />
    <Case_User2 />
    <Case_User3 />
    <Case_User4 />

    <Constituent_Abnormal>False</Constituent_Abnormal>

    <Constituent_AlwaysReport />

<Constituent_CategorySequence>2</Constituent_CategorySequence>

    <Constituent_Definition>OFBUPRENORPHINE-
C</Constituent_Definition>

    <Constituent_Description>Buprenorphine
</Constituent_Description>

    <Constituent_ImageName1 />
    <Constituent_ImageName2 />
    <Constituent_IsolateLabel />
    <Constituent_IsolateOrganism />
    <Constituent_IsolateOrganismDescription />
    <Constituent_IsolateOrganismModifier />
    <Constituent_KB />
    <Constituent_MIC />

```

```

<Constituent_ObjectID>15382509</Constituent_ObjectID>

<Constituent_PriorAbnormal>False</Constituent_PriorAbnormal>

    <Constituent_PriorR1>NEG</Constituent_PriorR1>

<Constituent_PriorR2>CONSISTENT</Constituent_PriorR2>

    <Constituent_PriorR3>1.5
ng/ml</Constituent_PriorR3>

    <Constituent_PriorR4>6 to 18
Hours</Constituent_PriorR4>

    <Constituent_PriorR5 />

    <Constituent_PriorR6 />

    <Constituent_PriorRange />

    <Constituent_PriorRangeComment />

<Constituent_PriorResult>0.1</Constituent_PriorResult>

    <Constituent_PriorResultDescription />

    <Constituent_PType>Immuno</Constituent_PType>

    <Constituent_R1>NEG</Constituent_R1>

    <Constituent_R1Desc>Detected</Constituent_R1Desc>

    <Constituent_R2>CONSISTENT</Constituent_R2>

<Constituent_R2Desc>Consistent</Constituent_R2Desc>

    <Constituent_R3>1.5 ng/ml</Constituent_R3>

    <Constituent_R3Desc>Cut Off</Constituent_R3Desc>

    <Constituent_R4>6 to 18 Hours</Constituent_R4>

    <Constituent_R4Desc>Detection
Window</Constituent_R4Desc>

    <Constituent_R5 />

```

```

    <Constituent_R5Desc>Result
Text</Constituent_R5Desc>

    <Constituent_R6 />
    <Constituent_R6Desc />
    <Constituent_Range />
    <Constituent_RangeComment />
    <Constituent_RangeCommentDesc />
    <Constituent_RangeDesc />
    <Constituent_ReportCategory />

<Constituent_ReportCategorySequence>283</Constituent_R
eportCategorySequence>

    <Constituent_Result>0.1</Constituent_Result>
    <Constituent_ResultDate />
    <Constituent_ResultDescription />
    <Constituent_ResultTime />
    <Constituent_Sensitivity />
    <Location_ConsultorCity />
    <Location_ConsultorCode />
    <Location_ConsultorName />
    <Location_ConsultorNPI />
    <Location_ConsultorPhone />
    <Location_ConsultorState />
    <Location_ConsultorStreet />
    <Location_ConsultorStreet2 />
    <Location_ConsultorStreet3 />
    <Location_ConsultorZip />
    <Location_GroupCity />

```

```

    <Location_GroupCode>HARTM</Location_GroupCode>

    <Location_GroupCounty />

    <Location_GroupFax />

    <Location_GroupName>Last,
First</Location_GroupName>

    <Location_GroupPhone />

    <Location_GroupState />

    <Location_GroupStreet />

    <Location_GroupStreet2 />

    <Location_GroupStreet3 />

    <Location_GroupUser1 />

    <Location_GroupZip />

    <Patient_City />

    <Patient_Code>13767197</Patient_Code>

    <Patient_County />

    <Patient_CreationDate>2018-05-07
12:42:08</Patient_CreationDate>

    <Patient_CREATOR>johnb</Patient_CREATOR>

    <Patient_Dob>1970-01-01</Patient_Dob>

    <Patient_Ethnicity />

    <Patient_Fax />

    <Patient_FirstName>J</Patient_FirstName>

    <Patient_LastName>Smith</Patient_LastName>

    <Patient_MedicareNumber />

    <Patient_MiddleName />

    <Patient_MODIFIER>johnb</Patient_MODIFIER>

    <Patient_Name>Smith, J</Patient_Name>

    <Patient_ObjectID>13767197</Patient_ObjectID>

```

```

    <Patient_Phone />
    <Patient_Race />
    <Patient_Sex />
    <Patient_Soundex>SMTHJ</Patient_Soundex>
    <Patient_Ssn>111-11-1111</Patient_Ssn>
    <Patient_State />
    <Patient_Street />
    <Patient_Street2 />
    <Patient_UpdateTime>2018-05-07
09:42:08</Patient_UpdateTime>
    <Patient_USER1 />
    <Patient_USER2 />
    <Patient_USER3 />
    <Patient_USER4 />
    <Patient_Zip />
    <Profile_Abnormal>False</Profile_Abnormal>
    <Profile_AlwaysReport />
    <Profile_ApproveDate />
    <Profile_Approver />

    <Profile_ApproverObjectID>0</Profile_ApproverObjectID>
    <Profile_ApproverTitle />
    <Profile_CompleteDate>2019-09-
24</Profile_CompleteDate>
    <Profile_CompletedBy>Vanger,
John</Profile_CompletedBy>

    <Profile_CompleterObjectID>2325137</Profile_CompleterO
bjectID>

```

```

<Profile_CompleteTime />
<Profile_Definition>OFBASIC</Profile_Definition>
<Profile_DefinitionUser1 />
<Profile_DefinitionUser2 />
<Profile_DefinitionUser3 />
<Profile_Description>OFBASIC (ORAL
FLUID)</Profile_Description>
<Profile_Designation />
<Profile_DNASequenceInterpretation />
<Profile_ImageName1 />
<Profile_ImageName2 />
<Profile_ImageName3 />
<Profile_ImageName4 />
<Profile_ImageName5 />
<Profile_ImageName6 />
<Profile_Images />
<Profile_Interpretation />
<Profile_KaryotypeInterpretation />
<Profile_Methodology />
<Profile_ObjectID>15382508</Profile_ObjectID>
<Profile_OrderDate />
<Profile_OrderTime />
<Profile_PDFs />

<Profile_PriorAbnormal>False</Profile_PriorAbnormal>
<Profile_PriorDNASequenceInterpretation />
<Profile_PriorInterpretation />
<Profile_PriorR1 />

```



```
<Profile_PriorR2 />
<Profile_PriorR3 />
<Profile_PriorR4 />
<Profile_PriorR5 />
<Profile_PriorR6 />
<Profile_PriorRange />
<Profile_PriorRangeComment />
<Profile_PriorResult />
<Profile_PriorResultDescription />
<Profile_PType>Immuno</Profile_PType>
<Profile_Questions />
<Profile_R1 />
<Profile_R1Desc />
<Profile_R2 />
<Profile_R2Desc />
<Profile_R3 />
<Profile_R3Desc />
<Profile_R4 />
<Profile_R4Desc />
<Profile_R5 />
<Profile_R5Desc />
<Profile_R6 />
<Profile_R6Desc />
<Profile_Range />
<Profile_RangeComment />
<Profile_RangeCommentDesc />
<Profile_RangeDesc />
```

```

<Profile_ReportCategory />
<Profile_Result />
<Profile_ResultDate />
<Profile_ResultDescription />
<Profile_ResultTime />
<Profile_Signer />
<Profile_SignerDate />
<Profile_SignerLocationAddress />
<Profile_SignerLocationAddress2 />
<Profile_SignerLocationCity />
<Profile_SignerLocationDescription />
<Profile_SignerLocationName />

<Profile_SignerLocationObjectID>0</Profile_SignerLocationObjectID>

<Profile_SignerLocationState />
<Profile_SignerLocationZip />
<Profile_SignerObjectID>0</Profile_SignerObjectID>
<Profile_SignerTitle />
<Profile_StartDate />
<Profile_StartTime />
<Profile_TestingLocationAddress />
<Profile_TestingLocationAddress2 />
<Profile_TestingLocationCity />
<Profile_TestingLocationDescription />
<Profile_TestingLocationName />

<Profile_TestingLocationObjectID>0</Profile_TestingLocationObjectID>

```

```

    <Profile_TestingLocationState />
    <Profile_TestingLocationZip />
    <Specimen_CollectionDate>2018-06-
29</Specimen_CollectionDate>
    <Specimen_CollectionTime />
    <Specimen_Comment />
    <Specimen_Description />
    <Specimen_Designation />
    <Specimen_DiagnosisFull />

<Specimen_DiagnosisMalignant>False</Specimen_Diagnosis
Malignant>
    <Specimen_FastingStatus />
    <SpecimenFrozenSection_CompleteDate />
    <SpecimenFrozenSection_CompletedBy />
    <SpecimenFrozenSection_CompleteTime />
    <SpecimenFrozenSection_Diagnosis />
    <SpecimenFrozenSection_Gross />
    <SpecimenFrozenSection_ImageName1 />
    <SpecimenFrozenSection_ImageName2 />
    <SpecimenFrozenSection_Images />

<SpecimenFrozenSection_ObjectID>0</SpecimenFrozenSecti
on_ObjectID>
    <SpecimenFrozenSection_PDFs />
    <SpecimenFrozenSection_Signer />
    <SpecimenFrozenSection_SignerDate />
    <SpecimenFrozenSection_SignerTime />
    <SpecimenFrozenSection_SignerTitle />

```

```

<SpecimenFrozenSection_Specimen>0</SpecimenFrozenSection_Specimen>

    <Specimen_Gross />
    <Specimen_ImageName1 />
    <Specimen_ImageName2 />
    <Specimen_Images />
    <Specimen_IndexNumber>1</Specimen_IndexNumber>
    <Specimen_Interpretation />
    <Specimen_Microscopic />
    <Specimen_Number>13770892</Specimen_Number>
    <Specimen_ObjectID>13770892</Specimen_ObjectID>
    <Specimen_OID>13770892</Specimen_OID>
    <Specimen_PatientAge>48</Specimen_PatientAge>
    <Specimen_PDFs />
    <Specimen_PendingResultMessage />
    <Specimen_Priority>ROUTINEx</Specimen_Priority>
    <Specimen_Questions />
    <Specimen_ReceivedDate>2018-06-29</Specimen_ReceivedDate>

    <Specimen_ReceivedTime>05:55:13</Specimen_ReceivedTime>

    <Specimen_Site />

    <Specimen_SourceCode>APPENDIX</Specimen_SourceCode>

    <Specimen_SourceDescription>APPENDIX</Specimen_SourceDescription>

    <Specimen_User1 />

```

```

    <Specimen_User2 />
    <Specimen_User3 />
    <Specimen_User4 />
    <Location_PhysicianCity />

    <Location_PhysicianCode>T124</Location_PhysicianCode>
    <Location_PhysicianCounty />
    <Location_PhysicianDegree />
    <Location_PhysicianFax />
    <Location_PhysicianGreetingDisplay>Earl Gentry
    M.D.</Location_PhysicianGreetingDisplay>
    <Location_PhysicianLocationCity />

    <Location_PhysicianLocationLocation>None</Location_PhysicianLocationLocation>
    <Location_PhysicianLocationPhone />
    <Location_PhysicianLocationState />
    <Location_PhysicianLocationStreet />
    <Location_PhysicianLocationStreet2 />
    <Location_PhysicianLocationStreet3 />
    <Location_PhysicianLocationZip />
    <Location_PhysicianName>Gentry M.D.,
    Earl</Location_PhysicianName>
    <Location_PhysicianNPI />
    <Location_PhysicianPhone />
    <Location_PhysicianState />
    <Location_PhysicianStreet />
    <Location_PhysicianStreet2 />
    <Location_PhysicianStreet3 />

```

```

    <Location_PhysicianUPIN />
    <Location_PhysicianUser1 />
    <Location_PhysicianZip />
    <Visit_ADMITTERFAX />
    <Visit_Code>13767202</Visit_Code>
    <Visit_Comment />
    <Visit_CONSULTORFAX />
    <Visit_CreationDate>2018-05-07
12:42:52</Visit_CreationDate>
    <Visit_DischargeDate />
    <Visit_Location />
    <Visit_MODIFIER>johnb</Visit_MODIFIER>
    <Visit_ObjectID>13767202</Visit_ObjectID>

<Visit_SuperObjectID>13767197</Visit_SuperObjectID>
    <Visit_Type />
    <Visit_UpdateTime>2018-05-07
09:42:52</Visit_UpdateTime>
    <Visit_USER1 />
    <Visit_USER2 />
    <Visit_USER3 />
    <Visit_USER4 />
    <Visit_VisitDate>2018-05-07</Visit_VisitDate>
</ResultReportTable>
</NewDataSet>

```

## GET REPORTS

<b>Endpoint</b>	api/N/GetReports
<b>Purpose</b>	Returns a collection of report records based on filter parameters passed.
<b>Request Type</b>	GET
<b>Parameters</b>	<pre>case_number={case_number} actives_only = {actives_only} report_date_from = {report_date_from} report_date_to = {report_date_to} submitting_physicians = {submitting_physicians} submitting_groups = {submitting_groups} filter_expression = {filter_expression} user_ID={user_ID} return_format={return_format}</pre>

Where:

**case\_number** is a record identifier which has the specified value in the case number. Wildcards can be embedded in the argument as "\*".

**actives\_only** if set to true, only returns active cases. If set to false, returns all cases.

**report\_date\_from** is the oldest inclusive report date in the format: (MMDDYYYY OR YYYY-MM-DD).

**report\_date\_to** is the newest inclusive report date in the format: (MMDDYYYY OR YYYY-MM-DD).

**submitting\_physicians** is a pipe delimited string specifying the physician identifying codes.

**submitting\_groups** is a pipe delimited string specifying the physician facility identifying codes.

**filter\_expression** is a dynamic filter expression which may pertain to any data element returned in the dataset.

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the reports.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

#### Example

[https://www.SomeServer.com/api/N/GetReports?case\\_number=18\\*&activesOnly=false&report\\_date\\_from=2000-01-01&report\\_date\\_to=2021-02-21&submitting\\_physicians=&submitting\\_groups=&filter\\_expression=&user\\_ID=15492125&return\\_format=0](https://www.SomeServer.com/api/N/GetReports?case_number=18*&activesOnly=false&report_date_from=2000-01-01&report_date_to=2021-02-21&submitting_physicians=&submitting_groups=&filter_expression=&user_ID=15492125&return_format=0)

#### Returns

If the request succeeds, returns a collection of Patient, Case and Report level data based on filtering parameters passed including Patient\_ObjectID, Patient\_Name, Patient\_Code, Patient\_Dob, Patient\_Sex, Case\_Number, Case\_OrderDate, Case\_CreationDate, Case\_ReportStatus, Case\_Status, Case\_ObjectID, Report\_ObjectID, Report\_Type, Report\_Date, Report\_Time and Report\_TitleModifier:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>

  <ResultReportTable>

    <Patient_ObjectID>13767197</Patient_ObjectID>

    <Patient_Name>Smith, J</Patient_Name>

    <Patient_Code>13767197</Patient_Code>

    <Patient_Dob>1970-01-01</Patient_Dob>

    <Patient_Sex />

    <Case_Number>183</Case_Number>
```



```

    <Case_OrderDate />
    <Case_CreationDate>2018-05-07T09:42:40.91-
04:00</Case_CreationDate>
    <Case_ReportStatus>FINAL</Case_ReportStatus>
    <Case_Status>I</Case_Status>
    <Case_ObjectID>13767201</Case_ObjectID>
    <Report_ObjectID>13781492</Report_ObjectID>
    <Report_Type>FINAL</Report_Type>
    <Report_Date>2018-09-24</Report_Date>
    <Report_Time>05:57</Report_Time>
</ResultReportTable>
<ResultReportTable>
    <Patient_ObjectID>13767197</Patient_ObjectID>
    <Patient_Name>Smith, J</Patient_Name>
    <Patient_Code>13767197</Patient_Code>
    <Patient_Dob>1970-01-01</Patient_Dob>
    <Patient_Sex />
    <Case_Number>184</Case_Number>
    <Case_OrderDate />
    <Case_CreationDate>2018-05-07T09:42:53.527-
04:00</Case_CreationDate>
    <Case_ReportStatus>FINAL</Case_ReportStatus>
    <Case_Status>I</Case_Status>
    <Case_ObjectID>13767203</Case_ObjectID>
    <Report_ObjectID>15386635</Report_ObjectID>
    <Report_Type>PRELIMINARY</Report_Type>
    <Report_Date>2019-09-24</Report_Date>
    <Report_Time>11:25</Report_Time>

```

```
</ResultReportTable>

</NewDataSet>
```

## GET RESULT REPORT

<b>Endpoint</b>	api/N/GetResultReport
<b>Purpose</b>	Return a PDF result report file based on case number and report Id.
<b>Request Type</b>	GET
<b>Parameters</b>	<pre>case_number={case_number} report_id={report_id} user_ID={user_ID} return_format={return_format}</pre>

Where:

**case\_number** is the case number of the report record to return.

**report\_id** is the report object id of the report record to return.

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the result report.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

**Example** [https://www.SomeServer.com/api/N/GetResultReport?](https://www.SomeServer.com/api/N/GetResultReport?case_number=FSH2020-000002&report_id=15498801&user_ID=15492125&return_format=0)  
case\_number=FSH2020-  
000002&report\_id=15498801&user\_ID=15492125&return\_form  
at=0

**Returns** If the request succeeds, returns a base64 encoded copy of the requested PDF result report:

```
<?xml version="1.0" encoding="UTF-  
8" ?><ReportRecords>  
  
  <Reports>  
  
    <report_id>15498801</report_id>  
  
    <case_number>FSH2020-000002</case_number>  
  
    <base_64>JVBERi0xLjQNCiWxsrO0DQolQ3Jl....o=</base_64>  
  
    <message />  
  
  </Reports>  
  
</ReportRecords>
```

## GET CASE TYPE DEFINITIONS

<b>Endpoint</b>	api/N/GetCaseTypeDefinitions
<b>Purpose</b>	Returns a collection of case type definitions in the system build based on filter parameters passed.
<b>Request Type</b>	GET
<b>Parameters</b>	user_ID={user_ID}  return_format={return_format}  with_retired = {with_retired}

Where:

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the case type definitions.

**return\_format** is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

**with\_retired** If false, will exclude retired entries in the system build. If true, will return all entries in the system build.

#### Example

[https://www.SomeServer.com/api/N/GetCaseTypeDefinitions?user\\_ID=15492125&return\\_format=0&with\\_retired=false](https://www.SomeServer.com/api/N/GetCaseTypeDefinitions?user_ID=15492125&return_format=0&with_retired=false)

#### Returns

If the request succeeds, returns the system unique identifier, code, description and retired state for each item defined in the system build:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
```

```
<CaseTypeDefinitionTable>
```

```
<CaseTypeDefinition_ObjectID>2327075</CaseTypeDefinition_ObjectID>
```

```
<CaseTypeDefinition_Code>GEN</CaseTypeDefinition_Code>
```

```
<CaseTypeDefinition_Description>GENERAL</CaseTypeDefinition_Description>
```

<CaseTypeDefinition\_Retired>False</CaseTypeDefinition\_Retired>

</CaseTypeDefinitionTable>

<CaseTypeDefinitionTable>

<CaseTypeDefinition\_ObjectID>2327641</CaseTypeDefinition\_ObjectID>

<CaseTypeDefinition\_Code>Prenatal</CaseTypeDefinition\_Code>

<CaseTypeDefinition\_Description>Prenatal Testing</CaseTypeDefinition\_Description>

<CaseTypeDefinition\_Retired>False</CaseTypeDefinition\_Retired>

</CaseTypeDefinitionTable>

<CaseTypeDefinitionTable>

<CaseTypeDefinition\_ObjectID>2327727</CaseTypeDefinition\_ObjectID>

<CaseTypeDefinition\_Code>Enzyme</CaseTypeDefinition\_Code>

<CaseTypeDefinition\_Description>Enzyme</CaseTypeDefinition\_Description>

<CaseTypeDefinition\_Retired>False</CaseTypeDefinition\_Retired>

</CaseTypeDefinitionTable>

<CaseTypeDefinitionTable>

<CaseTypeDefinition\_ObjectID>2327728</CaseTypeDefinition\_ObjectID>

```

        <CaseTypeDefinition_Code>DNA
Seq</CaseTypeDefinition_Code>

        <CaseTypeDefinition_Description>DNA
Sequence</CaseTypeDefinition_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_
Retired>

    </CaseTypeDefinitionTable>

    <CaseTypeDefinitionTable>

<CaseTypeDefinition_ObjectID>2327729</CaseTypeDefiniti
on_ObjectID>

<CaseTypeDefinition_Code>Disorder</CaseTypeDefinition_
Code>

        <CaseTypeDefinition_Description>Genetic
Disorders</CaseTypeDefinition_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_
Retired>

    </CaseTypeDefinitionTable>

    <CaseTypeDefinitionTable>

<CaseTypeDefinition_ObjectID>2327730</CaseTypeDefiniti
on_ObjectID>

<CaseTypeDefinition_Code>Cytogenetics</CaseTypeDefinit
ion_Code>

<CaseTypeDefinition_Description>Cytogenetics</CaseType
Definition_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_
Retired>

    </CaseTypeDefinitionTable>

```

```

<CaseTypeDefinitionTable>

<CaseTypeDefinition_ObjectID>2327732</CaseTypeDefiniti
on_ObjectID>

<CaseTypeDefinition_Code>Cytogenetics2</CaseTypeDefini
tion_Code>

<CaseTypeDefinition_Description>FISH</CaseTypeDefiniti
on_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_
Retired>

</CaseTypeDefinitionTable>

<CaseTypeDefinitionTable>

<CaseTypeDefinition_ObjectID>9758103</CaseTypeDefiniti
on_ObjectID>

<CaseTypeDefinition_Code>Cytogenetics_1</CaseTypeDefin
ition_Code>

<CaseTypeDefinition_Description>Karyotype_1</CaseTypeD
efinition_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_
Retired>

</CaseTypeDefinitionTable>

<CaseTypeDefinitionTable>

<CaseTypeDefinition_ObjectID>10866436</CaseTypeDefinit
ion_ObjectID>

<CaseTypeDefinition_Code>pGX</CaseTypeDefinition_Code>

```

```

        <CaseTypeDefinition_Description>TSI
Integration</CaseTypeDefinition_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_
Retired>

    </CaseTypeDefinitionTable>

    <CaseTypeDefinitionTable>

<CaseTypeDefinition_ObjectID>11363917</CaseTypeDefinit
ion_ObjectID>

<CaseTypeDefinition_Code>SURG</CaseTypeDefinition_Code
>

<CaseTypeDefinition_Description>SURGICAL</CaseTypeDefi
nition_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_
Retired>

    </CaseTypeDefinitionTable>

    <CaseTypeDefinitionTable>

<CaseTypeDefinition_ObjectID>13792557</CaseTypeDefinit
ion_ObjectID>

<CaseTypeDefinition_Code>Gastro</CaseTypeDefinition_Co
de>

<CaseTypeDefinition_Description>Gastro</CaseTypeDefini
tion_Description>

<CaseTypeDefinition_Retired>False</CaseTypeDefinition_
Retired>

    </CaseTypeDefinitionTable>

</NewDataSet>

```



## GET SOURCE DEFINITIONS

<b>Endpoint</b>	api/N/GetSourceDefinitions
<b>Purpose</b>	Returns a collection of source definitions in the system build based on filter parameters passed.
<b>Request Type</b>	GET
<b>Parameters</b>	<code>user_ID={user_ID}</code> <code>return_format={return_format}</code> <code>with_retired={with_retired}</code>

Where:

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the source definitions.

**return\_format** is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

**with\_retired** If false, will exclude retired entries in the system build. If true, will return all entries in the system build.

<b>Example</b>	<code>https://www.SomeServer.com/api/N/GetSourceDefinitions?user_ID=15492125&amp;return_format=0&amp;with_retired=false</code>
----------------	--

**Returns**

If the request succeeds, returns the system unique identifier, code, description and retired state for each item defined in the system build:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>

  <SourceDefinitionTable>

    <SourceDefinition_ObjectID>2326047</SourceDefinition_ObjectID>

    <SourceDefinition_Code>APPENDIX</SourceDefinition_Code>

    <SourceDefinition_Description>APPENDIX</SourceDefinition_Description>

    <SourceDefinition_Retired>False</SourceDefinition_Retired>

  </SourceDefinitionTable>

  <SourceDefinitionTable>

    <SourceDefinition_ObjectID>2326116</SourceDefinition_ObjectID>

    <SourceDefinition_Code>BLOOD</SourceDefinition_Code>

    <SourceDefinition_Description>BLOOD</SourceDefinition_Description>

    <SourceDefinition_Retired>False</SourceDefinition_Retired>

  </SourceDefinitionTable>

  <SourceDefinitionTable>
```

<SourceDefinition\_ObjectID>2326606</SourceDefinition\_ObjectID>

<SourceDefinition\_Code>DBlood</SourceDefinition\_Code>

<SourceDefinition\_Description>Dried Blood  
</SourceDefinition\_Description>

<SourceDefinition\_Retired>False</SourceDefinition\_Retired>

</SourceDefinitionTable>

<SourceDefinitionTable>

<SourceDefinition\_ObjectID>2326607</SourceDefinition\_ObjectID>

<SourceDefinition\_Code>Serum</SourceDefinition\_Code>

<SourceDefinition\_Description>Serum</SourceDefinition\_Description>

<SourceDefinition\_Retired>False</SourceDefinition\_Retired>

</SourceDefinitionTable>

<SourceDefinitionTable>

<SourceDefinition\_ObjectID>2326608</SourceDefinition\_ObjectID>

<SourceDefinition\_Code>Urine</SourceDefinition\_Code>

<SourceDefinition\_Description>Urine</SourceDefinition\_Description>

<SourceDefinition\_Retired>False</SourceDefinition\_Retired>

</SourceDefinitionTable>

<SourceDefinitionTable>

<SourceDefinition\_ObjectID>2326609</SourceDefinition\_ObjectID>

<SourceDefinition\_Code>Plasma</SourceDefinition\_Code>

<SourceDefinition\_Description>Plasma</SourceDefinition\_Description>

<SourceDefinition\_Retired>False</SourceDefinition\_Retired>

</SourceDefinitionTable>

<SourceDefinitionTable>

<SourceDefinition\_ObjectID>2326610</SourceDefinition\_ObjectID>

<SourceDefinition\_Code>Leukocytes</SourceDefinition\_Code>

<SourceDefinition\_Description>Leukocytes</SourceDefinition\_Description>

<SourceDefinition\_Retired>False</SourceDefinition\_Retired>

</SourceDefinitionTable>

<SourceDefinitionTable>

<SourceDefinition\_ObjectID>2326611</SourceDefinition\_ObjectID>

```

<SourceDefinition_Code>Fibroblasts</SourceDefinition_Code>

<SourceDefinition_Description>Fibroblasts</SourceDefinition_Description>

<SourceDefinition_Retired>False</SourceDefinition_Retired>

    </SourceDefinitionTable>

</NewDataSet>

```

## GET PROFILE DEFINITIONS

<b>Endpoint</b>	api/N/GetProfileDefinitions
<b>Purpose</b>	Returns a collection of profile definitions in the system build based on filter parameters passed.
<b>Request Type</b>	GET
<b>Parameters</b>	user_ID={user_ID} return_format={return_format} with_retired = {with_retired}

Where:

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the profile definitions.

**return\_format** is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

**with\_retired** If false, will exclude retired entries in the system build. If true, will return all entries in the system build.

**Example** `https://www.SomeServer.com/api/N/GetProfileDefinitions?user_ID=15492125&return_format=0&with_retired=false`

**Returns** If the request succeeds, returns the system unique identifier, code, description and retired state for each item defined in the system build.

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <ProfileDefinitionTable>

    <ProfileDefinition_ObjectID>2325220</ProfileDefinition_ObjectID>

    <ProfileDefinition_Code>BAESO</ProfileDefinition_Code>
      <ProfileDefinition_Description>Barrett's Esophagus
Progression: check for Esophageal
Adenocarcinoma</ProfileDefinition_Description>

    <ProfileDefinition_Retired>False</ProfileDefinition_Retired>

    <ProfileDefinition_LIBRARYTYPE>FISH</ProfileDefinition_LIBRARYTYPE>

  </ProfileDefinitionTable>

  <ProfileDefinitionTable>

    <ProfileDefinition_ObjectID>2335695</ProfileDefinition_ObjectID>
```

```

        <ProfileDefinition_Code>BCL-
6</ProfileDefinition_Code>

        <ProfileDefinition_Description>Break Apart
Probe</ProfileDefinition_Description>

<ProfileDefinition_Retired>False</ProfileDefinition_Re
tired>

<ProfileDefinition_LIBRARYTYPE>FISH</ProfileDefinition
_LIBRARYTYPE>

    </ProfileDefinitionTable>

    <ProfileDefinitionTable>

<ProfileDefinition_ObjectID>2326451</ProfileDefinition
_ObjectID>

<ProfileDefinition_Code>OVCA</ProfileDefinition_Code>

    <ProfileDefinition_Description>Ovarian Cancer
Kit</ProfileDefinition_Description>

<ProfileDefinition_Retired>False</ProfileDefinition_Re
tired>

<ProfileDefinition_LIBRARYTYPE>PCR</ProfileDefinition_
_LIBRARYTYPE>

    </ProfileDefinitionTable>

    <ProfileDefinitionTable>

<ProfileDefinition_ObjectID>2326532</ProfileDefinition
_ObjectID>

<ProfileDefinition_Code>RESPVIRUS</ProfileDefinition_C
ode>

    <ProfileDefinition_Description>Respiratory Virus
Panel</ProfileDefinition_Description>

```

<ProfileDefinition\_Retired>False</ProfileDefinition\_Retired>

<ProfileDefinition\_LIBRARYTYPE>PCR</ProfileDefinition\_LIBRARYTYPE>

</ProfileDefinitionTable>

<ProfileDefinitionTable>

<ProfileDefinition\_ObjectID>2326565</ProfileDefinition\_ObjectID>

<ProfileDefinition\_Code>BL</ProfileDefinition\_Code>

<ProfileDefinition\_Description>Beta-Lactamase  
Detection Comp Gram Neg  
Panel</ProfileDefinition\_Description>

<ProfileDefinition\_Retired>False</ProfileDefinition\_Retired>

<ProfileDefinition\_LIBRARYTYPE>PCR</ProfileDefinition\_LIBRARYTYPE>

</ProfileDefinitionTable>

<ProfileDefinitionTable>

<ProfileDefinition\_ObjectID>2326570</ProfileDefinition\_ObjectID>

<ProfileDefinition\_Code>HPV</ProfileDefinition\_Code>

<ProfileDefinition\_Description>HPV Genotyping-  
TMA</ProfileDefinition\_Description>

<ProfileDefinition\_Retired>False</ProfileDefinition\_Retired>



```

<ProfileDefinition_LIBRARYTYPE>PCR</ProfileDefinition_
LIBRARYTYPE>

  </ProfileDefinitionTable>

</NewDataSet>

```

## GET QUESTION DEFINITIONS

<b>Endpoint</b>	api/N/GetQuestionDefinitions
<b>Purpose</b>	Returns a collection of profile definitions in the system build based on filter parameters passed.
<b>Request Type</b>	GET
<b>Parameters</b>	user_ID={user_ID} return_format={return_format} with_retired={with_retired}

Where:

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the question definitions.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

**with\_retired** If false, will exclude retired entries in the system build. If true, will return all entries in the system build.

**Example** `https://www.SomeServer.com/api/N/GetQuestionDefinitions?user_ID=15492125&return_format=0&with_retired=false`

**Returns** If the request succeeds, returns the system unique identifier, code, description and retired state for each item defined in the system build for each question and response defined:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
<QuestionDefinitionTable>

<AAOEQuestionDefinition_ObjectID>1722188</AAOEQuestion
Definition_ObjectID>

<AAOEQuestionDefinition_Code>CATCH</AAOEQuestionDefini
tion_Code>

    <AAOEQuestionDefinition_Description>Enter catch
method</AAOEQuestionDefinition_Description>

<AAOEQuestionDefinition_RETIREDD>False</AAOEQuestionDef
inition_RETIREDD>

<AAOEQuestionDefinition_AllowOther>False</AAOEQuestion
Definition_AllowOther>

<AAOEQuestionDefinition_SingleSelect>False</AAOEQuesti
onDefinition_SingleSelect>

<AAOEQuestionDefinition_Required>False</AAOEQuestionDe
finition_Required>
```

<AAOEResponseDefinition\_Code>CLEAN</AAOEResponseDefinition\_Code>

<AAOEResponseDefinition\_Description>Clean  
Catch</AAOEResponseDefinition\_Description>

<AAOEResponseDefinition\_RETIREDD>False</AAOEResponseDefinition\_RETIREDD>

</QuestionDefinitionTable>

<QuestionDefinitionTable>

<AAOEQuestionDefinition\_ObjectID>1722188</AAOEQuestionDefinition\_ObjectID>

<AAOEQuestionDefinition\_Code>CATCH</AAOEQuestionDefinition\_Code>

<AAOEQuestionDefinition\_Description>Enter catch  
method</AAOEQuestionDefinition\_Description>

<AAOEQuestionDefinition\_RETIREDD>False</AAOEQuestionDefinition\_RETIREDD>

<AAOEQuestionDefinition\_AllowOther>False</AAOEQuestionDefinition\_AllowOther>

<AAOEQuestionDefinition\_SingleSelect>False</AAOEQuestionDefinition\_SingleSelect>

<AAOEQuestionDefinition\_Required>False</AAOEQuestionDefinition\_Required>

<AAOEResponseDefinition\_Code>RANDOM</AAOEResponseDefinition\_Code>

<AAOEResponseDefinition\_Description>Random</AAOEResponseDefinition\_Description>

<AAOEResponseDefinition\_RETIREDD>False</AAOEResponseDef  
inition\_RETIREDD>

</QuestionDefinitionTable>

<QuestionDefinitionTable>

<AAOEQuestionDefinition\_ObjectID>1722191</AAOEQuestion  
Definition\_ObjectID>

<AAOEQuestionDefinition\_Code>TEXT</AAOEQuestionDefinit  
ion\_Code>

<AAOEQuestionDefinition\_Description>Enter  
Text</AAOEQuestionDefinition\_Description>

<AAOEQuestionDefinition\_RETIREDD>False</AAOEQuestionDef  
inition\_RETIREDD>

</QuestionDefinitionTable>

<QuestionDefinitionTable>

<AAOEQuestionDefinition\_ObjectID>1722239</AAOEQuestion  
Definition\_ObjectID>

<AAOEQuestionDefinition\_Code>HOW</AAOEQuestionDefiniti  
on\_Code>

<AAOEQuestionDefinition\_Description>How are you  
today?</AAOEQuestionDefinition\_Description>

<AAOEQuestionDefinition\_RETIREDD>False</AAOEQuestionDef  
inition\_RETIREDD>

<AAOEQuestionDefinition\_AllowOther>False</AAOEQuestion  
Definition\_AllowOther>

<AAOEQuestionDefinition\_SingleSelect>True</AAOEQuestio  
nDefinition\_SingleSelect>

```

<AAOEQuestionDefinition_Required>True</AAOEQuestionDef
inition_Required>

<AAOEQuestionDefinition_DefaultResponse>Good</AAOEQues
tionDefinition_DefaultResponse>

<AAOEResponseDefinition_Code>GOOD</AAOEResponseDefinit
ion_Code>

<AAOEResponseDefinition_Description>Good</AAOEResponse
Definition_Description>

<AAOEResponseDefinition_RETIRED>False</AAOEResponseDef
inition_RETIRED>

  </QuestionDefinitionTable>

</NewDataSet>

```

## SET INSTRUMENT TEST RESULTS

<b>Endpoint</b>	api/N/SetInstrumentTestResults
<b>Purpose</b>	Sets result values for test orders.
<b>Request Type</b>	GET
<b>Parameters</b>	object_ids={object_ids} at_profile={at_profile} result_value={result_value} r1_value={r1_value} r2_value={r2_value}

```
r3_value={r3_value}

r4_value={r4_value}

r5_value={r5_value}

r6_value={r6_value}range_value is a (mandatory)
String, the a range value to set


rangeComment_value={rangeComment_value}

abnormal_value={abnormal_value}

sent_from={sent_from}

batch_id={batch_id}

user_ID={user_ID}

return_format={return_format}
```

Where:

**object\_ids** is a (mandatory) String, containing a pipe delimited list of unique identifiers of profile records.

**test\_code** is a (mandatory) String, containing profile or constituent code to result.

**at\_profile** is a (mandatory) Boolean, true for profile level result, false for constituent result.

**result\_value** is a (mandatory) String, a result value to set.

**r1\_value** is a (mandatory) String, the r1 result value to set.

**r2\_value** is a (mandatory) String, the r2 result value to set.

**r3\_value** is a (mandatory) String, the r3 result value to set.

**r4\_value** is a (mandatory) String, the r4 result value to set.

**r5\_value** is a (mandatory) String, the r5 result value to set.

**r6\_value** is a (mandatory) String, the r6 result value to set.

**range\_value** is a (mandatory) String, the a range value to set.

**rangeComment\_value** is a (mandatory) String, the result comment value to set.

**abnormal\_value** is a (mandatory) Boolean, true or false signifying if this value is abnormal.

**sent\_from** is (mandatory), the instrument or system in which the data originated

**batch\_id** is a (mandatory) String, the id of the batch associated with the results (if in use)

**user\_id** is a (mandatory) String, the authorized user

**return\_format** is the desired return format (mandatory)  
- Integer, 0 = XML, 1 = JSON

### Example

```
https://www.SomeServer.com/api/N/SetInstrumentTestResults?object_ids=15403525|15494579|15494594&test_code=Escherichia coli
0157&at_profile=false&result_value=Negative&r1_value=&r2_value=&r3_value=&r4_value=&r5_value=&r6_value=&range_value=&rangeComment_value=&abnormal_value=False&sent_from=Luminex&batch_id=Q123456&user_id=15492125&return_format=0
```

### Returns

If the request succeeds, returns the results of processing the result updates:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <InstrumentResultUpdate>
    <object_id>15403525</object_id>
```

```

        <object_type>PCR</object_type>
        <message>success</message>
    </InstrumentResultUpdate>
    <InstrumentResultUpdate>
        <object_id>15494579</object_id>
        <object_type>PCR</object_type>
        <message>success</message>
    </InstrumentResultUpdate>
    <InstrumentResultUpdate>
        <object_id>15494594</object_id>
        <object_type>PCR</object_type>
        <message>success</message>
    </InstrumentResultUpdate>
</NewDataSet>

```

## SET UPLOADED TEXT

<b>Endpoint</b>	api/N/SetUploadedText
<b>Purpose</b>	Sends inbound text and test processor identification for customized parsing and processing.
<b>Request Type</b>	GET
<b>Parameters</b>	processor_id={processor_id} file_text={file_text} user_ID={user_ID} return_format={return_format}



Where:

**processor\_id** is the (mandatory) the id of the text processor to use.

**file\_text** is the (mandatory) text to parse and process.

**user\_id** is the user code for the user in the NucleoLIS user table that is uploading the text.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

#### Example

`https://www.SomeServer.com/api/N/SetUploadedText?processor_id=MyProcessor&file_text=MSH|&user_ID=15492125&return_format=0`

#### Returns

If the request succeeds, returns the status of the text that was parsed and processed.:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <CustomTable>
    <object_type>CustomTable</object_type>
    <message>File Processed</message>
    <status>Success</status>
  </PatientTable>
</NewDataSet>
```

## ADD PATIENT

#### Endpoint

`api/N/Add Patient`

**Purpose** Adds a patient record

**Request Type** GET

**Parameters**

```
patient_lastname={patient_lastname}
patient_firstname={patient_firstname}
patient_middlename={patient_middlename}
patient_street={patient_street}
patient_street2={patient_street2}
patient_city={patient_city}
patient_state={patient_state}
patient_zip={patient_zip}
patient_sex={patient_sex}
patient_dob={patient_dob}
patient_ssn={patient_ssn}
user_ID={user_ID}
return_format={return_format}
```

Where:

**patient\_lastname** is the last name of the patient.

**patient\_firstname** is the first name of the patient.

**patient\_middleName** is the middle name of the patient.

**patient\_street** is the street of the patient.

**patient\_street2** is the street of the patient.

**patient\_city** is the city of the patient.

**patient\_state** is the state of the patient.

**patient\_zip** is the zip code of the patient.

**patient\_sex** is the sex of the patient.

**patient\_dob** is the date of birth of the patient in the format: yyyy-mm-dd.

**patient\_ssn** is the Social Security Number of the patient.

**additional\_fields** is not yet supported. Custom implementation, additional fields to update when in use. PipeDelimited Field^Value

**user\_id** is the user code for the user in the NucleoLIS user table that is adding the patient.

**return\_format** is the format for returning the data. Where:

0 = XML Format

1 = JSON Format

#### Example

```
https://www.SomeServer.com/api/N/
Addpatient?patient_lastname=Smith&patient_firstname=Francisco&patient_middleName=John&patient_street=126
Maple Street&patient_street2=Apt
9&patient_city=Milford&patient_state=MA&patient_zip=01
757&patient_Sex=M&patient_dob=1970-01-
01&patient_ssn=123-456-
7890&additional_fields=&user_ID=15492125&return_format
=0
```

#### Returns

If the request succeeds, returns the status of the patient add request:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <PatientTable>
    <object_type>PatientTable</object_type>
    <message>Patient Added</message>
    <status>Success</status>
  </PatientTable>
</NewDataSet>
```

## UPDATE PATIENT

<b>Endpoint</b>	api/N/UpdatePatient
<b>Purpose</b>	Updates an existing patient record.
<b>Request Type</b>	GET
<b>Parameters</b>	<pre> patient_id={patient_id} patient_lastname={patient_lastname} patient_firstname={patient_firstname} patient_middlename={patient_middlename} patient_street={patient_street} patient_street2={patient_street2} patient_city={patient_city} patient_state={patient_state} patient_zip={patient_zip} patient_sex={patient_sex} patient_dob={patient_dob} patient_ssn={patient_ssn} user_ID={user_ID} return_format={return_format} </pre>

Where:

**patient\_id** is the unique system identifier of the patient to update.

**patient\_lastname** is the last name of the patient.

**patient\_firstname** is the first name of the patient.

**patient\_middleName** is the middle name of the patient.

**patient\_street** is the street of the patient.

**patient\_street2** is the street of the patient.

**patient\_city** is the city of the patient.

**patient\_state** is the state of the patient.

**patient\_zip** is the zip code of the patient.

**patient\_sex** is the sex of the patient.

**patient\_dob** is the date of birth of the patient in the format: yyyy-mm-dd.

**patient\_ssn** is the Social Security Number of the patient.

**additional\_fields** is not yet supported. Custom implementation, additional fields to update when in use. PipeDelimited Field^Value

**user\_id** is the user code for the user in the NucleoLIS user table that is adding the patient.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

### Example

```
https://www.SomeServer.com/api/N/
UpdatePatient?patient_id=15498920&patient_lastname=Smith&patient_firstname=John&patient_middleName=J&patient_street=1 Pine Street&patient_street2=&patient_city=Milford&patient_state=MA&patient_zip=01757&patient_Sex=M&patient_dob=1976-07-01&patient_ssn=123-456-0000&additional_fields=&user_ID=15492125&return_format=0
```

### Returns

If the request succeeds, returns the updated patient identifier and status of processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <PatientTable>
    <object_type>PatientTable</object_type>
    <message>Patient Updated</message>
    <status>Success</status>
  </PatientTable>
```

</NewDataSet>

## ADD CASE

<b>Endpoint</b>	api/N/AddCase
<b>Purpose</b>	Adds a case record.
<b>Request Type</b>	GET
<b>Parameters</b>	<pre>patient_id={patient_id} caseType_code={caseType_code} submittingPhysician_code={submittingPhysician_code} icds={icds} clinical_information={clinical_information} aoe_questions={aoe_questions} additional_fields={additional_fields} user_ID={user_ID} return_format={return_format}</pre>

Where:

**patient\_id** is The unique system identifier of the patient in which to create a case.

**caseType\_code** is The unique code of the case type definition.

**submittingPhysician\_code** is The unique code of the submitting physician.

**submittingLocation** is The location of the submitting physician.

**icds** is Pipe Delimited icd codes.

**clinical\_information** is Clinical Information.

**aoe\_questions** is Not yet supported - Pipe Delimited AOE questions and answers Question^Response.

**additional\_fields** is Not yet supported - Custom implementation, Pipe delimited fields to update when in use. FieldName^FieldValue|

**user\_id** is the user code for the user in the NucleoLIS user table that is adding the case.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

#### Example

```
https://www.SomeServer.com/api/N/
AddCase?patient_id=15498920&caseType_code=GEN&submittingPhysician_code=000000&submittingLocation=OFFICE&icds=001.0|001.1|001.9|002.0&clinical_information=assorted ailments&aoe_questions=&additional_fields=&user_ID=15492125&return_format=0
```

#### Returns

If the request succeeds, returns the system generated case identifier and status of processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <CaseTable>
    <object_type>CaseTable</object_type>
    <message>Case Added</message>
    <status>Success</status>
  </CaseTable>
</NewDataSet>
```

## UPDATE CASE

<b>Endpoint</b>	api/N/UpdateCase
<b>Purpose</b>	Updates a case record.
<b>Request Type</b>	GET
<b>Parameters</b>	<pre> patient_id={patient_id} case_id={case_id} caseType_code={caseType_code} submittingPhysician_code={submittingPhysician_code} icds={icds} clinical_information={clinical_information} aoe_questions={aoe_questions} additional_fields={additional_fields} user_ID={user_ID} return_format={return_format} </pre>

Where:

**patient\_id** is the unique system identifier of the patient in which to create a case.

**case\_id** is the unique system identifier of the case to update.

**caseType\_code** is the unique code of the case type definition.

**submittingPhysician\_code** is the unique code of the submitting physician.

**submittingLocation** is the location of the submitting physician.

**icds** is a list of pipe delimited icd codes.

**clinical\_information** is Clinical Information.



**aoe\_questions** is not yet supported - Pipe Delimited AOE questions and answers Question^Response.

**additional\_fields** is not yet supported - Custom implementation, Pipe delimited fields to update when in use. FieldName^FieldValue|

**user\_id** is the user code for the user in the NucleoLIS user table that is updating the case.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

#### Example

```
https://www.SomeServer.com/api/N/
UpdateCase?patient_id=15498920&case_id=15498922&caseTy
pe_code=GEN&submittingPhysician_code=000000&submitting
Location=000000&icds=001.0&clinical_information=&aoe_q
uestions=&additional_fields=&user_ID=15492125&return_f
ormat=0
```

#### Returns

If the request succeeds, returns the system generated case identifier and status of processing. Note: Case updates may be disallowed depending on the state of case processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <CaseTable>
    <object_type>CaseTable</object_type>
    <message>Case Updated</message>
    <status>Success</status>
  </CaseTable>
</NewDataSet>
```

## ADD SPECIMEN

#### Endpoint

api/N/AddSpecimen

<b>Purpose</b>	Add a specimen record.
<b>Request Type</b>	GET
<b>Parameters</b>	<pre> patient_id={patient_id}  case_id={case_id}  sourceDefinition_code={sourceDefinition_code}  collection_date={collection_date}  collection_time={collection_time}  test_orders={test_orders}  aoe_questions={aoe_questions}  additional_fields={additional_fields}  user_ID={user_ID}  return_format={return_format} </pre>

Where:

**patient\_id** is the unique system identifier of the patient in which to add a specimen.

**case\_id** is the unique system identifier of the case in which to add a specimen.

**sourceDefinition\_code** is the unique system code of the source of the specimen.

**collection\_date** is the collection date in yyyy-mm-dd format.

**collection\_time** is the collection time in hh:mm format.

**test\_orders** is a list of pipe delimited Test Codes to order on the specimen.

**aoe\_questions** is Pipe Delimited AOE questions and answers Question^Response.

**additional\_fields** is Not yet supported - Custom implementation, Pipe delimited fields to update when in use. `FieldName^FieldValue|`

**user\_id** is the user code for the user in the NucleoLIS user table that is updating the case.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

**Example** `https://www.SomeServer.com/api/N/AddSpecimen?patient_id=15498920&case_id=15498922&sourceDefinition_code=BLOOD&collection_date=2021-04-19&collection_time=11:25&test_Orders=ZIKA|XTAG|RPR&aoe_questions=&additional_fields=&user_ID=15492125&return_format=0`

**Returns** If the request succeeds, returns the system generated case identifier and status of processing. Note: Case updates may be disallowed depending on the state of processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <SpecimenTable>
    <object_type>SpecimenTable</object_type>
    <message>Specimen Added</message>
    <status>Success</status>
  </SpecimenTable>
</NewDataSet>
```

## UPDATE SPECIMEN

**Endpoint** `api/N/UpdateSpecimen`

**Purpose** Updates a specimen record.

**Request Type**

GET

**Parameters**

```
patient_id={patient_id}
case_id={case_id}
specimen_id = {specimen_id}
sourceDefinition_code={sourceDefinition_code}
collection_date={collection_date}
collection_time={collection_time}
aoe_questions={aoe_questions}
additional_fields={additional_fields}
user_ID={user_ID}
return_format={return_format}
```

Where:

**patient\_id** is the unique system identifier of the patient in which to add a specimen.

**case\_id** is the unique system identifier of the case in which to add a specimen.

**specimen\_id** is the unique system identifier of the specimen in which to update.

**sourceDefinition\_code** is the unique system code of the source of the specimen.

**collection\_date** is the collection date in yyyy-mm-dd format.

**collection\_time** is the collection time in hh:mm format.

**aoe\_questions** is Pipe Delimited AOE questions and answers Question^Response.

**additional\_fields** is Not yet supported - Custom implementation, Pipe delimited fields to update when in use. FieldName^FieldValue|

**user\_id** is the user code for the user in the NucleoLIS user table that is updating the specimen.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

**Example**

```
https://www.SomeServer.com/api/N/
UpdateSpecimen?patient_id=15498920&case_id=15498922&specimen_id=15498927&sourceDefinition_code=BLOOD&collection_date=2021-04-20&collection_time=11:00&aoe_questions=Last Dose^15mg|Last Does Date^2021-03-19&additional_fields=&user_ID=15492125&return_format=0
```

**Returns**

If the request succeeds, returns the system generated specimen identifier and status of processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <SpecimenTable>
    <object_type>SpecimenTable</object_type>
    <message>Specimen Updated</message>
    <status>Success</status>
  </SpecimenTable>
</NewDataSet>
```

## ADD TEST ORDER

<b>Endpoint</b>	api/N/AddTestOrder
<b>Purpose</b>	Adds a test order.
<b>Request Type</b>	GET
<b>Parameters</b>	patient_id={patient_id}

```
case_id={case_id}

specimen_id = {specimen_id}

profileDefinition_code={profileDefinition_code}

user_ID={user_ID}

return_format={return_format}
```

Where:

**patient\_id** is the unique system identifier of the patient in which to add a specimen.

**case\_id** is the unique system identifier of the case in which to add a test.

**specimen\_id** is the unique system identifier of the specimen in which to add to a test.

**profileDefinition\_code** is the unique system code of the profile to order.

**user\_id** is the user code for the user in the NucleoLIS user table that is adding the test order.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

#### Example

```
https://www.SomeServer.com/api/N/
AddTestOrder?patient_id=15498920&case_id=15498948&specimen_id=15498949&profileDefinition_code=CLL&user_ID=15492125&return_format=0
```

#### Returns

If the request succeeds, returns the status of processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <SpecimenTable>
    <object_type>SpecimenTable</object_type>
    <message>Test Added</message>
    <status>Success</status>
```

```
</SpecimenTable>
</NewDataSet>
```

## CANCEL CASE

<b>Endpoint</b>	api/N/CancelCase
<b>Purpose</b>	Cancels a case record.
<b>Request Type</b>	GET
<b>Parameters</b>	patient_id={patient_id} case_id={case_id} user_ID={user_ID} return_format={return_format}

Where:

**patient\_id** is the unique system identifier of the patient in which to cancel the case.

**case\_id** is the unique system identifier of the case to cancel.

**user\_id** is the user code for the user in the NucleoLIS user table that is canceling the case.

**return\_format** is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

**Example** `https://www.SomeServer.com/api/N/CancelCase?patient_id=15498977&case_id=15498979&user_ID=15492125&return_format=0`

**Returns** If the request succeeds, returns the status of processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <CaseTable>
    <object_type>CaseTable</object_type>
    <message>Case Canceled</message>
    <status>Success</status>
  </CaseTable>
</NewDataSet>
```

## CANCEL SPECIMEN

**Endpoint** `api/N/CancelSpecimen`

**Purpose** Cancels a specimen record.

**Request Type** GET

**Parameters**

```
patient_id={patient_id}
case_id={case_id}
specimen_id={specimen_id}
user_ID={user_ID}
return_format={return_format}
```

Where:

**patient\_id** is the unique system identifier of the patient in which to cancel a specimen.



**case\_id** is the unique system identifier of the case of the specimen to cancel.

**specimen\_id** is the unique system identifier of the specimen to cancel.

**user\_id** is the user code for the user in the NucleoLIS user table that is canceling the specimen.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

**Example** `https://www.SomeServer.com/api/N/CancelSpecimen?patient_id=15498977&case_id=15498979&specimen_id=15498980&user_ID=15492125&return_format=0`

**Returns** If the request succeeds, returns the status of processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <SpecimenTable>
    <object_type>SpecimenTable</object_type>
    <message>Specimen Canceled</message>
    <status>Success</status>
  </SpecimenTable>
</NewDataSet>
```

## CANCEL TEST ORDER

<b>Endpoint</b>	api/N/CancelTestOrder
<b>Purpose</b>	Cancels a test order.
<b>Request Type</b>	GET
<b>Parameters</b>	patient_id={patient_id}

```
case_id={case_id}
specimen_id={specimen_id}
profileDefinition_Code={profileDefinition_Code}
user_ID={user_ID}
return_format={return_format}
```

Where:

**patient\_id** is the unique system identifier of the patient in which to cancel a test.

**case\_id** is the unique system identifier of the case in which to cancel a test.

**specimen\_id** is the unique system identifier of the specimen to cancel a test.

**profileDefinition\_Code** is the unique system code of the profile to cancel.

**user\_id** is the user code for the user in the NucleoLIS user table that is canceling the specimen.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

#### Example

```
https://www.SomeServer.com/api/N/
CancelTestOrder?patient_id=15498977&case_id=15498979&specimen_id=15498980&profileDefinition_Code=ZIKA&user_ID=15492125&return_format=0
```

#### Returns

If the request succeeds, returns the status of processing:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <SpecimenTable>
    <object_type>SpecimenTable</object_type>
    <message>Test Canceled</message>
    <status>Success</status>
```

```
</SpecimenTable>
</NewDataSet>
```

## GET PLATE MAP DEFINITIONS

**Endpoint**            `api/N/GetPlateMapDefinitions`

**Purpose**             Returns a collection of worksheet templates defined as plate maps.

**Request Type**       GET

**Parameters**        `user_ID={user_ID}`  
  
                      `return_format={return_format}`  
  
                      `with_retired={with_retired}`

Where:

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the plate map definitions.

**return\_format** is the format for returning the data.  
Where:

0 = XML Format

1 = JSON Format

**with\_retired** if false will exclude retired entries in the system build. If true, all entries will be returned.

**Example**            `https://www.SomeServer.com/api/N/GetPlateMapDefinitions?user_ID=15492125&return_format=0&with_retired=false`

## Returns

If the request succeeds, returns the system unique identifier, and other descriptive information:

```
<PlateMapDefinitionTable>

<WorksheetTemplate_ObjectID>15394236</WorksheetTemplate_ObjectID>
  <WorksheetTemplate_CreationDate>2019-12-03T06:05:40.123-05:00</WorksheetTemplate_CreationDate>
  <WorksheetTemplate_UpdateTime>2021-04-23T07:06:39.687-04:00</WorksheetTemplate_UpdateTime>
  <WorksheetTemplate_User1 />
  <WorksheetTemplate_User2 />
  <WorksheetTemplate_User3 />
  <WorksheetTemplate_User4 />

<WorksheetTemplate_CREATOR>johnb</WorksheetTemplate_CREATOR>

<WorksheetTemplate_MODIFIER>johnb</WorksheetTemplate_MODIFIER>

<WorksheetTemplate_Code>Extraction</WorksheetTemplate_Code>

<WorksheetTemplate_Description>Extraction</WorksheetTemplate_Description>

<WorksheetTemplate_RETIRED>False</WorksheetTemplate_RETIRED>

<WorksheetTemplate_CurrentVersion>8</WorksheetTemplate_CurrentVersion>

<WorksheetTemplate_Locked>False</WorksheetTemplate_Locked>
  <WorksheetTemplate_LockedBy />
  <WorksheetTemplate_LockedDate />
  <WorksheetTemplate_LockedTime />

<WorksheetTemplate_IsVariant>False</WorksheetTemplate_IsVariant>
  <WorksheetTemplate_VariantMap />
```

<WorksheetTemplate\_HasPlateMap>True</WorksheetTemplate\_HasPlateMap>

<WorksheetTemplate\_PlateMapTabID>Extraction</WorksheetTemplate\_PlateMapTabID>

<WorksheetTemplate\_PlateMapStartRow>2</WorksheetTemplate\_PlateMapStartRow>

<WorksheetTemplate\_PlateMapEndRow>25</WorksheetTemplate\_PlateMapEndRow>

<WorksheetTemplate\_PlateMapPositionColumn>A</WorksheetTemplate\_PlateMapPositionColumn>

<WorksheetTemplate\_PlateMapIDColumn>B</WorksheetTemplate\_PlateMapIDColumn>

<WorksheetTemplate\_PlateMapTypeColumn>C</WorksheetTemplate\_PlateMapTypeColumn>  
</PlateMapDefinitionTable>  
<PlateMapDefinitionTable>

<WorksheetTemplate\_ObjectID>15398875</WorksheetTemplate\_ObjectID>

<WorksheetTemplate\_CreationDate>2020-05-06T10:07:53.493-

04:00</WorksheetTemplate\_CreationDate>

<WorksheetTemplate\_UpdateTime>2021-04-

23T07:15:00.657-04:00</WorksheetTemplate\_UpdateTime>

<WorksheetTemplate\_User1 />

<WorksheetTemplate\_User2 />

<WorksheetTemplate\_User3 />

<WorksheetTemplate\_User4 />

<WorksheetTemplate\_CREATOR>johnb</WorksheetTemplate\_CREATOR>

<WorksheetTemplate\_MODIFIER>johnb</WorksheetTemplate\_MODIFIER>

<WorksheetTemplate\_Code>PCRFILL</WorksheetTemplate\_Code>

<WorksheetTemplate\_Description>PCRFILL</WorksheetTemplate\_Description>

```

<WorksheetTemplate_RETIREDD>False</WorksheetTemplate_RETIREDD>

<WorksheetTemplate_CurrentVersion>3</WorksheetTemplate_CurrentVersion>

<WorksheetTemplate_Locked>False</WorksheetTemplate_Locked>
    <WorksheetTemplate_LockedBy />
    <WorksheetTemplate_LockedDate />
    <WorksheetTemplate_LockedTime />

<WorksheetTemplate_IsVariant>False</WorksheetTemplate_IsVariant>
    <WorksheetTemplate_VariantMap />

<WorksheetTemplate_HasPlateMap>True</WorksheetTemplate_HasPlateMap>

<WorksheetTemplate_PlateMapTabID>Sheet1</WorksheetTemplate_PlateMapTabID>

<WorksheetTemplate_PlateMapStartRow>2</WorksheetTemplate_PlateMapStartRow>

<WorksheetTemplate_PlateMapEndRow>97</WorksheetTemplate_PlateMapEndRow>

<WorksheetTemplate_PlateMapPositionColumn>F</WorksheetTemplate_PlateMapPositionColumn>

<WorksheetTemplate_PlateMapIDColumn>C</WorksheetTemplate_PlateMapIDColumn>

<WorksheetTemplate_PlateMapTypeColumn>D</WorksheetTemplate_PlateMapTypeColumn>
    </PlateMapDefinitionTable>
</NewDataSet>

```

## GET PLATES

<b>Endpoint</b>	api/N/GetPlates
<b>Purpose</b>	Returns a collection of committed plate maps.
<b>Request Type</b>	GET

<b>Parameters</b>	<code>template_code={template_code}</code> <code>min_date={min_date}</code> <code>max_date={max_date}</code> <code>user_ID={user_ID}</code> <code>return_format={return_format}</code>
-------------------	--

Where:

**template\_code** is the system unique code value for the plate map query.

**min\_date** is the earliest commit date to include in the result set, inclusive in the format YYYY-MM-DD.

**max\_date** is the latest commit date to include in the result set, inclusive in the format YYYY-MM-DD.

**profileDefinition\_Code** is the unique system code of the profile to cancel.

**user\_id** is the user code for the user in the NucleoLIS user table that is making the plate map request.

**return\_format** is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

<b>Example</b>	<code>https://www.SomeServer.com/api/N/GetPlates?template_code=PCRFILL&amp;min_date=2021-04-</code>
----------------	---

```
23&max_date=2021-04-  
26&user_ID=15492125&return_format=0
```

## Returns

If the request succeeds, returns the system unique plate identifier, committed date and time, notes and template definition code for plate matching the criteria:

```
<PlateMapTable>  
  <Worksheet_ObjectID>15499389</Worksheet_ObjectID>  
  
  <Worksheet_WorksheetName>15499389</Worksheet_WorksheetName>  
  
  <Worksheet_CommittedBy>johnb</Worksheet_CommittedBy>  
    <Worksheet_CommittedDate>2021-04-  
23</Worksheet_CommittedDate>  
  
  <Worksheet_CommittedTime>07:24</Worksheet_CommittedTime>  
    <Worksheet_Notes />  
  
  <WorksheetTemplate_Code>PCRFILL</WorksheetTemplate_Code>  
  
  <WorksheetTemplate_ObjectID>15398875</WorksheetTemplate_ObjectID>  
    </PlateMapTable>  
</NewDataSet>
```

## GET PLATE MAP

<b>Endpoint</b>	api/N/GetPlateMap
<b>Purpose</b>	Returns the map for the unique identifier passed.
<b>Request Type</b>	GET
<b>Parameters</b>	<pre>platemap_id={platemap_id}  user_ID={user_ID}</pre>



```
return_format={return_format}
```

Where:

**platemap\_id** is the unique identifier of the plate.

**user\_id** is the user code for the user in the NucleoLIS user table that is requesting the specified plate map.

**return\_format** is the format for returning the data.

Where:

0 = XML Format

1 = JSON Format

#### Example

```
https://www.SomeServer.com/api/N/
GetPlateMap?platemad_id=15499389&user_ID=15492125&retu
rn_format=0
```

#### Returns

If the request succeeds, returns the plate map including plate position, specimen identifier, and additional configured plate map related information:

```
<?xml version="1.0" encoding="UTF-8" ?><NewDataSet>
  <PlateMapTable>
    <plate_map_position>A1</plate_map_position>
    <plate_map_id>GEN2020-1000185</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>A2</plate_map_position>
    <plate_map_id>GEN2020-1000185</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>A3</plate_map_position>
    <plate_map_id>GEN2020-1000188</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>A4</plate_map_position>
    <plate_map_id>GEN2020-1000188</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
```

```

<PlateMapTable>
  <plate_map_position>A5</plate_map_position>
  <plate_map_id>GEN2020-1000187</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>A6</plate_map_position>
  <plate_map_id>GEN2020-1000187</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>A7</plate_map_position>
  <plate_map_id>GEN2020-1000191</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>A8</plate_map_position>
  <plate_map_id>GEN2020-1000191</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>B1</plate_map_position>
  <plate_map_id>GEN2020-1000200</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>B2</plate_map_position>
  <plate_map_id>GEN2020-1000200</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>B3</plate_map_position>
  <plate_map_id>GEN2020-1000195</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>B4</plate_map_position>
  <plate_map_id>GEN2020-1000195</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>B5</plate_map_position>
  <plate_map_id>GEN2020-1000201</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
<PlateMapTable>
  <plate_map_position>B6</plate_map_position>

```

```

    <plate_map_id>GEN2020-1000201</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>B7</plate_map_position>
    <plate_map_id>GEN2020-1000202</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>B8</plate_map_position>
    <plate_map_id>GEN2020-1000202</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>C1</plate_map_position>
    <plate_map_id>GEN2020-1000193</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>C2</plate_map_position>
    <plate_map_id>GEN2020-1000193</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>C3</plate_map_position>
    <plate_map_id>GEN2020-1000203</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>C4</plate_map_position>
    <plate_map_id>GEN2020-1000203</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>C5</plate_map_position>
    <plate_map_id>GEN2020-1000196</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>C6</plate_map_position>
    <plate_map_id>GEN2020-1000196</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>
  <PlateMapTable>
    <plate_map_position>C7</plate_map_position>
    <plate_map_id>GEN2020-1000190</plate_map_id>
    <plate_map_type>sample</plate_map_type>
  </PlateMapTable>

```

```

</PlateMapTable>
<PlateMapTable>
  <plate_map_position>C8</plate_map_position>
  <plate_map_id>GEN2020-1000190</plate_map_id>
  <plate_map_type>sample</plate_map_type>
</PlateMapTable>
</NewDataSet>

```

## APPENDIX B – CODING EXAMPLES

### CODING EXAMPLES INTRODUCTION

Note: The coding examples in this appendix assume that the URL being connected to is <https://www.SomeServer.com>, and the username = LAB, and the password = 123456. These examples show how to connect to the Backbone API using a variety of programming languages and techniques. These examples are solely to help you get started communicating with the API, and are not intended to be used as production code.

For more information on the various frameworks mentioned in the code examples, please refer to the links below:

Language	Framework
C	<a href="https://curl.haxx.se/libcurl/c/">https://curl.haxx.se/libcurl/c/</a>
C#	<a href="http://restsharp.org/">http://restsharp.org/</a>
cURL	<a href="https://curl.haxx.se/">https://curl.haxx.se/</a>
Go	<a href="https://golang.org/pkg/net/http/">https://golang.org/pkg/net/http/</a>
HTTP	None (Raw HTTP request)
Java OkHTTP	<a href="https://github.com/square/okhttp">https://github.com/square/okhttp</a>
Java UniREST	<a href="https://github.com/Kong/unirest-java">https://github.com/Kong/unirest-java</a>
JavaScript Fetch	<a href="https://developer.mozilla.org/en-US/docs/Web/API/Fetch_API">https://developer.mozilla.org/en-US/docs/Web/API/Fetch_API</a>

JavaScript jQuery	<a href="http://api.jquery.com/jquery.ajax/">http://api.jquery.com/jquery.ajax/</a>
JavaScript XHR	<a href="https://developer.mozilla.org/en-US/docs/Web/API/XMLHttpRequest">https://developer.mozilla.org/en-US/docs/Web/API/XMLHttpRequest</a>
NodeJS Axios	<a href="https://github.com/axios/axios">https://github.com/axios/axios</a>
NodeJS Http	<a href="https://nodejs.org/api/http.html">https://nodejs.org/api/http.html</a>
NodeJS Request	<a href="https://github.com/request/request">https://github.com/request/request</a>
NodeJS UniRest	<a href="https://github.com/Kong/unirest-nodejs">https://github.com/Kong/unirest-nodejs</a>
Objective-C	<a href="https://developer.apple.com/library/ios/documentation/Foundation/Reference/NSURLSession_class/">https://developer.apple.com/library/ios/documentation/Foundation/Reference/NSURLSession_class/</a>
OCaml	<a href="https://github.com/mirage/ocaml-cohttp">https://github.com/mirage/ocaml-cohttp</a>
PHP CURL	<a href="http://php.net/manual/en/ref.curl.php">http://php.net/manual/en/ref.curl.php</a>
PHP HttpRequest	<a href="https://www.php.net/manual/en/reserved.variables.request.php">https://www.php.net/manual/en/reserved.variables.request.php</a>
PHP pecl_http	<a href="https://mdref.m6w6.name/http">https://mdref.m6w6.name/http</a>
PowerShell RestMethod	<a href="https://docs.microsoft.com/en-us/powershell/module/microsoft.powershell.utility/invoke-restmethod?view=powershell-7">https://docs.microsoft.com/en-us/powershell/module/microsoft.powershell.utility/invoke-restmethod?view=powershell-7</a>
Python Built-in http.client (Python 3)	<a href="https://docs.python.org/3/library/http.client.html">https://docs.python.org/3/library/http.client.html</a>
Python Requests	<a href="https://requests.readthedocs.io/en/master/">https://requests.readthedocs.io/en/master/</a>
Ruby Built-in NET::Http	<a href="http://docs.ruby-lang.org/en/2.0.0/Net/HTTP.html">http://docs.ruby-lang.org/en/2.0.0/Net/HTTP.html</a>
Shell HTTPie	<a href="https://github.com/jkbrzt/httpie">https://github.com/jkbrzt/httpie</a>
Shell wget	<a href="https://www.gnu.org/software/wget/">https://www.gnu.org/software/wget/</a>
Shell cURL	<a href="https://curl.haxx.se/">https://curl.haxx.se/</a>

Swift Built-in <https://developer.apple.com/documentation/foundation/urlsession>  
NSURLSession

VB.NET <http://restsharp.org/>  
RestSharp

## LOGON CODING EXAMPLES

Note: The Logon endpoint is a POST command. Any API call that uses POST would follow a similar pattern.

### C – LIBCURL

```
CURL *curl;
CURLcode res;
curl = curl_easy_init();
if(curl) {
    curl_easy_setopt(curl, CURLOPT_CUSTOMREQUEST, "POST");
    curl_easy_setopt(curl, CURLOPT_URL, "https://www.SomeServer.com/api/authenticate/logon");
    curl_easy_setopt(curl, CURLOPT_FOLLOWLOCATION, 1L);
    curl_easy_setopt(curl, CURLOPT_DEFAULT_PROTOCOL, "https");
    struct curl_slist *headers = NULL;
    headers = curl_slist_append(headers, "Content-Type: application/x-www-form-urlencoded");
    curl_easy_setopt(curl, CURLOPT_HTTPHEADER, headers);
    const char *data = "username=LAB&password=123456";
    curl_easy_setopt(curl, CURLOPT_POSTFIELDS, data);
    res = curl_easy_perform(curl);
}
curl_easy_cleanup(curl);
```

### C# - RESTSHARP

```
var client = new RestClient("https://www.SomeServer.com/api/authenticate/logon");
client.Timeout = -1;
var request = new RestRequest(Method.POST);
request.AddHeader("Content-Type", "application/x-www-form-urlencoded");
request.AddParameter("username", "LAB");
request.AddParameter("password", "123456");
IRestResponse response = client.Execute(request);
Console.WriteLine(response.Content);
```

## CURL

```
curl --location --request POST 'https://www.SomeServer.com/api/authenticate/logon' \  
--header 'Content-Type: application/x-www-form-urlencoded' \  
--data-urlencode 'username=LAB' \  
--data-urlencode 'password=123456'
```

## GO - NATIVE

```
package main  
  
import (  
    "fmt"  
    "strings"  
    "net/http"  
    "io/ioutil"  
)  
  
func main() {  
  
    url := "https://www.SomeServer.com/api/authenticate/logon"  
    method := "POST"  
  
    payload := strings.NewReader("username=LAB&password=123456")  
  
    client := &http.Client {  
    }  
    req, err := http.NewRequest(method, url, payload)  
  
    if err != nil {  
        fmt.Println(err)  
        return  
    }  
    req.Header.Add("Content-Type", "application/x-www-form-urlencoded")  
  
    res, err := client.Do(req)  
    if err != nil {  
        fmt.Println(err)  
        return  
    }  
    defer res.Body.Close()  
  
    body, err := ioutil.ReadAll(res.Body)  
    if err != nil {  
        fmt.Println(err)  
        return  
    }  
}
```

```

    }
    fmt.Println(string(body))
}

```

## HTTP

POST /WebApiBackboneQA/api/authenticate/logon? HTTP/1.1

Host: SomeServer.com

Content-Type: application/x-www-form-urlencoded

Content-Length: 28

username=LAB&password=123456

## JAVA - OKHTTP

```

OkHttpClient client = new OkHttpClient().newBuilder()
    .build();
MediaType mediaType = MediaType.parse("application/x-www-form-urlencoded");
RequestBody body = RequestBody.create(mediaType, "username=LAB&password=123456");
Request request = new Request.Builder()
    .url("https://www.SomeServer.com/api/authenticate/logon")
    .method("POST", body)
    .addHeader("Content-Type", "application/x-www-form-urlencoded")
    .build();
Response response = client.newCall(request).execute();

```

## JAVA – UNIREST

```

Unirest.setTimeouts(0, 0);
HttpResponse<String> response = Unirest.post("https://www.SomeServer.com/api/authenticate/logo")
    .header("Content-Type", "application/x-www-form-urlencoded")
    .field("username", "LAB")
    .field("password", "123456")
    .asString();

```

## JAVASCRIPT - FETCH

```

var myHeaders = new Headers();

```



```

myHeaders.append("Content-Type", "application/x-www-form-urlencoded");

var urlencoded = new URLSearchParams();
urlencoded.append("username", "LAB");
urlencoded.append("password", "123456");

var requestOptions = {
  method: 'POST',
  headers: myHeaders,
  body: urlencoded,
  redirect: 'follow'
};

fetch("https://www.SomeServer.com/api/authenticate/logon", requestOptions)
  .then(response => response.text())
  .then(result => console.log(result))
  .catch(error => console.log('error', error));

```

## JAVASCRIPT - JQUERY

```

var settings = {
  "url": "https://www.SomeServer.com/api/authenticate/logon",
  "method": "POST",
  "timeout": 0,
  "headers": {
    "Content-Type": "application/x-www-form-urlencoded",
  },
  "data": {
    "username": "LAB",
    "password": "123456"
  }
};

$.ajax(settings).done(function (response) {
  console.log(response);
});

```

## JAVASCRIPT – XHR

```

var data = "username=LAB&password=123456";

var xhr = new XMLHttpRequest();
xhr.withCredentials = true;

xhr.addEventListener("readystatechange", function() {
  if(this.readyState === 4) {

```

```

        console.log(this.responseText);
    }
});

xhr.open("POST", "https://www.SomeServer.com/api/authenticate/logon");
xhr.setRequestHeader("Content-Type", "application/x-www-form-urlencoded");

xhr.send(data);

```

## NODEJS – AXIOS

```

var axios = require('axios');
var qs = require('qs');
var data = qs.stringify({
    'username': 'LAB',
    'password': '123456'
});
var config = {
    method: 'post',
    url: 'https://www.SomeServer.com/api/authenticate/logon',
    headers: {
        'Content-Type': 'application/x-www-form-urlencoded',
    },
    data : data
};

axios(config)
    .then(function (response) {
        console.log(JSON.stringify(response.data));
    })
    .catch(function (error) {
        console.log(error);
    });

```

## NODEJS – NATIVE

```

var https = require('follow-redirects').https;
var fs = require('fs');

var qs = require('querystring');

var options = {
    'method': 'POST',
    'hostname': 'SomeServer.com',
    'path': '/WebApiBackboneQA/api/authenticate/logon?',

```

```

    'headers': {
      'Content-Type': 'application/x-www-form-urlencoded',
    },
    'maxRedirects': 20
  };

  var req = https.request(options, function (res) {
    var chunks = [];

    res.on("data", function (chunk) {
      chunks.push(chunk);
    });

    res.on("end", function (chunk) {
      var body = Buffer.concat(chunks);
      console.log(body.toString());
    });

    res.on("error", function (error) {
      console.error(error);
    });
  });

  var postData = qs.stringify({
    'username': 'LAB',
    'password': '123456'
  });

  req.write(postData);

  req.end();

```

## NODEJS – REQUEST

```

var request = require('request');
var options = {
  'method': 'POST',
  'url': 'https://www.SomeServer.com/api/authenticate/logon',
  'headers': {
    'Content-Type': 'application/x-www-form-urlencoded',
  },
  form: {
    'username': 'LAB',
    'password': '123456'
  }
};
request(options, function (error, response) {

```

```

    if (error) throw new Error(error);
    console.log(response.body);
  });

```

## NODEJS – UNIREST

```

var unirest = require('unirest');
var req = unirest('POST', 'https://www.SomeServer.com/api/authenticate/logon')
  .headers({
    'Content-Type': 'application/x-www-form-urlencoded',
  })
  .send('username=LAB')
  .send('password=123456')
  .end(function (res) {
    if (res.error) throw new Error(res.error);
    console.log(res.raw_body);
  });

```

## OBJECTIVE-C NSURLSESSION

```

#import <Foundation/Foundation.h>

```

```

dispatch_semaphore_t sema = dispatch_semaphore_create(0);

```

```

NSMutableURLRequest *request = [NSMutableURLRequest requestWithURL:[NSURL URLWithString:@"https://www.SomeServer.com/api/authenticate/logon"]
  cachePolicy:NSURLRequestUseProtocolCachePolicy
  timeoutInterval:10.0];
NSDictionary *headers = @{
  @"Content-Type": @"application/x-www-form-urlencoded",
};

```

```

[request setAllHTTPHeaderFields:headers];
NSMutableData *postData = [[NSMutableData alloc] initWithData:@"username=LAB" dataUsingEncoding:NSUTF8StringEncoding];
[postData appendData:@"&password=123456" dataUsingEncoding:NSUTF8StringEncoding];
[request setHTTPBody:postData];

```

```

[request setHTTPMethod:@"POST"];

```

```

NSURLSession *session = [NSURLSession sharedSession];
NSURLSessionDataTask *dataTask = [session dataTaskWithRequest:request
  completionHandler:^(NSData *data, NSURLResponse *response, NSError *error) {
  if (error) {
    NSLog(@"%@", error);
  }
}];

```

```

    } else {
        NSHTTPURLResponse *httpResponse = (NSHTTPURLResponse *) response;
        NSError *parseError = nil;
        NSDictionary *responseDictionary = [NSJSONSerialization JSONObjectWithData:data options:0
error:&parseError];
        NSLog(@"%@",responseDictionary);
        dispatch_semaphore_signal(sema);
    }
}];
[dataTask resume];
dispatch_semaphore_wait(sema, DISPATCH_TIME_FOREVER);

```

## OCAML – COHTTP

```

open Lwt
open Cohttp
open Cohttp_lwt_unix

let postData = ref "username=LAB&password=123456";;

let reqBody =
    let uri = Uri.of_string "https://www.SomeServer.com/api/authenticate/logon" in
    let headers = Header.init ()
    |> fun h -> Header.add h "Content-Type" "application/x-www-form-urlencoded"
    in
    let body = Cohttp_lwt.Body.of_string !postData in

    Client.call ~headers ~body `POST uri >>= fun (_resp, body) ->
    body |> Cohttp_lwt.Body.to_string >|= fun body -> body

let () =
    let respBody = Lwt_main.run reqBody in
    print_endline (respBody)

```

## PHP – CURL

```

<?php

$curl = curl_init();

curl_setopt_array($curl, array(
    CURLOPT_URL => "https://www.SomeServer.com/api/authenticate/logon",
    CURLOPT_RETURNTRANSFER => true,
    CURLOPT_ENCODING => "",
    CURLOPT_MAXREDIRS => 10,
    CURLOPT_TIMEOUT => 0,

```

```

CURLOPT_FOLLOWLOCATION => true,
CURLOPT_HTTP_VERSION => CURL_HTTP_VERSION_1_1,
CURLOPT_CUSTOMREQUEST => "POST",
CURLOPT_POSTFIELDS => "username=LAB&password=123456",
CURLOPT_HTTPHEADER => array(
    "Content-Type: application/x-www-form-urlencoded",
),
));

$response = curl_exec($curl);

curl_close($curl);
echo $response;

```

## PHP – HTTP\_REQUEST2

```

<?php
require_once 'HTTP/Request2.php';
$request = new HTTP_Request2();
$request->setUrl('https://www.SomeServer.com/api/authenticate/login');
$request->setMethod(HTTP_Request2::METHOD_POST);
$request->setConfig(array(
    'follow_redirects' => TRUE
));
$request->setHeader(array(
    'Content-Type' => 'application/x-www-form-urlencoded',
));
$request->addPostParameter(array(
    'username' => 'LAB',
    'password' => '123456'
));
try {
    $response = $request->send();
    if ($response->getStatus() == 200) {
        echo $response->getBody();
    }
    else {
        echo 'Unexpected HTTP status: ' . $response->getStatus() . ' ' .
            $response->getReasonPhrase();
    }
}
catch(HTTP_Request2_Exception $e) {
    echo 'Error: ' . $e->getMessage();
}

```

## PHP – PECL\_HTTP

```

<?php
$client = new http\Client;
$request = new http\Client\Request;
$request->setRequestUrl('https://www.SomeServer.com/api/authenticate/logon');
$request->setRequestMethod('POST');
$body = new http\Message\Body;
$body->append(new http\QueryString(array(
    'username' => 'LAB',
    'password' => '123456'))); $request->setBody($body);
$request->setOptions(array());
$request->setHeaders(array(
    'Content-Type' => 'application/x-www-form-urlencoded',
));
$client->enqueue($request)->send();
$response = $client->getResponse();
echo $response->getBody();

```

## POWERSHELL – RESTMETHOD

```

$headers = New-Object "System.Collections.Generic.Dictionary[[String],[String]]"
$headers.Add("Content-Type", "application/x-www-form-urlencoded")

$body = "username=LAB&password=123456"

$response = Invoke-RestMethod 'https://www.SomeServer.com/api/authenticate/logon' -
Method 'POST' -Headers $headers -Body $body

```

## PYTHON – HTTP.CLIENT

```

import http.client
import mimetypes
conn = http.client.HTTPSConnection("SomeServer.com")
payload = 'username=LAB&password=123456'
headers = {
    'Content-Type': 'application/x-www-form-urlencoded',
}
conn.request("POST", "/WebApiBackboneQA/api/authenticate/logon", payload, headers)
res = conn.getresponse()
data = res.read()
print(data.decode("utf-8"))

```

## PYTHON – REQUESTS

```

import requests

url = "https://www.SomeServer.com/api/authenticate/logon"

```

```

payload='username=LAB&password=123456'
headers = {
    'Content-Type': 'application/x-www-form-urlencoded',
}

response = requests.request("POST", url, headers=headers, data=payload)

print(response.text)

```

*Note: For an example Python integration with the API, please go to the following link: <https://github.com/marcuslehr/NucleoLIS-API>*

## RUBY – NET::HTTP

```

require "uri"
require "net/http"

url = URI("https://www.SomeServer.com/api/authenticate/logon")

https = Net::HTTP.new(url.host, url.port)
https.use_ssl = true

request = Net::HTTP::Post.new(url)
request["Content-Type"] = "application/x-www-form-urlencoded"
request.body = "username=LAB&password=123456"

response = https.request(request)
puts response.read_body

```

## SHELL – HTTPIE

```

http --ignore-stdin --form --follow --
timeout 3600 POST https://www.SomeServer.com/api/authenticate/logon \
'username'='LAB' \
'password'='123456' \
Content-Type: 'application/x-www-form-urlencoded' \

```

## SHELL – WGET

```

wget --no-check-certificate --quiet \
--method POST \
--timeout=0 \
--header 'Content-Type: application/x-www-form-urlencoded' \
--body-data 'username=LAB&password=123456' \
'https://www.SomeServer.com/api/authenticate/logon'

```



## SWIFT - URLSESSION

```
import Foundation

var semaphore = DispatchSemaphore (value: 0)

let parameters = "username=LAB&password=123456"
let postData = parameters.data(using: .utf8)

var request = URLRequest(url: URL(string: "https://Www.SomeServer.com/api/authenticate/logon")
!,timeoutInterval: Double.infinity)
request.addValue("application/x-www-form-urlencoded", forHTTPHeaderField: "Content-Type")

request.httpMethod = "POST"
request.httpBody = postData

let task = URLSession.shared.dataTask(with: request) { data, response, error in
    guard let data = data else {
        print(String(describing: error))
        semaphore.signal()
        return
    }
    print(String(data: data, encoding: .utf8)!)
    semaphore.signal()
}

task.resume()
semaphore.wait()
```

## VB.NET - RESTSHARP

```
Dim client = New
RestClient("https://Www.SomeServer.com/api/authenticate/logon")
client.Timeout = -1
Dim request = New RestRequest(Method.POST)
request.AddHeader("Content-Type", "application/x-www-form-urlencoded")
request.AddParameter("username", "LAB")
request.AddParameter("password", "123456")
Dim response As IRestResponse = client.Execute(request)
Console.WriteLine(response.Content)
```

## LOGOUT CODING EXAMPLES

Note: The Logout endpoint is a GET command. Any API call that uses GET would follow a similar pattern.

## C - LIBCURL

```

CURL *curl;
CURLcode res;
curl = curl_easy_init();
if(curl) {
    curl_easy_setopt(curl, CURLOPT_CUSTOMREQUEST, "GET");
    curl_easy_setopt(curl, CURLOPT_URL, "https://www.SomeServer.com/api/authenticate/logout");
    curl_easy_setopt(curl, CURLOPT_FOLLOWLOCATION, 1L);
    curl_easy_setopt(curl, CURLOPT_DEFAULT_PROTOCOL, "https");
    struct curl_slist *headers = NULL;
    headers = curl_slist_append(headers, "Cookie: session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900");
    curl_easy_setopt(curl, CURLOPT_HTTPHEADER, headers);
    res = curl_easy_perform(curl);
}
curl_easy_cleanup(curl);

```

## C# - RESTSHARP

```

var client = new RestClient("https://www.SomeServer.com/api/authenticate/logout");
client.Timeout = -1;
var request = new RestRequest(Method.GET);
request.AddHeader("Cookie", "session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900");
IRestResponse response = client.Execute(request);
Console.WriteLine(response.Content);

```

## CURL

```

curl --location --request GET 'https://www.SomeServer.com/api/authenticate/logout' \
--header 'Cookie: session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'

```

## GO – NATIVE

```

package main

import (
    "fmt"
    "net/http"
    "io/ioutil"
)

func main() {

```

```

url := "https://www.SomeServer.com/api/authenticate/logout"
method := "GET"

client := &http.Client {
}
req, err := http.NewRequest(method, url, nil)

if err != nil {
    fmt.Println(err)
    return
}
req.Header.Add("Cookie", "session=session-
id=kUn8axg5aDzfBOuVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900")

res, err := client.Do(req)
if err != nil {
    fmt.Println(err)
    return
}
defer res.Body.Close()

body, err := ioutil.ReadAll(res.Body)
if err != nil {
    fmt.Println(err)
    return
}
fmt.Println(string(body))
}

```

## HTTP

GET /WebApiBackboneQA/api/authenticate/logout HTTP/1.1

Host: SomeServer.com

Cookie: session=session-

id=kUn8axg5aDzfBOuVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelimit=900

## JAVA – OKHTTP

```

OkHttpClient client = new OkHttpClient().newBuilder()
    .build();
Request request = new Request.Builder()

```

```

        .url("https://www.SomeServer.com/api/authenticate/logout")
        .method("GET", null)
        .addHeader("Cookie", "session=session-id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSqt2mFq000%3d&username=346650&epoch=1604350665&timelimit=900")
        .build();
Response response = client.newCall(request).execute();

```

## JAVA – UNIREST

```

Unirest.setTimeouts(0, 0);
HttpResponse<String> response = Unirest.get("https://www.SomeServer.com/api/authenticate/logout")
    .header("Cookie", "session=session-id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSqt2mFq000%3d&username=346650&epoch=1604350665&timelimit=900")
    .asString();

```

## JAVASCRIPT – FETCH

```

var myHeaders = new Headers();
myHeaders.append("Cookie", "session=session-id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSqt2mFq000%3d&username=346650&epoch=1604350665&timelimit=900");

var requestOptions = {
  method: 'GET',
  headers: myHeaders,
  redirect: 'follow'
};

fetch("https://www.SomeServer.com/api/authenticate/logout", requestOptions)
  .then(response => response.text())
  .then(result => console.log(result))
  .catch(error => console.log('error', error));

```

## JAVASCRIPT – JQUERY

```

var settings = {
  "url": "https://www.SomeServer.com/api/authenticate/logout",
  "method": "GET",
  "timeout": 0,
  "headers": {

```

```

        "Cookie": "session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900"
    },
};

$.ajax(settings).done(function (response) {
    console.log(response);
});

```

## JAVASCRIPT – XHR

```

var xhr = new XMLHttpRequest();
xhr.withCredentials = true;

xhr.addEventListener("readystatechange", function() {
    if(this.readyState === 4) {
        console.log(this.responseText);
    }
});

xhr.open("GET", "https://www.SomeServer.com/api/authenticate/logout");
xhr.setRequestHeader("Cookie", "session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900");

xhr.send();

```

## NODEJS – AXIOS

```

var axios = require('axios');

var config = {
    method: 'get',
    url: 'https://www.SomeServer.com/api/authenticate/logout',
    headers: {
        'Cookie': 'session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'
    }
};

axios(config)
.then(function (response) {
    console.log(JSON.stringify(response.data));
});

```

```

})
.catch(function (error) {
  console.log(error);
});

```

## NODEJS – NATIVE

```

var https = require('follow-redirects').https;
var fs = require('fs');

var options = {
  'method': 'GET',
  'hostname': 'SomeServer.com',
  'path': '/WebApiBackboneQA/api/authenticate/logout',
  'headers': {
    'Cookie': 'session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSqt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'
  },
  'maxRedirects': 20
};

var req = https.request(options, function (res) {
  var chunks = [];

  res.on("data", function (chunk) {
    chunks.push(chunk);
  });

  res.on("end", function (chunk) {
    var body = Buffer.concat(chunks);
    console.log(body.toString());
  });

  res.on("error", function (error) {
    console.error(error);
  });
});

req.end();

```

## NODEJS – REQUEST

```

var request = require('request');
var options = {

```

```

    'method': 'GET',
    'url': 'https://www.SomeServer.com/api/authenticate/logout',
    'headers': {
        'Cookie': 'session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'
    }
};
request(options, function (error, response) {
    if (error) throw new Error(error);
    console.log(response.body);
});

```

## NODEJS – UNIREST

```

var unirest = require('unirest');
var req = unirest('GET', 'https://www.SomeServer.com/api/authenticate/logout')
    .headers({
        'Cookie': 'session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'
    })
    .end(function (res) {
        if (res.error) throw new Error(res.error);
        console.log(res.raw_body);
    });

```

## OBJECTIVE-C – NSURLSESSION

```
#import <Foundation/Foundation.h>
```

```
dispatch_semaphore_t sema = dispatch_semaphore_create(0);
```

```

NSMutableURLRequest *request = [NSMutableURLRequest requestWithURL:[NSURL URLWithString:@"http
s://www.SomeServer.com/api/authenticate/logout"]
    cachePolicy:NSURLRequestUseProtocolCachePolicy
    timeoutInterval:10.0];
NSDictionary *headers = @{
    @"Cookie": @"session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900"
};

```

```
[request setAllHTTPHeaderFields:headers];
```

```

[request setHTTPMethod:@"GET"];

NSURLSession *session = [NSURLSession sharedSession];
NSURLSessionDataTask *dataTask = [session dataTaskWithRequest:request
completionHandler:^(NSData *data, NSURLResponse *response, NSError *error) {
    if (error) {
        NSLog(@"%@", error);
    } else {
        NSHTTPURLResponse *httpResponse = (NSHTTPURLResponse *) response;
        NSError *parseError = nil;
        NSDictionary *responseDictionary = [NSJSONSerialization JSONObjectWithData:data options:0
error:&parseError];
        NSLog(@"%@", responseDictionary);
        dispatch_semaphore_signal(sema);
    }
}];
[dataTask resume];
dispatch_semaphore_wait(sema, DISPATCH_TIME_FOREVER);

```

## OCAML – COHTTP

```

open Lwt
open Cohttp
open Cohttp_lwt_unix

let reqBody =
    let uri = Uri.of_string "https://www.SomeServer.com/api/authenticate/logout" in
    let headers = Header.init ()
    |> fun h -> Header.add h "Cookie" "session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSqt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900"
    in
    Client.call ~headers `GET uri >>= fun (_resp, body) ->
    body |> Cohttp_lwt.Body.to_string >|= fun body -> body

let () =
    let respBody = Lwt_main.run reqBody in
    print_endline (respBody)

```

## PHP – CURL

```

<?php

$curl = curl_init();

curl_setopt_array($curl, array(

```



```

CURLOPT_URL => "https://www.SomeServer.com/api/authenticate/logout",
CURLOPT_RETURNTRANSFER => true,
CURLOPT_ENCODING => "",
CURLOPT_MAXREDIRS => 10,
CURLOPT_TIMEOUT => 0,
CURLOPT_FOLLOWLOCATION => true,
CURLOPT_HTTP_VERSION => CURL_HTTP_VERSION_1_1,
CURLOPT_CUSTOMREQUEST => "GET",
CURLOPT_HTTPHEADER => array(
    "Cookie: session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900"
),
));

$response = curl_exec($curl);

curl_close($curl);
echo $response;

```

## PHP – HTTP\_REQUEST2

```

<?php
require_once 'HTTP/Request2.php';
$request = new HTTP_Request2();
$request->setUrl('https://www.SomeServer.com/api/authenticate/logout');
$request->setMethod(HTTP_Request2::METHOD_GET);
$request->setConfig(array(
    'follow_redirects' => TRUE
));
$request->setHeader(array(
    'Cookie' => 'session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'
));
try {
    $response = $request->send();
    if ($response->getStatus() == 200) {
        echo $response->getBody();
    }
    else {
        echo 'Unexpected HTTP status: ' . $response->getStatus() . ' ' .
        $response->getReasonPhrase();
    }
}
catch(HTTP_Request2_Exception $e) {
    echo 'Error: ' . $e->getMessage();
}

```

```
}
```

## PHP – PECL\_HTTP

```
<?php
$client = new http\Client;
$request = new http\Client\Request;
$request->setRequestUrl('https://www.SomeServer.com/api/authenticate/logout');
$request->setRequestMethod('GET');
$request->setOptions(array());
$request->setHeaders(array(
    'Cookie' => 'session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSqt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'
));
$client->enqueue($request)->send();
$response = $client->getResponse();
echo $response->getBody();
```

## POWERSHELL – RESTMETHOD

```
$headers = New-Object "System.Collections.Generic.Dictionary[[String],[String]]"
$headers.Add("Cookie", "session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSqt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900")

$response = Invoke-RestMethod 'https://www.SomeServer.com/api/authenticate/logout' -
Method 'GET' -Headers $headers -Body $body
```

## PYTHON – HTTP.CLIENT

```
import http.client
import mimetypes
conn = http.client.HTTPSConnection("SomeServer.com")
payload = ''
headers = {
    'Cookie': 'session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSqt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'
}
conn.request("GET", "/WebApiBackboneQA/api/authenticate/logout", payload, headers)
res = conn.getresponse()
data = res.read()
print(data.decode("utf-8"))
```

## PYTHON – REQUESTS

```
import requests

url = "https://www.SomeServer.com/api/authenticate/logout"

payload={}
headers = {
    'Cookie': 'session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'
}

response = requests.request("GET", url, headers=headers, params=payload)

print(response.text)
```

## RUBY – NET::HTTP

```
require "uri"
require "net/http"

url = URI("https://www.SomeServer.com/api/authenticate/logout")

https = Net::HTTP.new(url.host, url.port)
https.use_ssl = true

request = Net::HTTP::Get.new(url)
request["Cookie"] = "session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900"

response = https.request(request)
puts response.read_body
```

## SHELL – HTTPIE

```
http --follow --timeout 3600 GET https://www.SomeServer.com/api/authenticate/logout \
  Cookie:'session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900'
```

## SHELL – WGET

```
wget --no-check-certificate --quiet \
  --method GET \
```

```
--timeout=0 \
--header 'Cookie: session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900' \
'https://Www.SomeServer.com/api/authenticate/logout'
```

## SWIFT – URLSESSION

```
import Foundation

var semaphore = DispatchSemaphore (value: 0)

var request = URLRequest(url: URL(string: "https://Www.SomeServer.com/api/authenticate/logout"
)!,timeoutInterval: Double.infinity)
request.addValue("session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1604350665&timelim
it=900", forHTTPHeaderField: "Cookie")

request.httpMethod = "GET"

let task = URLSession.shared.dataTask(with: request) { data, response, error in
    guard let data = data else {
        print(String(describing: error))
        semaphore.signal()
        return
    }
    print(String(data: data, encoding: .utf8)!)
    semaphore.signal()
}

task.resume()
semaphore.wait()
```

## VB.NET – RESTSHARP

```
Dim client = New
RestClient("https://Www.SomeServer.com/api/authenticate/logout")
client.Timeout = -1
Dim request = New RestRequest(Method.[GET])
request.AddHeader("Cookie", "session=session-
id=kUn8axg5aDzfB0uVqs2QwF2pZiaI%2fv%2fQSSQt2mFq000%3d&username=346650&epoch=1
604350665&timelimit=900")
Dim response As IRestResponse = client.Execute(request)
Console.WriteLine(response.Content)
```

## FOR MORE INFORMATION

Please contact Psyche Systems technical support.